

Kidney Bowl

105AA260

NPC Code FTA108


Quantity 260 per Carton



KM 650633



Tested to Comply with
PAS29:1999



Outer Case Dimensions 595x 245 x 400mm



OUTER BAR CODE

Vernacare

- SINGLE USE
- MOULDED FIBRE PULP
- MACERATABLE
- MADE FROM RECYCLED PRODUCT

- ☐ General purpose receptacle (not for use with detergents)
- ☐ Ideal for holding a variety of medical items including dressings, wipes or surgical instruments.
- ☐ Disposable single use product to reduce risk of cross infection
- ☐ Moulded products made from 100% recycled cellulose fibres
- ☐ Fully disposable in Vernacare Disposal Units

Regulatory Compliance

- ☐ Compliant with: Medical Device Regulation 2017/745 Class 1 (non-sterile)
- ☐ Product Specification Assessment PAS 29:1999
- ☐ Disposability: By Maceration
- ☐ ISO 9001:2015 / ISO 14001:2015

Description:

- ☐ Capacity: 700ml (max)
- ☐ Material: 100% Paper Pulp
- ☐ Dimensions: 240 x 130 x 45mm
- ☐ Weight: 19.2g (+/- 5g)
- ☐ Colour: Grey

Fluid Retention Testing

- ☐ Tested in Accordance with PAS 29:1999
- ☐ Guaranteed to hold fluid for 4 Hours
- ☐ Resistance Temperature: 35 C (+/- 3)



Due to the nature of the manufacturing process the products may vary in appearance. This will not affect the performance of the product.

+44 1204 529494 (UK) | (416) 661 5552 (N.America) | (800) 268 2422 (Toll Free)
Vernacare Ltd, Matrix Park, 1 Western Ave, Chorley Lancs UK, PR7 7NB

Artwork Code VC282 Revision 2 14/07/2022 Nathan Shields