

A simple 3 section couch, with an adjustable backrest and foot section for supporting the patient whilst being treated or examined. An optional breathing hole and plug can also be added for prone treatments


MODEL SPECIFICATION

- Safe working load and lifting capacity (from minimum height) of 200kgs (440lbs=31.5 stone)
- Power assisted head section adjustment from -25 degrees to +40 degrees
- Power assisted foot section adjustment from horizontal to +80 degrees
- Excellent height range from 46-99cm (18"-39")
- 100mm diameter individual locking castors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard
- Standard width 64cm (25")
- Choice of upholstery colours - Dark Grey, Dark Blue or Black

OPTIONAL COUCH ACCESSORIES

- Breathing hole and plug
- Single foot switch operation
- Paper towel holder
- Matching therapy chair or stool
- Removable head/knee support
- Stretch towelling covers

STERLING 3 SECTION COUCH - Model No. SX3008



Model **SX3008** - with optional breathing hole. Seated position. Colour Dark Grey



Model **SX3008** - with optional breathing hole. Head section in positive position. Colour Dark Grey



Model **SX3008** - with optional breathing hole. Low transfer height. Colour Dark Grey