3M[™] Defib-Pads

Bright orange color for safe and effective defibrillation and cardioversions.





$3M^{M}$ Defib-Pads save time when life is on the line.

When you're working to save a life, every second is critical. So take advantage of every second you've got with 3M Defib-Pads for defibrillation and cardioversion. The ready-to-use pads (electro-conductive gel supported by a porous non-woven fabric) are faster to apply than gels and creams in multiple cardiac arrest situations and their bright orange color makes them easy to see.

Quick and convenient

- Saves time just open package, place on patient and follow normal defib protocol
- Easy to remove from patient no slippery residue to interfere with CPR
- No cleanup

Safe

- High visibility permits proper pad positioning without overlap
- Self-contained pads eliminate "gel bridging," which can cause severe skin trauma or burns to the patient or operator
- No gel bridging also ensures that available energy is delivered to the patient
- No gel means no worry about improperly cleaned paddles

Effective

- Years of proven success
- Success rate equivalent to lowest resistivity gels
- No pitting or corrosion of paddles



For More Information

Contact your 3M Sales Representative, or call the 3M Health Care Customer Helpline at 1-800-228-3957. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.

3M Alaska 11151 Calaska Circle Anchorage, AK 99515 (907) 522-5200

3M Hawaii Box 30048 Honolulu, HI 96820 (808) 422-2721

3M Puerto Rico

Puerto Rico Industrial Park P.O. Box 100 Pueblo Station Carolina, Puerto Rico 00986-0100 (809) 750-3000

Ordering Information

3M™ Defib-Pads				
Catalog No.	Description	Pairs/Pkg	Pkg/Bx	Bx/Cs
2345N	4-1/2" x 4-1/2" (11,4cm x 11,4cm)	1	10	10
2346N	4-1/2" x 6" (11,4cm x 15,2cm)	1	10	10

2



Health Care

3M Center, Building 275-4W-01 St. Paul, MN 55144-1000 USA 1-800-228-3957 **3M Canada Inc.** Post Office Box 5757

Post Office Box 5757 London, Ontario N6A 4T1 Canada 1-800-563-2921



Recycled Paper 40% pre-consumer 10% post-consumer Printed in U.S.A. H.I. 4509 © 3M 2002 70-2008-9860-2 (22.23) DPI