

GAMMEX® Latex DermaShield

S U R G I C A L G L O V E S

Latex, powder-free surgical glove you can rely on for everyday performance

KEY FEATURES AND BENEFITS:

- Smooth with micro-textured finish for secure handling of instruments under dry and damp conditions
- Made using Ansell's unique anatomic formers to enhance comfort and fit to help reduce hand fatigue
- The DERMASHIELD™ Technology internal coating facilitates dry and damp donning as well as intra-operative changes
- Additional cleaning processes ensure low residues of both latex proteins and chemical accelerators to minimize the risk of developing latex Type I and chemical Type IV allergies and skin sensitivities

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- Brief, low-risk examination procedures
- Extra comfort during long, sustained procedures
- Sterile procedures outside of the O.R.
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES

 **A.R.T.**
Pinhole Barrier Technology

 **DERMASHIELD™**
Donning Aid Technology



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Latex DermaShield

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Straight with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	No, not tested in accordance with ASTM D6978
Freedom from Holes (Inspection level I)	0.65 AQL
Product Standards	AS/NZS 4179, EN 420:2003 + A1:2009, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthogen Polysulphides, Sodium di-n-butylthiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330063055	330063060	330063065	330063070	330063075	330063080	330063085	330063090

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: luvass.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com