

GAMMEX® Latex Textured

S U R G I C A L G L O V E S

Latex, powder-free surgical glove with a micro-textured finish for secure instrument handling and grip

KEY FEATURES AND BENEFITS:

- Micro-textured finish with special external surface treatment to ensure exceptional grip under dry and damp conditions
- Ideal thickness to optimize tactile sensitivity and protection
- DERMASHIELD™ Technology internal polymer coating facilitates dry and damp donning as well as intra-operative changes

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- Sterile procedures in the O.R.
- Enhanced grip in wet environments
- Dry and damp hand donning
- Extra comfort during long, sustained procedures
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES

 **A.R.T.**
Pinhole Barrier Technology

 **DERMASHIELD™**
Donning Aid Technology



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



GAMMEX[®] Latex Textured

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Straight with SUREFIT [™] Technology
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated with DERMASHIELD [™] Technology
Allergy Prevention	None
Grip Level	High
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Product Standards	EN 16523-1, EN 374-2 and -4, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 420:2003 + A1:2009, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	331300655	331300660	331300665	331300670	331300675	331300680	331300685	331300690

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking.
© 2020 Ansell Limited. All Rights Reserved.