

MICRO-TOUCH[®] AFFINITY[™]

Performance Series

Neoprene, powder-free examination glove with excellent elasticity and barrier protection

This non-sterile neoprene synthetic formulation provides superior protection, comfort, fit, elasticity and tactile sensitivity with excellent resistance to a broad range of chemicals.

KEY FEATURES AND BENEFITS:

- Offers comfort level similar to latex.
- Textured finger surfaces provide a secure grip in both dry and wet conditions for enhanced instrument handling.
- Excellent tear and puncture resistance
- Protects from latex Type I allergy
- Tested for use with chemotherapy drugs*

* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

WARNING: No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

Recommended For

- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients
- Enhanced grip in wet environments
- Sustained, high-risk examination procedures

FEATURED TECHNOLOGIES



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] AFFINITY[™]

Performance Series

PRODUCT INFORMATION

Material	Neoprene (Polychloroprene)	Not Made From Natural Rubber Latex	Yes
Color	Green	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.13 / 5.1
Finger Thickness (mm/mil)	0.16 / 6.3	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Performance Series

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	3770	3771	3772	3773	3774

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6977, CA Prop 65, Category III, EN 455 Part 1, EN 455 Part 3, FDA21 CFR 177-2600, ISO 13485, ISO 10993-10:2002

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking.
© 2020 Ansell Limited. All Rights Reserved.