



Easy to Apply



Contours to Skin



Cost Effective

AsGuard Flex + Non-Woven Island Dressings combine an absorbent pad with a soft and comfortable non-woven adhesive backing.

SOFT NON-WOVEN MATERIAL

Soft on skin and reduces edge lifting

LOW-ADHERENT PAD

Comfortable and absorbent with minimal pain on removal

FLEXIBLE

Ideal for use on areas difficult to dress (e.g. heel)



Suitable for:

Lightly exuding wounds

Post-operative and minor wounds where water resistance is not required

Reorder	Dressing Size	Padding Size	Sterility	Pcs/Unit	Units/Ctn
ISD001	7.5cm x 5cm	4cm x 2.5cm	Sterile	50	10
ISD002	8cm x 10cm	4cm x 5cm	Sterile	50	10
ISD030	8.25cm x 6cm	6cm x 3cm	Sterile	50	10
ISD006	10cm x 25cm	5cm x 20cm	Sterile	25	10
ISD008	10cm x 30cm	5cm x 25cm	Sterile	25	10
ISD003	15cm x 10cm	10.5cm x 5cm	Sterile	20	10
ISD004	20cm x 10cm	15.5cm x 5cm	Sterile	20	10

