

A Global Leader in Hypodermic Products



Helping all people
live healthy lives

Hypodermic Product Catalog

A Complete Line of
Medication Administration
and Safety-Engineered
Products For
The Healthcare Industry



BD™ Insulin Syringes (sterile, single use)

3/10 mL BD Lo-Dose™ Insulin Syringe with Permanently Attached Needle

Self-Contained Package. 100 per box, 500 per case

328438 31 G x 5/16 in. U-100 BD Ultra-Fine™ Short Needle
200 per box, 200 per case

324702 29 G x 1/2 in. U-100 BD Ultra-Fine



1/2 mL BD Lo-Dose™ Insulin Syringe with Permanently Attached Needle

Blister Package. 100 per box, 500 per case


329461 28 G x 1/2 in. U-100 BD Micro-Fine™ IV (Orange)

Self-Contained Package. 100 per box, 500 per case

329465 28 G x 1/2 in. U-100 BD Micro-Fine™ IV (Orange)
328468 31 G x 5/16 in. U-100 BD Ultra-Fine™ Short Needle

200 per box, 200 per case

324703 29 G x 1/2 in. U-100 BD Ultra-Fine



NOTE: All of the syringes in this section are in an Easy-Peel Blister Package unless otherwise noted as Self-Contained Package.

1 mL BD Insulin Syringe with Permanently Attached Needle

Self-Contained Package. 100 per box, 500 per case


329410 28 G x 1/2 in. U-100 BD Micro-Fine™ IV (Orange)
329412 27 G x 5/8 in. U-100 BD Micro-Fine™ IV (Orange)
328418 31 G x 5/16 in. U-100 BD Ultra-Fine™ Short Needle

Self-Contained Package. 200 per box, 200 per case

324704 29 G x 1/2 in. U-100 BD Ultra-Fine™
324705 27 G x 1/2 in. U-100 BD Ultra-Fine™

Blister Package. 100 per box, 500 per case

329424 28 G x 1/2 in. U-100 BD Micro-Fine™ IV (Orange)
329420 28 G x 1/2 in. U-100 BD Micro-Fine™ IV
Single Unit Graduation



1 mL BD Insulin Syringe with Detachable Needle

Blister Package. 100 per box, 800 per case - 2 unit scale

329651 25 G x 5/8 in. U-100 (Orange/Blue)
329652 26 G x 1/2 in. U-100 (Orange/Tan)
329622 25 G x 1 in. U-100 (Orange/Blue)



1 mL BD Insulin Syringe Only

Blister Package. 200 per box, 1600 per case- 2 unit scale

329654 U-100 Slip Tip (Orange)



100 per box, 800 per case - 1 unit scale

309629 BD Luer-Lok™, 1 mL

BD™ Allergy Syringes

BD SafetyGlide™ Allergy Syringe Tray with Permanently Attached Needle

25 per tray, 1,000 per case

305950 1 mL Allergy Tray with
27 G x 1/2 in. Regular Bevel
305951 1 mL Allergy Tray with
26 G x 3/8 in. Intradermal Bevel



1/2 mL BD Allergist Tray with Permanently Attached Needle

25 syringes per tray, 1,000 per case

305535 27 G x 1/2 in. Regular Bevel
305536 27 G x 3/8 in. Intradermal Bevel



1 mL BD Allergist Tray with Permanently Attached Needle

25 syringes per tray, 1,000 per case


305537 26 G x 1/2 in. Regular Bevel
305538 26 G x 1/2 in. Intradermal Bevel
305539 26 G x 3/8 in. Intradermal Bevel
305540 27 G x 1/2 in. Regular Bevel
305541 27 G x 3/8 in. Intradermal Bevel
305542 27 G x 3/8 in. Regular Bevel



1 mL BD Allergy Syringe with Permanently Attached Needle

100 per box, 500 per case

305500 28 G x 1/2 in. Regular Bevel



BD™ Tuberculin Syringes

1/2 mL BD Tuberculin Syringe with Permanently Attached Needle

100 per box, 500 per case 1/100 mL

305620 27 G x 1/2 in. Regular Bevel



1 mL BD Tuberculin Syringe Only

200 per box, 1600 per case 1/100 mL

309659 Slip Tip



1 mL BD Tuberculin Syringe with Detachable Needle, Slip Tip

100 per box, 800 per case 1/100 mL

309623 27 G x 1/2 in.
309625 26 G x 3/8 in. Intradermal Bevel
309626 25 G x 5/8 in.
309624 21 G x 1 in.



BD™ Syringes & Needles (sterile, single use)

1 mL BD Syringe/Needle Combination, Slip Tip

100 per box, 800 per case 1/100 mL

309597 26 G x 5/8 in. Sub-Q



1 mL BD Luer-Lok™ Syringe/Needle Combination

100 per box, 800 per case 1/100 mL

309637 20 G x 1 in.



1 mL BD Luer-Lok™ Syringe Only

100 per box, 800 per case 1/100 mL

309628 BD Luer-Lok™



3 mL BD Syringe/Needle Combination BD Luer-Lok™ Tip, Thin Wall

100 per box, 800 per case 1/10 mL

309587 26 G x 5/8 in. Sub-Q
309588 23 G x 1 in. IM
309589 23 G x 1 1/2 in. IM
309598 21 G x 1 in. IV
309599 21 G x 1 1/2 in. IV



3 mL BD Syringe/Needle Combination BD Luer-Lok™ Tip

100 per box, 800 per case 1/10 mL

309582 25 G x 1 1/2 in.	309577 21 G x 1 1/2 in.
309581 25 G x 1 in.	309575 21 G x 1 in.
309570 25 G x 5/8 in.	309579 20 G x 1 1/2 in.
309571 23 G x 1 in.	309578 20 G x 1 in.
309574 22 G x 1 1/2 in.	309580 18 G x 1 1/2 in.
309573 22 G x 1 1/4 in.	305060 18 G x 1 1/2 in.
309572 22 G x 1 in.	Blunt Fill Needle
309569 22 G x 3/4 in.	



3 mL BD Syringe Only

200 per box, 800 per case 1/10 mL

309657 BD Luer-Lok™ Tip
309656 Slip Tip



5 mL BD Syringe/Needle Combination BD Luer-Lok™ Tip

100 per box, 400 per case 1/5 mL


309630 22 G x 1 in.	309634 20 G x 1 in.
309631 22 G x 1 1/2 in.	309635 20 G x 1 1/2 in.
309632 21 G x 1 in.	305062 18 G x 1 1/2 in. Blunt Fill Needle
309633 21 G x 1 1/2 in.	



5 mL BD Syringe Only

125 per box, 500 per case 1/5 mL

309646 BD Luer-Lok™ Tip
309647 Slip Tip



BD™ Syringes & Needles (sterile, single use)

10 mL BD Syringe/Needle Combination BD Luer-Lok™ Tip

100 per box, 400 per case 1/5 mL

- 309640 22 G x 1 in.
- 309642 21 G x 1 in.
- 309643 21 G x 1 1/2 in.
- 309644 20 G x 1 in.
- 309645 20 G x 1 1/2 in.

100 per box, 400 per case

- 305064 18 G x 1 1/2 in. BD™ Blunt Fill Needle

BD Luer-Lok™ Tip, Thin Wall

100 per box, 400 per case 1/5 mL

- 309639 21 G x 1 1/2 in.



20 mL BD Syringe Only

48 per box, 192 per case 1 mL

- 302830 BD Luer-Lok™ Tip
- 302831 Slip Tip

120 per box, 480 per case

- 300613 Eccentric Tip



30 mL BD Syringe Only

56 per box, 224 per case 1 mL

- 302832 BD Luer-Lok™ Tip
- 302833 Slip Tip



10 mL BD Syringe Only

100 per box, 400 per case 1/5 mL

- 309604 BD Luer-Lok™ Tip
- 301604 Slip Tip
- 305482 Eccentric Tip

25 per box, 100 per case

- 309695 BD Luer-Lok™ Tip Control Syringe

100 per box, 400 per case

- 305462 Lab Syringe Eccentric Tip Non-Sterile, Clean, Ready-to-use Bulk-Pak




60 mL BD Syringe Only

40 per box, 160 per case 1 mL

- 309654 Slip Tip
- 309653 BD Luer-Lok™ Tip

60 per box, 240 per case


- 300866 Eccentric Tip



60 mL BD Catheter Tip Syringe Only (with Tip Shield)

40 per box, 160 per case 1 mL

- 309620 Catheter Tip



BD Pharmacy Convenience Paks

Tray Contents:

- | | | |
|--------|-----------------------|--|
| 309701 | 1 mL Luer Slip Tip | 25 syringes only per sterile tray, 12 trays per case of 300 |
| 309702 | 3 mL BD Luer-Lok™ Tip | 25 syringes only per sterile tray, 12 trays per case of 300 |
| 309703 | 5 mL BD Luer-Lok™ Tip | 25 syringes only per sterile tray, 12 trays per case of 300 |
| 309605 | 10 mLBD Luer-Lok™ Tip | 20 syringes only per sterile tray, 12 trays per case of 240 |
| 305617 | 20 mLBD Luer-Lok™ Tip | 10 each Syringes only per sterile tray, 12 trays per case of 120 |
| 305618 | 30 mLBD Luer-Lok™ Tip | 10 each Syringes only per sterile tray, 12 trays per case of 120 |
| 309680 | 60 mLBD Luer-Lok™ Tip | 20 each Syringes only per sterile tray, 6 trays per case of 120 |



BD Luer-Lok™ Sterile Syringe Tip Cap

200 per box, 2000 per case

- 305819 BD Syringe Tip Cap



BD Luer Slip™ Sterile Syringe Tip Cap Tray

50 per tray, 1,000 per box, 4,000 per case

- 305822 Syringe Tip Cap Tray



BD PrecisionGlide™ Needles (sterile, single use)

BD™ Regular Bevel Needles

100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305106	30 G x 1/2 in.	
305128	30 G x 1 in.	
305109	27 G x 1/2 in.	
305136	27 G x 1 1/4 in.	
301629	27 G x 1 1/2 in.	
305110	26 G x 3/8 in.	Intradermal Bevel
305111	26 G x 1/2 in.	
305122	25 G x 5/8 in.	
305124	25 G x 7/8 in.	
305125	25 G x 1 in.	
305127	25 G x 1 1/2 in.	
305143	23 G x 3/4 in.	
305145	23 G x 1 in.	
305120	23 G x 1 1/4 in.	
305155	22 G x 1 in.	
305156	22 G x 1 1/2 in.	
305165	21 G x 1 in.	
305167	21 G x 1 1/2 in.	
305129	21 G x 2 in.	
305175	20 G x 1 in.	
305176	20 G x 1 1/2 in.	
305186	19 G x 1 in.	Thin Wall
305187	19 G x 1 1/2 in.	Thin Wall
305195	18 G x 1 in.	
305196	18 G x 1 1/2 in.	
305197	16 G x 1 in.	
305198	16 G x 1 1/2 in.	
305115	26 G x 5/8 in.	Sub-Q
305193	23 G x 1 in.	Thin Wall IM
305194	23 G x 1 1/2 in.	Thin Wall IM
305185	18 G x 1 1/2 in.	Thin Wall Fill
305177	21 G x 1 in.	Thin Wall IV
305190	21 G x 1 1/2 in.	Thin Wall IV

BD™ Blunt Fill Needles

100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305180	18 G x 1 1/2 in.	
305181	18 G x 1 in.	
305183	20 G x 1 in.	
305211	18 G x 1 1/2 in.	Filter (5 Micron)

BD™ Short Bevel Needles

100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305159	22 G x 1 1/2 in.	
305178	20 G x 1 in.	
305179	20 G x 1 1/2 in.	
305188	19 G x 1 in.	Thin Wall
305189	19 G x 1 1/2 in.	Thin Wall
305199	18 G x 1 1/2 in.	

Types of needle walls and bevels

The Thin Wall needle is thinner than a Regular Wall needle. The inner diameter is bigger without any increase in the outside diameter. For example, a 21 G Thin Wall needle has the same outside diameter as a 21 G Regular Wall, but has the same inner diameter as a 20 G needle. Thin Wall allows the same flow volume as a Regular Wall needle one gauge larger. Thin Wall is available in 18, 19, 21, 22, and 23 gauge sizes.

Regular Wall




Thin Wall

Types of syringe tips


BD Luer-Lok™ Tip Sizes:

1 mL
3 mL
5 mL
10 mL
20 mL
30 mL
60 mL




Slip Tip Sizes:

1 mL
3 mL
5 mL
10 mL
20 mL
30 mL
60 mL



Eccentric Tip Sizes:
10 mL
20 mL
60 mL



Catheter Tip Sizes:
2 oz. only



Regular Bevel



Short Bevel



Intradermal Bevel



30 G

27 G

26 G

25 G

24 G

23 G

22 G

21 G

20 G

19 G

18 G

16 G

BD
Conventional
Needle &
Syringe
Color Chart

BD™ Disposal Solutions



BD Patient Room Collectors

Available in multiple colors

	Recykleen™	
305053	BD Recykleen™ Pearl, 3 gal.	10/Case
305054	BD Recykleen™ Red, 3 gal.	10/Case
	Horizontal Entry	
305422	Pearl, Nestable, 2 gal.	10/Case
	Side Entry	
305443	Red, Nestable, 5.4 qt.	20/Case
305444	Pearl, Nestable, 5.4 qt.	20/Case
305425	Pearl, Nestable, 5.4 qt.	12/Case

BD Phlebotomy Collectors

305487	1.5 qt.	36/Case
305635	1.0 qt.	60/Case

BD Multi-Use Nestable Collectors

305480	Clear Top, Hinge Cap with Petals, 14 qt.	20/Case
305464	Clear Top, Hinge Cap w/ Needle Ports, 14 qt.	20/Case
305456	Clear Top, Hinge Cap, 14 qt.	20/Case
305646	Clear Top, Tethered Cap, 14 qt.	20/Case
305457	Clear Top, Hinge Cap, 6 gal.	12/Case
305465	Clear Top, Hinge Cap with Petals, 6 gal.	12/Case

BD Large and Extra Large Collectors

305615	Red, Clear Hinged Top, 9 gal.	8/Case
305616	Red, Clear Slide Top, 9 gal.	8/Case
305601	Red, Clear Hinged Top, Gasketed, 9 gal.	8/Case
305602	Red, Clear Slide Top, Gasketed, 9 gal.	8/Case
305610	Red, Clear Hinged Top, 17 gal.	5/Case
305665	Red, Clear Hinged Top, Gasketed, 17 gal.	5/Case
305609	Red, Clear Slide Top, 19 gal.	5/Case
305666	Red, Clear Slide Top, Gasketed, 19 gal.	5/Case

BD Chemotherapy Collectors

305493	Yellow, Plug Cap, 5 gal.	8/Case
305578	Yellow, Screw Cap, 5 gal.	8/Case
305603	Yellow, Clear Hinged Top, 9 gal.	8/Case
305604	Yellow, Clear Slide Top, 9 gal.	8/Case
305613	Yellow, Clear Slide Top, 19 gal.	5/Case
305614	Yellow, Clear Hinged Top, 17 gal.	5/Case

BD Multi-Use One-Piece Collectors

305488	Small, 3.3 qt.	24/Case
305489	Medium, 6.9 qt.	12/Case
305490	Large, 8.2 qt.	12/Case
305491	Open Top, 5 gal.	8/Case
305477	Large Funnel, 5 gal.	8/Case
305577	Screw Cap, 5 gal.	8/Case
305557	Regular Funnel, 1.4 qt.	36/Case

BD Pharmaceutical Waste Sharps Collectors

305058	BD Recykleen™, 5.4 qt.	20/Case
305622	Counterbalanced door, 3 gal.	10/Case
305633	Blue Plug Cap, 3 gal.	12/Case
305634	Blue Hinged Top, 9 gal.	8/Case
305329	Slide Top, 9 gal.	8/Case
305028	Hinge Top, 17 gal.	5/Case
305328	Slide Top, 19 gal.	5/Case

BD RCRA Hazardous Waste Collectors

Non-DOT

305066	Hinge Top, 3 gal.	12/Case
305067	Hinge Top, 5 gal.	5/Case
305068	BD Recykleen™ Clear Hinge Top, 9 gal.	8/Case
305069	BD Recykleen™ Clear Slide Top, 9 gal.	8/Case
305071	BD Recykleen™ Clear Hinge Top, 17 gal.	5/Case
305072	BD Recykleen™ Clear Slide Top, 19 gal.	5/Case

BD Disposal Solutions Cabinets and Accessories

305096	BD Recykleen™ Brown Cabinet, 5.4 qt.	1/Case
305098	BD Recykleen™ Green Cabinet, 5.4 qt.	1/Case
305097	BD Recykleen™ Brown Cabinet, 2 and 3 gal.	1/Case
305099	BD Recykleen™ Green Cabinet, 2 and 3 gal.	1/Case
305017	Locking Wall Cabinet, 5.4 qt.	1/Case
305420	Locking Wall Cabinet, 2 and 3 gal.	1/Case
305475	Nestable, Vertical Entry Wall Cabinet, 3.2 qt.	12/Case
305409	Locking Wall Bracket, 2 and 3 gal.	10/Case
305447	Locking Wall Bracket, 5.4 qt. Collectors	10/Case
305451	One-Piece, Cup Bracket, 1.4 qt.	12/Case
305485	Nestable, Bracket, 8 qt., 14 qt. and 6 gal.	12/Case
305494	Security Lock	12/Case
305496	One-Piece, Bracket, 6.9 qt. and 8.2 qt.	12/Case
305497	One-Piece, Strap, 3.3 qt. through 5 gal.	12/Case
305552	Stabilizer, 2 and 3 gal.	5/Case
305570	Nestable, Stabilizer, 8 qt.	5/Case
305588	Wall/Cart, 2 and 3 gal.	6/Case
305021	Stabilizer, 5.4 qt.	5/Case
305494	Padlock	12/Case
305495	XL Security Bracket, 5 gal.	12/Case
305970	Cup Bracket, 1 qt.	12/Case

See BD Disposal Solutions Catalog for a complete listing of products available.

Interlink® System and BD™ Cannulas


BD™ Plastic Cannula Products for use with Baxter Interlink® System

100 per box, 1000 per case

303369 Threaded Lock Cannula



303370 Lever Lock Cannula



303365 Y Lock Cannula

BD Blood Collection Assembly for use with Baxter Interlink®, Abbott LifeShield®, or B. Braun SafeLine® Systems

25 per box, 200 per case

303380 BD Vacutainer® Specimen Collection Assembly with BD™ Blunt Plastic Cannula



Syringes with BD™ Cannula for use with Baxter Interlink®, Abbott LifeShield®, or B. Braun SafeLine® Systems

100 per box, 800 per case

303391 3 mL Syringe with BD Twinpak™ Cannula

303346 3 mL Syringe with BD Blunt Plastic Cannula

303401 3 mL Syringe with Vial Access Cannula

100 per box, 400 per case

303347 5 mL Syringe with BD Blunt Plastic Cannula

303348 10 mL Syringe with BD Blunt Plastic Cannula

303392 5 mL Syringe with BD Twinpak™ Cannula

303393 10 mL Syringe with BD Twinpak™ Cannula

303403 5 mL Syringe with Interlink® Vial Access Cannula

303405 10 mL Syringe with Interlink® Vial Access Cannula



BD Cannula Products for use with Baxter Interlink®, Abbott LifeShield®, or B. Braun SafeLine® Systems

100 per box, 1000 per case

303345 BD Blunt Plastic Cannula

303390 BD Twinpak™ Dual Cannula Device

303367 Interlink® Vial Access Cannula



BD Safety-Lok™ Syringes (sterile, single use)

1 mL BD Safety-Lok™ Insulin Syringe with Permanently Attached Needle

100 per box, 500 per case

329464 29 G x 1/2 in. U-100 BD Ultra-Fine™

1 mL BD Safety-Lok™ Tuberculin Syringe with Permanently Attached Needle

100 per box, 500 per case

305553 27 G x 1/2 in.

305554 25 G x 5/8 in.



3 mL BD Safety-Lok™ Syringe with Detachable Needle

100 per box, 800 per case

309592 25 G x 5/8 in.

309596 22 G x 1 in.

309594 23 G x 1 in.

309595 21 G x 1 1/2 in.

309593 22 G x 1 1/2 in.



5 mL BD Safety-Lok™ Syringe with Detachable Needle

50 per box, 400 per case

305561 21 G x 1 1/2 in.



10 mL BD Safety-Lok™ Syringe with Detachable Needle

50 per box, 400 per case

305564 21 G x 1 1/2 in.



BD Safety-Lok™ Syringes only, BD Luer-Lok™ Tip

309606 3 mL Syringe Only

100 box/800 case

305558 5 mL Syringe Only

50 box/400 case

305559 10 mL Syringe Only

50 box/400 case

BD SafetyGlide™ Syringes & Needles (sterile, single use)

BD SafetyGlide™ Hypodermic Needle

50 per box, 500 per case

305921	27 G x 5/8 in. Needle
305901	25 G x 5/8 in. Needle
305916	25 G x 1 in. Needle
305902	23 G x 1 in. Needle
305900	22 G x 1 1/2 in. Needle
305915	21 G x 1 in. Thin Wall
305917	21 G x 1 1/2 in. Thin Wall
305918	18 G x 1 1/2 in. Needle



BD SafetyGlide™ Insulin Syringe with Permanently Attached Needle

100 per box, 400 per case

305930	1 mL Insulin Syringe with 29 G x 1/2 in.
305932	1/2 mL Insulin Syringe with 29 G x 1/2 in.
305935	3/10 mL Insulin Syringe with 29 G x 1/2 in.
305934	1/2 mL Insulin Syringe with 30 G x 5/16 in. Thin Wall
305937	3/10 mL Insulin Syringe with 31 G x 5/16 in. Thin Wall



BD SafetyGlide™ Needle with Syringes

50 per box, 400 per case

1 mL BD Luer Slip Syringe w/ Detachable Needle

305903	25 G x 5/8 in.
305927	27 G x 5/8 in.

3 mL BD Luer-Lok™ Syringe with Detachable Needle

305904	25 G x 5/8 in.
305924	25 G x 1 in.
305905	23 G x 1 in.
305906	22 G x 1 1/2 in.
305909	21 G x 1 1/2 in. Thin Wall

5 mL BD Luer-Lok™ Syringe with Detachable Needle

305907	22 G x 1 1/2 in.
--------	------------------

10 mL BD Luer-Lok™ Syringe with Detachable Needle

305908	22 G x 1 1/2 in.
--------	------------------



BD SafetyGlide™ Tuberculin Syringe with Permanently Attached Needle

100 per box, 400 per case

305945	1 mL Tuberculin Syringe with 27 G x 1/2 in. Regular Bevel
305946	1 mL Tuberculin Syringe with 26 G x 3/8 in. Intradermal Bevel



BD SafetyGlide™ Allergy Syringe Tray with Permanently Attached Needle

25 per tray, 1000 per case

305950	1 mL Allergy Tray with 27 G x 1/2 in. Regular Bevel
305951	1 mL Allergy Tray with 26 G x 3/8 in. Intradermal Bevel



BD Eclipse™ Needles and Syringes

BD Eclipse™ Needles

100 per box, 1,200 per case

305757	30 G x 1/2 in.
305758	27 G x 1/2 in.
305759	25 G x 5/8 in.
305767	25 G x 1 1/2 in.
305761	25 G x 1 in.
305762	23 G x 1 in.
305768	22 G x 1 in.
305763	22 G x 1 1/2 in.
305765	21 G x 1 1/2 in. Thin Wall
305764	21 G x 1 in. Thin Wall
305766	18 G x 1 1/2 in.



BD Eclipse™ Needles with Syringes

50 per box, 300 per case

1 mL BD Luer-Lok™ Syringe with Detachable Needle

305778	30 G x 1/2 in.
305789	27 G x 1/2 in.
305780	25 G x 5/8 in.

3 mL BD Luer-Lok™ Syringe with Detachable Needle

305769	23 G x 1 1/4 in.
305781	25 G x 5/8 in.
305787	25 G x 1 in.
305782	23 G x 1 in.
305788	22 G x 1 in.
305783	22 G x 1 1/2 in.
305779	21 G x 1 in.
305784	21 G x 1 1/2 in.

5 mL BD Luer-Lok™ Syringe with Detachable Needle

305785	22 G x 1 1/2 in.
--------	------------------

10 mL BD Luer-Lok™ Syringe with Detachable Needle

305786	22 G x 1 1/2 in.
--------	------------------



BD Integra™ Needles and Syringes

BD Integra™ 3 mL Syringe with Detachable Needle

100 per box, 400 per case

305269	3 mL Syringe with 25 G x 5/8 in. Needle
305270	3 mL Syringe with 25 G x 1 in. Needle
305271	3 mL Syringe with 23 G x 1 in. Needle
305272	3 mL Syringe with 22 G x 1 1/2 in. Needle
305273	3 mL Syringe with 21 G x 1 in. Needle
305274	3 mL Syringe with 21 G x 1 1/2 in. Needle



BD Integra™ Needles

100 per box, 1,000 per case

305310	25 G x 5/8 in. Retracting Needle
305311	25 G x 1 in. Retracting Needle
305312	23 G x 1 in. Retracting Needle
305313	22 G x 1 1/2 in. Retracting Needle
305275 ¹	18 G x 1 1/2 in. Blunt Fill Needle
305833 ^{1,2}	18 G x 1 1/2 in. Blunt Filter Needle

1 – 18 G Blunt Fill and Blunt Filter Needles have a non-retracting needle
2 – Packaged 50 Box/500 Case



BD Integra™ 3 mL Syringe Only

100 per box, 400 per case

305283	3 mL Syringe Only
--------	-------------------

BD PosiFlush™ Prefilled Syringes

BD™ Prefilled Flush Syringes Are Available in a Wide Variety of Sizes and Solutions To Meet All Your Flush-Protocol Requirements.

BD Ref. #	Description	NHRIC Number	Box/Case
BD PosiFlush™ Normal Saline Syringes, all with 10 mL diameter			
306544	3 mL	8290-306544	30/Box 480/Case
306545	5 mL	8290-306545	30/Box 480/Case
306546	10 mL	8290-306546	30/Box 480/Case
306547	10 mL Pump Compatible Syringe*	8290-306547	30/Box 480/Case
306549	3 mL with Blunt Plastic Cannula	8290-306549	30/Box 120/Case
306550	5 mL with Blunt Plastic Cannula	8290-306550	30/Box 120/Case
306551	10 mL with Blunt Plastic Cannula	8290-306551	30/Box 120/Case
BD PosiFlush™ SF Saline Syringe			
306553	10mL Sterile Field Syringe	8290-306553	30/Box 240/Case
BD PosiFlush™ SureScrub™ Saline Syringe			
306559	10mL with Alcohol Pad	8290-306559	30/Box 480 Case
BD™ Pre-Filled Normal Saline Syringes			
306500	10 mL Fill in 10 mL Syringe	8290-306500	30/Box 120/Case
306502	5 mL Fill in 10 mL Syringe	8290-306502	30/Box 120/Case
306503	3 mL Fill in 10 mL Syringe	8290-306503	30/Box 120/Case
306504	5 mL Fill in 5 mL Syringe	8290-306504	30/Box 120/Case
306505	3 mL Fill in 5 mL Syringe	8290-306505	30/Box 120/Case
306507	3 mL Fill in 3 mL Syringe	8290-306507	30/Box 120/Case
306508	2 mL Fill in 3 mL Syringe	8290-306508	30/Box 120/Case
306518	10 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	8290-306518	30/Box 120/Case
306522	5 mL Fill in 5 mL Syringe with Blunt Plastic Cannula	8290-306522	30/Box 120/Case
306524	3 mL Fill in 3 mL Syringe with Blunt Plastic Cannula	8290-306524	30/Box 120/Case
BD PosiFlush™ Heparin Lock Flush Syringes, 10 USP units/mL			
306510	5 mL Fill in 10 mL Syringe	8290-306510	30/Box 120/Case
306511	5 mL Fill in 5 mL Syringe	8290-306511	30/Box 120/Case
306521	3 mL Fill in 10 mL Syringe	8290-306521	30/Box 120/Case
306512	3 mL Fill in 3 mL Syringe	8290-306512	30/Box 120/Case
306525	5 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	8290-306525	30/Box 120/Case
306527	2 mL Fill in 3 mL Syringe	8290-306527	30/Box 120/Case
BD PosiFlush™ Heparin Lock Flush Syringes, 100 USP units/mL			
306513	5 mL Fill in 10 mL Syringe	8290-306513	30/Box 120/Case
306514	3 mL Fill in 10 mL Syringe	8290-306514	30/Box 120/Case
306515	5 mL Fill in 5 mL Syringe	8290-306515	30/Box 120/Case
306516	3 mL Fill in 5 mL Syringe	8290-306516	30/Box 120/Case
306517	3 mL Fill in 3 mL Syringe	8290-306517	30/Box 120/Case
306528	2 mL Fill in 3 mL Syringe	8290-306528	30/Box 120/Case
306531	5 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	8290-306531	30/Box 120/Case
BD PosiFlush™ Heparin Lock Flush Syringes, 10 USP units/mL			
306413	3 mL	8290-306413	30/Box 480/Case
306414	5 mL	8290-306414	30/Box 480/Case
BD PosiFlush™ Heparin Lock Flush Syringes, 100 USP units/mL			
306423	3 mL	8290-306423	30/Box 480/Case
306424	5 mL	8290-306424	30/Box 480/Case

All products do not contain natural rubber latex or preservatives.

*Catalog No. 306547 10mL BD PosiFlush Normal Saline Flush Syringe is generally compatible for use with syringe pumps.



Brands are the property of their respective owners.
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD

Printed in. USA 2/15 MSS0158-1

1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com
2771 Bristol Circle
Oakville, Ontario L6H 6R5

