


BD SafetyGlide™ Needle

Safety without compromise



Efficiency through ease of activation

BD SafetyGlide™ Needles have Activation-Assist™ technology, making a safe injection easier and more efficient.

- Unique Activation-Assist™ technology feature means a fast and smooth activation while keeping finger well behind the needle.
- Fluid single-finger activation requiring minimal training and virtually no change in clinician technique.
- BD SafetyGlide™ Needles can help keep clinicians safe from accidental needlestick injuries.
- Detachable needle that can fit on any conventional luer syringe and facilitate standardization and SKU reduction across the facility.

Scan this code to watch a video on how to use the **BD SafetyGlide™ Needle**:



BD SafetyGlide™ Needle

BD SafetyGlide™ Needles		
Product number	Description	Packaging
305921	27 G x 5/8 in.	50/box, 500/case
305901	25 G x 5/8 in.	50/box, 500/case
305916	25 G x 1 in.	50/box, 500/case
305902	23 G x 1 in.	50/box, 500/case
304387	23 G x 1 1/2 in.	50/box, 500/case
305915	21 G x 1 in.	50/box, 500/case
305917	21 G x 1 1/2 in.	50/box, 500/case
305918	18 G x 1 1/2 in.	50/box, 500/case

BD SafetyGlide™ Needles with 1 mL luer slip syringe		
Product number	Description	Packaging
305927	27 G x 5/8 in.	50/box, 400/case
305903	25 G x 5/8 in.	50/box, 400/case

BD SafetyGlide™ Needles with 3 mL syringe		
Product number	Description	Packaging
305924	25 G x 1 in.	50/box, 400/case
305904	25 G x 5/8 in.	50/box, 400/case
305905	23 G x 1 in.	50/box, 400/case
303326	23 G x 1 1/12 in.	50/box, 400/case
305909	21 G x 1 1/2 in.	50/box, 400/case

Note: BD SafetyGlide™ Needles are compatible with both BD Luer-Lok™ Syringes and slip tip syringes.

For more information on our products, please contact Medical Support Services at **800.555.7422** or **medical.services@bd.com**

To speak with your BD sales representative, call **888.237.2762**

bd.com



BD, the BD Logo, Activation-Assist, Luer-Lok, PrecisionGlide and SafetyGlide are trademarks of Becton, Dickinson and Company or its affiliates. © 2023 BD. All rights reserved. BD-6145 (02/23)