

BD™ Sharps Collectors

Product Catalogue

A broad range of sharps collectors and bracketing devices designed for “point of use” disposal.

Delivered Clean & New every time



Helping all people
live healthy lives



Small Volume








Ideal for use in:

- Blood collection trays
- Laboratories
- Physician and Dental offices
- Community



Volume (litres) and Description	Catalogue No.	Case Quantity	Dimensions H x W x D (mm)
0.45 L	302434	100	170 x 105 x 50
0.87 L	305008	50	165 x 85 x 70
1.4 L High Top	303511	36	230 x 115 x 105
1.4 L One-piece	303205	36	190 x 90 x 90
1.4 L Temporary Closure	305011	36	210 x 80 x 80
1.5 L	305624	40	205 x 165 x 160

Compatible Brackets

Desktop Bracket for 302434	Small Cup Bracket for 305008 & 303511	Hanging Bracket for 303205 & 305011	Cup Bracket for 303205 & 305011	Lockable Cabinet for 303205 & 305011
				
305964	305970	305697	305451	305708
Hanging Bracket for 305624	Desktop Bracket for 305624			
				
305710	305645			

303511 and 305624 currently only available in New Zealand

Nestable

Available in Fluted top or Open top options. Fluted top is designed to limit accessibility to the collector contents. Open top features a larger diameter to accommodate disposal needs in non-patient or non-public areas.



Volume (litres) and Description	Catalogue No.	Case Quantity	Dimensions H x W x D (mm)
7.6 L Fluted top	303506	24	250 x 280 x 165
13.2 L Fluted top	303507	20	340 x 275 x 220
13.2 L Open top	303508	20	340 x 275 x 220
22.7 L Open top	303509	12	430 x 325 x 217
22.7 L Fluted top	300740	12	430 x 325 x 217

Compatible Brackets

Hanging Bracket for 303506	Wall Mount Bracket for 303507, 303508, 303509 & 300740	Mobile Trolley for 303507, 303508, 303509 & 300740
		
305706	305485	305702

Delivered Clean & New every time

One-Piece Range



- Single piece moulding, requires no assembly
- Fluted opening designed to limit accessibility to the collector contents
- Multiple bracketing options for stability

Volume (litres) and Description	Catalogue No.	Case Quantity	Dimensions H x W x D (mm)
1.4 L Temporary Closure	305011	36	210 x 80 x 80
1.4 L One-piece	303205	36	190 x 90 x 90
3 L One-piece	303207	24	139 x 212 x 330
8 L One-piece	303208	12	343 x 232 x 134
17 L One-piece	303209	8	450 x 266 x 190

Compatible Brackets

Hanging Bracket for 305011 & 303205 	Cup Bracket for 305011 & 303205 	Lockable Cabinet for 305011 & 303205 
305697	305451	305708
Hanging Bracket for 303207 & 303208 	Wall Mount Bracket for 303207 & 303208 	Quick Release Strap for 303207, 303208 & 303209 
305698	305699	305497
Hanging Bracket for 303209 	Wall Mount Bracket for 303209 	Mobile Trolley for 303209 
305700	305701	305702

Delivered Clean & New every time




Counter Balances & Next Generation

- Designed to improve tamper resistance by preventing hand access to the contents
- Compatible with locking cabinets to provide security and prevent unwanted access to collector contents



Volume (litres) and Description	Catalogue No.	Case Quantity	Dimensions H x W x D (mm)
5.1 L Next Generation	301272	20	272 x 297 x 108
7.6 L Next Generation	301273	10	289 x 333 x 150
7.6 L Counter Balance	305081	10	289 x 333 x 150
11.3 L Next Generation	301274	10	420 x 333 x 150
11.3 L Counter Balance	305082	10	420 x 333 x 150

Compatible Brackets

Wall Mount Bracket for 301273, 305081, 301274 & 305082	Hanging Bracket for 301273, 305081, 301274 & 305082	Mobile Trolley for 301273, 305081, 301274 & 305082
		
305409	305705	305702

Small Recykleen™ Locking Cabinet for 301272	Small Locking Cabinet for 301272	Wall Mount Bracket for 301272	Hanging Bracket for 301272
			
305096	305017	305447	305703

Mobile Point of Use Unit for 301272	Recykleen™ Locking Cabinet for 301273 & 301274	Locking Cabinet for 301273 & 301274
		
305704	305097	305420

BD Recykleen™ Cabinets are made from 100% recycled polycarbonate plastic.

Delivered Clean & New every time

Large Capacity

Ideal for use in:

- Operating rooms
- Laboratories
- Dialysis units



Volume (litres) and Description	Catalogue No.	Case Quantity	Dimensions H x W x D (mm)
34 L Slide top	305608	8	479 x 447 x 298
72 L Slide top	302406	5	660 x 488 x 357

Compatible Trolleys

Recykleen™ Foot-operated Trolley 34 L for 305608



305091

Recykleen™ Foot-operated Trolley 72 L for 302406



305093

BD Recykleen™ Trolleys



- Designed to help minimise the risk of collectors tipping over
- Wheels enable easy portability of the collector
- Locking wheels, ensuring the product is mobile only when desired
- The ability to open and close the collector for hands-free disposal




Chemotherapy



Volume (litres) and Description	Catalogue No.	Case Quantity	Dimensions H x W x D (mm)
8 L One-piece	303210	12	343 x 232 x 134
17 L One-piece	303510	8	450 x 266 x 190

Compatible Brackets

Hanging Bracket for 303210	Wall Mount Bracket for 303210	Quick Release Strap for 303210 & 303510
		
305698	305699	305497

Mobile Trolley for 303510	Hanging Bracket for 303510	Wall Mount Bracket for 303510
		
305702	305701	305700

Delivered Clean & New every time

BD™ Sharps Collector Dimensions

Style	Size	Catalogue No.	Sharps Collector Dimensions, H x W x D (mm)
Small Volume	0.45 L	302434	170 x 105 x 50
	0.87 L	305008	165 x 85 x 70
	1.4 L High Top	303511	230 x 115 x 105
	1.4 L	303205	190 x 90 x 90
	1.4 L Temp Closure	305011	210 x 80 x 80
	1.5 L	305624	205 x 165 x 160
Nestable	7.6 L	303506	250 x 280 x 165
	13.2 L Fluted top	303507	340 x 275 x 220
	13.2 L Open top	303508	340 x 275 x 220
	22.7 L Fluted top	300740	430 x 325 x 217
	22.7 L Open top	303509	430 x 325 x 217
One-piece	1.4 L	303205	190 x 90 x 90
	1.4 L Temp Closure	305011	210 x 80 x 80
	3 L	303207	139 x 212 x 330
	8 L	303208	343 x 232 x 134
	17 L	303209	450 x 266 x 190
Counter Balance	7.6 L	305081	289 x 333 x 150
	11.3 L	305082	420 x 333 x 150
Next Generation	5.1 L	301272	272 x 297 x 108
	7.6 L	301273	289 x 333 x 150
	11.3 L	301274	420 x 333 x 150
Large Capacity	34 L	305608	479 x 447 x 298
	72 L	302406	660 x 488 x 357
Chemotherapy	8 L	303210	343 x 232 x 134
	17 L	303510	450 x 266 x 190

Safe Sharps Disposal Tips

- Do not fill past the fill line
- Change the Sharps Collector as soon as the contents reach the fill line
- When removing a Sharps Collector always use the handles provided
- Secure the lid to the Sharps Collector base correctly, an audible click will be heard once the lid is in place
- Do not dispose of non-sharps waste in Sharps Collectors as this adds unnecessary cost to your waste disposal fees
- Always remember to aim for Point of Use positioning
- When placing a sharp into a Sharps Collector, place the sharp end in first. If using a Counter Balance or Next Generation Sharps Collector, you will still place the sharp end in first, but on a horizontal angle
- Do not recap needles
- Do not de-notch needles off syringes
- If you see that a collector has been positioned unsafely please alert the relevant person in your department

Delivered Clean & New every time



T 1300 111 250 | F 1300 111 251
E SALES@MARNEMEDICAL.COM.AU
W WWW.MARNEMEDICAL.COM.AU
A 109 MUNIBUNG RD, BOOLAROO NSW 2284



Australia: Becton Dickinson Pty Ltd.
4 Research Park Drive
North Ryde NSW 2113
Tel: 1800 656 100

New Zealand: Becton Dickinson Limited
14B George Bourke Drive
Mt Wellington, Auckland, 1060
Tel: 0800 572 468
www.bd.com/anz