



Light Weight



Conforms to Contours



Various Sizes

Conforming Gauze Bandages are constructed from a blend of polyamide yarn and rayon. These bandages are flexible and can be used for wounds in awkward places or to cover larger surface areas.

CONFORMS TO BODY CONTOURS

Allowing flexibility without restriction of movement

INDIVIDUALLY WRAPPED

Maintains cleanliness and hygiene

BREATHABLE

Porous non-woven material giving high comfort and air circulation



Suitable for:

Dressing retention

Wound covering

Reorder	Size	Pcs/Unit	Units/Ctn
CFB001	2.5cm x 1.5m	12	120
CFB002	5cm x 1.5m	12	60
CFB003	7.5cm x 1.5m	12	40
CFB004	10cm x 1.5m	12	30
CFB005	15cm x 1.5m	12	20

