

# eisco

I N D U S T R I A L



Glassware



Metalware



Plasticware



Woodware



Equipment



Rubberware



Ceramics

# LAB GLASS®

15 Years ago we decided to take our passion for manufacturing high quality products to lab supplies and glassware. Today we are proud to showcase our LABGLASS™ line.

Manufacturing of glassware is both an art and science that takes an incredible amount of skill. We employ the finest glass blowers and technicians to hand create glassware to the tightest tolerances while following state of the art manufacturing and calibration techniques. We use only the highest quality low expansion borosilicate 3.3 glass and follow manufacturing and calibration in accordance with ISO/DIN/ASTM standards. Each piece of volumetric glassware is calibrated individually in our ISO17025 compliant calibration lab.

The job is not done after blowing, calibrating and labeling a piece of glass, it then needs to be packaged in a way that delivers the piece unharmed to the end user. Our packaging is both innovative, protective and easy to use as we have unique package designs for each piece of glass.

We believe that you will enjoy the passion and attention to detail that can be seen when you use LABGLASS®.

## Thermal Properties

The coefficient of thermal expansion and temperature gradient of LabGlass® are low allowing the glass to withstand higher temperature gradients as well as sudden temperature changes (thermal shock).

|                                                    |                                               |
|----------------------------------------------------|-----------------------------------------------|
| Coefficient of Linear Expansion                    | $3.3 \times 10^{-6} \text{ }^{\circ}\text{C}$ |
| Strain Point                                       | 515 °C                                        |
| Annealing Point                                    | 565 °C                                        |
| Softening Point                                    | 820 °C                                        |
| Specific Heat (J/kg·°C)                            | 0.20                                          |
| Thermal Conductivity (Cal/cm <sup>3</sup> ·°C/Sec) | 0.0027                                        |

## Chemical Durability

LabGlass® glassware is highly resistant to water, neutral and acidic solutions, concentrated acids, chlorine, bromine iodine, and organic matter. It can withstand dry and wet sterilization without surface deterioration and subsequent contamination. Only hydrofluoric acid, very hot phosphoric acid and alkaline solutions increasingly attack the glass surface with rising concentration and temperature.

## Contact Chemical mg/m2

Water distilled at 100 °C  
Water Vapour Steam at 121°C  
Acid HCl  
80% H SO at 130°C  
Alkali -1N soln. of Na CO boiling  
Infusion Fluids Isotonic  
NaCl (0.85%) 121°C  
Glucose (5%) 121°C

## Duration in Hr.

6  
1  
6  
12  
6  
  
2 ½  
2 ½

## Loss in wt.

10  
75  
100  
140  
4000  
  
70  
50



## About Us...



Eisco Industrial is a division of Eisco Scientific, LLC. Over the past 50 years Eisco Scientific had grown into a multinational corporation supporting academic and industrial customers in over 100 countries through a growing network of distributors.

Eisco Scientific recently built a new manufacturing and distribution center at our corporate headquarters in India. We also have a facility in the United States that houses our research and development and customer support teams as well as our North American distribution center.

The Eisco Industrial Division has been created to focus on the specific needs of industrial, academic research, and botanical extraction laboratories. Eisco Industrial products have been carefully designed by laboratory professionals along-side experienced engineers with quality and durability in mind. They are manufactured with only the highest grade raw materials and according to international quality standards to ensure that they can withstand use in even the toughest laboratory environments.

## Certificates



WORLDIDAC MEMBER



ISO 9001:2008




CSA APPROVAL



CERTIFICATE OF  
EXPORT HOUSE



NABL APPROVED  
LABORATORY



USP

 **LAB** **GLASS<sup>®</sup>**



ISO

The ISO logo is a blue triangle pointing downwards, containing the text 'ISO' in white. The background image shows a close-up of a glass tube with an orange ring being heated by a flame.



ASTM

The ASTM logo is a blue triangle pointing downwards, containing the text 'ASTM' in white. The background image shows a rack of glass tubes with labels, including 'DIN 11812333-001'.




DIN

The DIN logo is a blue triangle pointing downwards, containing the text 'DIN' in white. The background image shows a glass tube with a small glass bulb being heated by a flame.



BORO 3.3

The BORO 3.3 logo is a blue triangle pointing downwards, containing the text 'BORO 3.3' in white. The background image shows a hand holding a glass tube with a small glass bulb, being heated by a flame.



NABL

The NABL logo is a blue triangle pointing downwards, containing the text 'NABL' in white. The background image shows a glass tube with a small glass bulb being heated by a flame, with a hand visible in the foreground.



INJECTION MOULDING



GLASS BLOWING




CALIBRATION




WAREHOUSE

# Understanding the Markings of Volumetric Glassware

## Burettes




## Cylinders




# Understanding the Markings of Volumetric Glassware

## Volumetric Flasks (ASTM)



## Volumetric Flasks (ISO)



### Beakers - Glass, Low Form with Spout

- Low form beaker with spout.
- Made of borosilicate glass.
- Graduated as per DIN 12231 & ISO 3819.

| Code     | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|----------|---------------|-----------|-------------|-----------|
| CH0126A  | 5             | 22        | 30          | 12        |
| CH0126B  | 10            | 26        | 35          | 12        |
| CH0126C  | 25            | 34        | 50          | 12        |
| CH0126D  | 50            | 42        | 60          | 12        |
| CH0126E  | 100           | 50        | 70          | 12        |
| CH0126F  | 150           | 60        | 80          | 12        |
| CH0126G  | 250           | 68        | 95          | 12        |
| CH0126H  | 400           | 80        | 110         | 6         |
| CH0126I  | 500           | 83        | 115         | 6         |
| CH0126J  | 600           | 90        | 125         | 6         |
| CH0126J8 | 800           | 100       | 135         | 6         |
| CH0126K  | 1000          | 105       | 145         | 6         |
| CH0126L  | 2000          | 130       | 185         | 4         |
| CH0126M  | 3000          | 150       | 210         | 1         |
| CH0126N  | 5000          | 170       | 270         | 1         |
| CH0126O  | 10,000        | 217       | 350         | 1         |



### Beakers - Glass, Tall Form with Spout

- Tall form beaker with spout.
- Made of borosilicate glass.
- Graduated as per DIN 12231 & ISO 3819.

| Code     | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|----------|---------------|-----------|-------------|-----------|
| CH0127A  | 25            | 30        | 55          | 12        |
| CH0127B  | 50            | 38        | 70          | 12        |
| CH0127C  | 100           | 48        | 80          | 12        |
| CH0127D  | 150           | 54        | 95          | 12        |
| CH0127E  | 250           | 60        | 120         | 12        |
| CH0127F  | 400           | 70        | 130         | 6         |
| CH0127G  | 500           | 75        | 137         | 6         |
| CH0127H  | 600           | 80        | 150         | 6         |
| CH0127H8 | 800           | 90        | 175         | 6         |
| CH0127I  | 1000          | 95        | 190         | 6         |
| CH0127J  | 2000          | 120       | 240         | 4         |
| CH0127K  | 3000          | 130       | 280         | 1         |





### Beakers - Glass, ASTM, Low Form

- Low form beaker with spout.
- Made of borosilicate glass.
- Graduated as per ASTM-E960.

| Code      | Capacity (ml) | O.D. (mm) | Height (mm) | Pack Size |
|-----------|---------------|-----------|-------------|-----------|
| CH0124A10 | 10            | 30        | 35          | 12        |
| CH0124A20 | 20            | 32        | 45          | 12        |
| CH0124A30 | 30            | 34        | 55          | 12        |
| CH0124B   | 50            | 42        | 60          | 12        |
| CH0124C   | 100           | 50        | 75          | 12        |
| CH0124D   | 150           | 60        | 90          | 12        |
| CH0124E   | 250           | 68        | 100         | 12        |
| CH0124F   | 400           | 80        | 120         | 6         |
| CH0124G   | 600           | 90        | 135         | 6         |
| CH0124H   | 800           | 100       | 145         | 6         |
| CH0124I   | 1000          | 105       | 160         | 6         |
| CH0124J   | 2000          | 130       | 200         | 4         |

### Beakers - Glass, ASTM, Tall Form

- Tall form beaker with spout.
- Made of borosilicate glass.
- Graduated as per ASTM-E960.

| Code    | Capacity (ml) | O.D. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0125A | 100           | 50        | 85          | 12        |
| CH0125B | 200           | 55        | 110         | 12        |
| CH0125C | 300           | 65        | 125         | 12        |
| CH0125D | 400           | 68        | 145         | 6         |
| CH0125E | 500           | 75        | 150         | 6         |
| CH0125F | 600           | 80        | 160         | 6         |
| CH0125G | 1000          | 90        | 205         | 6         |

### Beakers - Glass, Phillips Pattern

- Conical form beaker with spout.
- Made from borosilicate glass.

| Code    | Capacity (ml) | Height (mm) | Pack Size |
|---------|---------------|-------------|-----------|
| CH0128A | 250           | 105         | 12        |
| CH0128B | 500           | 142         | 6         |

### Crystallizing Dishes without Spout

- Flat bottom crystallizing dish with spout.
- Made of borosilicate glass.


| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0067A | 20            | 40        | 25          | 1         |
| CH0067B | 40            | 50        | 30          | 1         |
| CH0067C | 60            | 60        | 35          | 1         |
| CH0067D | 100           | 70        | 40          | 1         |
| CH0067E | 150           | 80        | 45          | 1         |
| CH0067F | 300           | 95        | 55          | 1         |
| CH0067G | 500           | 115       | 65          | 1         |



### Crystallizing Dishes with Spout

- Flat bottom crystallizing dish with spout.
- Made of borosilicate glass.

| Code    | Capacity (ml) | O.D. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0066A | 300           | 100       | 50          | 1         |
| CH0066B | 1000          | 150       | 75          | 1         |
| CH0066C | 2000          | 190       | 100         | 1         |
| CH0066D | 40            | 50        | 30          | 1         |
| CH0066E | 100           | 70        | 40          | 1         |
| CH0066F | 20            | 40        | 25          | 1         |
| CH0066G | 60            | 60        | 35          | 1         |
| CH0066H | 150           | 80        | 45          | 1         |
| CH0066I | 500           | 115       | 65          | 1         |
| CH0066J | 900           | 140       | 75          | 1         |
| CH0066K | 3500          | 230       | 100         | 1         |



### Evaporating Basins - Flat Bottom with Spout

- Flat bottom evaporating basin with spout.
- Made of borosilicate glass.

| Code      | Capacity (ml) | O.D. (mm) | Height (mm) | Pack Size |
|-----------|---------------|-----------|-------------|-----------|
| CH0068A01 | 15            | 50        | 25          | 1         |
| CH0068A02 | 45            | 60        | 30          | 1         |
| CH0068A   | 100           | 80        | 45          | 1         |
| CH0068B   | 250           | 105       | 55          | 1         |
| CH0068B01 | 500           | 131       | 62          | 1         |
| CH0068C   | 1000          | 150       | 80          | 1         |
| CH0068D   | 2000          | 200       | 100         | 1         |





### Reagent Bottles - Narrow Mouth, Polypropylene Stopper

- Narrow mouth reagent bottle with interchangeable acid proof polypropylene stopper.
- Made of clear borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0160A | 30            | 14/23       | 6         |
| CH0160B | 60            | 14/23       | 6         |
| CH0160C | 125           | 19/26       | 6         |
| CH0160D | 250           | 19/26       | 6         |
| CH0160E | 500           | 24/29       | 6         |
| CH0160F | 1000          | 29/32       | 6         |
| CH0160G | 2000          | 34/35       | 1         |



### Reagent Bottles - Narrow Mouth, Polypropylene Stopper, Amber

- Narrow mouth reagent bottle with interchangeable acid proof polypropylene stopper.
- Made of amber borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0161A | 30            | 14/23       | 6         |
| CH0161B | 60            | 14/23       | 6         |
| CH0161C | 125           | 19/26       | 6         |
| CH0161D | 250           | 19/26       | 6         |
| CH0161E | 500           | 24/29       | 6         |
| CH0161F | 1000          | 29/32       | 6         |
| CH0161G | 2000          | 34/35       | 1         |



GLASS  
BORO  
3.3



### Reagent Bottles - Narrow Mouth, Hexagonal Hollow Stopper

- Narrow mouth reagent bottle with interchangeable hollow glass stopper.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0162A | 30            | 14/23       | 6         |
| CH0162B | 60            | 14/23       | 6         |
| CH0162C | 125           | 19/26       | 6         |
| CH0162D | 250           | 19/26       | 6         |
| CH0162E | 500           | 24/29       | 6         |
| CH0162F | 1000          | 29/32       | 6         |
| CH0162G | 2000          | 34/35       | 1         |
| CH0162H | 5000          | 45/40       | 1         |
| CH0162I | 10,000        | 55/44       | 1         |

### Reagent Bottles - Wide Mouth, Hexagonal Hollow Stopper

- Wide mouth reagent bottle with interchangeable hexagonal hollow stopper.
- Made of borosilicate glass

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0163A | 60            | 24/29       | 6         |
| CH0163B | 125           | 29/32       | 6         |
| CH0163C | 250           | 34/35       | 6         |
| CH0163D | 500           | 45/40       | 6         |
| CH0163E | 1000          | 50/42       | 6         |
| CH0163F | 2000          | 55/44       | 1         |



GLASS  
BORO  
3.3



### Reagent Bottles - Graduated, Screw Cap

- Graduated reagent bottles with screw cap (GL 45).
- Pouring ring for drip free pouring.
- Bottle made of borosilicate glass.
- Screw cap & pouring ring autoclavable
- Pouring ring autoclavable.
- As per DIN ISO 4796-1 standard

| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0164A | 100           | 56        | 100         | 6         |
| CH0164B | 250           | 70        | 140         | 6         |
| CH0164C | 500           | 86        | 176         | 6         |
| CH0164D | 1000          | 101       | 225         | 6         |
| CH0164E | 2000          | 136       | 265         | 1         |



GLASS  
BORO  
3.3



### Reagent Bottles - Screw Cap, Amber

- Graduated reagent bottles with screw cap (GL 45).
- Pouring ring for drip free pouring.
- Bottle made of amber colored borosilicate glass.
- Screw cap & pouring ring autoclavable
- Pouring ring autoclavable.
- As per DIN ISO 4796-1 standard.

| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0165A | 100           | 56        | 100         | 6         |
| CH0165B | 250           | 70        | 140         | 6         |
| CH0165C | 500           | 81        | 176         | 6         |
| CH0165D | 1000          | 101       | 225         | 6         |
| CH0165E | 2000          | 136       | 265         | 1         |



GLASS  
BORO  
3.3





### BOD Bottles, Glass Stopper

- BOD bottle made of borosilicate glass with interchangeable glass stopper.
- As per ASTM E438 standard

| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0168A | 125           | 55        | 120         | 2         |
| CH0168B | 300           | 70        | 150         | 2         |



### Bottles - Pipette Dropping

- Clear pipette dropping bottles with three piece pipette section comprised of glass dropper.
- High density polythene stopper with dust-proof head & neoprene teat.
- Bottle made of borosilicate glass.
- Stopper and teat highly resistant to chemical attack by acids and alkalis.

| Code    | Capacity (ml) | Pack Size |
|---------|---------------|-----------|
| CH0194A | 30            | 6         |
| CH0194B | 60            | 6         |
| CH0194C | 125           | 6         |



### Bottles - Pipette Dropping, Amber

- Amber pipette dropping bottles with three piece pipette section comprised of glass dropper.
- High density polythene stopper with dust-proof head and neoprene teat.
- Bottle made of borosilicate glass.
- Stopper and teat highly resistant to chemical attack by acids and alkalis.

| Code    | Capacity (ml) | Pack Size |
|---------|---------------|-----------|
| CH0195A | 30            | 6         |
| CH0195B | 60            | 6         |
| CH0195C | 125           | 6         |



### Pycnometers

- Made of borosilicate glass.

| Code         |            | Capacity |
|--------------|------------|----------|
| Uncalibrated | Calibrated | (ml)     |
| BT1102       | BT1102-CL  | 10       |
| BT1104       | BT1104-CL  | 25       |
| BT1105       | BT1105-CL  | 50       |
| BT1106       | BT1106-CL  | 100      |

### Bottles - Gas Washing

- Dreschel gas washing bottle with interchangeable B-29 head joint.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Pack Size |
|---------|---------------|-----------|
| CH0196A | 125           | 1         |
| CH0196B | 250           | 1         |
| CH0196C | 500           | 1         |



GLASS BORO 3.3



### Bottles - Weighing, Low Form

- Low form weighing bottle with interchangeable ground stopper.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0198A | 20            | 40        | 30          | 10        |
| CH0198B | 20            | 50        | 25          | 10        |
| CH0198C | 35            | 50        | 35          | 10        |
| CH0198D | 50            | 50        | 50          | 10        |



GLASS BORO 3.3



### Bottles - Weighing, Tall Form

- Tall form weighing bottle with interchangeable ground stopper.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0199A | 5             | 20        | 40          | 10        |
| CH0199B | 15            | 25        | 50          | 10        |
| CH0199C | 25            | 30        | 60          | 10        |
| CH0199D | 60            | 40        | 80          | 10        |



GLASS BORO 3.3



### Bottles - Gas Wash Woulff, Two Necks

- Woulff gas wash bottles with two necks.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Socket Size | Pack Size |
|---------|---------------|-----------|-------------|-------------|-----------|
| CH0511C | 500           | 95        | 145         | 19/26       | 2         |
| CH0511D | 1000          | 110       | 180         | 24/29       | 2         |
| CH0511E | 2000          | 135       | 225         | 29/32       | 1         |



GLASS BORO 3.3



GLASS  
BORO  
3.3Chemical  
ResistantGLASS  
BORO  
3.3Chemical  
ResistantGLASS  
BORO  
3.3Chemical  
Resistant

### Bottles - Aspirator with Outlet for Tubing

- Aspirator bottle with outlet to connect to tubing.
- Made of borosilicate glass.

| Code    | Capacity (ml) | O. D. X Height (mm) | Joint (mm) | Pack Size |
|---------|---------------|---------------------|------------|-----------|
| CH0060A | 250           | 73 x 130            | 8          | 2         |
| CH0060B | 500           | 89 x 164            | 8          | 2         |
| CH0060C | 1000          | 110 x 200           | 10         | 2         |

### Bottles - Aspirator with Outlet for Stopcock

- Aspirator bottle with outlet to be used with stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | O. D. X Height (mm) | Joint | Pack Size |
|---------|---------------|---------------------|-------|-----------|
| CH0062A | 1000          | 107 x 200           | 19/26 | 1         |
| CH0062B | 2000          | 134 x 249           | 19/26 | 1         |
| CH0062C | 5000          | 182 x 320           | 29/32 | 1         |
| CH0062D | 10000         | 228 x 400           | 29/32 | 1         |

### Stopcock with standard ground joints

- For use with aspirator bottles CH0062A - D & above.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0063A | 1000 - 2000   | 19/26       | 1         |
| CH0063B | 5000 - 10000  | 29/32       | 1         |

### Petri Dishes

- Made of high quality soda lime glass with high clarity for better viewing.

| Code      | O. D. X Height (mm) | Pack Size |
|-----------|---------------------|-----------|
| CH0368C   | 90 x 15             | 10        |
| CH0368D   | 100 x 15            | 10        |
| CH0368D01 | 100 x 20            | 10        |
| CH0368E   | 120 x 20            | 6         |
| CH0368G   | 150 x 25            | 6         |
| CH0368G   | 180 x 30            | 6         |

### Petri Dishes


- Made of borosilicate glass, can with stand repeated wet or dry sterilization.
- Beaded edges provide greater mechanical strength.

| Code    | O. D. X Height (mm) | Pack Size |
|---------|---------------------|-----------|
| CH0370A | 50 x 17             | 10        |
| CH0370B | 80 x 17             | 10        |
| CH0370C | 100 x 17            | 10        |
| CH0370D | 150 x 20            | 6         |

## Flasks - Conical, Erlenmeyer, Narrow Neck

- Conical Erlenmeyer flask with narrow neck.
- Made of borosilicate glass.
- As per DIN ISO 1773 standard.


| Code      | Capacity (ml) | Dia. <sub>1</sub> (mm) | Dia. <sub>2</sub> (mm) | Height (mm) | Pack Size |
|-----------|---------------|------------------------|------------------------|-------------|-----------|
| CH0424B   | 25            | 42                     | 22                     | 75          | 12        |
| CH0424C   | 50            | 54                     | 22                     | 90          | 12        |
| CH0424D   | 100           | 65                     | 24                     | 105         | 12        |
| CH0424E   | 150           | 74                     | 28                     | 118         | 12        |
| CH0424E01 | 200           | 79                     | 34                     | 131         | 12        |
| CH0424F   | 250           | 85                     | 34                     | 145         | 12        |
| CH0424F01 | 300           | 87                     | 34                     | 156         | 12        |
| CH0424G   | 500           | 105                    | 34                     | 180         | 6         |
| CH0424H   | 1000          | 131                    | 42                     | 220         | 6         |
| CH0424I   | 2000          | 166                    | 50                     | 280         | 4         |
| CH0424J   | 3000          | 187                    | 50                     | 310         | 1         |
| CH0424K   | 5000          | 220                    | 50                     | 365         | 1         |



## Flasks - Conical, Erlenmeyer, Wide Neck

- Conical Erlenmeyer flask with wide neck.
- Made of borosilicate glass.
- As per DIN ISO 24450 standard.

| Code      | Capacity (ml) | Dia. <sub>1</sub> (mm) | Dia. <sub>2</sub> (mm) | Height (mm) | Pack Size |
|-----------|---------------|------------------------|------------------------|-------------|-----------|
| CH0426A02 | 25            | 42                     | 32                     | 70          | 12        |
| CH0426A03 | 50            | 51                     | 34                     | 85          | 12        |
| CH0426A   | 100           | 64                     | 34                     | 105         | 12        |
| CH0426B01 | 200           | 79                     | 50                     | 131         | 12        |
| CH0426B   | 250           | 85                     | 50                     | 140         | 12        |
| CH0426C01 | 300           | 87                     | 50                     | 156         | 12        |
| CH0426C   | 500           | 105                    | 50                     | 175         | 6         |
| CH0426D   | 1000          | 131                    | 50                     | 220         | 6         |
| CH0426E   | 2000          | 153                    | 72                     | 275         | 4         |





### Flasks - Conical, Erlenmeyer, Narrow Neck with Interchangeable Joint

- Conical Erlenmeyer flask with narrow neck and interchangeable joint.
- Made of borosilicate glass.
- As per DIN ISO 4797 standard.

| Code      | Capacity (ml) | Socket Size | Dia. (mm) | Height (mm) | Pack Size |
|-----------|---------------|-------------|-----------|-------------|-----------|
| CH0428A   | 25            | 14/23       | 42        | 75          | 12        |
| CH0428A01 | 25            | 19/26       | 42        | 75          | 12        |
| CH0428B   | 50            | 19/26       | 51        | 85          | 12        |
| CH0428B01 | 50            | 14/23       | 51        | 85          | 12        |
| CH0428B02 | 50            | 24/29       | 51        | 85          | 12        |
| CH0428B03 | 50            | 29/32       | 51        | 85          | 12        |
| CH0428C   | 100           | 19/26       | 64        | 105         | 12        |
| CH0428C01 | 100           | 14/23       | 64        | 105         | 12        |
| CH0428D   | 100           | 24/29       | 64        | 105         | 12        |
| CH0428D01 | 100           | 29/32       | 64        | 105         | 12        |
| CH0428D02 | 200           | 29/32       | 79        | 130         | 12        |
| CH0428E01 | 250           | 14/23       | 85        | 140         | 12        |
| CH0428E   | 250           | 19/26       | 85        | 140         | 12        |
| CH0428F   | 250           | 24/29       | 85        | 140         | 12        |
| CH0428G   | 250           | 29/32       | 85        | 140         | 12        |
| CH0428G01 | 300           | 29/32       | 87        | 155         | 6         |
| CH0428H   | 500           | 24/29       | 105       | 175         | 6         |
| CH0428H01 | 500           | 19/26       | 105       | 175         | 6         |
| CH0428I   | 500           | 29/32       | 105       | 175         | 6         |
| CH0428I01 | 500           | 45/40       | 105       | 175         | 6         |
| CH0428J   | 1000          | 24/29       | 131       | 220         | 6         |
| CH0428K   | 1000          | 29/32       | 131       | 220         | 6         |
| CH0428L   | 1000          | 34/35       | 131       | 220         | 6         |
| CH0428L01 | 1000          | 45/40       | 131       | 220         | 6         |
| CH0428M   | 2000          | 29/32       | 166       | 280         | 4         |
| CH0428N   | 2000          | 34/35       | 166       | 280         | 4         |
| CH0428N01 | 2000          | 45/40       | 166       | 280         | 4         |
| CH0428O   | 3000          | 45/40       | 187       | 310         | 1         |
| CH0428P   | 5000          | 45/40       | 220       | 365         | 1         |



### Flasks - Conical Erlenmeyer, Hollow Glass Stopper

- Conical Erlenmeyer flask with narrow neck and interchangeable hollow glass stopper.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0429A | 100           | 29/32       | 2         |
| CH0429B | 250           | 29/32       | 2         |
| CH0429C | 500           | 29/32       | 2         |
| CH0429D | 1000          | 29/32       | 2         |
| CH0429E | 2000          | 29/32       | 2         |

### Flasks - Conical, Erlenmeyer

- Conical Erlenmeyer flask.
- Made of borosilicate glass.
- Includes teflon lined screw cap.

| Code    | Capacity (ml) | Thread GL | Dia. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-----------|-------------|-----------|
| CH0430B | 50            | 25        | 54        | 94          | 12        |
| CH0430C | 100           | 25        | 65        | 109         | 12        |
| CH0430D | 250           | 32        | 85        | 149         | 12        |
| CH0430E | 500           | 32        | 105       | 180         | 6         |
| CH0430F | 1000          | 32        | 131       | 225         | 6         |



GLASS BORO 3.3



### Flasks - Iodine

- Conical shaped flask used for iodine value determination.
- Made of borosilicate glass with interchangeable glass stopper.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0434A | 250           | 24/29       | 2         |
| CH0434B | 250           | 29/32       | 2         |
| CH0434C | 500           | 24/29       | 2         |
| CH0434D | 500           | 29/32       | 2         |
| CH0434E | 1000          | 29/32       | 2         |



GLASS BORO 3.3



### Flasks - Filtering

- Conical shaped filtering flask with integral side arm.
- Made of borosilicate glass.
- Heavy wall for higher mechanical strength.

| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Pack Size |
|---------|---------------|-----------|-------------|-----------|
| CH0420A | 100           | 64        | 105         | 12        |
| CH0420B | 250           | 85        | 155         | 12        |
| CH0420C | 500           | 105       | 185         | 6         |
| CH0420D | 1000          | 131       | 230         | 6         |
| CH0420E | 2000          | 166       | 260         | 1         |
| CH0420F | 5000          | 220       | 360         | 1         |



GLASS BORO 3.3






### Flasks - Filtering, Buchner

- Buchner filtering flask with integral side arm.
- Made of borosilicate glass.

| Code      | Capacity (ml) | Socket Size | Pack Size |
|-----------|---------------|-------------|-----------|
| CH0422A   | 100           | 19/26       | 12        |
| CH0422B   | 100           | 24/29       | 12        |
| CH0422C   | 250           | 24/29       | 12        |
| CH0422C01 | 250           | 29/32       | 12        |
| CH0422D   | 500           | 24/29       | 6         |
| CH0422D01 | 500           | 29/32       | 6         |
| CH0422E   | 1000          | 24/29       | 6         |
| CH0422F   | 1000          | 29/32       | 6         |



### Flasks - Boiling, Flat Bottom, Narrow Neck with Rim

- Flat bottomed Florence boiling flask with beaded narrow neck.
- Made of borosilicate glass.
- As per DIN ISO 1773 standard.

| Code    | Capacity (ml) | Dia.1 (mm) | Dia.2 (mm) | Height (mm) | Pack Size |
|---------|---------------|------------|------------|-------------|-----------|
| CH0410B | 50            | 51         | 26         | 90          | 12        |
| CH0410C | 100           | 64         | 26         | 105         | 12        |
| CH0410E | 250           | 85         | 34         | 138         | 12        |
| CH0410F | 500           | 105        | 34         | 163         | 6         |
| CH0410G | 1000          | 131        | 42         | 190         | 6         |
| CH0410H | 2000          | 166        | 42         | 250         | 4         |
| CH0410I | 3000          | 185        | 50         | 250         | 1         |
| CH0410J | 5000          | 223        | 50         | 290         | 1         |
| CH0410K | 6000          | 237        | 65         | 315         | 1         |
| CH0410L | 10,000        | 280        | 65         | 360         | 1         |
| CH0410M | 20,000        | 345        | 76         | 500         | 1         |

## Flasks Florence - Boiling, Flat Bottom, Short Neck with Joint

- Flat bottomed boiling flask with short neck and interchangeable joint.
- Made of borosilicate glass.
- As per DIN ISO 4797 standard.


| Code      | Capacity (ml) | Socket Size | Dia. (mm) | Height (mm) | Pack Size |
|-----------|---------------|-------------|-----------|-------------|-----------|
| CH0412A   | 50            | 14/23       | 51        | 85          | 12        |
| CH0412B   | 50            | 19/26       | 51        | 85          | 12        |
| CH0412B01 | 50            | 24/29       | 51        | 85          | 12        |
| CH0412B02 | 50            | 29/32       | 51        | 85          | 12        |
| CH0412C   | 100           | 14/23       | 64        | 103         | 12        |
| CH0412D   | 100           | 19/26       | 64        | 103         | 12        |
| CH0412E   | 100           | 24/29       | 64        | 103         | 12        |
| CH0412E01 | 100           | 29/32       | 64        | 103         | 12        |
| CH0412F   | 250           | 19/26       | 85        | 130         | 12        |
| CH0412F01 | 250           | 14/23       | 85        | 130         | 12        |
| CH0412G   | 250           | 24/29       | 85        | 130         | 12        |
| CH0412H   | 250           | 29/32       | 85        | 130         | 12        |
| CH0412I   | 500           | 24/29       | 105       | 160         | 6         |
| CH0412I01 | 500           | 19/26       | 105       | 160         | 6         |
| CH0412J   | 500           | 29/32       | 105       | 160         | 6         |
| CH0412K   | 1000          | 24/29       | 131       | 187         | 6         |
| CH0412L   | 1000          | 29/32       | 131       | 187         | 6         |
| CH0412L01 | 1000          | 34/35       | 131       | 187         | 6         |
| CH0412L02 | 1000          | 45/50       | 131       | 187         | 6         |
| CH0412M   | 2000          | 29/32       | 166       | 230         | 4         |
| CH0412N   | 2000          | 34/35       | 166       | 230         | 4         |
| CH0412O   | 3000          | 29/32       | 185       | 235         | 1         |
| CH0412P   | 3000          | 34/35       | 185       | 235         | 1         |
| CH0412Q   | 3000          | 45/40       | 185       | 235         | 1         |
| CH0412R   | 5000          | 29/32       | 223       | 280         | 1         |
| CH0412S   | 5000          | 34/35       | 223       | 280         | 1         |
| CH0412T   | 5000          | 45/40       | 223       | 280         | 1         |
| CH0412U   | 10000         | 45/40       | 279       | 350         | 1         |




## Flasks - Boiling, Flat Bottom, Wide Neck with Rim

- Flat bottomed Florence boiling flask with beaded wide neck rim.
- Made of borosilicate glass.
- As per DIN ISO 24450 standard.

| Code    | Capacity (ml) | Dia. <sub>1</sub> (mm) | Dia. <sub>2</sub> (mm) | Height (mm) | Pack Size |
|---------|---------------|------------------------|------------------------|-------------|-----------|
| CH0411A | 50            | 51                     | 34                     | 90          | 12        |
| CH0411B | 100           | 64                     | 34                     | 105         | 12        |
| CH0411C | 250           | 85                     | 50                     | 140         | 12        |
| CH0411D | 500           | 105                    | 50                     | 165         | 6         |
| CH0411E | 1000          | 131                    | 50                     | 190         | 6         |
| CH0411F | 2000          | 166                    | 76                     | 230         | 4         |





## Flasks - Boiling, Round Bottom, Narrow Neck with Rim

- Round bottom Florence boiling flask.
- Narrow neck & beaded rim
- Made of borosilicate glass.
- As per DIN ISO 1773 standard.

| Code    | Capacity (ml) | Dia. <sub>1</sub> (mm) | Dia. <sub>2</sub> (mm) | Height (mm) | Pack Size |
|---------|---------------|------------------------|------------------------|-------------|-----------|
| CH0402B | 50            | 51                     | 26                     | 95          | 12        |
| CH0402C | 100           | 64                     | 26                     | 110         | 12        |
| CH0402E | 250           | 85                     | 34                     | 145         | 12        |
| CH0402F | 500           | 105                    | 34                     | 170         | 6         |
| CH0402G | 1000          | 131                    | 42                     | 200         | 6         |
| CH0402H | 2000          | 166                    | 50                     | 260         | 4         |
| CH0402I | 3000          | 185                    | 50                     | 260         | 1         |
| CH0402J | 5000          | 223                    | 50                     | 305         | 1         |
| CH0402K | 6000          | 236                    | 65                     | 350         | 1         |
| CH0402L | 10,000        | 279                    | 65                     | 380         | 1         |
| CH0402M | 20,000        | 345                    | 76                     | 515         | 1         |



## Flasks Florence - Boiling, Round Bottom with Joint

- Round bottom boiling flask with interchangeable joint.
- Made of borosilicate glass.
- As per DIN ISO 4797 standard.

| Code      | Capacity (ml) | Socket Size | Dia. (mm) | Height (mm) | Pack Size |
|-----------|---------------|-------------|-----------|-------------|-----------|
| CH0404A   | 50            | 14/23       | 51        | 90          | 12        |
| CH0404B   | 50            | 19/26       | 51        | 90          | 12        |
| CH0404B01 | 50            | 24/29       | 51        | 90          | 12        |
| CH0404B02 | 50            | 29/32       | 51        | 90          | 12        |
| CH0404C   | 100           | 14/23       | 64        | 105         | 12        |
| CH0404D   | 100           | 19/26       | 64        | 105         | 12        |
| CH0404E   | 100           | 24/29       | 64        | 105         | 12        |
| CH0404E01 | 100           | 29/32       | 64        | 105         | 12        |
| CH0404F   | 250           | 19/26       | 85        | 140         | 12        |
| CH0404F01 | 250           | 14/23       | 85        | 140         | 12        |
| CH0404G   | 250           | 24/29       | 85        | 140         | 12        |
| CH0404H   | 250           | 29/32       | 85        | 140         | 12        |
| CH0404I   | 500           | 24/29       | 105       | 163         | 6         |
| CH0404I01 | 500           | 19/26       | 105       | 163         | 6         |
| CH0404J   | 500           | 29/32       | 105       | 163         | 6         |
| CH0404K   | 1000          | 24/29       | 131       | 200         | 6         |
| CH0404L   | 1000          | 29/32       | 131       | 200         | 6         |
| CH0404L01 | 1000          | 34/35       | 131       | 200         | 6         |
| CH0404M   | 2000          | 29/32       | 166       | 240         | 4         |
| CH0404N   | 2000          | 34/35       | 166       | 240         | 4         |
| CH0404O   | 3000          | 29/32       | 185       | 250         | 1         |
| CH0404P   | 3000          | 34/35       | 185       | 250         | 1         |
| CH0404Q   | 3000          | 45/40       | 185       | 250         | 1         |
| CH0404R   | 5000          | 29/32       | 223       | 300         | 1         |
| CH0404S   | 5000          | 34/35       | 223       | 300         | 1         |
| CH0404T   | 5000          | 45/40       | 223       | 300         | 1         |
| CH0404U   | 10000         | 45/40       | 279       | 375         | 1         |

## Flasks - Boiling, Round Bottom, Wide Neck, Beaded Rim

- Round bottom Florence boiling flask with wide neck and beaded rim.
- Made of borosilicate glass.
- As per DIN ISO 24450 standard.

| Code    | Capacity (ml) | Dia. <sub>1</sub> (mm) | Dia. <sub>2</sub> (mm) | Height (mm) | Pack Size |
|---------|---------------|------------------------|------------------------|-------------|-----------|
| CH0403A | 50            | 51                     | 34                     | 105         | 12        |
| CH0403B | 100           | 64                     | 34                     | 110         | 12        |
| CH0403C | 250           | 85                     | 50                     | 145         | 12        |
| CH0403D | 500           | 105                    | 50                     | 168         | 6         |
| CH0403E | 1000          | 131                    | 50                     | 200         | 6         |
| CH0403F | 2000          | 165                    | 76                     | 240         | 4         |



## Flasks - Boiling, Pear Shape with Joint

- Pear shape boiling flask with short neck and interchangeable joint.
- Made from borosilicate glass.
- As per DIN 12383 standard.

| Code      | Capacity (ml) | Socket Size | Dia. (mm) | Height (mm) | Pack Size |
|-----------|---------------|-------------|-----------|-------------|-----------|
| CH0406A01 | 10            | 14/23       | 30        | 75          | 12        |
| CH0406A   | 25            | 14/23       | 40        | 90          | 12        |
| CH0406B   | 50            | 14/23       | 50        | 110         | 12        |
| CH0406B01 | 50            | 19/26       | 50        | 110         | 12        |
| CH0406C   | 100           | 14/23       | 60        | 125         | 12        |
| CH0406C01 | 100           | 19/26       | 60        | 125         | 12        |



## Flasks - Boiling, Pear Shape, Two Neck with Joint

- Pear shape two neck boiling flask with joint.
- Made from borosilicate glass.
- As per DIN 12383.

| Code    | Capacity (ml) | Center Size | Side Socket | Height (mm) | Pack Size |
|---------|---------------|-------------|-------------|-------------|-----------|
| CH0406H | 50            | 14/23       | 14/23       | 110         | 6         |
| CH0406I | 100           | 14/23       | 14/23       | 125         | 6         |



GLASS  
BORO  
3.3



### Flasks - Distilling, Pear Shape with Joint

- Pear shape distilling flask with short neck and interchangeable joint.
- Made from borosilicate glass.

| Code    | Capacity (ml) | Center Size | Side Socket | Height (mm) | Pack Size |
|---------|---------------|-------------|-------------|-------------|-----------|
| CH0407A | 10            | 14/23       | 14/23       | 120         | 1         |
| CH0407B | 25            | 14/23       | 14/23       | 135         | 1         |
| CH0407C | 50            | 14/23       | 14/23       | 150         | 1         |
| CH0407D | 100           | 14/23       | 14/23       | 165         | 1         |

### Flasks - Distilling, Round Bottom with Joint

- Round bottom distilling flask & one oblique neck with joint.
- Made from borosilicate glass.
- As per DIN 12392 standard.

| Code      | Capacity (ml) | Center Socket | Side Sockets | Height (mm) | Pack Size |
|-----------|---------------|---------------|--------------|-------------|-----------|
| CH0419A   | 50            | 24/29         | 14/23        | 90          | 1         |
| CH0419A01 | 50            | 29/32         | 14/23        | 90          | 1         |
| CH0419B   | 100           | 24/29         | 14/23        | 105         | 1         |
| CH0419C   | 100           | 24/29         | 19/26        | 105         | 1         |
| CH0419C01 | 100           | 29/32         | 14/23        | 105         | 1         |
| CH0419D   | 250           | 24/29         | 19/26        | 140         | 1         |
| CH0419E   | 250           | 24/29         | 14/23        | 140         | 1         |
| CH0419E01 | 250           | 29/32         | 14/23        | 140         | 1         |
| CH0419F   | 500           | 24/29         | 19/26        | 165         | 1         |
| CH0419F01 | 500           | 24/29         | 14/23        | 165         | 1         |
| CH0419F02 | 500           | 29/32         | 14/23        | 165         | 1         |
| CH0419G   | 1000          | 24/29         | 19/26        | 200         | 1         |
| CH0419G01 | 1000          | 24/29         | 14/23        | 200         | 1         |
| CH0419G02 | 1000          | 29/32         | 14/23        | 200         | 1         |
| CH0419H   | 2000          | 34/35         | 19/26        | 240         | 1         |
| CH0419I   | 2000          | 29/32         | 14/23        | 240         | 1         |

DIN  
12392

GLASS  
BORO  
3.3



### Flasks - Distilling, Round Bottom, 3 Neck

- Round bottom distilling flask with 3 parallel necks.
- Made of borosilicate glass.
- As per DIN 12392 standard.

| Code      | Capacity (ml) | Center Socket | Side Sockets | Height (mm) | Pack Size |
|-----------|---------------|---------------|--------------|-------------|-----------|
| CH0418A   | 250           | 24/29         | 19/26        | 140         | 1         |
| CH0418B   | 250           | 29/32         | 14/23        | 140         | 1         |
| CH0418B01 | 250           | 29/32         | 29/32        | 140         | 1         |
| CH0418C   | 500           | 24/29         | 19/26        | 165         | 1         |
| CH0418D   | 500           | 29/32         | 14/23        | 165         | 1         |
| CH0418E   | 500           | 29/32         | 19/26        | 165         | 1         |
| CH0418E01 | 500           | 29/32         | 29/32        | 165         | 1         |
| CH0418F   | 1000          | 24/29         | 19/26        | 200         | 1         |
| CH0418G   | 1000          | 29/32         | 14/23        | 200         | 1         |
| CH0418H   | 1000          | 29/32         | 19/26        | 200         | 1         |
| CH0418H01 | 1000          | 29/32         | 24/29        | 200         | 1         |
| CH0418H02 | 1000          | 29/32         | 29/32        | 200         | 1         |
| CH0418I   | 2000          | 24/29         | 19/26        | 240         | 1         |
| CH0418J   | 2000          | 29/32         | 19/26        | 240         | 1         |
| CH0418J01 | 2000          | 29/32         | 29/32        | 240         | 1         |
| CH0418J02 | 2000          | 34/35         | 19/26        | 240         | 1         |
| CH0418J03 | 2000          | 34/35         | 24/29        | 240         | 1         |

DIN  
12392

GLASS  
BORO  
3.3



### Flasks - Distilling, Round Bottom, 3 Neck

- Round bottom distilling flask with 3 angled necks.
- Made of borosilicate glass.
- As per DIN 12394.

| Code      | Capacity (ml) | Center Socket | Side Sockets | Height (mm) | Pack Size |
|-----------|---------------|---------------|--------------|-------------|-----------|
| CH0418K   | 100           | 19/26         | 14/23        | 105         | 1         |
| CH0418L   | 100           | 24/29         | 14/23        | 105         | 1         |
| CH0418L01 | 100           | 24/29         | 19/26        | 105         | 1         |
| CH0418L02 | 100           | 29/32         | 14/23        | 105         | 1         |
| CH0418M   | 250           | 19/26         | 19/26        | 140         | 1         |
| CH0418N   | 250           | 24/29         | 14/23        | 140         | 1         |
| CH0418O   | 250           | 24/29         | 19/26        | 140         | 1         |
| CH0418O01 | 250           | 29/32         | 14/23        | 140         | 1         |
| CH0418O02 | 250           | 29/32         | 19/26        | 140         | 1         |
| CH0418O03 | 250           | 29/32         | 29/32        | 140         | 1         |
| CH0418P   | 500           | 24/29         | 14/23        | 165         | 1         |
| CH0418Q   | 500           | 24/29         | 19/26        | 165         | 1         |
| CH0418R   | 500           | 29/32         | 19/26        | 165         | 1         |
| CH0418R01 | 500           | 29/32         | 14/23        | 165         | 1         |
| CH0418R02 | 500           | 29/32         | 29/32        | 165         | 1         |
| CH0418S   | 1000          | 24/29         | 14/23        | 200         | 1         |
| CH0418T   | 1000          | 24/29         | 19/26        | 200         | 1         |
| CH0418U   | 1000          | 29/32         | 19/26        | 200         | 1         |
| CH0418U01 | 1000          | 29/32         | 14/23        | 200         | 1         |
| CH0418U02 | 1000          | 29/32         | 29/32        | 200         | 1         |
| CH0418V   | 2000          | 24/29         | 19/26        | 240         | 1         |
| CH0418W   | 2000          | 29/32         | 24/29        | 240         | 1         |
| CH0418X   | 2000          | 29/32         | 29/32        | 240         | 1         |



**DIN**  
12394

**GLASS**  
**BORO**  
3.3

**Chemical**  
**Resistant**

### Flasks - Distilling, Round Bottom with Joint

- Round bottom distilling flask, one oblique neck with joint.
- Made of borosilicate glass.
- As per DIN 12394.

| Code      | Capacity (ml) | Center Socket | Side Sockets | Height (mm) | Pack Size |
|-----------|---------------|---------------|--------------|-------------|-----------|
| CH0419A   | 50            | 24/29         | 14/23        | 90          | 1         |
| CH0419A01 | 50            | 29/32         | 14/23        | 90          | 1         |
| CH0419B   | 100           | 24/29         | 14/23        | 105         | 1         |
| CH0419C   | 100           | 24/29         | 19/26        | 105         | 1         |
| CH0419C01 | 100           | 29/32         | 14/23        | 105         | 1         |
| CH0419D   | 250           | 24/29         | 19/26        | 140         | 1         |
| CH0419E   | 250           | 24/29         | 14/23        | 140         | 1         |
| CH0419E01 | 250           | 29/32         | 14/23        | 140         | 1         |
| CH0419F   | 500           | 24/29         | 19/26        | 165         | 1         |
| CH0419F01 | 500           | 24/29         | 14/23        | 165         | 1         |
| CH0419F02 | 500           | 29/32         | 14/23        | 165         | 1         |
| CH0419G   | 1000          | 24/29         | 19/26        | 200         | 1         |
| CH0419G01 | 1000          | 24/29         | 14/23        | 200         | 1         |
| CH0419G02 | 1000          | 29/32         | 14/23        | 200         | 1         |
| CH0419H   | 2000          | 34/35         | 19/26        | 240         | 1         |
| CH0419I   | 2000          | 29/32         | 14/23        | 240         | 1         |



**DIN**  
12394

**GLASS**  
**BORO**  
3.3


**Chemical**  
**Resistant**



## Flasks - Evaporating

- Evaporating flask made of borosilicate glass.

| Code       | Capacity (ml) | Socket Size | Dia. (mm) | Height (mm) | Pack Size |
|------------|---------------|-------------|-----------|-------------|-----------|
| CH0435A050 | 50            | 29/32       | 48        | 94          | 12        |
| CH0435A    | 100           | 24/29       | 60        | 110         | 12        |
| CH0435A01  | 100           | 29/32       | 60        | 110         | 12        |
| CH0435B    | 250           | 24/29       | 81        | 140         | 12        |
| CH0435C    | 250           | 29/32       | 81        | 140         | 12        |
| CH0435D    | 500           | 24/29       | 101       | 170         | 6         |
| CH0435E    | 500           | 29/32       | 101       | 170         | 6         |
| CH0435F    | 1000          | 24/29       | 126       | 210         | 6         |
| CH0435G    | 1000          | 29/32       | 126       | 210         | 6         |



## Flasks - Distilling with Side Arm

- Round bottom distilling flask with side arm.
- Made from borosilicate glass.

| Code      | Capacity (ml) | Dia. (mm) | Height (mm) | Length of Side Arm (mm) | Pack Size |
|-----------|---------------|-----------|-------------|-------------------------|-----------|
| CH0416A25 | 25            | 41        | 110         | 100                     | 2         |
| CH0416A50 | 50            | 51        | 120         | 100                     | 2         |
| CH0416A   | 100           | 64        | 215         | 100                     | 2         |
| CH0416B   | 250           | 85        | 225         | 125                     | 2         |
| CH0416C   | 500           | 105       | 240         | 140                     | 1         |
| CH0416D   | 1000          | 131       | 300         | 250                     | 1         |



## Flasks - Kjeldahl

- Round bottom
- Long neck with joint Kjeldahl flask designed for nitrogen determination as per Kjeldahl method.
- Made of borosilicate glass.

| Code      | Capacity (ml) | Socket Size | Dia. (mm) | Height (mm) | Pack Size |
|-----------|---------------|-------------|-----------|-------------|-----------|
| CH0436A   | 50            | 19/26       | 51        | 160         | 2         |
| CH0436B   | 100           | 19/26       | 60        | 165         | 2         |
| CH0436B01 | 100           | 24/29       | 60        | 260         | 2         |
| CH0436C   | 300           | 24/29       | 81        | 290         | 2         |
| CH0436D   | 500           | 24/29       | 101       | 305         | 2         |
| CH0436E   | 800           | 24/29       | 115       | 325         | 2         |
| CH0436F   | 800           | 29/32       | 115       | 325         | 2         |

### Flasks - Kjeldahl, Plain

- Round bottom
- Long neck Kjeldahl flask designed for nitrogen determination as per Kjeldahl method.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Dia. (mm) | Height (mm) | Neck O.D. (mm) | Pack Size |
|---------|---------------|-----------|-------------|----------------|-----------|
| CH0438A | 100           | 60        | 240         | 22             | 2         |
| CH0438B | 300           | 81        | 300         | 28             | 2         |
| CH0438C | 500           | 101       | 325         | 34             | 2         |
| CH0438D | 800           | 115       | 350         | 34             | 2         |



GLASS BORO 3.3



### Filter Funnels

- Borosilicate glass, plain & 60° angle.

| Code      | D <sub>1</sub> (mm) | D <sub>2</sub> (mm) | Pack Size |
|-----------|---------------------|---------------------|-----------|
| CH0456A01 | 25                  | 6                   | 10        |
| CH0456A02 | 35                  | 6                   | 10        |
| CH0456A03 | 45                  | 7                   | 10        |
| CH0456A   | 50                  | 7                   | 10        |
| CH0456B01 | 55                  | 8                   | 10        |
| CH0456B   | 65                  | 8                   | 10        |
| CH0456C   | 75                  | 8                   | 10        |
| CH0456C02 | 85                  | 9                   | 10        |
| CH0456C01 | 90                  | 9                   | 10        |
| CH0456D   | 100                 | 10                  | 10        |
| CH0456E   | 125                 | 16                  | 2         |
| CH0456F   | 150                 | 16                  | 2         |




GLASS BORO 3.3



### Filter Funnels - Heavy

- Made of borosilicate glass.
- Constructed with thick uniform walls for greater mechanical strength.


| Code    | D <sub>1</sub> (mm) | D <sub>2</sub> (mm) | Pack Size |
|---------|---------------------|---------------------|-----------|
| CH0458A | 35                  | 6                   | 10        |
| CH0458B | 55                  | 7                   | 10        |
| CH0458C | 75                  | 8                   | 10        |
| CH0458D | 100                 | 10                  | 10        |
| CH0458E | 125                 | 16                  | 2         |
| CH0458F | 150                 | 16                  | 2         |



GLASS BORO 3.3



GLASS  
BORO  
3.3



### Powder Funnel

- Made of borosilicate glass, plain & 60° angle.

| Code    | D <sub>1</sub><br>(mm) | D <sub>2</sub><br>(mm) | Height<br>(mm) | Pack<br>Size |
|---------|------------------------|------------------------|----------------|--------------|
| CH0459A | 45                     | 20                     | 30             | 10           |
| CH0459B | 55                     | 20                     | 30             | 10           |
| CH0459C | 75                     | 22                     | 30             | 10           |
| CH0459D | 85                     | 24                     | 30             | 10           |
| CH0459E | 100                    | 26                     | 30             | 10           |
| CH0459F | 120                    | 34                     | 30             | 2            |

GLASS  
BORO  
3.3



### Powder Funnels - Jointed

- Made of borosilicate glass & 60° angle with joint.

| Code    | Dia. <sub>1</sub><br>(mm) | Joint<br>Size | Pack<br>Size |
|---------|---------------------------|---------------|--------------|
| CH0460A | 45                        | 14/23         | 10           |
| CH0460B | 55                        | 14/23         | 10           |
| CH0460C | 70                        | 14/23         | 10           |
| CH0460D | 70                        | 29/32         | 10           |
| CH0460E | 100                       | 14/23         | 10           |
| CH0460F | 100                       | 29/32         | 10           |

GLASS  
BORO  
3.3



### Measuring Scoops

- Made of borosilicate glass with open tubular arm.
- Suitable for weighing small quantities of powder and dyes.

| Code    | Capacity<br>(ml) | Pack<br>Size |
|---------|------------------|--------------|
| CH0461J | 3                | 10           |
| CH0461K | 6                | 10           |
| CH0461L | 10               | 10           |

GLASS  
BORO  
3.3



### Gas Jars

- Made of borosilicate glass with heavy foot and grounded flange.

| Code    | Dia.<br>(mm) | Height<br>(mm) | Pack<br>Size |
|---------|--------------|----------------|--------------|
| CH0495A | 50           | 150            | 1            |
| CH0495B | 50           | 200            | 1            |
| CH0495C | 50           | 250            | 1            |

### Flasks, Volumetric Class-A with Stopper

- Flat bottom class- A volumetric flask with one graduation mark.
- Fitted with chemical resistant polypropylene stopper.
- Made of borosilicate glass.
- As per DIN ISO 1042 standard.

| Code          |                | Capacity | Tolerance | Socket | Pack |
|---------------|----------------|----------|-----------|--------|------|
| Blue Printing | White Printing | (ml)     | ± ml      | Size   | Size |
| CH0446A       | CH0446A-WT     | 5        | 0.025     | 10/19  | 10   |
| CH0446B       | CH0446B-WT     | 10       | 0.025     | 10/19  | 10   |
| CH0446B01     | CH0446B01-WT   | 20       | 0.040     | 10/19  | 10   |
| CH0446C       | CH0446C-WT     | 25       | 0.040     | 10/19  | 10   |
| CH0446D       | CH0446D-WT     | 50       | 0.060     | 12/21  | 2    |
| CH0446E       | CH0446E-WT     | 100      | 0.100     | 14/23  | 2    |
| CH0446F       | CH0446F-WT     | 200      | 0.150     | 14/23  | 2    |
| CH0446G       | CH0446G-WT     | 250      | 0.150     | 14/23  | 2    |
| CH0446H       | CH0446H-WT     | 500      | 0.250     | 19/26  | 2    |
| CH0446I       | CH0446I-WT     | 1000     | 0.400     | 24/29  | 2    |
| CH0446J       | CH0446J-WT     | 2000     | 0.600     | 29/32  | 1    |
| CH0446K       | CH0446K-WT     | 5000     | 1.200     | 34/35  | 1    |



### Flasks, Volumetric Class-B with Stopper

- Flat bottom class- B volumetric flask with one graduation mark.
- Fitted with chemical resistant polypropylene stopper.
- Made of borosilicate glass.
- As per DIN ISO 1042 standard.

| Code          |                | Capacity | Tolerance | Socket | Pack |
|---------------|----------------|----------|-----------|--------|------|
| Blue Printing | White Printing | (ml)     | ± ml      | Size   | Size |
| CH0448A       | CH0448A-WT     | 5        | 0.05      | 10/19  | 10   |
| CH0448B       | CH0448B-WT     | 10       | 0.05      | 10/19  | 10   |
| CH0448B01     | CH0448B01-WT   | 20       | 0.08      | 10/19  | 10   |
| CH0448C       | CH0448C-WT     | 25       | 0.08      | 10/19  | 10   |
| CH0448D       | CH0448D-WT     | 50       | 0.12      | 12/21  | 2    |
| CH0448E       | CH0448E-WT     | 100      | 0.20      | 14/23  | 2    |
| CH0448F       | CH0448F-WT     | 200      | 0.30      | 14/23  | 2    |
| CH0448G       | CH0448G-WT     | 250      | 0.30      | 14/23  | 2    |
| CH0448H       | CH0448H-WT     | 500      | 0.50      | 19/26  | 2    |
| CH0448I       | CH0448I-WT     | 1000     | 0.80      | 24/29  | 2    |
| CH0448J       | CH0448J-WT     | 2000     | 1.20      | 29/32  | 1    |
| CH0448K       | CH0448K-WT     | 5000     | 2.40      | 34/35  | 1    |



BATCH  
CERTIFICATES

GLASS  
BORO  
3.3


### Flasks, Volumetric Class-A

- Flat bottom class- A volumetric flask with one graduation mark and rimmed neck.
- Made of borosilicate glass.
- As per DIN ISO 1042 standard.

| Code      | Capacity (ml) | Tolerance $\pm$ ml | Pack Size |
|-----------|---------------|--------------------|-----------|
| CH0449A   | 5             | 0.025              | 10        |
| CH0449B   | 10            | 0.025              | 10        |
| CH0449B01 | 20            | 0.040              | 10        |
| CH0449C   | 25            | 0.040              | 10        |
| CH0449D   | 50            | 0.060              | 2         |
| CH0449E   | 100           | 0.100              | 2         |
| CH0449F   | 200           | 0.150              | 2         |
| CH0449G   | 250           | 0.150              | 2         |
| CH0449H   | 500           | 0.250              | 2         |
| CH0449I   | 1000          | 0.400              | 2         |
| CH0449J   | 2000          | 0.600              | 1         |


GLASS  
BORO  
3.3


### Flasks, Volumetric Class-B

- Flat bottom class- B volumetric flask with one graduation mark and rimmed neck.
- Made of borosilicate glass.
- As per DIN ISO 1042 standard.

| Code      | Capacity (ml) | Tolerance $\pm$ ml | Pack Size |
|-----------|---------------|--------------------|-----------|
| CH0450A   | 5             | 0.05               | 10        |
| CH0450B   | 10            | 0.05               | 10        |
| CH0450B01 | 20            | 0.08               | 10        |
| CH0450C   | 25            | 0.08               | 10        |
| CH0450D   | 50            | 0.12               | 2         |
| CH0450E   | 100           | 0.20               | 2         |
| CH0450F   | 200           | 0.30               | 2         |
| CH0450G   | 250           | 0.30               | 2         |
| CH0450H   | 500           | 0.50               | 2         |
| CH0450I   | 1000          | 0.80               | 2         |
| CH0450J   | 2000          | 1.20               | 1         |

BATCH  
CERTIFICATES

GLASS  
BORO  
3.3


### Flasks, Volumetric, Amber with Stopper, Class-A

- Flat bottom amber volumetric flask with one graduation mark.
- Fitted with chemical resistant polypropylene stopper.
- Made of amber borosilicate glass.
- As per DIN ISO 1042 standard.

| Code      | Capacity (ml) | Tolerance $\pm$ ml | Socket Size | Pack Size |
|-----------|---------------|--------------------|-------------|-----------|
| CH0451A   | 5             | 0.025              | 10/19       | 10        |
| CH0451B   | 10            | 0.025              | 10/19       | 10        |
| CH0451B01 | 20            | 0.040              | 10/19       | 10        |
| CH0451C   | 25            | 0.040              | 10/19       | 10        |
| CH0451D   | 50            | 0.060              | 12/21       | 2         |
| CH0451E   | 100           | 0.100              | 14/23       | 2         |
| CH0451F   | 200           | 0.150              | 14/23       | 2         |
| CH0451G   | 250           | 0.150              | 14/23       | 2         |
| CH0451H   | 500           | 0.250              | 19/26       | 2         |
| CH0451I   | 1000          | 0.400              | 24/29       | 2         |
| CH0451J   | 2000          | 0.600              | 29/32       | 1         |

## Flasks, Volumetric Class- A with Hollow Glass Stopper

- Flat bottom class- A volumetric flask with one graduation mark.
- Fitted with hollow hexagonal glass stopper.
- Made of borosilicate glass.
- As per DIN ISO 1042 standard.

| Code          |                | Capacity | Tolerance | Socket | Pack |
|---------------|----------------|----------|-----------|--------|------|
| Blue Printing | White Printing | (ml)     | ± ml      | Size   | Size |
| CH0439A       | CH0439A-WT     | 5        | 0.025     | 10/19  | 10   |
| CH0439B       | CH0439B-WT     | 10       | 0.025     | 10/19  | 10   |
| CH0439C       | CH0439C-WT     | 20       | 0.040     | 10/19  | 10   |
| CH0439D       | CH0439D-WT     | 25       | 0.040     | 10/19  | 10   |
| CH0439E       | CH0439E-WT     | 50       | 0.060     | 12/21  | 2    |
| CH0439F       | CH0439F-WT     | 100      | 0.100     | 14/23  | 2    |
| CH0439G       | CH0439G-WT     | 200      | 0.150     | 14/23  | 2    |
| CH0439H       | CH0439H-WT     | 250      | 0.150     | 14/23  | 2    |
| CH0439I       | CH0439I-WT     | 500      | 0.250     | 19/26  | 2    |
| CH0439J       | CH0439J-WT     | 1000     | 0.400     | 24/29  | 2    |
| CH0439K       | CH0439K-WT     | 2000     | 0.600     | 29/32  | 1    |
| CH0439L       | CH0439L-WT     | 5000     | 1.20      | 34/35  | 1    |



BATCH  
CERTIFICATES



GLASS  
BORO  
3.3



## Flasks, Volumetric Class- B with Hollow Glass Stopper

- Flat bottom class-B volumetric flask with one graduation mark.
- Fitted with hollow hexagonal glass stopper.
- Made of borosilicate glass.
- As per DIN ISO 1042 standard.

| Code          |                | Capacity | Tolerance | Socket | Pack |
|---------------|----------------|----------|-----------|--------|------|
| Blue Printing | White Printing | (ml)     | ± ml      | Size   | Size |
| CH0440A       | CH0440A-WT     | 5        | 0.05      | 10/19  | 10   |
| CH0440B       | CH0440B-WT     | 10       | 0.05      | 10/19  | 10   |
| CH0440C       | CH0440C-WT     | 20       | 0.08      | 10/19  | 10   |
| CH0440D       | CH0440D-WT     | 25       | 0.08      | 10/19  | 10   |
| CH0440E       | CH0440E-WT     | 50       | 0.12      | 12/21  | 2    |
| CH0440F       | CH0440F-WT     | 100      | 0.20      | 14/23  | 2    |
| CH0440G       | CH0440G-WT     | 200      | 0.30      | 14/23  | 2    |
| CH0440H       | CH0440H-WT     | 250      | 0.30      | 14/23  | 2    |
| CH0440I       | CH0440I-WT     | 500      | 0.50      | 19/26  | 2    |
| CH0440J       | CH0440J-WT     | 1000     | 0.80      | 24/29  | 2    |
| CH0440K       | CH0440K-WT     | 2000     | 1.20      | 29/32  | 1    |
| CH0440L       | CH0440L-WT     | 5000     | 2.40      | 34/35  | 1    |



GLASS  
BORO  
3.3



BATCH  
CERTIFICATES


ASTM

GLASS  
BORO  
3.3

Chemical  
Resistant


### Flasks, Volumetric Class- A ASTM with Solid Glass Stopper

- Flat bottom class- A ASTM volumetric flask with one graduation mark.
- Fitted with solid glass stopper.
- Made of borosilicate glass.
- As per ASTM-E288 standard.

| Code      | Capacity<br>(ml) | Tolerance<br>± ml | Stopper<br>No. | Pack<br>Size |
|-----------|------------------|-------------------|----------------|--------------|
| CH0441A01 | 5                | 0.020             | 9              | 10           |
| CH0441A02 | 10               | 0.020             | 9              | 10           |
| CH0441A   | 25               | 0.030             | 9              | 10           |
| CH0441B   | 50               | 0.050             | 13             | 2            |
| CH0441C   | 100              | 0.080             | 13             | 2            |
| CH0441D   | 200              | 0.100             | 16             | 2            |
| CH0441E   | 250              | 0.120             | 16             | 2            |
| CH0441F   | 500              | 0.200             | 19             | 2            |
| CH0441G   | 1000             | 0.300             | 22             | 2            |
| CH0441H   | 2000             | 0.500             | 27             | 1            |

### Flasks, Volumetric Class-B ASTM with Solid Glass Stopper

- Flat bottom class-B ASTM volumetric flask with one graduation mark.
- Fitted with solid glass stopper.
- Made of borosilicate glass.
- As per ASTM-E288 standard.

| Code      | Capacity<br>(ml) | Tolerance<br>± ml | Stopper<br>No. | Pack<br>Size |
|-----------|------------------|-------------------|----------------|--------------|
| CH0442A01 | 5                | 0.040             | 9              | 10           |
| CH0442A02 | 10               | 0.040             | 9              | 10           |
| CH0442A   | 25               | 0.060             | 9              | 10           |
| CH0442B   | 50               | 0.100             | 13             | 2            |
| CH0442C   | 100              | 0.160             | 13             | 2            |
| CH0442D   | 200              | 0.200             | 16             | 2            |
| CH0442E   | 250              | 0.240             | 16             | 2            |
| CH0442F   | 500              | 0.400             | 19             | 2            |
| CH0442G   | 1000             | 0.600             | 22             | 2            |
| CH0442H   | 2000             | 1.000             | 27             | 1            |

ASTM

GLASS  
BORO  
3.3

Chemical  
Resistant


### Flasks, Volumetric Class-A ASTM with PP Stopper

- Flat bottom class- A volumetric flask with one graduation mark.
- Fitted with chemical resistant polypropylene stopper.
- Made of borosilicate glass.
- As per ASTM-E288 standard.

| Code          |                | Capacity | Tolerance | Stopper | Pack |
|---------------|----------------|----------|-----------|---------|------|
| Blue Printing | White Printing | (ml)     | ± ml      | No.     | Size |
| CH0445A01     | CH0445A01-WT   | 5        | 0.020     | 9       | 10   |
| CH0445A02     | CH0445A02-WT   | 10       | 0.020     | 9       | 10   |
| CH0445A       | CH0445A-WT     | 25       | 0.030     | 9       | 10   |
| CH0445B       | CH0445B-WT     | 50       | 0.050     | 13      | 2    |
| CH0445C       | CH0445C-WT     | 100      | 0.080     | 13      | 2    |
| CH0445D       | CH0445D-WT     | 200      | 0.100     | 16      | 2    |
| CH0445E       | CH0445E-WT     | 250      | 0.120     | 16      | 2    |
| CH0445F       | CH0445F-WT     | 500      | 0.200     | 19      | 2    |
| CH0445G       | CH0445G-WT     | 1000     | 0.300     | 22      | 2    |
| CH0445H       | CH0445H-WT     | 2000     | 0.500     | 27      | 1    |



ASTM

GLASS  
BORO  
3.3



### Flasks, Volumetric Class-B ASTM with PP Stopper

- Flat bottom class-B volumetric flask with one graduation mark.
- Fitted with chemical resistant polypropylene stopper.
- Made of borosilicate glass.
- As per ASTM-E288 standard.

| Code          |                | Capacity | Tolerance | Stopper | Pack |
|---------------|----------------|----------|-----------|---------|------|
| Blue Printing | White Printing | (ml)     | ± ml      | No.     | Size |
| CH0445I01     | CH0445I01-WT   | 5        | 0.040     | 9       | 10   |
| CH0445I02     | CH0445I02-WT   | 10       | 0.040     | 9       | 10   |
| CH0445I       | CH0445I-WT     | 25       | 0.060     | 9       | 10   |
| CH0445J       | CH0445J-WT     | 50       | 0.100     | 13      | 2    |
| CH0445K       | CH0445K-WT     | 100      | 0.160     | 13      | 2    |
| CH0445L       | CH0445L-WT     | 200      | 0.200     | 16      | 2    |
| CH0445M       | CH0445M-WT     | 250      | 0.240     | 16      | 2    |
| CH0445N       | CH0445N-WT     | 500      | 0.400     | 19      | 2    |
| CH0445O       | CH0445O-WT     | 1000     | 0.600     | 22      | 2    |
| CH0445P       | CH0445P-WT     | 2000     | 1.000     | 27      | 1    |



ASTM

GLASS  
BORO  
3.3



BATCH  
CERTIFICATES


ASTM

GLASS  
BORO  
3.3

Chemical  
Resistant


### Flasks, Volumetric Class-A ASTM with Snap Cap

- Flat bottom class-A volumetric flask with one graduation mark.
- Fitted with snap cap.
- Made of borosilicate glass.
- As per ASTM-E288 standard.

| Code      | Capacity<br>(ml) | Tolerance<br>± ml | Snap<br>Cap No. | Pack<br>Size |
|-----------|------------------|-------------------|-----------------|--------------|
| CH0443A01 | 5                | 0.020             | 1               | 10           |
| CH0443A02 | 10               | 0.020             | 1               | 10           |
| CH0443A03 | 25               | 0.030             | 2               | 10           |
| CH0443A   | 50               | 0.050             | 3               | 2            |
| CH0443B   | 100              | 0.080             | 3               | 2            |
| CH0443C   | 200              | 0.100             | 4               | 2            |
| CH0443D   | 250              | 0.120             | 4               | 2            |
| CH0443E   | 500              | 0.200             | 5               | 2            |
| CH0443F   | 1000             | 0.300             | 6               | 2            |
| CH0443G   | 2000             | 0.500             | 7               | 1            |

ASTM

GLASS  
BORO  
3.3

Chemical  
Resistant


### Flasks, Volumetric Class-B ASTM with Snap Cap

- Flat bottom class-B volumetric flask with one graduation mark.
- Fitted with snap cap.
- Made of borosilicate glass.
- As per ASTM-E288 standard.

| Code      | Capacity<br>(ml) | Tolerance<br>± ml | Stopper<br>No. | Pack<br>Size |
|-----------|------------------|-------------------|----------------|--------------|
| CH0444A01 | 5                | 0.020             | 1              | 10           |
| CH0444A02 | 10               | 0.020             | 1              | 10           |
| CH0444A03 | 25               | 0.030             | 2              | 10           |
| CH0444A   | 50               | 0.100             | 3              | 2            |
| CH0444B   | 100              | 0.160             | 3              | 2            |
| CH0444C   | 200              | 0.200             | 4              | 2            |
| CH0444D   | 250              | 0.240             | 4              | 2            |
| CH0444E   | 500              | 0.400             | 5              | 2            |
| CH0444F   | 1000             | 0.600             | 6              | 2            |
| CH0444G   | 2000             | 1.000             | 7              | 1            |

### Burettes - Glass with Pinch Clip

- Single bore burettes with detachable rubber tube with pinch clip & class-B.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH0234A | 10            | 0.05           | 0.05               | 10        |
| CH0234B | 25            | 0.10           | 0.10               | 10        |
| CH0234C | 50            | 0.10           | 0.10               | 10        |
| CH0234D | 100           | 0.20           | 0.20               | 10        |



Check out  
our new website  
[www.eiscoindustrial.com](http://www.eiscoindustrial.com)



GLASS  
BORO  
3.3



### Burettes - Glass with Glass Key Stopcock

- Single straight bore burettes with glass key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385 standard.

#### Class-A

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH0236A | 10            | 0.05           | 0.03               | 1         |
| CH0236B | 25            | 0.10           | 0.05               | 1         |
| CH0236C | 50            | 0.10           | 0.05               | 1         |
| CH0236D | 100           | 0.20           | 0.10               | 1         |

#### Class-B

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH0236E | 10            | 0.05           | 0.05               | 1         |
| CH0236F | 25            | 0.10           | 0.10               | 1         |
| CH0236G | 50            | 0.10           | 0.10               | 1         |
| CH0236H | 100           | 0.20           | 0.20               | 1         |



BATCH  
CERTIFICATES



GLASS  
BORO  
3.3



BATCH  
CERTIFICATESGLASS  
BORO  
3.3Chemical  
Resistant

### Burettes - Glass with Needle Valve Stopcock

- Burettes with PTFE needle valve stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385, standard.

#### Class-A

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0238A        | CH0238A-WT      | 10       | 0.05      | 0.03      | 1    |
| CH0238B        | CH0238B-WT      | 25       | 0.10      | 0.05      | 1    |
| CH0238C        | CH0238C-WT      | 50       | 0.10      | 0.05      | 1    |
| CH0238D        | CH0238D-WT      | 100      | 0.20      | 0.10      | 1    |

#### Class-B

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0238E        | CH0238E-WT      | 10       | 0.05      | 0.05      | 1    |
| CH0238F        | CH0238F-WT      | 25       | 0.10      | 0.10      | 1    |
| CH0238G        | CH0238G-WT      | 50       | 0.10      | 0.10      | 1    |
| CH0238H        | CH0238H-WT      | 100      | 0.20      | 0.20      | 1    |

BATCH  
CERTIFICATESGLASS  
BORO  
3.3Chemical  
Resistant

### Burettes - Glass with PTFE Key Stopcock

- Single straight bore burettes with interchangeable PTFE key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385 standard.

#### Class-A

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0240A        | CH0240A-WT      | 10       | 0.05      | 0.03      | 1    |
| CH0240B        | CH0240B-WT      | 25       | 0.10      | 0.05      | 1    |
| CH0240C        | CH0240C-WT      | 50       | 0.10      | 0.05      | 1    |
| CH0240D        | CH0240D-WT      | 100      | 0.20      | 0.10      | 1    |

#### Class-B

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0240E        | CH0240E-WT      | 10       | 0.05      | 0.05      | 1    |
| CH0240F        | CH0240F-WT      | 25       | 0.10      | 0.10      | 1    |
| CH0240G        | CH0240G-WT      | 50       | 0.10      | 0.10      | 1    |
| CH0240H        | CH0240H-WT      | 100      | 0.20      | 0.20      | 1    |

## Burettes - Schellbach, Glass with Glass Key Stopcock

- Single bore Schellbach burettes with glass key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385, standard.

### Class-A

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0241D01      | CH0241D01-WT    | 10       | 0.05      | 0.03      | 1    |
| CH0241D        | CH0241D-WT      | 25       | 0.10      | 0.05      | 1    |
| CH0241E        | CH0241E-WT      | 50       | 0.10      | 0.05      | 1    |
| CH0241F        | CH0241F-WT      | 100      | 0.20      | 0.10      | 1    |

### Class-B

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0241A01      | CH0241A01-WT    | 10       | 0.05      | 0.05      | 1    |
| CH0241A        | CH0241A-WT      | 25       | 0.10      | 0.10      | 1    |
| CH0241B        | CH0241B-WT      | 50       | 0.10      | 0.10      | 1    |
| CH0241C        | CH0241C-WT      | 100      | 0.20      | 0.20      | 1    |

## Burettes - Schellbach, Glass with PTFE Key Stopcock

- Single bore Schellbach burettes with PTFE key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385 standard.

### Class-A

| Code      | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|-----------|---------------|----------------|--------------------|-----------|
| CH0241J01 | 10            | 0.05           | 0.03               | 1         |
| CH0241J   | 25            | 0.10           | 0.05               | 1         |
| CH0241K   | 50            | 0.10           | 0.05               | 1         |
| CH0241L   | 100           | 0.20           | 0.10               | 1         |

### Class-B

| Code      | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|-----------|---------------|----------------|--------------------|-----------|
| CH0241M01 | 10            | 0.05           | 0.05               | 1         |
| CH0241M   | 25            | 0.10           | 0.10               | 1         |
| CH0241N   | 50            | 0.10           | 0.10               | 1         |
| CH0241O   | 100           | 0.20           | 0.20               | 1         |

## Burettes - Modular, Glass with PTFE Key Stopcock

- Modular burettes with PTFE key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385 standard.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH0243A | 10            | 0.05           | 0.03               | 1         |
| CH0243B | 25            | 0.10           | 0.05               | 1         |
| CH0243C | 50            | 0.10           | 0.05               | 1         |
| CH0243D | 100           | 0.20           | 0.10               | 1         |



BATCH CERTIFICATES



BATCH CERTIFICATES



BATCH CERTIFICATES



BATCH  
CERTIFICATESGLASS  
BORO  
3.3Chemical  
Resistant

### Micro Burettes

- Class-A micro burettes designed for precise analytical methods in microchemistry.
- Manufactured from precision bore tubing to assure accuracy between any two points throughout the body.
- Funnel top facilitates addition of reagents by pouring or by suction from below through a rubber-stoppered tube.
- Tips are tapered to assure smooth flow.
- The capacity tolerance on these burettes is established by ASTM E-1189 or ISO 385.

Class-A

| Code    | Capacity<br>(ml) | Sub. Div.<br>(ml) | Tolerance<br>± ml | Pack<br>Size |
|---------|------------------|-------------------|-------------------|--------------|
| CH0245A | 5                | 0.02              | 0.01              | 1            |
| CH0245B | 10               | 0.02              | 0.02              | 1            |

ASTM

GLASS  
BORO  
3.3Chemical  
Resistant

### Burettes - Glass, Class-B ASTM with PTFE Key Stopcock

- Single bore burettes with PTFE key stopcock,
- Made of borosilicate glass.
- Accuracy as per. ASTM-E287 compliance standards.

Class-A

| Code    | Capacity<br>(ml) | Sub. Div.<br>(ml) | Tolerance<br>± ml | Pack<br>Size |
|---------|------------------|-------------------|-------------------|--------------|
| CH0245E | 25               | 0.10              | 0.06              | 1            |
| CH0245F | 30               | 0.10              | 0.10              | 1            |

### Automatic Zero Burettes - Glass with PTFE Key Stopcock

- Automatic zero burettes mounted on glass reservoir with rubber below.
- PTFE key stopcock
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385.

#### Class-A

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH0242A | 10            | 0.05           | 0.03               | 1         |
| CH0242B | 25            | 0.10           | 0.05               | 1         |
| CH0242C | 50            | 0.10           | 0.05               | 1         |
| CH0242D | 100           | 0.20           | 0.10               | 1         |

#### Class-B

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH0242E | 10            | 0.05           | 0.05               | 1         |
| CH0242F | 25            | 0.10           | 0.10               | 1         |
| CH0242G | 50            | 0.10           | 0.10               | 1         |
| CH0242H | 100           | 0.20           | 0.20               | 1         |



BATCH CERTIFICATES



### Automatic Zero Burettes - Glass with Intermediate PTFE Key Stopcock

- Automatic zero burettes mounted on glass reservoir with rubber below.
- Intermediate PTFE key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385.

#### Class-A

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH0242P | 10            | 0.05           | 0.03               | 1         |
| CH0242Q | 25            | 0.10           | 0.05               | 1         |
| CH0242R | 50            | 0.10           | 0.05               | 1         |
| CH0242S | 100           | 0.20           | 0.10               | 1         |

#### Class-B

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH0242T | 10            | 0.05           | 0.05               | 1         |
| CH0242U | 25            | 0.10           | 0.10               | 1         |
| CH0242V | 50            | 0.10           | 0.10               | 1         |
| CH0242W | 100           | 0.20           | 0.20               | 1         |



BATCH CERTIFICATES



GLASS  
BORO  
3.3BATCH  
CERTIFICATES

### Automatic Zero Burettes - Schellbach, Glass with PTFE Key Stopcock

- Automatic zero burettes & blue Schellbach markings.
- PE reservoir has high chemical resistance and is extremely break resistant.
- For precise dosing press button.
- Fine titration with micro screw and exact spout tip.
- Made of borosilicate glass.

Class-A

| Code    | Capacity<br>(ml) | Sub. Div.<br>(ml) | Tolerance<br>± ml | Volume<br>Bottle(ml) | Pack<br>Size |
|---------|------------------|-------------------|-------------------|----------------------|--------------|
| CH0242J | 25               | 0.10              | 0.05              | 1000                 | 1            |
| CH0242K | 50               | 0.10              | 0.05              | 1000                 | 1            |



Check out  
our new website  
[www.eiscoindustrial.com](http://www.eiscoindustrial.com)

## Graduated Cylinders - Glass, Class-A with Round Base

- Class-A graduated cylinder with round base and spout.
- Made of borosilicate glass.
- As per DIN EN ISO 4788 standard.

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0344H        | CH0344H-WT      | 5        | 0.1       | 0.05      | 10   |
| CH0344I        | CH0344I-WT      | 10       | 0.2       | 0.10      | 2    |
| CH0344J        | CH0344J-WT      | 25       | 0.5       | 0.25      | 2    |
| CH0344K        | CH0344K-WT      | 50       | 1.0       | 0.50      | 2    |
| CH0344L        | CH0344L-WT      | 100      | 1.0       | 0.50      | 2    |
| CH0344M        | CH0344M-WT      | 250      | 2.0       | 1.00      | 2    |
| CH0344N        | CH0344N-WT      | 500      | 5.0       | 2.50      | 2    |
| CH0344O        | CH0344O-WT      | 1000     | 10.0      | 5.00      | 2    |
| CH0344P        | CH0344P-WT      | 2000     | 20.0      | 10.00     | 1    |



BATCH CERTIFICATES



## Graduated Cylinders - Glass, Class-B with Round Base

- Class-B graduated cylinder with round base and spout.
- Made of borosilicate glass.
- As per DIN EN ISO 4788 standard.

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0344A5       | CH0344A5-WT     | 5        | 0.1       | 0.10      | 10   |
| CH0344A10      | CH0344A10-WT    | 10       | 0.2       | 0.20      | 2    |
| CH0344A        | CH0344A-WT      | 25       | 0.5       | 0.50      | 2    |
| CH0344B        | CH0344B-WT      | 50       | 1.0       | 1.00      | 2    |
| CH0344C        | CH0344C-WT      | 100      | 1.0       | 1.00      | 2    |
| CH0344D        | CH0344D-WT      | 250      | 2.0       | 2.00      | 2    |
| CH0344E        | CH0344E-WT      | 500      | 5.0       | 5.00      | 2    |
| CH0344F        | CH0344F-WT      | 1000     | 10.0      | 10.00     | 2    |
| CH0344G        | CH0344G-WT      | 2000     | 20.0      | 20.00     | 1    |





### Graduated Cylinders - Glass, Class-A with Hexagonal Base

- Class-A graduated cylinder with hexagonal base and spout.
- Made of borosilicate glass.
- As per DIN EN ISO 4788 standard.

| Code      | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|-----------|---------------|----------------|--------------------|-----------|
| CH0345H05 | 5             | 0.1            | 0.05               | 10        |
| CH0345H   | 10            | 0.2            | 0.10               | 2         |
| CH0345I   | 25            | 0.5            | 0.25               | 2         |
| CH0345J   | 50            | 1.0            | 0.50               | 2         |
| CH0345K   | 100           | 1.0            | 0.50               | 2         |
| CH0345L   | 250           | 2.0            | 1.00               | 2         |
| CH0345M   | 500           | 5.0            | 2.50               | 2         |
| CH0345N   | 1000          | 10.0           | 5.00               | 2         |
| CH0345O   | 2000          | 20.0           | 10.00              | 1         |



### Graduated Cylinders - Glass, Class-B with Hexagonal Base

- Class-B graduated cylinder with hexagonal base and spout.
- Made of borosilicate glass.
- As per DIN EN ISO 4788 standard.

| Code      | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|-----------|---------------|----------------|--------------------|-----------|
| CH0345A05 | 5             | 0.1            | 0.10               | 10        |
| CH0345A10 | 10            | 0.2            | 0.20               | 2         |
| CH0345A   | 25            | 0.5            | 0.50               | 2         |
| CH0345B   | 50            | 1.0            | 1.00               | 2         |
| CH0345C   | 100           | 1.0            | 1.00               | 2         |
| CH0345D   | 250           | 2.0            | 2.00               | 2         |
| CH0345E   | 500           | 5.0            | 5.00               | 2         |
| CH0345F   | 1000          | 10.0           | 10.00              | 2         |
| CH0345G   | 2000          | 20.0           | 20.00              | 1         |

### Graduated Cylinders - Glass, Class-B with Hexagonal Base and Guard

- Class-B graduated cylinder with spout, plastic hexagonal base and protective guard.
- Made of borosilicate glass.
- As per DIN EN ISO 4788 standard.

| Code           |                 | Capacity | Sub. Div. | Tolerance | Pack |
|----------------|-----------------|----------|-----------|-----------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      | Size |
| CH0346A        | CH0346A-WT      | 10       | 0.2       | 0.20      | 2    |
| CH0346B        | CH0346B-WT      | 25       | 0.5       | 0.50      | 2    |
| CH0346C        | CH0346C-WT      | 100      | 1.0       | 1.00      | 2    |
| CH0346D        | CH0346D-WT      | 50       | 0.5       | 1.00      | 2    |
| CH0346E        | CH0346E-WT      | 250      | 2.0       | 2.00      | 2    |
| CH0346F        | CH0346F-WT      | 500      | 5.0       | 5.00      | 2    |
| CH0346G        | CH0346G-WT      | 1000     | 10.0      | 10.00     | 2    |



### Graduated Cylinders - Glass, Class-A with Round Base and Stoppers

- Class-A graduated cylinder with round base and spout.
- Made of borosilicate glass.
- Each graduated cylinder comes with interchangeable hexagonal glass or polypropylene stoppers.
- As per DIN EN ISO 4788 standard.

| Code           | Capacity (ml)   | Sub. Div. (ml) | Tolerance $\pm$ ml | Joint | Pack Size |
|----------------|-----------------|----------------|--------------------|-------|-----------|
| Blue Graduated | White Graduated |                |                    |       |           |
| CH0348I        | CH0348I-WT      | 10             | 0.2                | 0.10  | 10/19 2   |
| CH0348J        | CH0348J-WT      | 25             | 0.5                | 0.25  | 14/23 2   |
| CH0348K        | CH0348K-WT      | 50             | 1.0                | 0.50  | 14/23 2   |
| CH0348L        | CH0348L-WT      | 100            | 1.0                | 0.50  | 19/26 2   |
| CH0348M        | CH0348M-WT      | 250            | 2.0                | 1.00  | 24/29 2   |
| CH0348N        | CH0348N-WT      | 500            | 5.0                | 2.50  | 24/29 2   |
| CH0348O        | CH0348O-WT      | 1000           | 10.0               | 5.00  | 29/32 2   |
| CH0348P        | CH0348P-WT      | 2000           | 20.0               | 10.00 | 34/35 1   |



### Graduated Cylinders - Glass, Class-B with Round Base and Stoppers

- Class-B graduated cylinder with round base and spout.
- Made of borosilicate glass.
- Each graduated cylinder comes with interchangeable hexagonal glass or polypropylene stoppers.
- As per DIN EN ISO 4788 standard.

| Code           | Capacity (ml)   | Sub. Div. (ml) | Tolerance $\pm$ ml | Joint | Pack Size |
|----------------|-----------------|----------------|--------------------|-------|-----------|
| Blue Graduated | White Graduated |                |                    |       |           |
| CH0348A        | CH0348A-WT      | 10             | 0.2                | 0.20  | 10/19 2   |
| CH0348B        | CH0348B-WT      | 25             | 0.5                | 0.50  | 14/23 2   |
| CH0348C        | CH0348C-WT      | 50             | 1.0                | 1.00  | 14/23 2   |
| CH0348D        | CH0348D-WT      | 100            | 1.0                | 1.00  | 19/26 2   |
| CH0348E        | CH0348E-WT      | 250            | 2.0                | 2.00  | 24/29 2   |
| CH0348F        | CH0348F-WT      | 500            | 5.0                | 5.00  | 24/29 2   |
| CH0348G        | CH0348G-WT      | 1000           | 10.0               | 10.00 | 29/32 2   |
| CH0348H        | CH0348H-WT      | 2000           | 20.0               | 20.00 | 34/35 1   |



BATCH  
CERTIFICATES



GLASS  
BORO  
3.3



### Graduated Cylinders - Glass, Class-A with Hexagonal Base and Stoppers

- Class- A graduated cylinder with hexagonal base and spout.
- Made of borosilicate glass.
- Each graduated cylinder comes with interchangeable hexagonal glass or polypropylene stoppers.
- As per DIN EN ISO 4788 standard.

| Code           |                 | Capacity | Sub. Div. | Tolerance | Joint | Pack |
|----------------|-----------------|----------|-----------|-----------|-------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      |       | Size |
| CH0349I        | CH0349I-WT      | 10       | 0.2       | 0.10      | 10/19 | 2    |
| CH0349J        | CH0349J-WT      | 25       | 0.5       | 0.25      | 14/23 | 2    |
| CH0349K        | CH0349K-WT      | 50       | 1.0       | 0.50      | 14/23 | 2    |
| CH0349L        | CH0349L-WT      | 100      | 1.0       | 0.50      | 19/26 | 2    |
| CH0349M        | CH0349M-WT      | 250      | 2.0       | 1.00      | 24/29 | 2    |
| CH0349N        | CH0349N-WT      | 500      | 5.0       | 2.50      | 24/29 | 2    |
| CH0349O        | CH0349O-WT      | 1000     | 10.0      | 5.00      | 29/32 | 2    |
| CH0349P        | CH0349P-WT      | 2000     | 20.0      | 10.00     | 34/35 | 1    |



GLASS  
BORO  
3.3



### Graduated Cylinders - Glass, Class-B with Hexagonal Base and Stopper

- Class-B graduated cylinder with hexagonal base and spout.
- Made of borosilicate glass.
- Each graduated cylinder comes with interchangeable hexagonal glass or polypropylene stoppers.
- As per DIN EN ISO 4788 standard.

| Code           |                 | Capacity | Sub. Div. | Tolerance | Joint | Pack |
|----------------|-----------------|----------|-----------|-----------|-------|------|
| Blue Graduated | White Graduated | (ml)     | (ml)      | ± ml      |       | Size |
| CH0349A        | CH0349A-WT      | 10       | 0.2       | 0.20      | 10/19 | 2    |
| CH0349B        | CH0349B-WT      | 25       | 0.5       | 0.50      | 14/23 | 2    |
| CH0349C        | CH0349C-WT      | 50       | 1.0       | 1.00      | 14/23 | 2    |
| CH0349D        | CH0349D-WT      | 100      | 1.0       | 1.00      | 19/26 | 2    |
| CH0349E        | CH0349E-WT      | 250      | 2.0       | 2.00      | 24/29 | 2    |
| CH0349F        | CH0349F-WT      | 500      | 5.0       | 5.00      | 24/29 | 2    |
| CH0349G        | CH0349G-WT      | 1000     | 10.0      | 10.0      | 29/32 | 2    |
| CH0349H        | CH0349H-WT      | 2000     | 20.0      | 20.0      | 34/35 | 1    |

### Graduated Cylinders - Glass, Class-A ASTM with Hexagonal Base and Guard

- Class-A graduated cylinder with hexagonal base.
- Protective guard and spout.
- Made of borosilicate glass.
- As per ASTM E1272 standard.

| Code         | Blue Graduated | White Graduated | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|--------------|----------------|-----------------|---------------|----------------|--------------------|-----------|
| CH0341A05-BL |                | CH0341A05       | 5             | 0.1            | 0.05               | 2         |
| CH0341A-BL   |                | CH0341A         | 10            | 0.2            | 0.10               | 2         |
| CH0341B-BL   |                | CH0341B         | 25            | 0.5            | 0.17               | 2         |
| CH0341C-BL   |                | CH0341C         | 50            | 1.0            | 0.25               | 2         |
| CH0341D-BL   |                | CH0341D         | 100           | 1.0            | 0.50               | 2         |
| CH0341E-BL   |                | CH0341E         | 250           | 2.0            | 1.00               | 2         |
| CH0341F-BL   |                | CH0341F         | 500           | 5.0            | 2.00               | 2         |
| CH0341G-BL   |                | CH0341G         | 1000          | 10.0           | 3.0                | 2         |
| CH0341H-BL   |                | CH0341H         | 2000          | 20.0           | 6.0                | 1         |



BATCH CERTIFICATES

ASTM

GLASS BORO 3.3

Chemical Resistant

### Graduated Cylinders - Glass, Class-B ASTM with Hexagonal Base and Guard

- Class-B graduated cylinder with hexagonal base.
- Protective guard and spout
- Made of borosilicate glass.
- As per ASTM E1272 standard.

| Code         | Blue Graduated | White Graduated | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|--------------|----------------|-----------------|---------------|----------------|--------------------|-----------|
| CH0341I05-BL |                | CH0341I05       | 5             | 0.1            | 0.10               | 2         |
| CH0341I-BL   |                | CH0341I         | 10            | 0.2            | 0.20               | 2         |
| CH0341J-BL   |                | CH0341J         | 25            | 0.5            | 0.34               | 2         |
| CH0341K-BL   |                | CH0341K         | 50            | 1.0            | 0.50               | 2         |
| CH0341L-BL   |                | CH0341L         | 100           | 1.0            | 1.00               | 2         |
| CH0341M-BL   |                | CH0341M         | 250           | 2.0            | 2.00               | 2         |
| CH0341N-BL   |                | CH0341N         | 500           | 5.0            | 4.00               | 2         |
| CH0341O-BL   |                | CH0341O         | 1000          | 10.0           | 6.00               | 2         |
| CH0341P-BL   |                | CH0341P         | 2000          | 20.0           | 12.00              | 1         |



ASTM

GLASS BORO 3.3

Chemical Resistant

GLASS  
BORO  
3.3GLASS  
BORO  
3.3GLASS  
BORO  
3.3

### Nessler Cylinders

- Class- A Nessler cylinder with blue ceramic graduations.
- Used for color comparisons.
- Made of borosilicate glass.

| Code    | Capacity<br>(ml) | Tolerance<br>± ml | Pack<br>Size |
|---------|------------------|-------------------|--------------|
| CH0350A | 50               | 0.40              | 10           |
| CH0350B | 100              | 0.80              | 10           |

### Crow Receivers

- Class-B graduated Crow Receiver with base and spout.
- Made of borosilicate glass.

| Code    | Capacity<br>(ml) | Sub<br>Division | Pack<br>Size |
|---------|------------------|-----------------|--------------|
| CH0351A | 25               | 0.5             | 1            |
| CH0351B | 50               | 1.0             | 1            |
| CH0351C | 100              | 1.0             | 1            |

### Conical Measures

- Ungraduated conical measures with base and spout.
- Made of borosilicate glass.

| Code    | Capacity<br>(ml) | Pack<br>Size |
|---------|------------------|--------------|
| CH0355F | 125              | 1            |
| CH0355G | 250              | 1            |

### Pipettes - Glass, Class-A

- Bulb form pipettes, color coded with one mark accuracy & class-A.
- Made of borosilicate glass.
- As per DIN ISO 648 standard.

| Code    | Blue Graduated | White Graduated | Capacity (ml) | Tolerance $\pm$ ml | Colour Code | Pack Size |
|---------|----------------|-----------------|---------------|--------------------|-------------|-----------|
| CH0559A |                | CH0559A-WT      | 1             | 0.007              | Blue        | 1         |
| CH0559B |                | CH0559B-WT      | 2             | 0.010              | Orange      | 1         |
| CH0559C |                | CH0559C-WT      | 5             | 0.015              | White       | 1         |
| CH0559D |                | CH0559D-WT      | 10            | 0.020              | Red         | 1         |
| CH0559E |                | CH0559E-WT      | 15            | 0.025              | Green       | 1         |
| CH0559F |                | CH0559F-WT      | 20            | 0.030              | Yellow      | 1         |
| CH0559G |                | CH0559G-WT      | 25            | 0.030              | Blue        | 1         |
| CH0559H |                | CH0559H-WT      | 50            | 0.050              | Red         | 1         |
| CH0559I |                | CH0559I-WT      | 100           | 0.080              | Yellow      | 1         |



### Pipettes, Glass, Class-B

- Bulb form pipettes, color coded with one mark accuracy & class-B.
- Made of borosilicate glass.
- As per DIN ISO 648 standard.

| Code    | Blue Graduated | White Graduated | Capacity (ml) | Tolerance $\pm$ ml | Color Code | Pack Size |
|---------|----------------|-----------------|---------------|--------------------|------------|-----------|
| CH0558A |                | CH0558A-WT      | 1             | 0.015              | Blue       | 1         |
| CH0558B |                | CH0558B-WT      | 2             | 0.020              | Orange     | 1         |
| CH0558C |                | CH0558C-WT      | 5             | 0.030              | White      | 1         |
| CH0558D |                | CH0558D-WT      | 10            | 0.040              | Red        | 1         |
| CH0558E |                | CH0558E-WT      | 15            | 0.050              | Green      | 1         |
| CH0558F |                | CH0558F-WT      | 20            | 0.060              | Yellow     | 1         |
| CH0558G |                | CH0558G-WT      | 25            | 0.060              | Blue       | 1         |
| CH0558H |                | CH0558H-WT      | 50            | 0.100              | Red        | 1         |
| CH0558I |                | CH0558I-WT      | 100           | 0.160              | Yellow     | 1         |





BATCH  
CERTIFICATES



GLASS  
BORO  
3.3



Chemical  
Resistant



GLASS  
BORO  
3.3



Chemical  
Resistant

## Serological Pipettes - Glass, Class-A

- Color coded serological pipettes calibrated for delivery & class-A.
- Made of borosilicate glass.
- As per DIN ISO 835 standard.

| Code<br>Blue Graduated | White Graduated | Capacity<br>(ml) | Tolerance<br>± ml | Colour<br>Code | Pack<br>Size |
|------------------------|-----------------|------------------|-------------------|----------------|--------------|
| CH0573A                | CH0573A-WT      | 1.0              | 0.006             | Yellow         | 10           |
| CH0573B                | CH0573B-WT      | 2.0              | 0.010             | Black          | 10           |
| CH0573C                | CH0573C-WT      | 5.0              | 0.030             | Red            | 10           |
| CH0573D                | CH0573D-WT      | 10.0             | 0.050             | Orange         | 10           |
| CH0573E                | CH0573E-WT      | 25.0             | 0.100             | White          | 5            |

Ask us about OEM and  
Private Label Opportunities

## Serological Pipettes - Glass, Class-B

- Color coded serological pipettes calibrated for delivery & class-B.
- Made of borosilicate glass.
- As per DIN ISO 835 standard.

| Code<br>Blue Graduated | White Graduated | Capacity<br>(ml) | Tolerance<br>± ml | Colour<br>Code | Pack<br>Size |
|------------------------|-----------------|------------------|-------------------|----------------|--------------|
| CH0572D                | CH0572D-WT      | 1.0              | 0.010             | Yellow         | 10           |
| CH0572E                | CH0572E-WT      | 2.0              | 0.020             | Black          | 10           |
| CH0572F                | CH0572F-WT      | 5.0              | 0.050             | Red            | 10           |
| CH0572G                | CH0572G-WT      | 10.0             | 0.100             | Orange         | 10           |
| CH0572H                | CH0572H-WT      | 25.0             | 0.200             | White          | 5            |

### Crucibles - Gooch Type, Sintered Disc

- Gooch type crucible made of borosilicate glass.
- Sintered disc with G1, G2, G3 or G4 porosity.

| Code    | Capacity (ml) | Pack Size |
|---------|---------------|-----------|
| CH0465F | 15            | 10        |
| CH0465G | 30            | 10        |
| CH0465H | 50            | 10        |



GLASS  
BORO  
3.3

Chemical  
Resistant



Come chat with us  
Live on [www.eiscoindustrial.com](http://www.eiscoindustrial.com)

### Buchner Funnels - Glass with Sintered Disc

- Buchner funnel with plain stem, made of borosilicate glass.
- Sintered disc with G0, G1, G2, G3 or G4 porosity.

| Code    | Capacity (ml) | Disc Dia (mm) | Pack Size |
|---------|---------------|---------------|-----------|
| CH0466A | 35            | 30            | 2         |
| CH0466B | 80            | 40            | 2         |
| CH0466C | 200           | 65            | 2         |
| CH0466D | 500           | 90            | 2         |
| CH0466E | 1000          | 120           | 2         |



GLASS  
BORO  
3.3

Chemical  
Resistant

GLASS  
BORO  
3.3

Chemical  
Resistant



### Buchner Funnels - Glass with Sintered Disk

- Buchner funnel with cone at bottom of stem.
- Made of borosilicate glass.
- Sintered disc with G0, G1, G2, G3 or G4 porosity

| Code      | Capacity (ml) | Disc Dia. (mm) | Cone Size | Pack Size |
|-----------|---------------|----------------|-----------|-----------|
| CH0466F   | 35            | 30             | 14/23     | 2         |
| CH0466F01 | 35            | 30             | 19/26     | 2         |
| CH0466G   | 80            | 40             | 14/23     | 2         |
| CH0466G01 | 80            | 40             | 19/26     | 2         |
| CH0466H   | 200           | 65             | 19/26     | 2         |
| CH0466H01 | 200           | 65             | 24/29     | 2         |
| CH0466I   | 500           | 90             | 19/26     | 1         |
| CH0466I01 | 500           | 90             | 24/29     | 1         |
| CH0466I02 | 500           | 90             | 29/32     | 1         |
| CH0466J   | 1000          | 120            | 19/26     | 1         |
| CH0466J01 | 1000          | 120            | 24/29     | 1         |
| CH0466J02 | 1000          | 120            | 29/32     | 1         |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Filtration Assembly - 47 mm

- This filtration assembly can be used for vacuum filtration of tissue culture media, biological fluids, HPLC and other aqueous solutions.
- Set includes 300 ml borosilicate glass funnel.
- Highly integrated sealed 47 mm fitted glass support base, one litre flask
- An aluminium clamp to attach the filter funnel to the support base properly.
- The unit is designed to maximize flow and reduce foaming.
- The 47 mm filter provides an effective filtration area of 9.6 cm<sup>2</sup>.

| Code   | Item Description            |
|--------|-----------------------------|
| FSAS16 | Filtration Assembly - 47 mm |

- Spares for above

| Includes  | Item Description                        |
|-----------|-----------------------------------------|
| FSAS16-AA | 1 Funnel - 300 ml                       |
| FSAS16-BB | 2 Vacuum Base with Sintered Disc- 47 mm |
| FSAS16-CC | 3 Clamp - 47 mm                         |



1



2




3

### Dropping Funnels - Open with PTFE Key Stopcock

- Ungraduated cylindrical dropping funnel with open top and PTFE key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Pack Size |
|---------|---------------|-----------|
| CH0470A | 50            | 1         |
| CH0470B | 100           | 1         |
| CH0470C | 250           | 1         |
| CH0470D | 500           | 1         |
| CH0470E | 1000          | 1         |



### Dropping Funnel - Open with Glass Key Stopcock

- Ungraduated cylindrical dropping funnel with open top and glass key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Pack Size |
|---------|---------------|-----------|
| CH0470F | 50            | 1         |
| CH0470G | 100           | 1         |
| CH0470H | 250           | 1         |
| CH0470I | 500           | 1         |
| CH0470J | 1000          | 1         |



### Dropping Funnels - Glass Key Stopcock, Plastic Stopper

- Cylindrical dropping funnel with interchangeable plastic stopper and glass key stopcock.
- Made of borosilicate glass.
- Graduated and ungraduated versions available.

| Code        | Capacity (ml) | Socket Size | Pack Size |
|-------------|---------------|-------------|-----------|
| Ungraduated | Graduated     |             |           |
| CH0472A     | CH0472A-GR    | 50          | 19/26     |
| CH0472B     | CH0472B-GR    | 100         | 19/26     |
| CH0472C     | CH0472C-GR    | 250         | 29/32     |
| CH0472D     | CH0472D-GR    | 500         | 29/32     |
| CH0472E     | CH0472E-GR    | 1000        | 29/32     |



GLASS  
BORO  
3.3GLASS  
BORO  
3.3GLASS  
BORO  
3.3

### Dropping Funnel - Needle Valve Stopcock, Plastic Stopper

- Cylindrical dropping funnel with interchangeable plastic stopper and needle valve stopcock.
- Made of borosilicate glass.
- Graduated and ungraduated versions available.

| Code        | Capacity (ml) | Socket Size | Pack Size |
|-------------|---------------|-------------|-----------|
| UnGraduated | Graduated     |             |           |
| CH0473A     | CH0473A-GR    | 50          | 19/26     |
| CH0473B     | CH0473B-GR    | 100         | 19/26     |
| CH0473C     | CH0473C-GR    | 250         | 29/32     |
| CH0473D     | CH0473D-GR    | 500         | 29/32     |
| CH0473E     | CH0473E-GR    | 1000        | 29/32     |

### Dropping Funnel - PTFE Key Stopcock, Plastic Stopper

- Cylindrical dropping funnel with interchangeable plastic stopper and PTFE key stopcock.
- Made of borosilicate glass.
- Graduated and ungraduated versions available.

| Code        | Capacity (ml) | Socket Size | Pack Size |
|-------------|---------------|-------------|-----------|
| UnGraduated | Graduated     |             |           |
| CH0474A     | CH0474A-GR    | 50          | 19/26     |
| CH0474B     | CH0474B-GR    | 100         | 19/26     |
| CH0474C     | CH0474C-GR    | 250         | 29/32     |
| CH0474D     | CH0474D-GR    | 500         | 29/32     |
| CH0474E     | CH0474E-GR    | 1000        | 29/32     |

### Dropping Funnel - Glass Key Stopcock, Stem with Cone, Plastic Stopper

- Ungraduated cylindrical dropping funnel with cone on stem.
- Interchangeable plastic stopper and glass key stopcock.
- Made of borosilicate glass.

| Code      | Capacity (ml) | Socket/Cone Size | Pack Size |
|-----------|---------------|------------------|-----------|
| CH0469A   | 50            | 14/23            | 1         |
| CH0469A01 | 50            | 19/26            | 1         |
| CH0469A02 | 50            | 29/32            | 1         |
| CH0469B   | 100           | 14/23            | 1         |
| CH0469B01 | 100           | 19/26            | 1         |
| CH0469B02 | 100           | 29/32            | 1         |
| CH0469C   | 250           | 19/26            | 1         |
| CH0469D   | 250           | 29/32            | 1         |
| CH0469E   | 500           | 29/32            | 1         |
| CH0469F   | 1000          | 29/32            | 1         |

### Dropping Funnel - PTFE Key Stopcock, Stem with Cone, Plastic Stopper

- Ungraduated cylindrical dropping funnel with cone on stem.
- Interchangeable plastic stopper and PTFE key stopcock.
- Made of borosilicate glass.

| Code      | Capacity (ml) | Socket/Cone Size | Pack Size |
|-----------|---------------|------------------|-----------|
| CH0469G   | 50            | 14/23            | 1         |
| CH0469G01 | 50            | 19/26            | 1         |
| CH0469G02 | 50            | 29/32            | 1         |
| CH0469H   | 100           | 14/23            | 1         |
| CH0469H01 | 100           | 19/26            | 1         |
| CH0469H02 | 100           | 29/32            | 1         |
| CH0469I   | 250           | 19/26            | 1         |
| CH0469J   | 250           | 29/32            | 1         |
| CH0469K   | 500           | 29/32            | 1         |
| CH0469L   | 1000          | 29/32            | 1         |



### Dropping Funnel - Needle Valve Stopcock, Stem with Cone, Plastic Stopper

- Ungraduated cylindrical dropping funnel with cone on stem.
- Interchangeable plastic stopper and needle valve stopcock.
- Made of borosilicate glass.

| Code      | Capacity (ml) | Socket/Cone Size | Pack Size |
|-----------|---------------|------------------|-----------|
| CH0469M   | 50            | 14/23            | 1         |
| CH0469M01 | 50            | 19/26            | 1         |
| CH0469M02 | 50            | 29/32            | 1         |
| CH0469N   | 100           | 14/23            | 1         |
| CH0469N01 | 100           | 19/26            | 1         |
| CH0469N02 | 100           | 29/32            | 1         |
| CH0469O   | 250           | 19/26            | 1         |
| CH0469P   | 250           | 29/32            | 1         |
| CH0469Q   | 500           | 29/32            | 1         |
| CH0469R   | 1000          | 29/32            | 1         |



GLASS  
BORO  
3.3Chemical  
Resistant

### Dropping Funnel - Pressure Equalizing with PTFE Key Stopcock and Plastic Stopper

- Pressure equalizing ungraduated cylindrical dropping funnel with interchangeable plastic stopper and PTFE key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket/Cone Size | Pack Size |
|---------|---------------|------------------|-----------|
| CH0471A | 50            | 14/23            | 1         |
| CH0471B | 50            | 19/26            | 1         |
| CH0471C | 50            | 29/32            | 1         |
| CH0471D | 100           | 14/23            | 1         |
| CH0471E | 100           | 19/26            | 1         |
| CH0471F | 100           | 29/32            | 1         |
| CH0471G | 250           | 29/32            | 1         |
| CH0471H | 500           | 29/32            | 1         |

GLASS  
BORO  
3.3Chemical  
Resistant

### Dropping Funnel - Pressure Equalizing with Glass Key Stopcock and Plastic Stopper

- Pressure equalizing ungraduated cylindrical dropping funnel with interchangeable plastic stopper and glass key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket/Cone Size | Pack Size |
|---------|---------------|------------------|-----------|
| CH0471K | 50            | 14/23            | 1         |
| CH0471L | 50            | 19/26            | 1         |
| CH0471M | 50            | 29/32            | 1         |
| CH0471N | 100           | 14/23            | 1         |
| CH0471O | 100           | 19/26            | 1         |
| CH0471P | 100           | 29/32            | 1         |
| CH0471Q | 250           | 29/32            | 1         |
| CH0471R | 500           | 29/32            | 1         |

GLASS  
BORO  
3.3Chemical  
Resistant

### Squibb Separating Funnel - Glass Key Stopcock, Plastic Stopper

- Squibb separating funnel with interchangeable plastic stopper and glass key stopcock.
- Made of borosilicate glass.
- Graduated and ungraduated versions available.

| Code        |              | Capacity | Socket | Pack |
|-------------|--------------|----------|--------|------|
| UnGraduated | Graduated    | (ml)     | Size   | Size |
| CH0476A50   | CH0476A50-GR | 50       | 14/23  | 2    |
| CH0476A     | CH0476A-GR   | 100      | 19/26  | 2    |
| CH0476B     | CH0476B-GR   | 250      | 19/26  | 2    |
| CH0476C     | CH0476C-GR   | 500      | 24/29  | 2    |
| CH0476D     | CH0476D-GR   | 1000     | 29/32  | 2    |
| CH0476E     | CH0476E-GR   | 2000     | 29/32  | 1    |

### Squibb Separating Funnel - Rotaflow Stopcock, Plastic Stopper

- Squibb separating funnel with interchangeable plastic stopper and Rotaflow stopcock.
- Made of borosilicate glass.
- Graduated and ungraduated versions available.

| Code        | Capacity (ml) | Socket Size | Pack Size |
|-------------|---------------|-------------|-----------|
| UnGraduated | Graduated     |             |           |
| CH0478A     | CH0478A-GR    | 50          | 14/23     |
| CH0478B     | CH0478B-GR    | 100         | 19/26     |
| CH0478C     | CH0478C-GR    | 250         | 19/26     |
| CH0478D     | CH0478D-GR    | 500         | 24/29     |
| CH0478E     | CH0478E-GR    | 1000        | 29/32     |
| CH0478F     | CH0478F-GR    | 2000        | 29/32     |



GLASS BORO 3.3



### Squibb Separating Funnel - PTFE Key Stopcock, Plastic Stopper

- Squibb separating funnel with interchangeable plastic stopper and PTFE key stopcock.
- Made of borosilicate glass.
- Graduated and ungraduated versions available.

| Code        | Capacity (ml) | Socket Size | Pack Size |
|-------------|---------------|-------------|-----------|
| UnGraduated | Graduated     |             |           |
| CH0479A     | CH0479A-GR    | 50          | 14/23     |
| CH0479B     | CH0479B-GR    | 100         | 19/26     |
| CH0479C     | CH0479C-GR    | 250         | 19/26     |
| CH0479D     | CH0479D-GR    | 500         | 24/29     |
| CH0479E     | CH0479E-GR    | 1000        | 29/32     |
| CH0479F     | CH0479F-GR    | 2000        | 29/32     |



GLASS BORO 3.3



### Separating Funnel - Pear Shaped, Glass Key Stopcock, Plastic Stopper

- Ungraduated pear shaped separating funnel with interchangeable plastic stopper and glass key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0482A | 50            | 14/23       | 2         |
| CH0482B | 100           | 19/26       | 2         |
| CH0482C | 250           | 19/26       | 2         |
| CH0482D | 500           | 24/29       | 2         |
| CH0482E | 1000          | 29/32       | 2         |
| CH0482F | 2000          | 29/32       | 1         |



GLASS BORO 3.3



GLASS  
BORO  
3.3



## Separating Funnel - Pear Shaped, Rotaflow Stopcock, Plastic Stopper

- Ungraduated pear shaped separating funnel with interchangeable plastic stopper and Rotaflow stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0484A | 50            | 14/23       | 2         |
| CH0484B | 100           | 19/26       | 2         |
| CH0484C | 250           | 19/26       | 2         |
| CH0484D | 500           | 24/29       | 2         |
| CH0484E | 1000          | 29/32       | 2         |
| CH0484F | 2000          | 29/32       | 1         |

GLASS  
BORO  
3.3



## Separating Funnel - Pear Shaped, PTFE Key Stopcock, Plastic Stopper

- Ungraduated pear shaped separating funnel with interchangeable plastic stopper and PTFE key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0485A | 50            | 14/23       | 2         |
| CH0485B | 100           | 19/26       | 2         |
| CH0485C | 250           | 19/26       | 2         |
| CH0485D | 500           | 24/29       | 2         |
| CH0485E | 1000          | 29/32       | 2         |
| CH0485F | 2000          | 29/32       | 1         |

GLASS  
BORO  
3.3



## Separating Funnel - Pear Shaped, PTFE Key Stopcock, Plastic Stopper, Stem with Cone

- Ungraduated pear shaped separating funnel with interchangeable plastic stopper.
- Stem with cone and PTFE key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket/Cone Size | Pack Size |
|---------|---------------|------------------|-----------|
| CH0481A | 50            | 14/23            | 2         |
| CH0481B | 50            | 19/26            | 2         |
| CH0481C | 100           | 19/26            | 2         |
| CH0481D | 250           | 19/26            | 2         |
| CH0481E | 250           | 29/32            | 2         |
| CH0481F | 500           | 29/32            | 2         |
| CH0481G | 1000          | 29/32            | 2         |

### Separating Funnel - Pear Shaped, Glass Key Stopcock, Plastic Stopper, Stem with Cone

- Ungraduated pear shaped separating funnel with interchangeable plastic stopper.
- Stem with cone & glass key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0481H | 50            | 14/23       | 2         |
| CH0481I | 50            | 19/26       | 2         |
| CH0481J | 100           | 19/26       | 2         |
| CH0481K | 250           | 19/26       | 2         |
| CH0481L | 250           | 29/32       | 2         |
| CH0481M | 500           | 29/32       | 2         |
| CH0481N | 1000          | 29/32       | 2         |



### Pressure Equalizing Separating Funnel - Pear Shaped, PTFE Key Stopcock, Plastic Stopper

- Ungraduated pressure equalizing pear shaped separating funnel with interchangeable plastic stopper and PTFE key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket/Cone Size | Pack Size |
|---------|---------------|------------------|-----------|
| CH0483A | 50            | 14/23            | 1         |
| CH0483B | 100           | 14/23            | 1         |
| CH0483C | 100           | 19/26            | 1         |
| CH0483D | 250           | 19/26            | 1         |
| CH0483E | 250           | 29/32            | 1         |
| CH0483F | 500           | 29/32            | 1         |
| CH0483G | 1000          | 29/32            | 1         |



### Pressure Equalizing Separating Funnel - Pear Shaped, Glass Key Stopcock, Plastic Stopper

- Ungraduated pressure equalizing pear shaped separating funnel with interchangeable plastic stopper and glass key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket/Cone Size | Pack Size |
|---------|---------------|------------------|-----------|
| CH0483J | 50            | 14/23            | 1         |
| CH0483K | 100           | 14/23            | 1         |
| CH0483L | 100           | 19/26            | 1         |
| CH0483M | 250           | 19/26            | 1         |
| CH0483N | 250           | 29/32            | 1         |
| CH0483O | 500           | 29/32            | 1         |
| CH0483P | 1000          | 29/32            | 1         |





## Gilson Separating Funnel - Gas Key Stopcock, Plastic Stopper

- Ungraduated Gilson separating funnel with interchangeable plastic stopper and glass key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0487A | 50            | 14/23       | 1         |
| CH0487B | 100           | 19/26       | 1         |
| CH0487C | 250           | 29/32       | 1         |
| CH0487D | 500           | 29/32       | 1         |
| CH0487E | 1000          | 29/32       | 1         |

Ask us about **OEM** and **Private Label Opportunities**



## Gilson Separating Funnel - PTFE Key Stopcock, Plastic Stopper

- Ungraduated Gilson separating funnel with interchangeable plastic stopper and PTFE key stopcock.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Pack Size |
|---------|---------------|-------------|-----------|
| CH0487G | 50            | 14/23       | 1         |
| CH0487H | 100           | 19/26       | 1         |
| CH0487I | 250           | 29/32       | 1         |
| CH0487J | 500           | 29/32       | 1         |
| CH0487K | 1000          | 29/32       | 1         |

## Liebig Condensers

1

- Inner tube integral with body.
- Narrow jacket for efficient cooling.
- Made of borosilicate glass.

| Code    | Length of Jacket (cm) | Total Length (cm) | Pack Size |
|---------|-----------------------|-------------------|-----------|
| CH0301A | 30                    | 48                | 1         |
| CH0301B | 40                    | 59                | 1         |
| CH0301C | 50                    | 70                | 1         |
| CH0301D | 20                    | 35                | 1         |

## Liebig Condensers - Jointed

2

- Body and integral tube made from borosilicate glass.
- Plastic or glass connectors available in several socket & cone sizes.

| Glass Connectors | Plastic Connectors | Socket/Cone Size | Effective Length (cm) | Pack Size |
|------------------|--------------------|------------------|-----------------------|-----------|
| CH0304AN         | CH0304AN-PST       | 14/23            | 16                    | 1         |
| CH0304BN         | CH0304BN-PST       | 19/26            | 16                    | 1         |
| CH0304CN         | CH0304CN-PST       | 19/26            | 25                    | 1         |
| CH0304DN         | CH0304DN-PST       | 24/29            | 25                    | 1         |
| CH0304EN         | CH0304EN-PST       | 19/26            | 30                    | 1         |
| CH0304FN         | CH0304FN-PST       | 24/29            | 30                    | 1         |
| CH0304GN         | CH0304GN-PST       | 29/32            | 30                    | 1         |
| CH0304HN         | CH0304HN-PST       | 24/29            | 40                    | 1         |
| CH0304IN         | CH0304IN-PST       | 29/32            | 40                    | 1         |

## Allihn Bulb Condensers

3

- Body and integral tube made from borosilicate glass.
- Plastic or glass connectors available in several socket & cone sizes.

| Glass Connectors | Plastic Connectors | Socket/Cone Size | Effective Length (cm) | Pack Size |
|------------------|--------------------|------------------|-----------------------|-----------|
| CH0306AN         | CH0306AN-PST       | 14/23            | 16                    | 1         |
| CH0306BN         | CH0306BN-PST       | 19/26            | 16                    | 1         |
| CH0306CN         | CH0306CN-PST       | 14/23            | 25                    | 1         |
| CH0306DN         | CH0306DN-PST       | 19/26            | 25                    | 1         |
| CH0306EN         | CH0306EN-PST       | 24/29            | 25                    | 1         |
| CH0306FN         | CH0306FN-PST       | 24/29            | 30                    | 1         |
| CH0306GN         | CH0306GN-PST       | 29/32            | 30                    | 1         |
| CH0306HN         | CH0306HN-PST       | 24/29            | 40                    | 1         |
| CH0306IN         | CH0306IN-PST       | 29/32            | 40                    | 1         |
| CH0306IN01       | CH0306IN01-PST     | 19/26            | 30                    | 1         |





GLASS  
BORO  
3.3




## Dimroth Condenser <sup>1</sup>

- Body and integral tube made from borosilicate glass.
- Glass or plastic connectors available in several socket & cone sizes.

| Glass Connectors | Plastic Connectors | Socket/Cone Size | Effective Length (cm) | Pack Size |
|------------------|--------------------|------------------|-----------------------|-----------|
| CH0307A          | CH0307A-PST        | 14/23            | 16                    | 1         |
| CH0307B          | CH0307B-PST        | 19/26            | 25                    | 1         |
| CH0307C          | CH0307C-PST        | 19/26            | 30                    | 1         |
| CH0307D          | CH0307D-PST        | 24/29            | 30                    | 1         |
| CH0307E          | CH0307E-PST        | 29/32            | 30                    | 1         |



Come chat with us  
Live on [www.eiscoindustrial.com](http://www.eiscoindustrial.com)



GLASS  
BORO  
3.3



## Graham Condenser - Coiled Distillate <sup>2</sup>

- Body and integral tube made from borosilicate glass.
- Glass or plastic connectors available in several socket & cone sizes.

| Glass Connectors | Plastic Connectors | Socket/Cone Size | Effective Length (cm) | Pack Size |
|------------------|--------------------|------------------|-----------------------|-----------|
| CH0308AN         | CH0308AN-PST       | 14/23            | 16                    | 1         |
| CH0308BN         | CH0308BN-PST       | 19/26            | 16                    | 1         |
| CH0308DN         | CH0308DN-PST       | 19/26            | 25                    | 1         |
| CH0308EN         | CH0308EN-PST       | 24/29            | 25                    | 1         |
| CH0308GN         | CH0308GN-PST       | 24/29            | 30                    | 1         |
| CH0308HN         | CH0308HN-PST       | 29/32            | 30                    | 1         |
| CH0308IN         | CH0308IN-PST       | 29/32            | 40                    | 1         |
| CH0308JN         | CH0308JN-PST       | 24/29            | 40                    | 1         |

## Expansion Adapters

- Made of borosilicate glass as per DIN 12257 standard.

| Code    | Socket Size | Cone Size | Pack Size |
|---------|-------------|-----------|-----------|
| CH0820A | 19/26       | 14/23     | 10        |
| CH0820B | 24/29       | 14/23     | 10        |
| CH0820C | 24/29       | 19/26     | 10        |
| CH0820D | 29/32       | 19/26     | 10        |
| CH0820E | 29/32       | 24/29     | 10        |
| CH0820F | 34/35       | 19/26     | 10        |
| CH0820G | 34/35       | 24/29     | 10        |
| CH0820H | 34/35       | 29/32     | 10        |
| CH0820I | 40/38       | 24/29     | 10        |



**DIN**  
12257

**GLASS**  
**BORO**  
3.3



## Reduction Adapters

- Made of borosilicate glass as per DIN 12257 standard.

| Code    | Socket Size | Cone Size | Pack Size |
|---------|-------------|-----------|-----------|
| CH0821A | 14/23       | 19/26     | 10        |
| CH0821B | 14/23       | 24/29     | 10        |
| CH0821C | 19/26       | 24/29     | 10        |
| CH0821D | 19/26       | 29/32     | 10        |
| CH0821E | 19/26       | 40/38     | 10        |
| CH0821F | 24/29       | 29/32     | 10        |
| CH0821G | 24/29       | 34/35     | 10        |
| CH0821H | 24/29       | 40/38     | 10        |
| CH0821I | 29/32       | 45/40     | 10        |



**DIN**  
12257

**GLASS**  
**BORO**  
3.3



## Single Sockets

- Made of borosilicate glass as per DIN 12249 standard.

| Code    | Socket Size | Total Length (mm) | Pack Size |
|---------|-------------|-------------------|-----------|
| CH0823A | 10/19       | 120               | 10        |
| CH0823B | 12/21       | 120               | 10        |
| CH0823C | 14/23       | 120               | 10        |
| CH0823D | 19/26       | 120               | 10        |
| CH0823E | 24/29       | 135               | 10        |
| CH0823F | 29/32       | 135               | 10        |
| CH0823G | 34/35       | 150               | 10        |



**DIN**  
12249

**GLASS**  
**BORO**  
3.3



**DIN**  
12249



### Single Cones

- Made of borosilicate glass as per DIN 12249 standard.

| Code    | Cone Size | Total Length (mm) | Pack Size |
|---------|-----------|-------------------|-----------|
| CH0825A | 10/19     | 120               | 10        |
| CH0825B | 12/21     | 120               | 10        |
| CH0825C | 14/23     | 120               | 10        |
| CH0825D | 19/26     | 120               | 10        |
| CH0825E | 24/29     | 150               | 10        |
| CH0825F | 29/32     | 150               | 10        |
| CH0825G | 34/35     | 150               | 10        |

**DIN**  
12594

**GLASS BORO 3.3**



### Multiple Adapters - Two Neck, Parallel

- Designed with two parallel necks.
- Made of borosilicate glass as per DIN 12594 standard.

| Code    | Socket Size | Cone Size | Pack Size |
|---------|-------------|-----------|-----------|
| CH0822A | 14/23       | 14/23     | 1         |
| CH0822B | 19/26       | 19/26     | 1         |
| CH0822C | 19/26       | 24/29     | 1         |
| CH0822D | 24/29       | 24/29     | 1         |
| CH0822E | 19/26       | 34/35     | 1         |
| CH0822F | 29/32       | 29/32     | 1         |

**GLASS BORO 3.3**



### Multiple Adapters - Two Neck, Angled

- Designed with one vertical neck and one at 45°.
- Made of borosilicate glass.

| Code    | Socket Size | Cone Size | Pack Size |
|---------|-------------|-----------|-----------|
| CH0824A | 19/26       | 19/26     | 1         |
| CH0824B | 19/26       | 24/29     | 1         |
| CH0824C | 19/26       | 34/35     | 1         |



**GLASS BORO 3.3**



### Multiple Adapters - Three Neck

- Designed with two vertical necks and one at 45°.
- Made of borosilicate glass.

| Code    | Socket Size | Cone Size | Pack Size |
|---------|-------------|-----------|-----------|
| CH0826A | 19/26       | 24/29     | 1         |
| CH0826B | 19/26       | 29/32     | 1         |
| CH0826C | 19/26       | 34/35     | 1         |
| CH0826D | 19/26       | 45/40     | 1         |

### Receiver Delivery Adapters

- Short stem and long stem available.
- Made of borosilicate glass.

| Code<br>Short Stem | Socket<br>Size | Length<br>(mm) | Pack<br>Size |
|--------------------|----------------|----------------|--------------|
| CH0830A            | 14/23          | 65             | 10           |
| CH0830B            | 19/26          | 65             | 10           |
| CH0830C            | 24/29          | 65             | 10           |
| CH0830D            | 29/32          | 65             | 5            |

| Code<br>Long Stem | Socket<br>Size | Length<br>(mm) | Pack<br>Size |
|-------------------|----------------|----------------|--------------|
| CH0830E           | 10/19          | 120            | 10           |
| CH0830F           | 14/23          | 190            | 10           |
| CH0830G           | 19/26          | 200            | 10           |
| CH0830H           | 24/29          | 200            | 5            |
| CH0830I           | 29/32          | 200            | 5            |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Bend Receiver Adapters with Vent

- Made of borosilicate glass.

| Code    | Socket<br>Size | Cone<br>Size | Pack<br>Size |
|---------|----------------|--------------|--------------|
| CH0834A | 14/23          | 14/23        | 1            |
| CH0834B | 19/26          | 19/26        | 1            |
| CH0834C | 19/26          | 24/29        | 1            |
| CH0834D | 24/29          | 24/29        | 1            |
| CH0834E | 29/32          | 29/32        | 1            |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Cone/Screw Thread Adapters - Straight

- Thread size 13 & diameter through screw cap 6-7 mm.
- Cone made of borosilicate glass.

| Code    | Cone<br>Size | Pack<br>Size |
|---------|--------------|--------------|
| CH0835A | 14/23        | 10           |
| CH0835B | 19/26        | 10           |
| CH0835C | 24/29        | 10           |
| CH0835D | 29/32        | 10           |



GLASS  
BORO  
3.3

Chemical  
Resistant



GLASS  
BORO  
3.3

Chemical  
Resistant

### Receiver Adapters - Bend with Vacuum

- Connects the outlet of a condenser to a receiving flask.
- Made of borosilicate glass.

| Code    | Socket Size | Cone Size | Pack Size |
|---------|-------------|-----------|-----------|
| CH0836A | 14/23       | 14/23     | 1         |
| CH0836B | 19/26       | 19/26     | 1         |
| CH0836C | 19/26       | 24/29     | 1         |
| CH0836D | 24/29       | 24/29     | 1         |
| CH0836E | 24/29       | 29/32     | 1         |
| CH0836F | 29/32       | 29/32     | 1         |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Still Heads - Plain with Thermometer Socket

- Plain still head with B-14 thermometer socket.
- Made of borosilicate glass.

| Code    | Cone Size to fit Flask | Cone Size to fit Condenser | Pack Size |
|---------|------------------------|----------------------------|-----------|
| CH0848A | 14/23                  | 14/23                      | 1         |
| CH0848B | 19/26                  | 19/26                      | 1         |
| CH0848C | 24/29                  | 19/26                      | 1         |
| CH0848D | 29/32                  | 19/26                      | 1         |
| CH0848E | 34/35                  | 19/26                      | 1         |
| CH0848F | 24/29                  | 24/29                      | 1         |
| CH0848G | 34/35                  | 24/29                      | 1         |
| CH0848H | 34/35                  | 34/35                      | 1         |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Splash Heads Adapter - Pear Shape, Straight

- Used between a flask and condenser to prevent carry over of raw material.
- Made of borosilicate glass.

| Code    | Socket Size | Cone Size | Pack Size |
|---------|-------------|-----------|-----------|
| CH0852A | 14/23       | 14/23     | 1         |
| CH0852B | 19/26       | 24/29     | 1         |
| CH0852C | 24/29       | 29/32     | 1         |

### Splash Head Adapters - Pear Shape, Verticle

- Used between a flask and condenser to prevent carry over of raw material.
- Made of borosilicate glass.

| Code    | Cone Size<br>to fit Flask | Cone Size<br>to fit Condenser | Pack<br>Size |
|---------|---------------------------|-------------------------------|--------------|
| CH0854A | 19/26                     | 19/26                         | 1            |
| CH0854B | 24/29                     | 19/26                         | 1            |
| CH0854C | 24/29                     | 24/29                         | 1            |
| CH0854D | 29/32                     | 29/32                         | 1            |



### Splash Head Adapters - Pear Shape, Sloping

- Used between a flask and condenser to prevent carry over of raw material.
- Made of borosilicate glass.

| Code    | Cone Size<br>to fit Flask | Cone Size<br>to fit Condenser | Pack<br>Size |
|---------|---------------------------|-------------------------------|--------------|
| CH0856A | 19/26                     | 19/26                         | 1            |
| CH0856B | 24/29                     | 19/26                         | 1            |
| CH0856C | 24/29                     | 24/29                         | 1            |
| CH0856D | 29/32                     | 29/32                         | 1            |



### Stoppers - Hollow, Hexagonal with Flat End

- Hollow hexagonal stopper with flat end.
- Made of borosilicate glass.

| Code    | Cone<br>Size | Pack<br>Size |
|---------|--------------|--------------|
| CH0868A | 7/16         | 10           |
| CH0868B | 10/19        | 10           |
| CH0868C | 12/21        | 10           |
| CH0868D | 14/23        | 10           |
| CH0868E | 19/26        | 10           |
| CH0868F | 24/29        | 10           |
| CH0868G | 29/32        | 10           |



GLASS  
BORO  
3.3



### Stoppers - Penny Head

- Solid penny head stopper with flat end.
- Made of borosilicate glass.

| Code    | Cone Size | Pack Size |
|---------|-----------|-----------|
| CH0870A | 10/19     | 10        |
| CH0870B | 14/23     | 10        |
| CH0870C | 19/26     | 10        |
| CH0870D | 24/29     | 10        |
| CH0870E | 29/32     | 10        |



GLASS  
BORO  
3.3



### Stoppers - Flat Head

- Solid flat head stopper with flat end.
- Made of borosilicate glass.

| Code    | Cone Size | Pack Size |
|---------|-----------|-----------|
| CH0872A | 10/19     | 10        |
| CH0872B | 14/23     | 10        |
| CH0872C | 19/26     | 10        |
| CH0872D | 24/29     | 10        |
| CH0872E | 29/32     | 10        |



GLASS  
BORO  
3.3



### Stoppers - Hexagonal

- Hollow hexagonal stopper with round end.
- Made of borosilicate glass.

| Code    | Cone Size | Pack Size |
|---------|-----------|-----------|
| CH0872F | 10/19     | 10        |
| CH0872G | 12/21     | 10        |
| CH0872H | 14/23     | 10        |
| CH0872I | 19/26     | 10        |
| CH0872J | 24/29     | 10        |
| CH0872K | 29/32     | 10        |

## Kipp's Gas Generator

- Rapid and effective means of producing controllable amounts of common gases.
- Comprised of three separate pieces.
- Made of borosilicate glass

| Code    | Capacity (ml) | Pack Size |
|---------|---------------|-----------|
| CH0516A | 250           | 1         |
| CH0516B | 500           | 1         |
| CH0516C | 1000          | 1         |



Custom Packaging for Glass Products to ensure safe transit



GLASS BORO 3.3

Chemical Resistant

## Gas Preparation Kit

This kit contains all the essential glassware needed for production of gases in a laboratory including Ammonia, Bromine, Carbon Dioxide, Carbon Monoxide, Chlorine, Hydrogen, Hydrogen Bromide, Hydrogen Iodide, Hydrogen Chloride, Hydrogen Sulphide, Hydrogen Peroxide, Nitric Oxide, Nitrogen, Nitrogen Dioxide, Nitrous Oxide, Oxygen, Sulphur Dioxide, Sulphur Trioxide etc. Complete kit is packaged in 2 Grattnels® trays with lids, individual glassware pieces are secured in die-cut foam.

Includes:

|                                    |                                          |
|------------------------------------|------------------------------------------|
| Recovery bend 24/29                | Bend tube                                |
| Adjustable wash bottle head 24/29  | 'T' adapter                              |
| Test tube 24/19                    | Combustion tube                          |
| Dropping funnel 50 cm <sup>3</sup> | Gas mixing tube                          |
| Test tube with side arm            | 2x Cone/screw cap adapter                |
| Delivery adapter 24/29             | 2x Gas cover jars 75 mm dia              |
| Dropping funnel plug 10/19         | Dreshel bottle 125 cm <sup>3</sup> 24/29 |
| Beehive shelf 75 mm dia. 24/29     | Flask r.B. 150 Cm <sup>3</sup> 24/29     |
| 5X Glass tuber 7 mm o.D.           | Plastic trough                           |
| 2X Gas jar 200 x 50 mm             | Instruction manual                       |

| Code   | Item Name           |
|--------|---------------------|
| CH0889 | Gas Preparation Kit |



GLASS  
BORO  
3.3



GLASS  
BORO  
3.3



GLASS  
BORO  
3.3



### Soxhlet Extraction Apparatus

- Includes Soxhlet extractor bulb, condenser.
- Flask with interchangeable standard joints
- Made of borosilicate glass.

| Code    | Capacity (ml) | Flask Size (ml) | Pack Size |
|---------|---------------|-----------------|-----------|
| CH0888A | 60            | 150             | 1         |
| CH0888B | 100           | 250             | 1         |
| CH0888C | 200           | 500             | 1         |
| CH0888D | 400           | 1000            | 1         |
| CH0888E | 600           | 2000            | 1         |
| CH0888F | 1000          | 3000            | 1         |
| CH0888G | 2000          | 5000            | 1         |

### Allihn Condenser for Soxhlet Apparatus

- Made of borosilicate glass.
- Various cone sizes available.

| Code     | Cone Size | Extractor (ml) | Pack Size |
|----------|-----------|----------------|-----------|
| CH0306JN | 45/40     | 100-250        | 1         |
| CH0306KN | 60/46     | 500            | 1         |
| CH0306LN | 71/51     | 1000           | 1         |

### Extractors for Soxhlet Apparatus

- For use with Soxhlet Apparatus.
- Made of borosilicate glass.

| Code    | Capacity (ml) | Socket Size | Cone Size | Pack Size |
|---------|---------------|-------------|-----------|-----------|
| CH0888H | 60            | 34/35       | 24/29     | 1         |
| CH0888I | 100           | 34/35       | 24/29     | 1         |
| CH0888J | 200           | 50/42       | 24/29     | 1         |
| CH0888K | 400           | 50/42       | 24/29     | 1         |
| CH0888L | 600           | 50/42       | 24/29     | 1         |
| CH0888M | 1000          | 55/44       | 34/35     | 1         |
| CH0888N | 2000          | 55/44       | 34/35     | 1         |

### Test Tubes - Light Wall with Rim

- Light walled test tubes with rim.
- Made of borosilicate glass with strong resistance to heat and chemicals.

| Code           | O.D. x Length (mm) | Pack Size |
|----------------|--------------------|-----------|
| CH0718AA       | 10 x 75            | 48        |
| CH0718A/PK48   | 12 x 75            | 48        |
| CH0718B/PK48   | 12 x 100           | 48        |
| CH0718C/PK48   | 16 x 100           | 48        |
| CH0718D01/PK48 | 15 x 130           | 48        |
| CH0718E/PK48   | 16 x 150           | 48        |
| CH0718F/PK24   | 18 x 150           | 24        |
| CH0718F01/PK24 | 20 x 150           | 24        |
| CH0718G/PK24   | 24 x 150           | 24        |
| CH0718H        | 25 x 250           | 24        |

### Test Tubes - Rimless

- Test tubes with fire polished ends to allow use with plugs and easier rack stacking.
- Made of borosilicate glass.

| Code    | O.D. x Length (mm) | Pack Size |
|---------|--------------------|-----------|
| CH0718K | 12 x 75            | 48        |
| CH0718L | 12 x 100           | 48        |
| CH0718M | 16 x 100           | 48        |
| CH0718N | 16 x 125           | 48        |
| CH0718O | 16 x 150           | 48        |
| CH0718P | 18 x 150           | 24        |
| CH0718Q | 24 x 150           | 24        |

### Test Tubes - Round Bottom with Screw Cap

- Test tubes with round bottom and bakelite screw cap with rubber liner.
- Tube made of borosilicate glass.

| Code           | O.D. x Length (mm) | Pack Size |
|----------------|--------------------|-----------|
| CH0724A01/PK24 | 12 x 100           | 24        |
| CH0724A/PK24   | 16 x 75            | 24        |
| CH0724B01/PK24 | 16 x 100           | 24        |
| CH0724B/PK24   | 16 x 125           | 24        |
| CH0724C        | 16 x 150           | 24        |
| CH0724D/PK24   | 18 x 150           | 24        |
| CH0724E/PK12   | 25 x 100           | 12        |
| CH0724F/PK12   | 25 x 150           | 12        |

### Test Tubes - Flat Bottom with Screw Cap

- Test tubes with flat bottom and bakelite screw cap with rubber liner.
- Tube made of borosilicate glass.

| Code    | O.D. x Length (mm) | Pack Size |
|---------|--------------------|-----------|
| CH0724K | 12 x 100           | 24        |
| CH0724L | 16 x 75            | 24        |
| CH0724M | 16 x 100           | 24        |
| CH0724N | 16 x 125           | 24        |





### Flasks - Boiling, Flat Bottom with Screw Thread Neck

- Round bottom boiling flask with screw thread neck.
- Made of borosilicate glass.

| Code     | Capacity (ml) | Socket Size | Pack Size |
|----------|---------------|-------------|-----------|
| CH01002A | 50            | 14/23       | 12        |
| CH01002B | 50            | 19/26       | 12        |
| CH01002C | 100           | 14/23       | 12        |
| CH01002D | 100           | 19/26       | 12        |
| CH01002E | 100           | 24/29       | 12        |
| CH01002F | 250           | 19/26       | 12        |
| CH01002G | 250           | 24/29       | 12        |
| CH01002H | 500           | 19/26       | 6         |
| CH01002I | 500           | 24/29       | 6         |
| CH01002J | 1000          | 24/29       | 6         |



### Flasks - Boiling, Flat Bottom with Screw Thread Neck

- Flat bottom boiling flask with screw thread neck.
- Made of borosilicate glass.

| Code     | Capacity (ml) | Socket Size | Pack Size |
|----------|---------------|-------------|-----------|
| CH01004A | 50            | 14/23       | 12        |
| CH01004B | 50            | 19/26       | 12        |
| CH01004C | 100           | 14/23       | 12        |
| CH01004D | 100           | 19/26       | 12        |
| CH01004E | 100           | 24/29       | 12        |
| CH01004F | 250           | 19/26       | 12        |
| CH01004G | 250           | 24/29       | 12        |
| CH01004H | 500           | 19/26       | 6         |
| CH01004I | 500           | 24/29       | 6         |
| CH01004J | 1000          | 24/29       | 6         |

## Flasks - Conical, Erlenmeyer with Screw Thread Neck

- Conical Erlenmeyer flask with screw thread neck.
- Made of borosilicate glass.

| Code     | Capacity (ml) | Socket Size | Pack Size |
|----------|---------------|-------------|-----------|
| CH01006A | 50            | 14/23       | 12        |
| CH01006B | 50            | 19/26       | 12        |
| CH01006C | 100           | 14/23       | 12        |
| CH01006D | 100           | 19/26       | 12        |
| CH01006E | 100           | 24/29       | 12        |
| CH01006F | 250           | 19/26       | 12        |
| CH01006G | 250           | 24/29       | 12        |
| CH01006H | 500           | 19/26       | 6         |
| CH01006I | 500           | 24/29       | 6         |
| CH01006J | 1000          | 24/29       | 6         |



GLASS  
BORO  
3.3



## Distillation Flasks - Round Bottom, Oblique Neck with Screw Thread

- Round bottom distillation flask with oblique screw thread necks.
- Made of borosilicate glass.

| Code     | Capacity (ml) | Socket Size | Pack Size |
|----------|---------------|-------------|-----------|
| CH01008A | 100           | 14/23       | 1         |
| CH01008B | 100           | 19/26       | 1         |
| CH01008C | 250           | 19/26       | 1         |
| CH01008D | 250           | 24/29       | 1         |
| CH01008E | 500           | 24/29       | 1         |
| CH01008F | 1000          | 24/29       | 1         |



GLASS  
BORO  
3.3



## Distillation Flasks - Round Bottom, Tow Oblique Necks with Screw Thread

- Round bottom distillation flask with two oblique necks with screw threads.
- Made of borosilicate glass.

| Code     | Capacity (ml) | Socket Size | Pack Size |
|----------|---------------|-------------|-----------|
| CH01010A | 100           | 14/23       | 1         |
| CH01010B | 100           | 19/26       | 1         |
| CH01010C | 250           | 19/26       | 1         |
| CH01010D | 250           | 24/29       | 1         |
| CH01010E | 500           | 19/26       | 1         |
| CH01010F | 500           | 24/29       | 1         |



GLASS  
BORO  
3.3



GLASS  
BORO  
3.3GLASS  
BORO  
3.3GLASS  
BORO  
3.3

### Liebig Condensers - Screw Thread

- Liebig condenser with screw thread joints.
- Made of borosilicate glass exclusive of bottom screw thread.
- Tubing connector for CH01012D is 12 mm, for all other sizes 8 mm.

| Code     | Socket Size | Cone Size | Effective Length (mm) | Pack Size |
|----------|-------------|-----------|-----------------------|-----------|
| CH01012A | 14/23       | 14/23     | 200                   | 1         |
| CH01012B | 19/26       | 19/26     | 200                   | 1         |
| CH01012C | 19/26       | 19/26     | 300                   | 1         |
| CH01012D | 24/29       | 24/29     | 300                   | 1         |

### Bulb Condensers - Screw Thread

- Bulb condenser with screw thread bottom joint.
- Made of borosilicate glass exclusive of screw thread.
- Tubing connector for CH01012D is 12 mm, for all other sizes 8 mm.

| Code     | Socket Size | Cone Size | Effective Length (mm) | Pack Size |
|----------|-------------|-----------|-----------------------|-----------|
| CH01014A | 14/23       | 14/23     | 200                   | 1         |
| CH01014B | 19/26       | 19/26     | 200                   | 1         |
| CH01014C | 19/26       | 19/26     | 300                   | 1         |
| CH01014D | 24/29       | 24/29     | 300                   | 1         |

### Vigrex Columns - Screw Thread

- Vigrex column with screw thread joints.
- Made of borosilicate glass exclusive of bottom screw thread.

| Code     | Socket Size | Cone Size | Effective Length (mm) | Pack Size |
|----------|-------------|-----------|-----------------------|-----------|
| CH01016A | 14/23       | 14/23     | 100                   | 1         |
| CH01016B | 19/26       | 19/26     | 150                   | 1         |
| CH01016C | 24/29       | 24/29     | 200                   | 1         |

### Receiver Adapters - Screw Thread

- Receiver adapter with screw thread joint.
- Made of borosilicate glass.

| Code     | Length (mm) | Socket Size | Pack Size |
|----------|-------------|-------------|-----------|
| CH01018A | 100         | 14/23       | 1         |
| CH01018B | 100         | 19/26       | 1         |
| CH01018C | 100         | 24/29       | 1         |



GLASS  
BORO  
3.3



### Receiver Adapters with Vacuum Connection - Screw Thread

- Receiver adapter with vacuum connection and screw thread joints.
- Made of borosilicate glass.
- Exclusive of bottom joint screw thread.

| Code     | Cone Size | Socket Size | Pack Size |
|----------|-----------|-------------|-----------|
| CH01020A | 14/23     | 14/23       | 1         |
| CH01020B | 19/26     | 19/26       | 1         |
| CH01020C | 24/29     | 24/29       | 1         |



GLASS  
BORO  
3.3



### Still Heads - Screw Thread

- Still head with screw thread joints and thermometer socket that accepts 6 - 7 mm dia. thermometers.
- Made of borosilicate glass exclusive of screw threads.

| Code     | Cone Size | Pack Size |
|----------|-----------|-----------|
| CH01022A | 14/23     | 1         |
| CH01022B | 19/26     | 1         |
| CH01022C | 24/29     | 1         |




GLASS  
BORO  
3.3



### Dropping Funnels - Screw Thread

- Ungraduated cylindrical dropping funnel with open top, PTFE key stopcock.
- Screw thread joints
- Made of borosilicate glass exclusive of stop cock and bottom screw thread.

| Code     | Capacity (ml) | Joint Size | Pack Size |
|----------|---------------|------------|-----------|
| CH01024A | 50            | 14/23      | 1         |
| CH01024B | 50            | 19/26      | 1         |
| CH01024C | 50            | 24/29      | 1         |
| CH01024D | 100           | 19/26      | 1         |
| CH01024E | 100           | 24/29      | 1         |



GLASS  
BORO  
3.3



GLASS  
BORO  
3.3

### Simple Distillation Apparatus - Secure Joint

4 piece simple distillation apparatus with size 24/29 screw thread joints, made of borosilicate glass exclusive of screw threads.

Includes:

- Round bottom flask - 250 ml
- Still head
- Liebig condenser
- Receiver adaptor

| Code   | Item Name                                    |
|--------|----------------------------------------------|
| CH1026 | Simple Distillation Apparatus - Secure Joint |

GLASS  
BORO  
3.3

### Simple Vacuum Distillation Apparatus - Secure Joint

5 piece simple vacuum distillation apparatus with size 24/29 screw thread joints, made of borosilicate glass exclusive of screw threads.

Includes:

- Round bottom flask - 250 ml
- Still head
- Liebig condenser
- Adapter receiver with vacuum
- Erlenmeyer flask - 100 ml

| Code   | Item Name                                           |
|--------|-----------------------------------------------------|
| CH1028 | Simple Vacuum Distillation Apparatus - Secure Joint |

GLASS  
BORO  
3.3

### Fractional Distillation Apparatus - Secure Joint

5 piece fractional distillation apparatus with size 24/29 screw thread, made of borosilicate glass exclusive of screw threads.

Includes:

- Round bottom flask - 250 ml
- Vigrex column
- Still head
- Liebig condenser
- Receiver adapter with vacuum

| Code   | Item Name                                        |
|--------|--------------------------------------------------|
| CH1030 | Fractional Distillation Apparatus - Secure Joint |

GLASS  
BORO  
3.3

### Fractional Vacuum Distillation Apparatus - Secure Joint

5 piece fractional vacuum distillation apparatus with size 24/29 screw thread joints, made of borosilicate glass exclusive of screw threads.

Includes:

- Vigrex column
- Still head
- Liebig condenser
- Receiver adapter vacuum
- Erlenmeyer flask - 100 ml

| Code   | Item Name                                               |
|--------|---------------------------------------------------------|
| CH1032 | Fractional Vacuum Distillation Apparatus - Secure Joint |

### Hydro Distillation Apparatus - Secure Joint

5 piece hydro distillation apparatus with size 24/29 screw thread joints, made of borosilicate glass exclusive of screw threads. Includes:

- Round bottom boiling flask - two neck - 250 ml
- Dropping funnel - 50 ml
- Still head
- Liebig condenser
- Receiver adapter

| Code   | Item Name          |
|--------|--------------------|
| CH1034 | Hydro Distillation |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Dean and Stark

- Graduated Dean and Stark with interchangeable stopper and open top to allow weighing with steel balls.
- Used with measuring expansion of edible oils and fat.
- Made of borosilicate glass.

| Code      | Capacity | Pack Size |
|-----------|----------|-----------|
| DN-STRK02 | 2 ml     | 1         |
| DN-STRK10 | 10 ml    | 1         |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Plastic Screw Caps

| Code    | Joint Size | Pack Size |
|---------|------------|-----------|
| CH1036A | 14/23      | 10        |
| CH1036B | 19/26      | 10        |
| CH1036C | 24/29      | 10        |



### Rubber O Rings

| Code    | Joint Size | Pack Size |
|---------|------------|-----------|
| CH1038A | 14/23      | 10        |
| CH1038B | 19/26      | 10        |
| CH1038C | 24/29      | 10        |





### Set 10 BU Organic Chemistry Kit

This set contains 22 items with 10/19 ground joint size for maximum versatility. Designed for up to 15 g starting material. Packaged in a die-cut foam lined box. Glassware is made of borosilicate glass.

Includes:

- 1 Liebig condenser
- 1 Dropping funnel, 20 ml
- 1 Flask, pear shape, 25 ml
- 1 Flask, three neck, 25 ml
- 3 Flask, round bottom, 5 ml
- 1 Air leak / steam inlet
- 1 Thermometer -10°C to 250°C
- 1 Dropping pipette
- 1 Rubber teat
- 1 Receiver adapter
- 2 Stoppers
- 1 Still head
- 1 Air condenser
- 1 Filter funnel
- 1 Receiver tube
- 1 Willstätter nail
- 1 Gas inlet tube
- 1 Drying tube
- 1 Stirrer

| Code   | Item Name                       |
|--------|---------------------------------|
| CH0874 | Set 10 BU Organic Chemistry Kit |



### Set 29 BU Organic Chemistry Kit

This set contains 6 items with 10/19 ground joint size for maximum versatility, packaged in a die-cut foam lined box. Glassware is made of borosilicate glass.

Includes:

- 1 Pear shaped flask, 50 ml
- 1 Stillhead
- 1 Liebig condenser
- 1 Screwcap adapter
- 1 Receiver adapter
- 1 Stopper

| Code   | Item Name                       |
|--------|---------------------------------|
| CH0875 | Set 29 BU Organic Chemistry Kit |

### Set 27 BU/M Organic Chemistry Kit

This set contains 9 items with 14/23 ground joint size for maximum versatility. Packaged in a die-cut foam lined box. Glassware made of borosilicate glass.

Includes:

- 1 Flask, pear shaped, 50 ml
- 1 Still head
- 1 Liebig condenser
- 1 Screwcap adapter
- 1 Receiver adapter
- 1 Air leak / steam inlet tube
- 1 Dropping funnel, 50 ml with GP rotaflow tap
- 1 Stopper
- 1 Thermometer

| Code   | Item Name                         |
|--------|-----------------------------------|
| CH0876 | Set 27 BU/M Organic Chemistry Kit |



### Set 27 X/M - Organic Chemistry Kit

This extension kit contains the 8 items required to convert Set 27 BU/M to Set 23 BU/M. All glassware has size 14/23 ground joints. Packaged in a die-cut foam lined box. Glassware made of borosilicate glass.

Includes:

- 1 Round bottom flask, 25 ml
- 1 Air condenser
- 1 Sintered glass funnel
- 1 Drying tube
- 1 Pear shaped 2-neck flask, 50 ml
- 1 Steam inlet tube
- 1 Adapter with T - connector
- 1 Screwcap adapter

| Code   | Item Name                          |
|--------|------------------------------------|
| CH0877 | Set 27 X/M - Organic Chemistry Kit |





### Set 27 BU Organic Chemistry Kit - Secure Joint

This set contains 9 items with size 10/19 screw thread secure joints for maximum versatility and safety. Packaged in a die-cut foam lined Grattells® tray with lid. Glassware made of borosilicate glass.

Includes:

- 1 Flask, pear shaped, 50 ml
- 1 Still head
- 1 Liebig condenser
- 1 Screwcap adapter
- 1 Receiver adapter
- 1 Air leak / steam inlet tube
- 1 Dropping funnel, 50 ml with GP rotaflow tap
- 1 Stopper
- 1 Thermometer

| Code      | Item Name                                      |
|-----------|------------------------------------------------|
| CH0876-SJ | Set 27 BU Organic Chemistry Kit - Secure Joint |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Set 34 BU Organic Chemistry Kit

- This set contains 16 items with various size joints.
- Designed for up to 150 g starting material.
- Packaged in a die-cut foam lined box.
- Glassware is made of borosilicate glass.

Includes:

| Item Name                     | Socket Size | Joint Size |
|-------------------------------|-------------|------------|
| 1 Still head                  | 14/23       | 19/26      |
| 1 Receiver adapter            | 19/26       | 24/29      |
| 1 Air leak / steam inlet      | ----        | 19/26      |
| 1 Liebig condenser            | 19/26       | 19/26      |
| 1 Dropping funnel, 100 ml     | 19/26       | 19/26      |
| 1 Reduction adapter           | 19/26       | 24/29      |
| 1 Multiple adapter            | 19/26       | 24/29      |
| 1 Flask, round bottom, 50 ml  | 24/29       | ----       |
| 1 Flask, Erlenmeyer, 250 ml   | 24/29       | ----       |
| 1 Flask, round bottom, 250 ml | 24/29       | ----       |
| 2 Stoppers                    | ----        | 19/26      |
| 1 Thermometer packet          | ----        | 14/23      |
| 1 Receiver adapter            | 19/26       | ----       |
| 1 Flask, round bottom, 100 ml | 24/29       | ----       |
| 1 Stopper                     | ----        | 24/29      |

| Code   | Item Name                       |
|--------|---------------------------------|
| CH0878 | Set 34 BU Organic Chemistry Kit |



GLASS  
BORO  
3.3

Chemical  
Resistant

### Set 34 BU Organic Chemistry Kit in Grattells® Tray with Lid

- This set contains 16 items with various size joints.
- Designed for up to 150 g starting material.
- Packaged in a die-cut foam lined Grattells® tray with lid.
- Glassware is made of borosilicate glass.

Includes:

| Item Name                     | Socket Size | Joint Size |
|-------------------------------|-------------|------------|
| 1 Still head                  | 14/23       | 19/26      |
| 1 Receiver adapter            | 19/26       | 24/29      |
| 1 Air leak / steam inlet      | ----        | 19/26      |
| 1 Liebig condenser            | 19/26       | 19/26      |
| 1 Dropping funnel, 100 ml     | 19/26       | 19/26      |
| 1 Reduction adapter           | 19/26       | 24/29      |
| 1 Multiple adapter            | 19/26       | 24/29      |
| 1 Flask, round bottom, 50 ml  | 24/29       | ----       |
| 1 Flask, Erlenmeyer, 250 ml   | 24/29       | ----       |
| 1 Flask, round bottom, 250 ml | 24/29       | ----       |
| 2 Stoppers                    | ----        | 19/26      |
| 1 Thermometer packet          | ----        | 14/23      |
| 1 Receiver adapter            | 19/26       | ----       |
| 1 Flask, round bottom, 100 ml | 24/29       | ----       |
| 1 Stopper                     | ----        | 24/29      |

| Code      | Item Name                                      |
|-----------|------------------------------------------------|
| CH0878-GT | Set 34 BU Organic Chemistry Kit in Garter Tray |



### Set 46 BU Organic Chemistry Kit

This set contains 11 items with size 19/26 ground joints, packaged in a die-cut foam lined box. Glassware made of borosilicate glass.

Includes:

- 1 Dropping funnel, 100 ml
- 1 Air leak / steam inlet tube
- 1 Liebig condenser
- 1 Receiver adapter
- 1 Stopper
- 1 Flask, round bottom, 250 ml
- 1 Cone / screwcap adapter
- 1 Still head
- 1 Flask, pear shaped 100 ml
- 1 Filter funnel, 24 mm. dia. perforated plate
- 1 Flask, buchner, 100 ml

| Code   | Item Name                       |
|--------|---------------------------------|
| CH0880 | Set 46 BU Organic Chemistry Kit |



GLASS  
BORO  
3.3GLASS  
BORO  
3.3

### Simple Organic Chemistry Set - Micro

This set contains 8 items with size 14/23 ground joints, packaged in a die-cut foam lined box. Glassware made of borosilicate glass.

Includes:

- Beaker, 25 ml
- Vigrex column, 100 mm
- Liebig condenser, 100 mm
- Still head
- Screw cap adapter
- Receiver adapter
- Pear shaped flask, 25 ml
- Thermometer, -10 to 110°C

| Code   | Item Name                            |
|--------|--------------------------------------|
| CH0881 | Simple Organic Chemistry Set - Micro |

### Simple Organic Chemistry Set - in Gratnells® Tray with Lid

This set contains 8 items with size 14/23 ground joints, packaged in a die-cut foam lined Gratnells® tray with lid. Glassware made of borosilicate glass.

Includes:

- Beaker, 25 ml
- Vigrex column, 100 mm
- Liebig condenser, 100 mm
- Still head
- Screw cap adapter
- Receiver adapter
- Pear shaped flask, 25 ml
- Thermometer, -10 to 110°C

| Code      | Item Name                                                         |
|-----------|-------------------------------------------------------------------|
| CH0881-GT | Simple Organic Chemistry Set - Micro, in Gratnells® Tray with Lid |

### Micro Distillation Apparatus

This set contains 9 items allowing for small scale distillation. Everything is packed in a die-cut foam lined box, glassware made of borosilicate glass.

Includes:

- Round bottom flask, 25 ml
- Conical flask, 50 ml
- Liebig condenser, 100 mm
- Still head
- Receiver adapter
- Funnel, 35 mm
- Beaker, 50 ml
- Burette, 10 ml
- Spirit burner

| Code   | Item Name                     |
|--------|-------------------------------|
| CH0879 | Micro Glass Distillation Unit |

### Kjeldahl Distillation Apparatus

This set contains 6 items with size 19/26 ground joints. Everything is packed in a die-cut foam lined box, glassware made of borosilicate glass.

Includes:

- Kjeldahl Flask, 500 ml
- Multiple adapter
- Dropping funnel, 50 ml
- Splash head
- Leibig condenser
- 4 Joint clips

| Code   | Item Name                       |
|--------|---------------------------------|
| CH0885 | Kjeldahl Distillation Apparatus |



### Simple Distillation Apparatus

This set contains 5 items with size 19/26 ground joints. Everything is packed in a die-cut foam lined box, glassware made of borosilicate glass.

Includes:

- Round bottom flask, 500 ml
- Still head with thermometer support
- Liebig condenser, 300 mm
- Receiver adapter
- 3 Joint clips, 19/26

| Code   | Item Name                     |
|--------|-------------------------------|
| CH0886 | Simple Distillation Apparatus |



### Standard Distillation Apparatus

This set contains 5 items with size 19/26 ground joints. Everything is packed in a die-cut foam lined box, glassware made of borosilicate glass.

Includes:

- Round bottom boiling flask, 500 ml
- Still head
- Liebig condenser, 300 mm
- Receiver adapter
- Separating funnel with stopper
- Supplied without glass beaker and metal stands.

| Code   | Item Name                       |
|--------|---------------------------------|
| CH0890 | Standard Distillation Apparatus |



GLASS  
BORO  
3.3Chemical  
Resistant


### Fractional Distillation Apparatus

This set contains 6 items with size 19/26 ground joints. Everything is packed in a die-cut foam lined box, glassware made of borosilicate glass.

Includes:

- Round-bottom boiling flask, 250ml
- Vigrex column
- Liebig condenser, 300 mm
- Vacuum receiver adapter
- Erlenmeyer flask, 250ml
- 4 Joint clips, 19/26

| Code    | Item Name                         |
|---------|-----------------------------------|
| CH0890B | Fractional Distillation Apparatus |

GLASS  
BORO  
3.3Chemical  
Resistant

### Hydro Distillation Apparatus

This set contains 6 items with size 19/26 ground joints. Everything is packed in a die-cut foam lined box, glassware made of borosilicate glass.

Includes:

- Round bottom, two neck boiling flask, 250ml
- Still head, plain
- Liebig condenser, 300 mm
- Delivery receiver adapter
- Dropping funnel, 100 ml with stopper
- 4 Joint clips, 19/26

| Code    | Item Name                    |
|---------|------------------------------|
| CH0890C | Hydro Distillation Apparatus |

GLASS  
BORO  
3.3Chemical  
Resistant

### Markham Distillation Apparatus

This apparatus is fitted with a spherical ground glass joint between the distillation vessel and condenser for easy assembly, includes screw thread outlets and plastic connectors. The apparatus automatically sucks back when the steam supply to the distillation is cut off.

Glassware made of borosilicate glass.

Includes:

- Funnel plug with cone, size 10/19
- Spherical Joint clip
- Condenser, 1.0 x 10-2m2
- Markham still body, 25 ml capacity

| Code     | Item Name                      |
|----------|--------------------------------|
| MRKH-STL | Markham Distillation Apparatus |



# CLASS - A INDIVIDUALLY CERTIFIED GLASSWARE

 **EISCO**  
CALIBRATION  
LABORATORY

 **CLASS**  
LAB

**CLASS "A" TEST CERTIFICATE**

Certificate No. : BU50-09027/EIS-014 Certificate Issue Date : 17-09-2018  
Calibration Method : Gravimetric Method Calibration Date : 14-09-2018

**PRODUCT DETAIL**

Description of the Item: Burette 50 ml Cat No. : CH0240C Batch No. : BU50-09-027  
Nominal Volume: 50 ml Tolerance: ± 0.05 ml Serial No. : EIS-014

**STANDARD EQUIPMENTS USED FOR CALIBRATION**

| Name        | Weighting Balance | Thermometer    |
|-------------|-------------------|----------------|
| Make        | OHAUS             | CE             |
| I.D. No.    | WB-18             | DTM-01         |
| Range       | 0.210 gm.         | -50°C to 200°C |
| Least Count | 0.1 mg            | 0.1°C          |

**CALIBRATION RESULT**

| Nominal Vol. | Certified Vol. | Delivery Time (sec.) |
|--------------|----------------|----------------------|
| 50 ml        | 50.02 ml       | 78 Sec.              |

**TERMS**

- The above product has been examined for compliance to DIN ISO 385, Class A.
- This certificate shall not be reproduced, except in full, without the written permission from EISCO.
- Results reported are valid under lab conditions at the time of calibration.
- Laboratory Standards are traceable to NPL India, traceable to International standards.
- It has been found to satisfy the specified requirements. The above item has been calibrated at 24.51°C at RH 42%.

Calibrated by : JNY Approved by :   
Quality Manager

**Eisco Scientific LLC**

INDIVIDUAL  
WORK  
CERTIFICATESGLASS  
BORO  
3.3**Burettes - Glass Key Stopcock Class A with individual work certificate**

- Single straight bore burettes with glass key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1900A | 10            | 0.05           | 0.03               | 1         |
| CH1900B | 25            | 0.10           | 0.05               | 1         |
| CH1900C | 50            | 0.10           | 0.05               | 1         |
| CH1900D | 100           | 0.20           | 0.10               | 1         |

INDIVIDUAL  
WORK  
CERTIFICATESGLASS  
BORO  
3.3**Burettes - Schellbach, Glass Key Stopcock Class A with individual work certificate**

- Single bore Schellback burettes with glass key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385, standard.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1904A | 10            | 0.05           | 0.03               | 1         |
| CH1904B | 25            | 0.10           | 0.05               | 1         |
| CH1904C | 50            | 0.10           | 0.05               | 1         |
| CH1904D | 100           | 0.20           | 0.10               | 1         |

INDIVIDUAL  
WORK  
CERTIFICATESGLASS  
BORO  
3.3**Burettes - Glass with PTFE Key Stopcock Class A & individual work certificate**

- Single straight bore burettes with interchangeable PTFE key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1908A | 10            | 0.05           | 0.03               | 1         |
| CH1908B | 25            | 0.10           | 0.05               | 1         |
| CH1908C | 50            | 0.10           | 0.05               | 1         |
| CH1908D | 100           | 0.20           | 0.10               | 1         |

### Burettes - Schellbach, Glass with PTFE Key Stopcock Class A & individual work certificate

- Single bore Schellbach burettes with PTFE key stopcock.
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1912A | 10            | 0.05           | 0.03               | 1         |
| CH1912B | 25            | 0.10           | 0.05               | 1         |
| CH1912C | 50            | 0.10           | 0.05               | 1         |
| CH1912D | 100           | 0.20           | 0.10               | 1         |



INDIVIDUAL WORK CERTIFICATES



GLASS BORO 3.3



### Automatic Zero Burettes - Glass with PTFE Key Stopcock Class A & individual work certificate

- Automatic zero burettes mounted on glass reservoir with rubber below.
- PTFE key stopcock
- Made of borosilicate glass.
- Accuracy as per DIN ISO 385.
- Printed in Blue.

Class-A

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1916A | 10            | 0.05           | 0.03               | 1         |
| CH1916B | 25            | 0.10           | 0.05               | 1         |
| CH1916C | 50            | 0.10           | 0.05               | 1         |
| CH1916D | 100           | 0.20           | 0.10               | 1         |



INDIVIDUAL WORK CERTIFICATES



GLASS BORO 3.3



### Graduated Cylinders - Glass, Class-A with Round Base & individual work certificate

- Class-A graduated cylinder with round base and spout.
- Made of borosilicate glass.
- As per DIN EN ISO 4788 standard.
- Accuracy as per DIN ISO 385.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1922A | 5             | 0.1            | 0.05               | 2         |
| CH1922B | 10            | 0.2            | 0.10               | 2         |
| CH1922C | 25            | 0.5            | 0.25               | 2         |
| CH1922D | 50            | 1.0            | 0.50               | 2         |
| CH1922E | 100           | 1.0            | 0.50               | 2         |
| CH1922F | 250           | 2.0            | 1.00               | 2         |
| CH1922G | 500           | 5.0            | 2.50               | 2         |
| CH1922H | 1000          | 10.0           | 5.00               | 2         |
| CH1922I | 2000          | 20.0           | 10.00              | 1         |



INDIVIDUAL WORK CERTIFICATES



GLASS BORO 3.3





### Graduated Cylinders - Glass, Class-A with Hexagonal Base & individual work certificate

- Class-A graduated cylinder with hexagonal base and spout.
- Made of borosilicate glass.
- As per DIN EN ISO 4788 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1926A | 5             | 0.1            | 0.05               | 2         |
| CH1926B | 10            | 0.2            | 0.10               | 2         |
| CH1926C | 25            | 0.5            | 0.25               | 2         |
| CH1926D | 50            | 1.0            | 0.50               | 2         |
| CH1926E | 100           | 1.0            | 0.50               | 2         |
| CH1926F | 250           | 2.0            | 1.00               | 2         |
| CH1926G | 500           | 5.0            | 2.50               | 2         |
| CH1926H | 1000          | 10.0           | 5.00               | 2         |
| CH1926I | 2000          | 20.0           | 10.00              | 1         |



### Graduated Cylinders - Glass, Class-A with Hexagonal Base, USP & individual work certificate

- Class-A graduated cylinder with hexagonal base and spout.
- Made of borosilicate glass.
- USP <31> Complaint.
- As per DIN EN ISO 4788 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1930A | 5             | 0.1            | 0.05               | 2         |
| CH1930B | 10            | 0.2            | 0.10               | 2         |
| CH1930C | 25            | 0.5            | 0.17               | 2         |
| CH1930D | 50            | 1.0            | 0.25               | 2         |
| CH1930E | 100           | 1.0            | 0.50               | 2         |
| CH1930F | 250           | 2.0            | 1.00               | 2         |
| CH1930G | 500           | 5.0            | 2.00               | 2         |
| CH1930H | 1000          | 10.0           | 3.00               | 2         |
| CH1930I | 2000          | 20.0           | 6.00               | 1         |

### Graduated Cylinders - Glass, Class-A ASTM with Hexagonal Base, Guard & individual work certificate

- Class-A graduated cylinder with hexagonal base.
- Protective guard and spout.
- Made of borosilicate glass.
- As per ASTM E1272 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Pack Size |
|---------|---------------|----------------|--------------------|-----------|
| CH1932A | 5             | 0.1            | 0.05               | 2         |
| CH1932B | 10            | 0.2            | 0.10               | 2         |
| CH1932C | 25            | 0.5            | 0.17               | 2         |
| CH1932D | 50            | 1.0            | 0.25               | 2         |
| CH1932E | 100           | 1.0            | 0.50               | 2         |
| CH1932F | 250           | 2.0            | 1.00               | 2         |
| CH1932G | 500           | 5.0            | 2.00               | 2         |
| CH1932H | 1000          | 10.0           | 3.0                | 2         |
| CH1932I | 2000          | 20.0           | 6.0                | 1         |



### Cylinders - Glass, Class-A with Hexagonal Base, Stoppers & individual work certificate

- Class-A graduated cylinder with hexagonal base and spout.
- Made of borosilicate glass.
- Each graduated cylinder comes with interchangeable hexagonal glass or polypropylene stoppers.
- As per DIN EN ISO 4788 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Sub. Div. (ml) | Tolerance $\pm$ ml | Joint | Pack Size |
|---------|---------------|----------------|--------------------|-------|-----------|
| CH1936A | 10            | 0.2            | 0.10               | 10/19 | 2         |
| CH1936B | 25            | 0.5            | 0.25               | 14/23 | 2         |
| CH1936C | 50            | 1.0            | 0.50               | 14/23 | 2         |
| CH1936D | 100           | 1.0            | 0.50               | 19/26 | 2         |
| CH1936E | 250           | 2.0            | 1.00               | 24/29 | 2         |
| CH1936F | 500           | 5.0            | 2.50               | 24/29 | 2         |
| CH1936G | 1000          | 10.0           | 5.00               | 29/32 | 2         |
| CH1936I | 2000          | 20.0           | 10.00              | 34/35 | 1         |



INDIVIDUAL  
WORK  
CERTIFICATES



### Flasks, Volumetric Class-A with Stopper & individual work certificate

- Flat bottom class- A volumetric flask with one graduation mark.
- Fitted with chemical resistant polypropylene stopper.
- Made of borosilicate glass.
- As per DIN ISO 1042 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Tolerance $\pm$ ml | Socket Size | Pack Size |
|---------|---------------|--------------------|-------------|-----------|
| CH1952A | 5             | 0.025              | 10/19       | 2         |
| CH1952B | 10            | 0.025              | 10/19       | 2         |
| CH1952C | 20            | 0.040              | 10/19       | 2         |
| CH1952D | 25            | 0.040              | 10/19       | 2         |
| CH1952E | 50            | 0.060              | 12/21       | 2         |
| CH1952F | 100           | 0.100              | 14/23       | 2         |
| CH1952G | 200           | 0.150              | 14/23       | 2         |
| CH1952H | 250           | 0.150              | 14/23       | 2         |
| CH1952I | 500           | 0.250              | 19/26       | 2         |
| CH1952J | 1000          | 0.400              | 24/29       | 2         |
| CH1952K | 2000          | 0.600              | 29/32       | 1         |

INDIVIDUAL  
WORK  
CERTIFICATES



### Flasks, Volumetric, Amber with Stopper, Class-A & individual work certificate

- Flat bottom amber volumetric flask with one graduation mark.
- Fitted with chemical resistant polypropylene stopper.
- Made of amber borosilicate glass.
- As per DIN ISO 1042 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Tolerance $\pm$ ml | Socket Size | Pack Size |
|---------|---------------|--------------------|-------------|-----------|
| CH1960A | 5             | 0.025              | 10/19       | 2         |
| CH1960B | 10            | 0.025              | 10/19       | 2         |
| CH1960C | 20            | 0.040              | 10/19       | 2         |
| CH1960D | 25            | 0.040              | 10/19       | 2         |
| CH1960E | 50            | 0.060              | 12/21       | 2         |
| CH1960F | 100           | 0.100              | 14/23       | 2         |
| CH1960G | 200           | 0.150              | 14/23       | 2         |
| CH1960H | 250           | 0.150              | 14/23       | 2         |
| CH1960I | 500           | 0.250              | 19/26       | 2         |
| CH1960J | 1000          | 0.400              | 24/29       | 2         |
| CH1960K | 2000          | 0.600              | 29/32       | 1         |

### Flasks, Volumetric Class- A ASTM with Solid Glass, Stopper & individual work certificate

- Flat bottom class- A ASTM volumetric flask with one graduation mark.
- Fitted with solid glass stopper.
- Made of borosilicate glass.
- As per ASTM-E288 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Tolerance $\pm$ ml | Stopper No. | Pack Size |
|---------|---------------|--------------------|-------------|-----------|
| CH1964A | 5             | 0.020              | 9           | 2         |
| CH1964B | 10            | 0.020              | 9           | 2         |
| CH1964C | 25            | 0.030              | 9           | 2         |
| CH1964D | 50            | 0.050              | 13          | 2         |
| CH1964E | 100           | 0.080              | 13          | 2         |
| CH1964F | 200           | 0.100              | 16          | 2         |
| CH1964G | 250           | 0.120              | 16          | 2         |
| CH1964H | 500           | 0.200              | 19          | 2         |
| CH1964I | 1000          | 0.300              | 22          | 2         |
| CH1964J | 2000          | 0.500              | 27          | 1         |



INDIVIDUAL WORK CERTIFICATES

ASTM

GLASS BORO 3.3

Chemical Resistant

### Flasks, Volumetric Class- A, USP with PP Stopper & individual work certificate

- Flat bottom class- A volumetric flask with one graduation mark.
- Fitted with chemical resistant PP stopper.
- USP <31> Compliant.
- As per ISO 1042 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Tolerance $\pm$ ml | Socket Size | Pack Size |
|---------|---------------|--------------------|-------------|-----------|
| CH1972A | 5             | 0.020              | 10/19       | 2         |
| CH1972B | 10            | 0.020              | 10/19       | 2         |
| CH1972C | 25            | 0.030              | 10/19       | 2         |
| CH1972D | 50            | 0.050              | 10/19       | 2         |
| CH1972E | 100           | 0.080              | 12/21       | 2         |
| CH1972F | 200           | 0.100              | 14/23       | 2         |
| CH1972G | 250           | 0.120              | 14/23       | 2         |
| CH1972H | 500           | 0.200              | 14/23       | 2         |
| CH1972I | 1000          | 0.300              | 19/26       | 2         |
| CH1972J | 2000          | 0.500              | 24/29       | 1         |



INDIVIDUAL WORK CERTIFICATES

ISO 1042

USP

GLASS BORO 3.3

Chemical Resistant

INDIVIDUAL  
WORK  
CERTIFICATESGLASS  
BORO  
3.3

### Pipettes - Glass, Class-A with individual work certificate

- Bulb form pipettes, color coded with one mark accuracy & class-A.
- Made of borosilicate glass.
- As per DIN ISO 648 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Tolerance $\pm$ ml | Colour Code | Pack Size |
|---------|---------------|--------------------|-------------|-----------|
| CH1990A | 1             | 0.007              | Blue        | 1         |
| CH1990B | 2             | 0.010              | Orange      | 1         |
| CH1990C | 5             | 0.015              | White       | 1         |
| CH1990D | 10            | 0.020              | Red         | 1         |
| CH1990E | 15            | 0.025              | Green       | 1         |
| CH1990F | 20            | 0.030              | Yellow      | 1         |
| CH1990G | 25            | 0.030              | Blue        | 1         |
| CH1990H | 50            | 0.050              | Red         | 1         |
| CH1990I | 100           | 0.080              | Yellow      | 1         |

INDIVIDUAL  
WORK  
CERTIFICATESGLASS  
BORO  
3.3

### Serological Pipettes - Glass, Class-A with individual work certificate

- Color coded serological pipettes calibrated for delivery & class-A.
- Made of borosilicate glass.
- As per DIN ISO 835 standard.
- Printed in Blue.

| Code    | Capacity (ml) | Tolerance $\pm$ ml | Colour Code | Pack Size |
|---------|---------------|--------------------|-------------|-----------|
| CH1996A | 1.0           | 0.006              | Yellow      | 10        |
| CH1996B | 2.0           | 0.010              | Black       | 10        |
| CH1996C | 5.0           | 0.030              | Red         | 10        |
| CH1996D | 10.0          | 0.050              | Orange      | 10        |
| CH1996E | 25.0          | 0.100              | White       | 5         |

**eisco**  
INDUSTRIAL

# METALWARE



## EISCO PREMIUM METAL CLAMPS

- Research grade quality for industrial and academic labs
- Solid zinc die-cast construction made from Zamac 3 Alloy + aluminium alloy
- Complete non-ferrous construction, rust resistant and durable in even the toughest scientific testing environments
- Wide range of options for various laboratory applications
- High grade non-toxic powdered coating



### Premium Clamp - Three Prong, Vinyl Coated

- 3 prong vinyl coated clamp attached to stainless steel rod
- Clamp constructed of powder coated high grade aluminium alloy
- Maximum clamp opening is 8.5 cm
- 14 cm stainless steel connecting rod
- Complete metal construction, thread fitted actuating joints
- 10 mm diameter stainless steel connecting rod

| Code    | Specs                                       |
|---------|---------------------------------------------|
| PRCMP01 | 1 Premium Clamp - Three Prong, Vinyl Coated |



### Premium Clamp - Three Prong, Vinyl Coated, with Bosshead

- 3 prong vinyl coated clamp attached to swivel bosshead
- Clamp constructed of powder coated high grade aluminium alloy
- Maximum clamp opening is 8.5 cm
- Swivel between clamp and bosshead
- Bosshead accommodates up to 18 mm rod
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      | Specs                                                      |
|-----------|------------------------------------------------------------|
| PRCMP01BH | 2 Premium Clamp - Three Prong, Vinyl Coated, with Bosshead |

**Premium Clamp - Four Prong, Cork Lined**

- 4 prong cork lined clamp attached to stainless steel rod
- Clamp constructed of powder coated high grade aluminium alloy
- Maximum clamp opening is 10.5 cm
- 14 cm stainless steel connecting rod
- Complete metal construction, thread fitted actuating joints
- 10 mm diameter stainless steel connecting rod

| Code    | Specs                                    |
|---------|------------------------------------------|
| PRCMP02 | 1 Premium Clamp - Four Prong, Cork Lined |

**Premium Clamp - Four Prong, Cork Lined, with Bosshead**

- 4 prong cork lined clamp attached to swivel bosshead
- Clamp constructed of powder coated high grade aluminium alloy
- Maximum clamp opening is 10.5 cm
- Swivel between clamp and bosshead
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      | Specs                                                   |
|-----------|---------------------------------------------------------|
| PRCMP02BH | 2 Premium Clamp - Four Prong, Cork Lined, with Bosshead |

**Premium Clamp - Four Prong, Vinyl Coated**

- 4 prong vinyl coated clamp attached to stainless steel rod
- Clamp constructed of powder coated high grade zinc alloy
- Maximum clamp opening is 10.5 cm
- 14 cm stainless steel connecting rod
- Complete metal construction, thread fitted actuating joints
- 10 mm diameter stainless steel connecting rod

| Code    | Specs                                      |
|---------|--------------------------------------------|
| PRCMP03 | 3 Premium Clamp - Four Prong, Vinyl Coated |

**Premium Clamp - Four Prong, Vinyl Coated, with Bosshead**

- 4 prong vinyl coated clamp attached to swivel bosshead
- Clamp constructed of powder coated high grade zinc alloy
- Maximum clamp opening is 10.5 cm
- Swivel between clamp and bosshead
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      | Specs                                                     |
|-----------|-----------------------------------------------------------|
| PRCMP03BH | 4 Premium Clamp - Four Prong, Vinyl Coated, with Bosshead |





### Premium Clamp - Two Prong, Vinyl Coated, Dual Adjustment

- 2 prong vinyl coated clamp attached to stainless steel rod
- Heavy duty spring loaded clamp jaw with dual screw adjustment to fine tune any fitting
- Made of zinc and powder coated alloy with all metal adjustment screws
- Maximum clamp opening is 7 cm
- 12.5 mm diameter stainless steel connecting rod

| Code    | Specs                                                      |
|---------|------------------------------------------------------------|
| PRCMP04 | 1 Premium Clamp - Two Prong, Vinyl Coated, Dual Adjustment |



### Premium Clamp - Two Prong, Vinyl Coated, with Bosshead and Dual Adjustment

- 2 prong vinyl coated clamp attached to swivel bosshead
- Heavy duty spring loaded clamp jaw with dual screw adjustment to fine tune any fitting
- Clamp constructed of powder coated high grade zinc alloy with all metal screws
- Maximum clamp opening is 7 cm
- Swivel between clamp and bosshead
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      | Specs                                                                        |
|-----------|------------------------------------------------------------------------------|
| PRCMP04BH | 2 Premium Clamp - Two Prong, Vinyl Coated, with Bosshead and Dual Adjustment |



### Premium Clamp - Two Prong, Vinyl Coated

- 2 prong vinyl coated clamp attached to stainless steel rod
- Clamp constructed of powder coated high grade zinc alloy
- Maximum clamp opening is 11 cm
- 14 cm stainless steel connecting rod
- Complete metal construction, thread fitted actuating joints
- 10 mm diameter stainless steel connecting rod

| Code    | Specs                                     |
|---------|-------------------------------------------|
| PRCMP05 | 3 Premium Clamp - Two Prong, Vinyl Coated |



### Premium Clamp - Two Prong, Vinyl Coated, with Bosshead

- 2 prong vinyl coated clamp attached to swivel bosshead
- Clamp constructed of powder coated high grade zinc alloy
- Maximum clamp opening is 11 cm
- Swivel between clamp and bosshead
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      | Specs                                                    |
|-----------|----------------------------------------------------------|
| PRCMP05BH | 4 Premium Clamp - Two Prong, Vinyl Coated, with Bosshead |



### Premium Clamp - Three Prong, Vinyl Coated, Dual Adjustment

- 3 prong vinyl coated clamp attached to stainless steel rod
- Heavy duty spring loaded clamp jaw with dual screw adjustment to fine tune any fitting
- Made of zinc and powder coated alloy with all metal adjustment screws
- Maximum clamp opening is 6 cm
- 12.5 mm diameter stainless steel connecting rod

| Code    |   | Specs                                                      |
|---------|---|------------------------------------------------------------|
| PRCMP06 | 1 | Premium Clamp - Three Prong, Vinyl Coated, Dual Adjustment |



### Premium Clamp - Three Prong, Vinyl Coated, with Bosshead and Dual Adjustment

- 3 prong vinyl coated clamp attached to swivel bosshead
- Heavy duty spring loaded clamp jaw with dual screw adjustment to fine tune any fitting
- Clamp constructed of powder coated high grade zinc alloy with all metal screws
- Maximum clamp opening is 6 cm
- Swivel between clamp and bosshead
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      |   | Specs                                                                        |
|-----------|---|------------------------------------------------------------------------------|
| PRCMP06BH | 2 | Premium Clamp - Three Prong, Vinyl Coated, with Bosshead and Dual Adjustment |



### Premium Clamp - Three Prong, Vinyl Coated

- 3 prong vinyl coated clamp attached to stainless steel rod
- Clamp constructed of powder coated high grade aluminium alloy
- Maximum clamp opening is 8.5 cm
- 14 cm stainless steel connecting rod
- Complete metal constructions, thread fitted actuating joints
- 10 mm diameter stainless steel connecting rod

| Code    |   | Specs                                     |
|---------|---|-------------------------------------------|
| PRCMP07 | 3 | Premium Clamp - Three Prong, Vinyl Coated |



### Premium Clamp - Three Prong, Vinyl Coated, with Bosshead

- 3 prong vinyl coated clamp attached to swivel bosshead
- Clamp constructed of powder coated high grade aluminium alloy
- Maximum clamp opening is 8.5 cm
- Swivel between clamp and bosshead
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      |   | Specs                                                    |
|-----------|---|----------------------------------------------------------|
| PRCMP07BH | 4 | Premium Clamp - Three Prong, Vinyl Coated, with Bosshead |





### Premium Clamp - Two Prong, Cork Lined

- 2 prong cork lined clamp attached to stainless steel rod
- Clamp constructed of powder coated high grade zinc alloy.
- Maximum clamp opening is 12 cm
- 14 cm stainless steel connecting rod
- Complete metal construction, thread fitted actuating joints
- 10 mm diameter stainless steel connecting rod

| Code    | Specs                                   |
|---------|-----------------------------------------|
| PRCMP08 | 1 Premium Clamp - Two Prong, Cork Lined |



### Premium Clamp - Two Prong, Cork Lined, with Bosshead

- 2 prong cork lined clamp attached to swivel bosshead
- Clamp constructed of powder coated high grade zinc alloy
- Maximum clamp opening is 11 cm
- Swivel between clamp and bosshead
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      | Specs                                                  |
|-----------|--------------------------------------------------------|
| PRCMP08BH | 2 Premium Clamp - Two Prong, Cork Lined, with Bosshead |



### Premium Clamp - Three Prong, Cork Lined

- 3 prong cork lined clamp attached to stainless steel rod
- Clamp constructed of powder coated high grade aluminium alloy.
- Maximum clamp opening is 8.5 cm
- 14 cm stainless steel connecting rod
- Complete metal construction, thread fitted actuating joints
- 10 mm diameter stainless steel connecting rod

| Code    | Specs                                     |
|---------|-------------------------------------------|
| PRCMP09 | 3 Premium Clamp - Three Prong, Cork Lined |



### Premium Clamp - Three Prong, Cork Lined, with Bosshead

- 3 prong cork lined clamp attached to swivel bosshead
- Clamp constructed of powder coated high grade aluminium alloy
- Maximum clamp opening is 8.5 cm
- Swivel between clamp and bosshead
- High strength alloy bosshead
- Bosshead accommodates up to 18 mm rod

| Code      | Specs                                                    |
|-----------|----------------------------------------------------------|
| PRCMP09BH | 4 Premium Clamp - Three Prong, Cork Lined, with Bosshead |

### Premium Bosshead - Die Cast

- Made from high grade zinc alloy, finished in non-toxic powder coating
- Threaded metal one piece clamping screws
- Metal thumb screws
- Able to hold rods up to 16 mm in diameter
- Holds rods orthogonal (90 degrees) to one another
- All metal construction supports heavy loads and prevents torque damage

| Code   | Specs                         |
|--------|-------------------------------|
| PRBH01 | 1 Premium Bosshead - Die Cast |



1

### Premium Bosshead - Metal

- Threaded metal one piece clamping screws
- Metal thumb screws
- Able to hold rods up to 15 mm in diameter
- Holds rods orthogonal (90 degrees) to one another
- All metal construction supports heavy loads and prevents torque damage

| Code   | Specs                      |
|--------|----------------------------|
| PRBH02 | 2 Premium Bosshead - Metal |



2

### Premium Bosshead - Heavy Duty

- Made from high grade zinc alloy, finished in non-toxic powder coating
- Thicker than normal 7 mm threaded clamping screws
- Metal thumb screws
- Able to hold rods up to 16 mm in diameter
- Holds rods orthogonal (90 degrees) to one another
- All metal construction with thicker casting to support more load and prevent torque damage

| Code   | Specs                           |
|--------|---------------------------------|
| PRBH03 | 3 Premium Bosshead - Heavy Duty |



3

### Premium Bosshead - Universal Swivel

- Made from high grade zinc alloy, finished in non-toxic powder coating
- Threaded metal one piece clamping screws
- Metal thumb screws
- Able to hold rods up to 16 mm in diameter
- Adjusts to any angle, both bossheads can swivel 360 degrees
- All metal construction supports heavy loads and prevents torque damage

| Code   | Specs                                 |
|--------|---------------------------------------|
| PRBH04 | 4 Premium Bosshead - Universal Swivel |



4

### Premium Bosshead - Square Head

- Made from high grade zinc alloy, finished in non-toxic powder coating
- Threaded metal one piece clamping screws
- Metal thumb screws
- Able to hold rods up to 15 mm in diameter
- Holds rods orthoangle (90 degrees) to one another
- All metal construction supports heavy loads and prevents torque damage

| Code   | Specs                            |
|--------|----------------------------------|
| PRBH05 | 5 Premium Bosshead - Square Head |



5

# Next Gen™ Burners

**eisco**  
INDUSTRIAL

## Teachers and Lab Tech Testimonials

*"I love the simplicity of the safety features to make my class safer."*

*- High school chemistry teacher*

*"This burner does not tip and still fits under all our stands."*


*- Community college teacher*

*"Should have been thought of 100 years ago. Great design, nice work."*

*- 8th grade physical science teacher*



# Burner Features





1

Tube Dia.- 13 mm  
Ht.- 145 mm  
Wt.- 300 g

### Bunsen Burner - Simple

- Nickel plated brass tube with adjustable air regulator.
- Fitted with integrated flame stabilizer and mounted on tip resistant die-cast base.
- Height 145 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                       |
|------------|-----------------------------|
| CH0971LP 1 | Bunsen Burner - LPG         |
| CH0971NG   | Bunsen Burner - Natural Gas |



2

Tube Dia.- 11 mm  
Ht.- 125 mm  
Wt.- 290 g

### Bunsen Burner

- Nickel plated brass tube with adjustable air regulator.
- Fitted with integrated flame stabilizer and mounted on tip resistant die-cast base.
- Height 125 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                       |
|------------|-----------------------------|
| CH0972LP 2 | Bunsen Burner - LPG         |
| CH0972NG   | Bunsen Burner - Natural Gas |



3

Tube Dia.- 8 mm  
Ht.- 100 mm  
Wt.- 280 g

### Bunsen Burner - Semi Micro

- Nickel plated brass tube with adjustable air regulator.
- Fitted with integrated flame stabilizer and mounted on tip resistant die-cast base.
- Height 100 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                       |
|------------|-----------------------------|
| CH0973LP 3 | Bunsen Burner - LPG         |
| CH0973NG   | Bunsen Burner - Natural Gas |



4

Tube Dia.- 13 mm  
Ht.- 155 mm  
Wt.- 370 g

### Bunsen Burner - Teclu

- This premium burner features a nickel plated brass tube with an adjustable air regulator and fuel adjustment needle valve.
- It is fitted with a integrated flame stabilizer and mounted on a tip resistant die-cast base.
- Height 155 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                       |
|------------|-----------------------------|
| CH0974LP 4 | Bunsen Burner - LPG         |
| CH0974NG   | Bunsen Burner - Natural Gas |

### Bunsen Burner - Basic

- This burner is our most economical burner.
- Includes a nickel plated brass tube with adjustable air regulator but not a flame stabilizer.
- It is mounted on a tip resistant die-cast base.
- Height 145 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                       |
|------------|-----------------------------|
| CH0991LP ① | Bunsen Burner - LPG         |
| CH0991NG   | Bunsen Burner - Natural Gas |



### Bunsen Burner - Flame Stabilizer

- This burner includes a nickel plated brass tube with adjustable air regulator and flame stabilizer mounted on a tip resistant die-cast base.
- Height 160 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                                                                                  |
|------------|----------------------------------------------------------------------------------------|
| CH0992LP ② | Bunsen Burner - Flame Stabilizer - LPG<br>BTU: 2000 - 3200; Max. Temp.: 1550 °C        |
| CH0992NG   | Bunsen Burner - Flame Stabilizer - Natural Gas<br>BTU: 800 - 1200; Max. Temp.: 1550 °C |



### Bunsen Burner - Flame Stabilizer, Gas adjustment

- This burner includes a nickel plated brass tube with adjustable air regulator.
- Gas adjustment needle and flame stabilizer mounted on a tip resistant die-cast base.
- Height 160 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                                                                                  |
|------------|----------------------------------------------------------------------------------------|
| CH0993LP ③ | Bunsen Burner - Flame Stabilizer - LPG<br>BTU: 2000 - 3200; Max. Temp.: 1550 °C        |
| CH0993NG   | Bunsen Burner - Flame Stabilizer - Natural Gas<br>BTU: 800 - 1200; Max. Temp.: 1550 °C |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS



1

**Bunsen Burner - Micro**

- This burner is perfect for tight spaces.
- It includes a nickel plated brass tube with adjustable air regulator mounted on a tip resistant die-cast base.
- Height 90 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                       |
|------------|-----------------------------|
| CH0994LP 1 | Burner Bunsen - LPG         |
| CH0994NG   | Burner Bunsen - Natural Gas |



2

**Bunsen Burner - Tirrill**

- The Tirrill burner is a type of Bunsen burner that features a wing nut on the bottom of the tube for gas adjustment.
- The nickel plated brass tube screws for fine adjustment of the air to fuel ratio.
- Mounted on a tip resistant die-cast base.
- Height 150 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                                                                                  |
|------------|----------------------------------------------------------------------------------------|
| CH0995LP 2 | Bunsen Burner - Flame Stabilizer - LPG<br>BTU: 2000 - 3200; Max. Temp.: 1550 °C        |
| CH0995NG   | Bunsen Burner - Flame Stabilizer - Natural Gas<br>BTU: 800 - 1200; Max. Temp.: 1550 °C |



3

**Bunsen Burner - Meker**

- The Meker burner is a relatively radical type of Bunsen burner that is designed to produce an extremely hot and stable flame.
- The Meker burner was designed around the principle that the hottest point of a flame is the outer tip of the inner flame.
- It features a flame stabilizer with a grid that separates the flame into an array of much smaller flames creating several "outer tips" or zones of highest temperature.
- This allows the average temperature of the flame to rise from 935°C to 1180°C.
- This burner also includes a wing nut on the bottom of the nickel plated brass tube for gas adjustment as well as an airflow regulator.
- It is mounted on a tip resistant die-cast base.
- Height 220 mm. Liquid propane gas (LPG) or Natural Gas options available.

| Code       | Specs                                                                                  |
|------------|----------------------------------------------------------------------------------------|
| CH0996LP 3 | Bunsen Burner - Flame Stabilizer - LPG<br>BTU: 1900 - 3100; Max. Temp.: 1780 °C        |
| CH0996NG   | Bunsen Burner - Flame Stabilizer - Natural Gas<br>BTU: 600 - 1100; Max. Temp.: 1780 °C |

**Bunsen Burner**

- Same as above but mounted on more safer, easy to handle diecasted base.

| Code          | Specs                       |
|---------------|-----------------------------|
| CH0093QP-LP ① | Bunsen Burner - LPG         |
| CH0093QP-NG ② | Bunsen Burner - Natural Gas |



Check out  
our new website  
[www.eiscoindustrial.com](http://www.eiscoindustrial.com)

GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Clips

**Mohr's Clips**

- Nickel plated
- Maximum opening 18 mm.

| Code      | Specs                |
|-----------|----------------------|
| CH0290A 1 | Mohr's Clips - Iron  |
| CH0290B   | Mohr's Clips - Brass |

**Hoffman's Clips**

- Nickel plated, hinged bottom bar for direct fitting and removal without access to ends of tube.
- Protecting plate to prevent damage to tubing by end of screw.
- Maximum opening 15 mm.

| Code      | Specs                                     |
|-----------|-------------------------------------------|
| CH0298B   | Hoffman's Clips - Brass, Closed-Closed    |
| CH0298C 2 | Hoffman's Clips - Brass, Open-Closed-Open |
| CH0298D 3 | Hoffman's Clips - Brass, Open-Closed      |



Come chat with us  
Live on [www.eiscoindustrial.com](http://www.eiscoindustrial.com)

**Hoffman's Clips - Heavy Duty**

- Die-cast
- Nickel plated
- Maximum opening 20 mm

| Code      | Specs                            |
|-----------|----------------------------------|
| CH0299A 4 | Hoffman's Clips - Die-Cast, Open |
| CH0299B 5 | Clips Hoffman Die-Cast, Closed   |

**Wire Gauze - Iron with Ceramic Center**

- Square wire gauze made of iron with ceramic center.
- Pack of 10.

| Code                 | Size    |
|----------------------|---------|
| CH0525A <sup>1</sup> | 12.5 cm |
| CH0525B <sup>2</sup> | 15 cm   |

**Wire Gauze - Iron Plain**


- Square wire gauze.
- Pack of 10.

| Code                 | Size    |
|----------------------|---------|
| CH0526A <sup>3</sup> | 12.5 cm |
| CH0526B              | 15 cm   |

**Wire Gauze - Stainless Steel with Ceramic Center**


- Square wire gauze made of stainless steel with ceramic center.
- Pack of 10.

| Code                 | Size    |
|----------------------|---------|
| CH0527A <sup>4</sup> | 12.5 cm |
| CH0527B <sup>5</sup> | 15 cm   |

**Wire Gauze - Stainless Steel Plain**

- Square wire gauze made of stainless steel.
- Pack of 10.

| Code                 | Size    |
|----------------------|---------|
| CH0527C <sup>6</sup> | 12.5 cm |
| CH0527D              | 15 cm   |

**Clay Pipe Triangle**

- Clay pipes on thick iron wire.

| Code                 | Size   |
|----------------------|--------|
| CH0530A <sup>7</sup> | 5 cm   |
| CH0530B              | 6 cm   |
| CH0530C <sup>8</sup> | 7.5 cm |

**Spirit Lamp - Metal**

- Die-pressed sheet metal with woven wick and cap.

| Code                 | Specs                     |
|----------------------|---------------------------|
| CH0535B <sup>5</sup> | Spirit Lamp Glass, 80 ml  |
| CH0535C              | Spirit Lamp Glass, 150 ml |



## Spatulas &amp; Scoops

**Spatula Spoon - Stainless Steel, Wide Handle**

- Flat handle with one concave spoon end and one flat spatula end.
- Made of polished stainless steel.

| Code    | Specs              |
|---------|--------------------|
| CH0632A | 1 Length - 12.5 cm |
| CH0632B | Length - 15 cm     |
| CH0632C | Length - 20 cm     |

**Spatula - Stainless Steel, Wide Handle, Dual End**

- Flat handle with two flat spatula ends.
- Made of polished stainless steel.

| Code    | Specs            |
|---------|------------------|
| CH0634C | 3 Length - 20 cm |

**Scoop - Stainless Steel**

- Semicircular cross section with one blunt end and one tapered end, 16 cm long.
- Made of stainless steel.
- Lacquered hardwood handle available separately (not pictured).

| Code    | Specs                  |
|---------|------------------------|
| CH0635A | 4 Spatula - Scoop Type |
| CH0635B | Spatula Handle         |

**Micro Spatula Spoon - Stainless Steel**

- 15 cm long with one concave spoon end and one flat spatula end (30 x 8 mm).
- Handle diameter 2 mm
- Made of stainless steel.

| Code    | Specs          |
|---------|----------------|
| CH0635C | Length - 15 cm |

**Micro Spatula Scoop - Stainless Steel**

- 15 cm long with one concave scoop end and one flat spatula end (30 x 8 mm)
- Handle diameter 2 mm
- Made of stainless steel.

| Code    | Specs            |
|---------|------------------|
| CH0635D | 6 Length - 15 cm |



### Semi Micro Arrowhead Spatula - Stainless Steel Dual End

- 20 cm long with one flat blunt end and one flat tapered end.
- Handle diameter 4 mm
- Made of stainless steel.

| Code                   | Specs          |
|------------------------|----------------|
| CH0635E <span>1</span> | Length - 20 cm |

### Chattaway Spatula

- Traditional Chattaway spatula with one flat blunt end
- One rounded end with bent tip.
- Made of polished stainless steel.

| Code                   | Specs          |
|------------------------|----------------|
| CH0636A                | Length - 15 cm |
| CH0636B                | Length - 20 cm |
| CH0636C <span>2</span> | Length - 10 cm |

### Trowel Scoop - Stainless Steel, Dual End

- 17.5 cm long with two trowel shaped ends.
- Larger is 22 mm wide and the smaller is 11 mm wide.
- Made of stainless steel.

| Code                  | Specs            |
|-----------------------|------------------|
| CH0639 <span>3</span> | Length - 17.5 cm |

### Palette Knife Spatula

- Flat, flexible metal blade with wooden handle

| Code                  | Specs                         |
|-----------------------|-------------------------------|
| CH0640 <span>4</span> | Palette Knife Spatula - 15 cm |
| CH0640B               | Palette Knife Spatula - 10 cm |
| CH0640C               | Palette Knife Spatula - 20 cm |

### Spatula & Scoops - Teflon Coated

- The Teflon coating is non-stick and resistant to most chemicals.
- Items sold separately.

| Code                       | Specs                                               |
|----------------------------|-----------------------------------------------------|
| CH0632A-TFC <span>5</span> | Spatula spoon, wide handle - 15 cm long             |
| CH0634-TFC <span>6</span>  | Spatula, wide handle, dual end - 15 cm long         |
| CH0635-TFC <span>7</span>  | Micro spatula spoon - 23 cm long                    |
| CH0635A-TFC <span>8</span> | Scoop                                               |
| CH0635E-TFC <span>9</span> | Semi micro arrowhead spatula, dual end - 20 cm long |
| CH0636A-TFC                | Chattaway spatula - 15 cm long                      |
| CH0636B-TFC                | Chattaway spatula - 20 cm long                      |
| CH0637-TFC <span>10</span> | Nuffield spatula                                    |



## EISCO PREMIUM SPATULAS & SPOONS

- Research grade quality for industrial and academic lab.
- Made of high quality stainless steel 304 grade.
- Rust and corrosion free.
- Perfect for every industry.



1

### Premium Double Sided Spoons - Wide Handle, Stainless Steel

- Wide handle, 15 cm long with a spoon on both ends.
- Made of grade SS-304 stainless steel.

| Code          | Length |
|---------------|--------|
| CH0632F-SUP 1 | 15 cm  |



2

### Premium Double Sided Spatula - Wide Handle, Stainless Steel

- Wide handle with flat, blunt blades on both ends.
- Made of polished grade 304 stainless steel.

| Code          | Length |
|---------------|--------|
| CH0634A-SUP 2 | 12 cm  |
| CH0634B-SUP 3 | 15 cm  |



3

### Premium Nuffield Spatula - Wide Handle, Stainless Steel

- Wide, raised handle
- 14 cm long with flat rounded spatula at one end.
- A short angled rounded spatula at the other end.
- Made of polished grade 304 stainless steel.

| Code         | Length |
|--------------|--------|
| CH0637-SUP 4 | 14 cm  |




4

### Premium Double Sided Micro Spatulas - Stainless Steel

- Double sided flat spatulas with one rounded and one blunt end.
- Made of grade 304 stainless steel.

| Code          | Length x Dia. Handle |
|---------------|----------------------|
| CH0635H-SUP 5 | 13 x 0.2 cm          |
| CH0635I-SUP 6 | 10 x 0.3 cm          |
| CH0635J-SUP 7 | 15 x 0.3 cm          |
| CH0635K-SUP 8 | 15 x 0.4 cm          |




5



6



7



8

**Premium Double Sided Micro Spatula - Curved, Stainless Steel**


- Double sided spatula with round curved blade at both ends.
- Made of grade 304 stainless steel.

| Code          | Length x Dia. Handle |
|---------------|----------------------|
| CH0635O-SUP 1 | 13 x 0.5 cm          |
| CH0635P-SUP 2 | 21 x 0.5 cm          |

**Premium Double Sided Micro Spatula Scoops - Stainless Steel**

- Double sided spatula with flat blunt blade at one end
- Rounded scoop at the other end.
- Made of grade SS-304 stainless steel.

| Code          | Length x Dia. Handle |
|---------------|----------------------|
| CH0635Q-SUP 3 | 13 x 0.4 cm          |
| CH0635R-SUP 4 | 150 x 0.4 cm         |
| CH0635S-SUP 5 | 18.5 x 0.4 cm        |

**Premium Combustion Spoon - Stainless Steel**

- Designed for use in test tubes.
- Made of grade 304 stainless steel.

| Code         | Length |
|--------------|--------|
| CH0521-SUP 6 | 21 cm  |

**Premium Micro Spoon Spatula - Stainless Steel**

- Double sided spatula spoon with flat blunt blade at one end
- Micro spoon at the other end.
- Made of grade 304 stainless steel.

| Code          | Length |
|---------------|--------|
| CH0635T-SUP 7 | 10 cm  |

**Premium Fusion Spoon - Round Bottom, Stainless Steel with Wooden Handle**

- Round bottom fusion spoon with wooden handle & 100 ml capacity.
- Made of grade 304 stainless steel.

| Code          | Length |
|---------------|--------|
| CH0736E-SUP 8 | 7.5 cm |

**Premium Sodium Scoop - Brass with Wooden Handle**

- Sodium spoon with 2 cm diameter brass basket.
- Wooden handle.

| Code          | Specs          |
|---------------|----------------|
| CH0737B-SUP 9 | Length - 19 cm |





### Laboratory Jacks - Stainless Steel

- Made with stainless steel plates.
- Slip proof rubber feet and a screw operated mechanism that allows smooth, easy height adjustment.

| Code     | Dimensions          | Max. Open Height | Top Plate    | Weight   |
|----------|---------------------|------------------|--------------|----------|
| SSJCT-10 | 14.5 x 7.5 x 4.5 cm | 15 cm            | 8 x 7.3 cm   | 0.58 kg  |
| SSJCT-20 | 23.5 x 13 x 6.5 cm  | 28 cm            | 16.5 x 13 cm | 1.86 kg  |
| SSJCT-30 | 25 x 20 x 7 cm      | 28 cm            | 20 x 20 cm   | 2.58 kg  |
| SSJCT-40 | 25 x 30 x 6 cm      | 33 cm            | 25 x 25 cm   | 6.60 kg  |
| SSJCT-50 | 35 x 30 x 8 cm      | 48 cm            | 30 x 30 cm   | 6.60 kg  |
| SSJCT-60 | 45 x 30 x 8 cm      | 50 cm            | 40 x 40 cm   | 10.00 kg |



### Laboratory Jacks - Standard

- Light weight steel construction with a screw operated mechanism that allows smooth, easy height adjustment.
- Minimum height 60 mm & maximum height 25 cm.
- Maximum load 25 kg.

| Code    |   | Specs                   |
|---------|---|-------------------------|
| CH0642A | 1 | Plate size 15 x 13.5 cm |
| CH0642B | 2 | Plate size 20 x 20 cm   |



### Burette Clamps

- Spring loaded jaws made from high strength zinc die-cast alloy with non-toxic powdered coating.
- Suitable for all burettes up to 100 ml
- Bosshead accommodates rods up to 16 mm dia.


| Code    |   | Specs                  |
|---------|---|------------------------|
| CH0258A | 3 | Burette Clamp - Single |
| CH0258B | 4 | Burette Clamp - Double |




### Burette Clamps - Die Casted - Economy

- 3 prong vinyl coated clamp with in built bosshead suitable for all burettes up to 100 ml
- Bosshead accommodates rods up to 16 mm dia.

| Code    |   | Specs                  |
|---------|---|------------------------|
| CH0256A | 5 | Burette Clamp - Single |
| CH0256B | 6 | Burette Clamp - Double |



## Burner Stands



### Stands for Alcohol Burner

- Chrome plated steel wire gauze burner.
- Designed to accommodate most standard spirit lamps, not included.

| Code    | Specs                                                |
|---------|------------------------------------------------------|
| CH0643A | 1 Stand for Alcohol Burner - Height 90 mm            |
| CH0643B | 2 Stand for Alcohol Burner - Two way - Height 125 mm |



### Stands for Alcohol Burner - Stainless Steel

- Stainless steel legs and iron mesh gauze with ceramic center.
- Designed to accommodate most standard spirit lamps, not included.

| Code    | Specs                                                 |
|---------|-------------------------------------------------------|
| CH0643C | 3 Wire gauze burner stand with ceramic center - 9 cm  |
| CH0643D | Wire gauze burner stand with ceramic center - 12.5 cm |
| CH0645E | Wire gauze burner stand with ceramic center - 11 cm   |



### Tripod Stand - Circular

- Circular tripod stand made of black sheet metal.
- For use with spirit lamps, not included.

| Code   | I.D.     | Height |
|--------|----------|--------|
| CH0643 | 5 9.5 cm | 11 cm  |



### Tripod Stands - Circular

- Circular tripod stand with zinc plated cast iron top.
- Mild steel legs.

| Code    | I.D.    | Height |
|---------|---------|--------|
| CH0644A | 5 12 cm | 21 cm  |
| CH0644B | 6 15 cm | 21 cm  |

**Tripod Stand - Circular with Removable Legs**

- Cast iron circular tripod stand, with heat resistant painted top.
- Removable legs for easy storage.
- For use with spirit lamps, not included.

| Code         | I.D.    | Height |
|--------------|---------|--------|
| CH0644-SPL 1 | 10.5 cm | 21 cm  |

**Tripod Stand - Triangular with Removable Legs**

- Cast iron triangular tripod stand, with heat resistant painted top.
- Removable legs for easy storage.
- For use with spirit lamps, not included.

| Code          | Length of side | Height |
|---------------|----------------|--------|
| CH0645A-SPL 2 | 12 cm          | 21 cm  |
| CH0645A       | 12 cm          | 21 cm  |
| CH0645B       | 12 cm          | 21 cm  |

**Tripod Stand - Triangular**

- Triangular tripod stands are 20 cm tall.
- Length of each side is 12.5 cm.
- For use with spirit lamps, not included.

| Code      | Specs                          |
|-----------|--------------------------------|
| CH0646A 3 | Tripod Stand - Stainless Steel |
| CH0646B 4 | Tripod Stand - Iron            |



## Retort Base Stands

**Retort Stand Bases - Triangular** ①

- Black cast iron base with three points for stability.
- Tapered hole with 10 x 1.5 mm thread.

| Code    | Specs       | Radius | Weight |
|---------|-------------|--------|--------|
| CH0648A | Medium      | 100 mm | 800 g  |
| CH0648B | Large       | 130 mm | 1000 g |
| CH0648C | Extra Large | 160 mm | 1500 g |

**Deluxe Retort Stand Bases - Sheet Metal** ②

- Steel base, tapered hole with 10 x 1.5 mm thread in center of one short side.
- Finished with chemical resistant blue paint & fitted with rubber feet.

| Code    | Length | Width   | Weight |
|---------|--------|---------|--------|
| CH0649A | 16 cm  | 10 cm   | 750 g  |
| CH0649B | 20 cm  | 12.5 cm | 1200 g |
| CH0649C | 25 cm  | 16 cm   | 1800 g |

**Premium Retort Stand Bases - Cast Iron** ③

- Heavy duty cast iron base designed with expanded feet for excellent stability.
- Tapered hole with 10 x 1.5 mm thread in center of one short side.
- Finished with chemical resistant blue paint.

| Code    | Length  | Width   | Weight |
|---------|---------|---------|--------|
| CH0654F | 16 cm   | 10 cm   | 900 g  |
| CH0654G | 21 cm   | 13.5 cm | 1600 g |
| CH0654H | 23.5 cm | 16 cm   | 2400 g |
| CH0654I | 28.5 cm | 16 cm   | 2900 g |

**'A' Shape Retort Stand Bases - Cast Iron with Leveling Screws**

- 'A' shape cast iron base, tapered hole with 10 x 1.5 mm thread in center.
- Finished with chemical resistant paint and fit with leveling screws.

| Code    | Side  | Weight |
|---------|-------|--------|
| CH0655A | 25 cm | 3.5 kg |
| CH0655B | 25 cm | 4.5 kg |
| CH0655C | 25 cm | 3.8 kg |

**'A' Shape Retort Stand Bases - Cast Iron Nesting**

- 'A' shape cast iron base, tapered hole with 10 x 1.5 mm thread in center.
- Finished with chemical resistant powder coat.
- Bases can nest in one another for compact storage.

| Code    | Side  | Weight  |
|---------|-------|---------|
| CH0655E | 18 cm | 0.95 kg |
| CH0655F | 22 cm | 1.50 kg |
| CH0655G | 28 cm | 2.50 kg |

**'H' Shape Retort Stand Base - Metal**

- 'H' shape metal retort stand base for use with rods 8-12 mm dia (not included).
- Finished in chemical resistant blue paint and fit with plastic ends.

| Code   | Length | Width | Weight |
|--------|--------|-------|--------|
| CH0656 | 20 cm  | 20 cm | 3 kg   |

**Burette Support Stand - Cast Iron**

- Heavy duty cast iron base designed for use with a double burette holder.
- Tapered hole with 10 x 1.5 mm thread in center of one long side.
- Finished with chemical resistant blue paint.

| Code    | Length | Width | Weight |
|---------|--------|-------|--------|
| CH0656A | 13 cm  | 7 cm  | 3 kg   |

**Steel Rods 3**

- A range of steel rods with threading for use with Eisco stand bases.

| Code    | Length | Dia.  |
|---------|--------|-------|
| CH0657A | 50 cm  | 12 mm |
| CH0657B | 60 cm  | 12 mm |
| CH0657C | 75 cm  | 12 mm |
| CH0657D | 100 cm | 12 mm |

**Aluminium Rods 4**

- A range of aluminium rods with threading for use with Eisco stand bases.

| Code    | Length | Dia.  |
|---------|--------|-------|
| CH0658A | 50 cm  | 12 mm |
| CH0658B | 60 cm  | 12 mm |
| CH0658C | 75 cm  | 12 mm |
| CH0658D | 100 cm | 12 mm |

**Stainless Steel Rods 5**

- A range of stainless steel rods with threading for use with Eisco stand bases.
- Made of grade 202 stainless steel.

| Code    | Length | Dia.  |
|---------|--------|-------|
| CH0659A | 50 cm  | 12 mm |
| CH0659B | 60 cm  | 12 mm |
| CH0659C | 75 cm  | 12 mm |
| CH0659D | 100 cm | 12 mm |



## Clamps Holders and Bossheads



### Multipurpose Clamp Holder

- Square section aluminium alloy clamp holder with metal clamping screw.
- Two rectangular slots can accept square or circular section rods and tubing up to 16 mm diameter.
- Can also attach to table edges and boards up to 16 mm thick.

| Code      | Specs                     |
|-----------|---------------------------|
| CH0663D 1 | Multipurpose Clamp Holder |



### Lattice Clamps

- For lattice building and rod coupling with 12.5 mm dia. rods.
- Fitted with screws for Allen wrench, Pack of 12.

| Code      | Specs                       |
|-----------|-----------------------------|
| CH0664C 2 | Lattice Clamps - Pack of 12 |



### Universal Clamp Holder - Jumbo

- This clamp holder accepts rods up to 20 mm in diameter or flat boards.
- Tables up to 18 mm wide.

| Code       | Specs                          |
|------------|--------------------------------|
| CH0664JM 3 | Universal Clamp Holder - Jumbo |



### Double Bossheads - Large

- Die cast double bosshead, accepts rods 5 - 20 mm in diameter.
- Fitted with plastic cap adjustment screws.

| Code      | Specs                    |
|-----------|--------------------------|
| CH0665A 4 | Bosshead - Chrome Plated |
| CH0665B 5 | Bosshead - Black Coated  |



### Double Bosshead - Swivel

- Zinc die-cast double swivel bosshead.
- Accepts rods up to 16 mm in diameter.
- Allows two rods to be fitted at any angle.

| Code     | Specs                  |
|----------|------------------------|
| CH0666 6 | Bosshead Double Swivel |

### Suspension Clamps

- Use to suspend a thermometer over a hot plate.
- To suspend specific heat specimens.
- Clamp holds items at a distance of 110 mm from the support rod.
- The hook is adjustable and locks at any angle.


| Code                | Specs             |
|---------------------|-------------------|
| CH0667 <sup>1</sup> | Suspension Clamps |



### Thermometer Clamp with Bosshead

- Thermometer clamp with thumbscrews and 360° jaw rotation.
- Accommodates support rods up to 16 mm in diameter.
- Thermometers 6 to 12 mm in diameter.

| Code                 | Specs                           |
|----------------------|---------------------------------|
| CH0667A <sup>2</sup> | Thermometer Clamp with Bosshead |



### Burette Clamp with Bosshead - PVC Coated Round Jaws

- Clamp features round PVC dipped steel jaws with a 15 mm to 45 mm opening to hold burettes and test tubes.
- Accepts rods up to 13 mm in diameter
- Threaded rod and wing nut adjust clamp angle.
- Overall length from 13 to 15 cm.

| Code                 | Specs                                               |
|----------------------|-----------------------------------------------------|
| CH0688A <sup>3</sup> | Burette Clamp with Bosshead - PVC Coated Round Jaws |



### Universal Clamp with Bosshead - PVC Coated Jaws

- Clamp features 1 flat angled and 1 'V' steel prongs dipped in PVC ensuring slip resistant grip of cylindrical objects.
- Accepts rods up to 13 mm in diameter.
- Threaded rod and wing nut adjust clamp angle.
- Overall length from 13 to 15 cm.

| Code                 | Specs                                           |
|----------------------|-------------------------------------------------|
| CH0688B <sup>4</sup> | Universal Clamp with Bosshead - PVC Coated Jaws |



## Support Rings | Table Clamp



1

**Closed Support Rings** 1

- Closed ring made of plated mild steel.

| Code    | I.D. of Ring | Length of Stem |
|---------|--------------|----------------|
| CH0668A | 5 cm         | 14 cm          |
| CH0668B | 7 cm         | 14 cm          |
| CH0668C | 10 cm        | 14 cm          |
| CH0668D | 15 cm        | 14 cm          |



2

**Closed Support Rings with Bossheads** 2

- Stainless steel ring with integrated stem.
- Bosshead fitted with thumb screw.
- Accepts rods from 8 mm to 12.5 mm in diameter.
- 14 cm from center of support to center of ring.

| Code       | Inner Dia. |
|------------|------------|
| CH0672A/SS | 5 cm       |
| CH0672B/SS | 7.5 cm     |
| CH0672C/SS | 10 cm      |



3

**Closed Support Rings with Bossheads - Zinc Plated** 3

- Zinc plated cast iron ring with integrated stem.
- Distance is 12 cm from center of support to center of the ring.
- Fitted with thumb screws.
- Accepts rods from 8 mm to 12.5 mm in diameter.

| Code    | Inner Dia. |
|---------|------------|
| CH0674A | 5 cm       |
| CH0674B | 7.5 cm     |
| CH0674C | 10 cm      |
| CH0674D | 12.5 cm    |
| CH0674E | 15 cm      |



4

**Table Clamp - Premium**

- Versatile multipurpose clamp for securing a wide variety of laboratory apparatus and accessories.
- Made from heavy duty die-cast metal.
- Non-ferrous metal construction that will not bend or break under heavy loads.
- Features a vinyl sleeve and swivel pad to ensure a strong, slip proof grip on a variety of surfaces.

| Code      | Specs         |
|-----------|---------------|
| CH0688FBG | 4 Table Clamp |

### Retort Stand Sets - Value

- Includes a pressed steel base with 10 x 1.5 mm threaded hole in the center of one short side.
- A steel rod with 10 x 1.5 mm thread.

| Code       | Base         | Rod length | Rod Dia. |
|------------|--------------|------------|----------|
| STD-RD ①   | 12.5 x 20 cm | 60 cm      | 12 mm    |
| STD-RDVL02 | 10 x 16 cm   | 50 cm      | 12 mm    |
| STD-RDVL03 | 16 x 25 cm   | 75 cm      | 12 mm    |

### Retort Stand Sets - Superior

- Includes a steel base with rubber feet and 10 x 1.5 mm threaded hole in the center of one short side.
- A steel rod with 10 x 1.5 mm thread.

| Code       | Base         | Rod length | Rod Dia. |
|------------|--------------|------------|----------|
| STD-RD1 ②  | 12.5 x 20 cm | 60 cm      | 12 mm    |
| STD-RDSP02 | 10 x 16 cm   | 50 cm      | 12 mm    |
| STD-RDSP03 | 16 x 25 cm   | 75 cm      | 12 mm    |

### Retort Stand Sets - Cast Iron

- Includes a cast iron base, 10 x 1.5 mm threaded hole in the center of one short side.
- A steel rod with 10 x 1.5 mm thread.

| Code         | Base         | Rod length | Rod Dia. |
|--------------|--------------|------------|----------|
| CH0653ARD1 ② | 10 x 16 cm   | 50 cm      | 11 mm    |
| CH0653BRD2   | 12.5 x 20 cm | 50 cm      | 11 mm    |
| CH0653D1RD3  | 16 x 23 cm   | 60 cm      | 11 mm    |
| CH0653E1RD4  | 16 x 28 cm   | 90 cm      | 12.5 mm  |

### 'A' Shaped Retort Stand Sets

- These sets include an 'A' shape cast iron base.
- Steel rod.

This set comprises of the following:

| Code        | Item                                                           | Qty |
|-------------|----------------------------------------------------------------|-----|
| STDRDAB01   | 'A' base, side 18 cm, 0.95 kg(CH0655E)<br>rod 500mm (CH0657A)  | 1   |
| STDRDAB02   | 'A' base, side 22 cm, 1.50 kg (CH0655F)<br>rod 600mm (CH0657B) | 1   |
| STDRDAB03 ⑤ | 'A' base, side 28 cm, 2.50 kg (CH0655G)<br>rod 750mm (CH0657C) | 1   |



# METALWARE SETS

## Practical For All Laboratories

- Perfect for a wide variety of experiments and lab procedures
- Suitable for every lab
- Heavy base provide excellent stability
- Each set has its own unique purpose



### Metalware Set - Premium

#### MT-ST1

This set comprises of the following:

| Code    | Item                            | Qty |
|---------|---------------------------------|-----|
| CH0652B | Steel base, 20 x 12.5 cm, 800 g | 1   |
| CH0657B | Steel Rod, 60 cm                | 1   |
| CH0663A | Bosshead                        | 1   |
| CH0682  | Cork lined 4 prong clamp        | 1   |
| CH0672D | Support ring, 75 mm outer dia.  | 1   |



### Metalware Set - Deluxe

#### MT-ST5

This set comprises of the following:

| Code    | Item                                     | Qty |
|---------|------------------------------------------|-----|
| CH0649B | Steel base, 20 x 12.5 cm, 1200 g         | 1   |
| CH0657B | Steel Rod, 60 cm                         | 1   |
| CH0672D | Support ring, 75 mm outer dia.           | 1   |
| CH0688  | Vinyl coated 3 prong clamp with bosshead | 1   |

**Metalware Set - Superior****MT-ST6**

This set comprises of the following:

| Code    | Item                            | Qty |
|---------|---------------------------------|-----|
| CH0654B | Steel base, 25 x 16 cm, 2000 g  | 1   |
| CH0659B | Stainless steel rod, 60 cm      | 1   |
| CH0663D | Bosshead                        | 1   |
| CH0676A | Vinyl coated 3 prong clamp      | 1   |
| CH0674A | Support ring, 75 mm inner dia.  | 1   |
| CH0674C | Support ring, 100 mm inner dia. | 1   |

**Metalware Set****CH0652B-R8**

This set comprises of the following:

| Code    | Item                            | Qty |
|---------|---------------------------------|-----|
| CH0652B | Steel base, 20 x 12.5 cm, 800 g | 1   |
| CH0659A | Stainless steel rod, 50 cm.     | 1   |
| CH0674B | Support ring, 75 mm inner dia.  | 1   |
| CH0674C | Support ring, 100 mm inner dia. | 1   |
| CH0674D | Support ring, 125 mm inner dia. | 1   |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Chemistry Lab Starter Set

## MT-ST7

This set comprises of the following:

| Code       | Item                                     | Qty |
|------------|------------------------------------------|-----|
| CH0649C    | Steel base, 25 x 16 cm, 1800 g           | 1   |
| CH0659C    | Stainless steel rod, 75 cm               | 1   |
| CH0258B    | Double burette clamp                     | 1   |
| CH0674B    | Support ring 75 mm inner dia.            | 1   |
| CH0674C    | Support ring 100 mm inner dia.           | 1   |
| CH0675B    | Vinyl coated 2 prong clamp with bosshead | 1   |
| CH0676A    | Vinyl coated 3 prong clamp               | 1   |
| CH0665A    | Bosshead                                 | 1   |
| CH0645-SPL | Tripod stand                             | 1   |
| CH0993     | NextGen bunsen burner, Natural Gas       | 1   |
| CH0100E    | Rubber tubing                            | 1   |
| CH0525B    | Wire gauze with ceramic center, 150 mm   | 1   |



## Lab Essential Kit - Superior

## LAB-KIT

This set comprises of the following:

| Code       | Item                                   | Qty |
|------------|----------------------------------------|-----|
| CH0649B    | Steel retort base, 20 x 12.5 cm, 1200g | 1   |
| CH0657B    | Steel rod, 60 cm                       | 1   |
| BKTG       | Beaker tongs                           | 1   |
| CH0729B    | Crucible tongs                         | 1   |
| CH0644-SPL | Tripod stand                           | 1   |
| CH0525A    | Wire gauze with ceramic center         | 1   |
| CH0672D    | Support ring, 75 mm inner dia          | 1   |
| CH0258B    | Double burette clamp                   | 1   |
| CH0730     | Test tube holder                       | 1   |
| CH0713A    | Test tube rack, 'Z' shape              | 1   |
| CH0204G    | Test tube brush                        | 1   |
| CH0632B    | Wide handle spatula/spoon              | 1   |
| CH0663A    | Bosshead                               | 1   |
| CH0682     | Cork lined 4 finger clamp              | 1   |



**Laboratory Tool Set - Basic****Lab-TS**

This set comprises of the following:

| Code    | Item                  | Qty |
|---------|-----------------------|-----|
| CH0730  | Test tube holder      | 1   |
| CH0632B | Spatula               | 1   |
| CH0729B | Crucible tongs        | 1   |
| CH0290A | Mohr clip             | 1   |
| CH0204G | Test tube brush       | 1   |
| BI0168A | Forceps               | 1   |
|         | Acrylic scale, 150 mm | 1   |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS



Check out  
our new website  
[www.eiscoindustrial.com](http://www.eiscoindustrial.com)



## Ultimate Lab Starter Kit

### ULTLABKIT

Includes:

- Steel retort base with rubber feet, 20 x 12.5 cm
- Stainless steel rod, 60 cm
- Beaker tongs
- Crucible tongs
- Tripod stand
- Bosshead
- Test tube clamp
- Sheet metal test tube rack
- Cork lined test tube holder
- Test tube brush
- Stainless steel wide handle spatula/spoon
- 2 x Polypropylene graduated cylinders, 100 ml
- 2 x Polypropylene beakers, 500 ml
- Acrylic burette, 50 ml
- Erlenmeyer flask, glass, 250 ml
- Cleaning/drying basket 14 x 21 x 11 cm



## EISCO PREMIUM TONGS & HOLDERS

- Research grade quality for industrial and academic lab.
- Made of high quality stainless steel 304 grade.
- Rust and corrosion free.
- Perfect for research and industrial applications.

### Premium Beaker Tongs - Stainless Steel

- SS 304 grade beaker tongs.
- Length 26 cm long.

| Code          | Specs                                         |
|---------------|-----------------------------------------------|
| CH0738D-SUP 1 | Premium Beaker Tongs - Stainless Steel, 26 cm |



### Premium Mercury Tongs - Stainless Steel

- SS 304 grade mercury tongs for safe handling globules of mercury.
- Length 22.5 cm.

| Code          | Specs                                            |
|---------------|--------------------------------------------------|
| CH0736B-SUP 2 | Premium Mercury Tongs - Stainless Steel, 22.5 cm |



### Crucible Tongs - Straight, Stainless Steel

- SS 304 grade straight crucible tongs.
- Flat hinge.


| Code          | Specs |
|---------------|-------|
| CH0734E-SUP   | 15 cm |
| CH0734F-SUP 3 | 20 cm |



### Premium Crucible Tongs - Stainless Steel, with Bow

- SS 304 grade crucible tongs with bow.

| Code          | Specs |
|---------------|-------|
| CH0735E-SUP   | 15 cm |
| CH0735F-SUP 4 | 20 cm |





### Premium Flask Tongs - Stainless Steel, Angular

- SS 304 grade angular jaws provide an extremely firm grip for flasks.
- 22 cm length.

| Code                       | Specs                                        |
|----------------------------|----------------------------------------------|
| CH07336-SUP <span>1</span> | Premium Flask Tongs - Stainless Steel, 22 cm |



### Premium Flask Tongs - Stainless Steel

- SS 304 grade accepts flasks.
- 31 cm length.

| Code                       | Specs                                        |
|----------------------------|----------------------------------------------|
| CH0738C-SUP <span>2</span> | Premium Flask Tongs - Stainless Steel, 31 cm |



### Premium Test Tube Clamp - Stainless Steel

- SS 304 grade test tube clamp.
- Accepts test tubes & 13 cm length.

| Code                       | Specs                                            |
|----------------------------|--------------------------------------------------|
| CH0730C-SUP <span>3</span> | Premium Test Tube Clamp - Stainless Steel, 13 cm |



### Premium Test Tube Clamp - Stainless Steel

- SS 304 grade test tube holder.
- 15 cm length.

| Code                       | Specs                                            |
|----------------------------|--------------------------------------------------|
| CH0730D-SUP <span>4</span> | Premium Test Tube Clamp - Stainless Steel, 15 cm |



### Premium Test Tube Rack - Stainless Steel, 24 tubes

- SS 304 grade test tube rack.
- Holds 24 tubes (2 rows x 12 holes) up to 20 mm diameter.

| Code                       | Specs                                              |
|----------------------------|----------------------------------------------------|
| CH0713G-SUP <span>5</span> | Premium Test Tube Rack - Stainless Steel, 24 Tubes |

## Test Tube Holder - Vinyl coated

- Steel test tube holder with strong spring clip and wooden handle.
- Accepts 15-25 mm test tubes.

| Code                 | Specs                                               |
|----------------------|-----------------------------------------------------|
| CH0728A              | Test Tube Holder - Steel, (add length cm)           |
| CH0728B <sup>1</sup> | Test Tube Holder - Stainless Steel (add length, cm) |



## Crucible Tongs - General Use

- Nickel plated steel crucible tongs with bow.
- 20 cm long

| Code                 | Specs                        |
|----------------------|------------------------------|
| CH0729A <sup>2</sup> | Crucible Tongs - General Use |



## Beaker Tongs - Vinyl Coated, 3 Prong

- Special 3 prong tong for holding odd shaped objects.
- Jaws are PVC coated.
- Overall length 30 cm.

| Code                | Specs                                |
|---------------------|--------------------------------------|
| CH0732 <sup>3</sup> | Beaker Tongs - Vinyl Coated, 3 Prong |



## Flask Tongs - Cork Lined, Stainless Steel

- Cork lined stainless steel angular jaws provide an extremely firm grip for flasks from 100 - 2000 ml.
- 30 cm length.

| Code                | Specs                                            |
|---------------------|--------------------------------------------------|
| CH0740 <sup>4</sup> | Flask Tongs - Cork Lined, Stainless Steel, 30 cm |



## Flask Tongs - Cork Lined, Heavy Duty

- Cork lined & heavy duty angular jaws for miscellaneous applications.
- 25 cm long.

| Code                | Specs                   |
|---------------------|-------------------------|
| CH0733 <sup>5</sup> | Tong Flask - Heavy Duty |



## Test Tube Racks



### Test Tube Rack - Epoxy Coated Steel Wire, 20 Tubes

- Epoxy-coated, steel wire test tube rack.
- Resists most solvents weak acids base with stands temperatures up to 121°C.
- Holds 20 test tubes (4x5 rows) & with 18-20 mm diameter.

| Code      | Specs                                                |
|-----------|------------------------------------------------------|
| CH182015A | 1 Test Tube Rack - Epoxy Coated Steel Wire, 20 Tubes |



### Test Tube Rack - Epoxy Coated Steel Wire, 24 Tubes

- Epoxy-coated, steel wire test tube rack.
- Resists most solvents weak acids base with stands temperatures up to 121°C.
- Holds 24 test tubes (4x6 rows) & with 15-16 mm diameter.

| Code      | Specs                                                |
|-----------|------------------------------------------------------|
| CH182015B | 2 Test Tube Rack - Epoxy Coated Steel Wire, 24 Tubes |



### Test Tube Rack - Epoxy Coated Steel Wire, 40 Tubes

- Epoxy-coated, steel wire test tube rack.
- Resists most solvents weak acids base with stands temperatures up to 121°C.
- Holds 40 test tubes (4x10 rows) & with 18-20 mm diameter.

| Code      | Specs                                                |
|-----------|------------------------------------------------------|
| CH182015C | 3 Test Tube Rack - Epoxy Coated Steel Wire, 40 Tubes |



### Test Tube Rack - Epoxy Coated Steel Wire, 40 Tubes


- Epoxy-coated, steel wire test tube rack.
- Resists most solvents weak acids base with stands temperatures up to 121°C.
- Holds 40 test tubes (4x10 rows) & with 20-22 mm diameter.

| Code      | Specs                                                |
|-----------|------------------------------------------------------|
| CH182015D | 4 Test Tube Rack - Epoxy Coated Steel Wire, 40 Tubes |

**Test Tube Rack - Epoxy Coated Steel Wire, 36 Tubes**

- Epoxy-coated, steel wire test tube rack.
- Resists most solvents weak acids base with stands temperatures up to 121°C.
- Holds 36 test tubes (6x6 rows) & with 10-13 mm diameter.

| Code      |   | Specs                                              |
|-----------|---|----------------------------------------------------|
| CH182015E | 1 | Test Tube Rack - Epoxy Coated Steel Wire, 36 Tubes |




1

**Test Tube Rack - Epoxy Coated Steel Wire, 72 Tubes**

- Epoxy-coated, steel wire test tube rack.
- Resists most solvents weak acids base with stands temperatures up to 121°C.
- Holds 72 test tubes (6x12 rows) & with 15-16 mm diameter.

| Code      |   | Specs                                              |
|-----------|---|----------------------------------------------------|
| CH182015F | 2 | Test Tube Rack - Epoxy Coated Steel Wire, 72 Tubes |



2

**Test Tube Rack - Epoxy Coated Steel Wire, 80 Tubes**

- Epoxy-coated, steel wire test tube rack.
- Resists most solvents weak acids base with stands temperatures up to 121°C.
- Holds 80 test tubes (8x10 rows) & with 18-20 mm diameter.

| Code      |   | Specs                                              |
|-----------|---|----------------------------------------------------|
| CH182015G | 3 | Test Tube Rack - Epoxy Coated Steel Wire, 80 Tubes |




3

**Test Tube Rack - Epoxy Coated Steel Wire, 80 Tubes**

- Epoxy-coated, steel wire test tube rack.
- Resists most solvents weak acids base with stands temperatures up to 121°C.
- Holds 80 test tubes (8x10 rows) & with 22-25 mm diameter.

| Code      |   | Specs                                              |
|-----------|---|----------------------------------------------------|
| CH182015H | 4 | Test Tube Rack - Epoxy Coated Steel Wire, 80 Tubes |




4



# EISCO PREMIUM PLASTICWARE

- Manufactured with high quality polymers
- Great chemical and heat resistance, autoclavable at 121°C
- Class-A and Class-B certified options available



Crystal clear graduations using newly developed intellectual property, materials and moulds.



# Technical guide for use of different plasticware products

## Physical properties & chemical resistance of plastics

### Polypropylene, PP

- Translucent rigid polymer
- Temperature range -20 to +135°C
- Autoclavable at 121°C
- Good to excellent chemical resistance
- Resistant to fatigue, making it tough
- Typically used for beakers, bottles, cylinders, funnels, jugs, etc.



### Polytetrafluoroethylene, PTFE

- Opaque rigid polymer
- Wide temperature range -200 to +260°C
- Autoclavable at 121°C
- Unrivalled resistance to almost all chemicals
- Extremely low friction coefficient
- Typically used for EISCO stopcocks in burettes



### Low Density Polyethylene, LDPE

- Translucent flexible polymer
- Narrow temperature range of -50 to +80°C
- Not autoclavable at 121°C
- Good to excellent chemical resistance
- Robust and virtually unbreakable
- Typically used for wash bottles



### Polymethylmethacrylate, Acrylic (PMMA)

- Transparent rigid polymer
- Narrow temperature range -60 to +50°C
- Not autoclavable at 121°C
- Moderate chemical resistance
- Very tough and high clarity
- Typically used for radiation shields



### High Density Polyethylene, HDPE

- Translucent rigid polymer
- Broad temperature range of -100 to +120°C
- Not autoclavable at 121°C
- Good to excellent chemical resistance
- High tensile strength making it very tough
- Typically used for bottles



### Polystyrene, PS

- Transparent rigid polymer
- Narrow temperature range -40 to +90°C
- Not autoclavable at 121°C
- Moderate chemical resistance
- Brittle yet has excellent clarity
- Typically used for container ware



### Polymethylpentene, PMP (TPX)

- Transparent rigid polymer
- Broad temperature range -180 to +145°C
- Autoclavable at 121°C
- Good to excellent chemical resistance
- Has a low density and a high clarity
- Typically used for beakers and cylinders



### Polycarbonate, PC

- Transparent rigid polymer
- Broad temperature range -135 to +135°C
- Autoclavable at 121°C
- Moderate chemical resistance
- High impact strength
- Typically used for safety shields



**Private label and OEM options with bulk supplies available,  
please contact us for more details.**

## Beakers



### Beaker - Polypropylene, Low Form <sup>1</sup>

- Low form polypropylene beaker with spout.
- Moulded graduations & autoclavable.

| Code    | Capacity                                          | Graduations |
|---------|---------------------------------------------------|-------------|
| CH0137  | 25 ml                                             | 5 ml        |
| CH0137A | 50 ml                                             | 5 ml        |
| CH0137B | 100 ml                                            | 5 ml        |
| CH0137C | 250 ml                                            | 10 ml       |
| CH0137D | 500 ml                                            | 10 ml       |
| CH0137E | 1000 ml                                           | 20 ml       |
| CH0137F | 2000 ml                                           | 50 ml       |
| CH0137G | 5000 ml                                           | 200 ml      |
| CH0137H | Set of 5 - 50 ml, 100 ml, 250 ml, 500 ml, 1000 ml |             |



### Beaker - Polypropylene, Low Form, Printed <sup>2</sup>

- Low form polypropylene beaker with spout.
- Easy to read chemical resistant printed graduations & autoclavable.

| Code        | Capacity | Graduations |
|-------------|----------|-------------|
| ECH0139B/PR | 100 ml   | 10 ml       |
| ECH0139C/PR | 250 ml   | 25 ml       |
| ECH0139D/PR | 600 ml   | 50 ml       |
| ECH0139E/PR | 1000 ml  | 50 ml       |
| ECH0139F/PR | 2000 ml  | 100 ml      |
| ECH0139F/PR | 5000 ml  | 500 ml      |



### Beaker - TPX, Low Form, Printed <sup>3</sup>

- Low form TPX (PMP) beaker with spout.
- Easy to read chemical resistant printed graduations & autoclavable.

| Code       | Capacity | Graduations |
|------------|----------|-------------|
| CH0138B/PR | 100 ml   | 10 ml       |
| CH0138C/PR | 250 ml   | 25 ml       |
| CH0138D/PR | 600 ml   | 50 ml       |
| CH0138E/PR | 1000 ml  | 50 ml       |
| CH0138F/PR | 2000 ml  | 100 ml      |
| CH0138F/PR | 5000 ml  | 500 ml      |

### Reagent Bottles - Polypropylene

- Polypropylene reagent bottles with leak proof screw top.
- Narrow and wide neck options available & autoclavable.

| Narrow Neck | Wide Neck   | Capacity |
|-------------|-------------|----------|
| CH0170A-N   |             | 4 ml     |
| CH0170B-N   |             | 8 ml     |
| CH0170C-N   |             | 15 ml    |
| CH0170D-N   | CH0170D-W   | 30 ml    |
| CH0172A-N ① | CH0172A-W   | 60 ml    |
| CH0172B-N   | CH0172B-W   | 125 ml   |
| CH0172C-N   | CH0172C-W ② | 250 ml   |
| CH0172D-N   | CH0172D-W   | 500 ml   |
| CH0172E-N   | CH0172E-W ③ | 1000 ml  |



### Wash bottles - Premium, Narrow Neck

- LDPE wash bottles with leak proof narrow neck.
- Screw cap
- Non flexible delivery tube & translucent.

| Code       | Capacity |
|------------|----------|
| ECH0181A ⑤ | 125 ml   |
| ECH0181B ⑥ | 250 ml   |
| ECH0181C   | 500 ml   |
| ECH0181D   | 1000 ml  |



### Wash bottles - Premium, Narrow Neck, Colored Top

- LDPE wash bottles with leak proof narrow neck.
- Blue screw cap
- Non flexible delivery tube, translucent.

| Code       | Capacity |
|------------|----------|
| ECH0171A   | 125 ml   |
| ECH0171B   | 250 ml   |
| ECH0171C ⑦ | 500 ml   |
| ECH0171D   | 1000 ml  |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Bottles



## Wash bottles - Premium, Labeled, Wide Neck

- LDPE wash bottles with leak proof wide neck.
- Screw cap and non flexible delivery tube.
- Translucent & preprinted with chemical resistant ink.

| Code                  | Chemical            | Capacity |
|-----------------------|---------------------|----------|
| ECH0179A <sup>1</sup> | Acetone             | 500 ml   |
| ECH0179B              | Distilled Water     | 500 ml   |
| ECH0179C              | Isopropanol         | 500 ml   |
| ECH0179D <sup>2</sup> | Ethanol             | 500 ml   |
| ECH0179E <sup>3</sup> | Methanol            | 500 ml   |
| ECH0179F <sup>4</sup> | Sodium Hypochlorite | 500 ml   |
| ECH0179G              | Acetone             | 1000 ml  |
| ECH0179H              | Distilled Water     | 1000 ml  |
| ECH0179I              | Isopropanol         | 1000 ml  |
| ECH0179J              | Ethanol             | 1000 ml  |
| ECH0179K              | Methanol            | 1000 ml  |
| ECH0179L              | Sodium Hypochlorite | 1000 ml  |



## Wash bottles - Premium, Wide Neck, Colored Top

- LDPE wash bottles with leak proof wide neck.
- Blue screw cap, non flexible delivery tube, translucent.

| Code                  | Capacity |
|-----------------------|----------|
| ECH0173A <sup>5</sup> | 250 ml   |
| ECH0173B              | 500 ml   |
| ECH0173C              | 1000 ml  |



## Dropping Bottles - LDPE, Euro Design

- LDPE dropping bottles with screw cap.
- Attached cap, captive spout, translucent.

| Code                 | Capacity |
|----------------------|----------|
| CH0182A <sup>5</sup> | 15 ml    |
| CH0182B              | 30 ml    |
| CH0182C              | 60 ml    |
| CH0182D              | 125 ml   |
| CH0182E <sup>6</sup> | 250 ml   |

**Burettes - Acrylic 1**

- Transparent acrylic body & chemical resistant blue printed graduations.
- Leakproof PTFE key stopcock.

| Code    | Capacity                            | Sub. Div. | Tolerance $\pm$ |
|---------|-------------------------------------|-----------|-----------------|
| CH0233A | 25 ml                               | 0.10 ml   | 0.10 ml         |
| CH0233B | 50 ml                               | 0.10 ml   | 0.10 ml         |
| CH0233C | Spare stopcock for above (Pk. of 5) |           |                 |

**Burette Clamps - Polypropylene**

- Polypropylene burette clamp with a spring loaded jack in the middle to provide secure grip on burette.
- Two tapered wedges are also provided to ensure firm grip on rods.
- For use with 10 mm or 12.5 mm rods.

| Code      | Specs                  |
|-----------|------------------------|
| CH0262A 2 | Burette Clamp - Single |
| CH0262B 3 | Burette Clamp - Double |



Check out  
our new website  
[www.eiscoindustrial.com](http://www.eiscoindustrial.com)

### Conical Centrifuge Tubes - Polypropylene

- Polypropylene conical centrifuge tubes.
- Leak proof with rubber seal screw cap.
- Printed graduations & autoclavable.



| Code    | Capacity          | Type  | Height | O.D.  |
|---------|-------------------|-------|--------|-------|
| CH0270A | 15 ml             | Clear | 120 mm | 17 mm |
| CH0270B | 50 ml             | Clear | 116 mm | 29 mm |
| CH0270C | 15 ml             | Amber | 120 mm | 17 mm |
| CH0270D | 50 ml             | Amber | 116 mm | 29 mm |
| CH0270G | 15 ml- sterilised | Clear | 120 mm | 17 mm |
| CH0270H | 50 ml- sterilised | Clear | 116 mm | 29 mm |
| CH0270I | 15 ml- sterilised | Amber | 120 mm | 17 mm |
| CH0270J | 50 ml- sterilised | Amber | 116 mm | 29 mm |

### Conical Centrifuge Tubes - Polypropylene, Self Standing


- Polypropylene self standing conical centrifuge tubes.
- Leak proof with rubber seal screw cap.
- Printed graduations & autoclavable.



| Code    | Capacity           | Height | O.D.  |
|---------|--------------------|--------|-------|
| CH0271A | 50 ml              | 118 mm | 29 mm |
| CH0271B | 50 ml - Sterilized | 118 mm | 29 mm |

### Centrifuge Tube Storage Box with Lid

- Indexed grid, transparent lid.
- Temperature range is -90°C to 120°C & autoclavable.



| Code    | Capacity | Specs    |
|---------|----------|----------|
| CH0272A | 15 ml    | 36 Tubes |
| CH0272B | 50 ml    | 16 Tubes |

### Microcentrifuge Tubes

- Polypropylene & attached leak proof snap lid.
- Conical bottom & autoclavable.



| Code     | Capacity                      |
|----------|-------------------------------|
| CH0284A  | 0.5 ml capacity, Pack of 1000 |
| CH0284B  | 1.5 ml capacity, Pack of 500  |
| CH0284BB | 2 ml capacity, Pack of 500    |

### PCR Tubes

- Polypropylene & ultra thin walls.
- Attached leak proof snap lid.
- Conical bottom & autoclavable.

| Code                 | Specs                                 |
|----------------------|---------------------------------------|
| CH0289A <sup>1</sup> | 0.2 ml - Clear Flat Cap, Pack of 1000 |



### PCR Tubes

- Polypropylene, ultra thin walls.
- Strip of 8 connected tubes attached leak proof snap lid.
- Conicle bottom & autoclavable.

| Code                 | Specs                      |
|----------------------|----------------------------|
| CH0289B <sup>2</sup> | 0.2 ml, Pack of 120 strips |



### Floating Microcentrifuge Tube Rack

- Polypropylene rack with raised alphanumeric index on top.
- Can accomodate 16 microcentrifuge tubes & floats
- Central handle for easy lifting & autoclavable.

| Code                | Capacity | Holes |
|---------------------|----------|-------|
| CH0285 <sup>4</sup> | 1.5 ml   | 16    |



### Dual Sided Microcentrifuge Tube Rack with Lid

- Dual sided polypropylene rack with alpha-numeric index and transparent polycarbonate lid.
- One side can accomodate up to 108 1.5 ml tubes.
- The other side can accomodate up to 108 0.5 ml tubes & fully autoclavable.

| Code                | Specs                                    |
|---------------------|------------------------------------------|
| CH0286 <sup>5</sup> | 108 x 1.5 ml tubes or 108 x 0.5 ml tubes |



### PCR Tube Rack with Lid

- Single sided polypropylene rack with alpha-numeric index and transparent polycarbonate base and lid.
- Can accomodate up to 96 0.2 ml PCR tubes.
- Fully autoclavable & compatible with multipipettes.


| Code                | Capacity | Holes |
|---------------------|----------|-------|
| CH0287 <sup>6</sup> | 0.2 ml   | 96    |



### Microcentrifuge Tube Storage Box with Lid

- Polycarbonate box and lid with numbered grid.
- Temperature range -190°C to 121°C & autoclavable.

| Code                 | Capacity | Holes |
|----------------------|----------|-------|
| CH0288A              | 0.5 ml   | 100   |
| CH0288B <sup>7</sup> | 1.5 ml   | 64    |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Graduated Cylinders

# EISCO PREMIUM CYLINDERS

- Research grade quality for academic and industrial labs
- Made from premium quality polypropylene & TPX
- Clear, easy to read moulded graduations
- High clarity
- Certified class-A and class-B for accuracy
- Superb quality and durability
- Autoclavable



### Graduated Cylinders - Polypropylene, Class-B, Octagonal Base <sup>1</sup>

- Single piece polypropylene graduated cylinder with easy to read moulded graduations.
- Octagonal base for stability, class-B & autoclavable.

| Code       | Capacity | Graduation |
|------------|----------|------------|
| CH0354A/ES | 10 ml    | 0.2        |
| CH0354B/ES | 25 ml    | 0.5        |
| CH0354C/ES | 50 ml    | 1.0        |
| CH0354D/ES | 100 ml   | 1.0        |
| CH0354E/ES | 250 ml   | 2.0        |
| CH0354F/ES | 500 ml   | 5.0        |
| CH0354G/ES | 1000 ml  | 10.0       |



### Graduated Cylinders - TPX, Class-B, Octagonal Base <sup>2</sup>


- Single piece TPX graduated cylinder with easy to read moulded graduations.
- Octagonal base for stability, Class-B & autoclavable.

| Code    | Capacity | Graduation |
|---------|----------|------------|
| CH0354I | 10 ml    | 0.2 ml     |
| CH0354J | 25 ml    | 0.5 ml     |
| CH0354K | 50 ml    | 1.0 ml     |
| CH0354L | 100 ml   | 1.0 ml     |
| CH0354M | 250 ml   | 2.0 ml     |
| CH0354N | 500 ml   | 5.0 ml     |
| CH0354O | 1000 ml  | 10.0 ml    |

## Graduated Cylinders - TPX, Class-A, Octagonal Base 1

- Single piece TPX graduated cylinder with easy to read moulded graduations.
- Octagonal base for stability, Class-A & autoclavable.

| Code    | Capacity | Graduation |
|---------|----------|------------|
| CH0354P | 10 ml    | 0.2 ml     |
| CH0354Q | 25 ml    | 0.5 ml     |
| CH0354R | 50 ml    | 1.0 ml     |
| CH0354S | 100 ml   | 1.0 ml     |
| CH0354T | 250 ml   | 2.0 ml     |
| CH0354U | 500 ml   | 5.0 ml     |
| CH0354V | 1000 ml  | 10.0 ml    |



Come chat with us  
Live on [www.eiscoindustrial.com](http://www.eiscoindustrial.com)

## Graduated Cylinder Protectors 2

- Hexagonal ABS shape to prevent an overturned cylinder from rolling off the bench.

| Code    | Fits    | Pack Size |
|---------|---------|-----------|
| CH0352A | 10 ml   | Pk/10     |
| CH0352B | 25 ml   | Pk/10     |
| CH0352C | 50 ml   | Pk/10     |
| CH0352D | 100 ml  | Pk/10     |
| CH0352E | 250 ml  | Pk/10     |
| CH0352F | 500 ml  | Pk/5      |
| CH0352G | 1000 ml | Pk/5      |



## Flasks



1

### Conical Measures - Polypropylene, Hexagonal Base

- Polypropylene conical measures with spout.
- Moulded graduations and hexagonal base for stability
- Autoclavable.

| Code      | Capacity | Graduation |
|-----------|----------|------------|
| CH0355A   | 12 ml    | 1 ml       |
| CH0355B   | 25 ml    | 1 ml       |
| CH0355C   | 50 ml    | 2 ml       |
| CH0355D 1 | 125 ml   | 5 ml       |
| CH0355E   | 200 ml   | 5 ml       |



2

### Filtrating Flasks - Polypropylene

- Polypropylene filtering flasks with integrated side arm.
- Autoclavable.

| Code      | Capacity |
|-----------|----------|
| CH0422A 2 | 500 ml   |
| CH0422B   | 1000 ml  |



3

### Conical Flasks - Polypropylene

- Polypropylene conical flasks.
- Autoclavable.

| Code      | Capacity |
|-----------|----------|
| CH0431F 3 | 100 ml   |
| CH0431G   | 250 ml   |
| CH0431H   | 500 ml   |

### Measuring Jugs - Polypropylene

- Polypropylene measuring jugs with thumb grip handle.
- Spout, easy to read moulded graduations.
- Autoclavable.

| Code                 | Capacity | Graduation |
|----------------------|----------|------------|
| CH0356A <sup>1</sup> | 250 ml   | 10 ml      |
| CH0356B              | 500 ml   | 10 ml      |
| CH0356C              | 1000 ml  | 10 ml      |
| CH0356D <sup>2</sup> | 2000 ml  | 20 ml      |



### Measuring Jugs - Polypropylene, EURO Design

- Polypropylene short form measuring jugs with thumb grip handle.
- Long spout, easy to read moulded graduations.
- Autoclavable.

| Code                 | Capacity | Graduation |
|----------------------|----------|------------|
| CH0356E              | 500 ml   | 10 ml      |
| CH0356F <sup>3</sup> | 1000 ml  | 10 ml      |
| CH0356G <sup>4</sup> | 2000 ml  | 20 ml      |
| CH0356H              | 3000 ml  | 20 ml      |
| CH0356I              | 5000 ml  | 50 ml      |



### Measuring Jug - Polypropylene, Printed Graduations

- Polypropylene short form transparent measuring jugs with thumb grip handle.
- Spout & chemical resistant printed graduations.
- Autoclavable.

| Code                 | Capacity | Graduation |
|----------------------|----------|------------|
| CH0357A              | 100 ml   | 10 ml      |
| CH0357B              | 250 ml   | 25 ml      |
| CH0357C <sup>5</sup> | 600 ml   | 50 ml      |
| CH0357D <sup>6</sup> | 1000 ml  | 50 ml      |
| CH0357E              | 2000 ml  | 100 ml     |
| CH0357E              | 5000 ml  | 500 ml     |



### Measuring Jugs - TPX, Printed Graduations

- TPX short form transparent measuring jugs with thumb grip handle.
- Spout & chemical resistant printed graduations.
- Autoclavable.

| Code                 | Capacity | Graduation |
|----------------------|----------|------------|
| CH0357I <sup>7</sup> | 100 ml   | 10 ml      |
| CH0357J              | 250 ml   | 25 ml      |
| CH0357K <sup>8</sup> | 600 ml   | 50 ml      |
| CH0357L              | 1000 ml  | 50 ml      |
| CH0357M              | 2000 ml  | 100 ml     |
| CH0357M              | 5000 ml  | 500 ml     |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS



1

### Desiccators - Plain

- Desiccator with polypropylene base.
- Transparent polycarbonate flange without stopcock.
- Grease must be applied to the flange to obtain an airtight seal.

| Code      | Base I.D. | Flange O.D. |
|-----------|-----------|-------------|
| CH0364A 1 | 150 mm    | 171 mm      |
| CH0364B   | 200 mm    | 232 mm      |
| CH0364C   | 240 mm    | 276 mm      |
| CH0364D   | 300 mm    | 340 mm      |



2

### Desiccators - Vacuum

- Desiccator with polypropylene base.
- Transparent polycarbonate flange with an internal groove holding a silicon rubber O-ring.
- Polypropylene stopcock fitted with a self lubricating PTFE plug.
- The plug can be used for creating, shutting off and releasing the vacuum.
- Can hold a vacuum up to 740 mm Hg for 24 hrs. Autoclavable.

| Code      | Base I.D. | Flange O.D. |
|-----------|-----------|-------------|
| CH0366A   | 150 mm    | 171 mm      |
| CH0366B   | 200 mm    | 232 mm      |
| CH0366C   | 240 mm    | 276 mm      |
| CH0366D 2 | 300 mm    | 340 mm      |



3

### Desiccator - Vacuum, Transparent

- Desiccator with transparent polypropylene base and flange.
- Flange includes an internal groove holding a silicon rubber O-ring.
- Polypropylene stopcock fitted with a self lubricating PTFE plug.
- The plug can be used for creating, shutting off and releasing the vacuum.
- Can hold a vacuum up to 740 mm Hg for 24 hrs.
- Autoclavable.

| Code      | I.D.   | Flange O.D. |
|-----------|--------|-------------|
| CH0367A   | 150 mm | 171 mm      |
| CH0367B   | 200 mm | 232 mm      |
| CH0367C 3 | 240 mm | 276 mm      |
| CH0367D   | 300 mm | 340 mm      |

### Petri Dishes - Reusable, Non Sterile

- Clear, shatter resistant polypropylene petri dishes.
- Autoclavable.


| Code                 | Outer Dia. |
|----------------------|------------|
| CH0372A <sup>1</sup> | 50 mm      |
| CH0372B              | 75 mm      |
| CH0372C              | 100 mm     |
| CH0372D              | 125 mm     |
| CH0372E              | 150 mm     |



### Petri Dishes - Disposable, Sterile

- Clear, polystyrene petri dishes.
- Tripple vented & nontoxic.
- Individually packaged in polyethylene bags & radiation sterilized.

| Code                 | I.D.  | Flange O.D. |
|----------------------|-------|-------------|
| CH0374A              | 90 mm | 14 mm       |
| CH0374B <sup>2</sup> | 60 mm | 15 mm       |



### Petri Dish Racks

- Clear acrylic rack with white polycarbonate posts, can hold 60 mm and 90 mm petri dishes.
- Can be used during inoculation, incubation and storage.
- Easy plate accessibility and visibility enables the cultures to be seen and checked during incubation.

| Code                 | Specs                       |
|----------------------|-----------------------------|
| CH0376A <sup>4</sup> | Rack for 90 mm Petri Dishes |
| CH0376B              | Rack for 60 mm Petri Dishes |





### Volumetric Flasks - Polypropylene

- Polypropylene volumetric flasks with screw caps.
- Autoclavable

| Code    | Capacity |
|---------|----------|
| CH0452A | 25 ml    |
| CH0452B | 50 ml    |
| CH0452C | 100 ml   |
| CH0452D | 250 ml   |
| CH0452E | 500 ml   |
| CH0452F | 1000 ml  |



### Filter Funnels - Polyethylene

| Code    | Dia.   | Stem Length |
|---------|--------|-------------|
| CH0454A | 50 mm  | 50 mm       |
| CH0454B | 62 mm  | 62 mm       |
| CH0454C | 75 mm  | 76 mm       |
| CH0454D | 100 mm | 80 mm       |
| CH0454E | 150 mm | 110 mm      |



### Funnels - Industrial

- PP funnel with handle for easy handling.

| Code    | Capacity |
|---------|----------|
| CH0455A | 20 cm    |
| CH0455B | 25 cm    |

### Buchner Funnels - Polypropylene

- Two piece polypropylene Buchner funnel.
- Light weight and autoclavable.
- Top portion can be fitted on the lower portion forming a vacuum seal which does not requires greasing.
- Base of top portion has a perforated filter plate.

| Code                 | Dia.   | Stem Length |
|----------------------|--------|-------------|
| CH0464A              | 70 mm  | 60 mm       |
| CH0464B <sup>1</sup> | 110 mm | 85 mm       |
| CH0464C              | 90 mm  | 65 mm       |



### Funnel Holders - Polypropylene

- Can hold funnels with 3" - 6" dia. or 1" - 3" with the extra plate provided.
- Two tapered wedges are also provided to ensure firm grip on rods having 1/2"-3/8" dia.

| Code                | Specs                  |
|---------------------|------------------------|
| CH0467 <sup>2</sup> | Funnel Holder - Single |
| CH0467B             | Funnel Holder - Double |



### Separating Funnel - Polypropylene

- Polypropylene separating funnel with excellent clarity and strength.
- Leak proof threaded screw cap.
- Fitted polypropylene stopcock with self lubricating PTFE plug.
- Autoclavable.

| Code                 | Capacity |
|----------------------|----------|
| CH0480A              | 100 ml   |
| CH0480B <sup>3</sup> | 250 ml   |
| CH0480C              | 500 ml   |



### Separatory Funnel Holder - Polypropylene

- Polypropylene separatory funnel holder with a front opening that allows unobstructed placement of separating funnel and a clear view of its contents.
- Two tapered wedges are provided to ensure firm grip on rods having 1/2" or 3/8" dia.

| Code                | Specs                                    |
|---------------------|------------------------------------------|
| CH0468 <sup>4</sup> | Separatory Funnel Holder - Polypropylene |





### Weighing Scoops - Polypropylene <sup>1</sup>

- Polypropylene weighing scoops
- Autoclavable.
- Economical alternative to metal weighing scoops.

| Code    | Dia.    |
|---------|---------|
| CH0461A | 2 ml    |
| CH0461B | 5 ml    |
| CH0461C | 10 ml   |
| CH0461D | 25 ml   |
| CH0461E | 50 ml   |
| CH0461F | 100 ml  |
| CH0461G | 250 ml  |
| CH0461H | 500 ml  |
| CH0461I | 1000 ml |



### Pneumatic Troughs - Polypropylene

- Polypropylene pneumatic can be used for storage
- Aterilizing & drying labware.

| Code                 | Diameter | Height |
|----------------------|----------|--------|
| CH0505A              | 6 cm     | 18 cm  |
| CH0505B              | 10 cm    | 20 cm  |
| CH0505C <sup>2</sup> | 10 cm    | 25 cm  |



### Kipp's Gas Generators - Polypropylene

- Rapid and effective means of producing controllable amounts of common gases.
- Comprised of three separate pieces.

| Code                 | Capacity |
|----------------------|----------|
| CH0518A              | 250 ml   |
| CH0518B <sup>3</sup> | 500 ml   |
| CH0518C              | 1000 ml  |

### Pipette Stand - Vertical

- Vertical Pipette Stand with adjustable disks to accommodate a variety of pipette lengths.
- The top disk secures the pipette while the second disk has pilot holes where the pipette rests also allowing excess liquid to drain down to the bottom disk or base.

| Code   | Specs                      |
|--------|----------------------------|
| CH0594 | 1 Pipette Stand - Vertical |



Ask us about **OEM** and  
Private Label Opportunities

### Pipette Stand - Rotary

- Polypropylene pipette stand can hold up to 94 pipettes & rotates on a vertical axis.
- The lower rotating disc has four concentric troughs with tapered sides.
- Open bottoms to cradle the pipette tips while allowing excess water to drain.
- Stand can be easily assembled & disassembled.

| Code   | Specs                    |
|--------|--------------------------|
| CH0595 | 2 Pipette Stand - Rotary |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY


RUBBERWARE

CERAMICS

## Test Tube Racks



1



2


### Test Tube Racks - Polypropylene, Blue

- 6 place blue polypropylene test tube racks with drying pins.
- Autoclavable.

| Code    | Specs                              |
|---------|------------------------------------|
| CH0694A | 1 2 x 25 mm tubes, 4 x 16 mm tubes |
| CH0694B | 2 6 x 25 mm tubes                  |



Check out  
our new website  
[www.eiscoindustrial.com](http://www.eiscoindustrial.com)



3

### Test Tube Rack - Polypropylene

- 12 place polypropylene test tube rack.
- Accommodates up to 6 x 30 mm tubes and 6 x 15 mm tubes
- Autoclavable.

| Code   | Specs                            |
|--------|----------------------------------|
| CH0696 | 3 Test Tube Rack - Polypropylene |



4

### Test Tube Rack - Round, Polypropylene

- 12 place round polypropylene test tube rack.
- Can hold up to 4 x 25 mm tubes and 8 x 19 mm tubes.
- Includes 16 vertical pins for drying tubes and a trough to catch drippings.
- Autoclavable.

| Code   | Specs                                   |
|--------|-----------------------------------------|
| CH0697 | 4 Test Tube Rack - Round, Polypropylene |

### Test Tube Racks - 3-Tier, White

- White polypropylene test tube racks with 3-tiers to ensure secure tube fit.
- Autoclavable.

| Code                 | Specs            |
|----------------------|------------------|
| CH0708A              | 31 x 13 mm Tubes |
| CH0708B              | 31 x 16 mm Tubes |
| CH0708C              | 20 x 20 mm Tubes |
| CH0708D              | 12 x 25 mm Tubes |
| CH0708E <sup>1</sup> | 18 x 25 mm Tubes |
| CH0708F <sup>2</sup> | 12 x 32 mm Tubes |
| CH0708G              | 36 x 25 mm Tubes |
| CH0708H              | 18 x 13 mm Tubes |
| CH0708I              | 15 x 18 mm Tubes |
| CH0708J              | 18 x 12 mm Tubes |
| CH0708K              | 13 x 62 mm Tubes |
| CH0708L              | 16 x 62 mm Tubes |
| CH0708M              | 20 x 40 mm Tubes |



### Test Tube Racks

- Clear polypropylene test tube racks provide clear view of tube contents with 3-tiers to ensure secure tube fit.
- Autoclavable.

| Code                 | Specs            |
|----------------------|------------------|
| CH0709A              | 31 x 13 mm Tubes |
| CH0709B <sup>3</sup> | 16 x 13 mm Tubes |
| CH0709C              | 20 mm x 20 Tubes |
| CH0709D              | 25 mm x 12 Tubes |
| CH0709E              | 13 mm x 62 Tubes |
| CH0709F              | 16 mm x 62 Tubes |
| CH0709G <sup>4</sup> | 20 mm x 40 Tubes |
| CH0709H              | 25 mm x 24 Tubes |
| CH0709I              | 25 mm x 18 Tubes |
| CH0709J              | 25 mm x 36 Tubes |
| CH0709K              | 32 mm x 12 Tubes |



## Test Tube Racks



### Test Tube Racks - Polypropylene, Stackable

- White polypropylene test tube racks with 3-tiers provides secure tube fit, stackable when empty, fold up sides and end bar lock into top tier.
- Autoclavable.


| Code    | Specs            |
|---------|------------------|
| CH0710A | 90 x 13 mm Tubes |
| CH0710B | 60 x 16 mm Tubes |
| CH0710C | 40 x 20 mm Tubes |
| CH0710D | 24 x 25 mm Tubes |
| CH0710E | 21 x 30 mm Tubes |



### Test Tube Racks - Polypropylene, Moulded

- One piece moulded test tube racks with three tier grid design.
- Good chemical and thermal resistance.

| Code    | Specs            |
|---------|------------------|
| CH0710F | 36 x 13 mm Tubes |
| CH0710G | 36 x 16 mm Tubes |



### Test Tube Peg Racks

- Place tubes upright between pegs or inverted for drying purposes.
- Can also be used to hold electrophoresis and chromatography plates.

| Code     | Specs                      |
|----------|----------------------------|
| CH0711B1 | 96 x 13 mm tubes (upright) |
| CH0711B2 | 50 x 16 mm tubes (upright) |



### Test Tube Racks - Slant Racks

- Slant racks allow 5° or 20° slant when rack is set down on either side.
- Perfect for autoclaving tubes of agar while upright, then set on bench at desired angle to cool so Agar will solidify with a consistent slant.
- Also useful for incubating liquid cultures on a slant.

| Code    | Specs            |
|---------|------------------|
| CH0712A | 40 x 20 mm tubes |
| CH0712B | 36 x 25 mm tubes |



### Universal Multi-Tube Rack

- Polypropylene tube rack
- Each of the four sides accommodate different sized tubes: 30 mm, 20 mm, 17 mm, 12 mm.

| Code    | Specs                     |
|---------|---------------------------|
| CH0711A | Universal Multi-Tube Rack |

### Rotatable Multi Tube Rack

- Polypropylene multi rack is a combination of four modules that can be rotated in different directions to accommodate different types and sizes of tubes at the same time.
- Each module can accommodate 2 x 50 ml, 6 x 20 ml, 8 x 15 ml, or 10 x 12 mm tubes.
- Autoclavable.

| Code                   | Specs                     |
|------------------------|---------------------------|
| CH0712C <span>1</span> | Rotatable Multi Tube Rack |

### Microcentrifuge Tube Racks - Polycarbonate, Stackable

- Polycarbonate microcentrifuge tube rack for use with 1.5 ml microcentrifuge tubes.
- Stackable
- Alphanumeric index & Autoclavable.

| Code                   | Specs                                        |
|------------------------|----------------------------------------------|
| CH0711C <span>2</span> | Microcentrifuge Tube Rack, 24 x 1.5 ml tubes |
| CH0711D                | Microcentrifuge Tube Rack, 48 x 1.5 ml tubes |

### Microcentrifuge Tube Rack - Polypropylene, Folding

- Polypropylene microcentrifuge tube rack for use with 1.5 ml microcentrifuge tubes.
- Stackable
- Alphanumeric grid
- Folds flat when empty & autoclavable.

| Code                   | Specs                                         |
|------------------------|-----------------------------------------------|
| CH0711E <span>3</span> | Microcentrifuge Tube Rack, 100 x 1.5 ml tubes |

### Test Tube Baskets

- Polypropylene, break resistant and autoclavable for easy handling.
- Compact storage of test tubes and general labware.

| Code                   | Specs                      |
|------------------------|----------------------------|
| CH0707A <span>4</span> | 14 x 12 x 11 cm            |
| CH0707B                | 16 x 16 x 16 cm            |
| CH0707F <span>5</span> | 11 x 12 x 15 cm with cover |
| CH0707G                | 18 x 17 x 16 cm with cover |
| CH0707H                | 23 x 23 x 23 cm with cover |

### Draining Basket

- Polypropylene basket with wide net like structure on all sides perfect for draining labware after washing or autoclaving.
- Size 400 x 400 x 110 mm.

| Code                   | Specs                                |
|------------------------|--------------------------------------|
| CH0707C <span>6</span> | Draining Basket - 400 x 400 x 110 mm |



## Tubing Connectors | Stir Bars



### Tubing Connectors - Polypropylene, T & Y Shaped

- Polypropylene tubing connectors.
- Three way & ruffled end for rubber tubing.

| Code    | Specs               |
|---------|---------------------|
| CH0765A | T Connector - 8 mm  |
| CH0765B | Y Connector - 8 mm  |
| CH0765C | T Connector - 6 mm  |
| CH0765D | Y Connector - 6 mm  |
| CH0765E | T Connector - 10 mm |
| CH0765F | Y Connector - 10 mm |



### Tubing Connectors - Polypropylene, L Shaped

- Polypropylene tubing connectors.
- Two way & provide leak-proof connections.

| Code    | Specs               |
|---------|---------------------|
| CH0765H | L Connector - 6 mm  |
| CH0765I | L Connector - 8 mm  |
| CH0765J | L Connector - 10 mm |



### Tubing Connectors - Polypropylene, Cross

- Polypropylene 4-way cross tubing connectors.
- Provide leak-proof connections.

| Code     | Specs                   |
|----------|-------------------------|
| CH0765EC | Cross Connector - 6 mm  |
| CH0765FC | Cross Connector - 8 mm  |
| CH0765GC | Cross Connector - 10 mm |



### Magnetic Stir Bars

- PTFE coated Alnico magnets for mixing aqueous solutions.

| Code      | Size (Length x Dia.) | Shape       |
|-----------|----------------------|-------------|
| CH0803A01 | 12 x 6 mm            | Centre Ring |
| CH0803B   | 20 x 8 mm            | Centre Ring |
| CH0803C   | 25 x 8 mm            | Centre Ring |
| CH0803D01 | 30 x 8 mm            | Centre Ring |
| CH0803E01 | 35 x 9 mm            | Centre Ring |
| CH0803F   | 50 x 8 mm            | Centre Ring |

## Magnet Retriever

- PTFE coated magnet retriever
- Used to remove stirring bars from most vessels
- 30 cm length.

| Code    | Specs                |
|---------|----------------------|
| CH0803I | 1 Magnetic Retriever |

## Stoppers - Polypropylene

- Assorted size available.

| Code    | Stopper Size | Pack Size |
|---------|--------------|-----------|
| CH0873A | 10/19        | Pk/10     |
| CH0873B | 12/21        | Pk/10     |
| CH0873C | 14/23        | Pk/10     |
| CH0873D | 19/26        | Pk/10     |
| CH0873E | 24/29        | Pk/10     |
| CH0873F | 29/32        | Pk/5      |

## Joint Clips

- Chemical and temperature resistant plastic clips for use with jointed glassware reduces the need for clamping.
- Clips are color coded by joint size for easy identification.
- Recommended maximum temperature 80°C.

| Code    | Joint Size | Color  |
|---------|------------|--------|
| CH0869A | 10/19      | Mint   |
| CH0869B | 14/23      | Yellow |
| CH0869C | 19/26      | Blue   |
| CH0869D | 24/29      | Green  |
| CH0869E | 29/32      | Red    |
| CH0869F | 34/35      | Orange |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

# eisco

INDUSTRIAL

---

## WOODWARE

- Burette Stands
- Pipette Stands
- Test Tube Stands
- Separating Funnel Stands
- Funnel Stands



## Bottles Stands

- Polished wooden bottle stands with 4.5 cm drilled holes.

| Code                 | Specs           |
|----------------------|-----------------|
| CH0167               | 6 Bottle Stand  |
| CH0167B <sup>1</sup> | 12 Bottle Stand |



## Burette Stands

- Burette stands made of polished and seasoned hardwood.
- Includes clamping screws
- Rods are approximately 46 cm high with a 16 mm diameter.

| Code                 | Specs                  | Base Size    |
|----------------------|------------------------|--------------|
| CH0264A <sup>2</sup> | Burette Stand - Single | 200 x 100 mm |
| CH0264B <sup>3</sup> | Burette Stand - Double | 290 x 100 mm |



## Burette Rack - Vertical

- Natural wood burette rack with metal support rods
- Accepts maximum of 12 burettes up to 50 ml capacity and up to 2 cm diameter.
- Special design for easy racking and removing of burettes
- Wall mounting brackets included.

| Code                    | Specs                   |
|-------------------------|-------------------------|
| BUR-RCK 12 <sup>4</sup> | Burette Rack - 12 place |



## Pipette Rack - Horizontal

- Polished wood pipette rack hold up to 12 pipettes horizontally.
- Stand is 23 x 26 x 12 cm (HxWxD).

| Code                | Specs                  |
|---------------------|------------------------|
| CH0589 <sup>5</sup> | Pipette Rack, 12 Place |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Stands and Racks



### Pipette Stands - Verticle

- Polished wood pipette stands hold assorted size pipettes vertically.
- Stand is 41 x 36 x 11 cm (HxWxD).


| Code      | Specs                         |
|-----------|-------------------------------|
| CH0590A 1 | Pipette Stand for 24 Pipettes |
| CH0590B   | Pipette Stand for 10 Pipettes |



### Funnel Stands

- Polished wood funnel stands with clamping screws.
- Rods are approximately 46 cm tall with 16 mm diameter.
- Single base size 200 x 110 mm, double base size 290 x 110 mm.

| Code      | Specs                         |
|-----------|-------------------------------|
| CH0689A 2 | Single Funnel Stand           |
| CH0689B 3 | Double Funnel Stand           |
| CH0689C   | Funnel Stand for four funnels |



### Separating Funnel Stand

- Height adjustable polished wood separating funnel stands.
- Can accommodate up to 6 (place volume range for funnels here) separating funnels.

| Code      | Specs                   |
|-----------|-------------------------|
| CH0689E 4 | Separating Funnel Stand |



### Funnel Stand Clamp

- Hardwood rack with clamp holds up to 4 funnels in beveled holes that are 11 cm apart.
- Attached iron clamp accommodates support rods up to 16 mm.

| Code      | Specs              |
|-----------|--------------------|
| CH0689D 5 | Funnel Stand Clamp |

## Combination Funnel and Burette Stand

- Polished wood stand can accommodate 1 burette and 1 funnel of any size.
- Clamping screws included.
- Rod is approx.
- 46 tall with a 16 mm diameter.
- Base is 290 x 110 mm.

| Code   | Specs                                  |
|--------|----------------------------------------|
| CH0690 | 1 Combination Funnel and Burette Stand |



## Test Tube Racks - 4 Piece

- Wooden test tube racks clip together quickly
- Easily with two plastic supports.


| Code    | Specs                                         |
|---------|-----------------------------------------------|
| CH0702  | 2 Test Tube Rack - 2 x 30 mm, 4 x 20 mm Tubes |
| CH0702B | Test Tube Rack - 6 x 20 mm Tubes              |
| CH0702C | Test Tube Rack - 8 x 22 mm, 4 x 32 mm tubes   |



## Test Tube Racks

- Polished wood test tube racks with 22 mm diameter drilled holes.

| Code    | Specs                       |
|---------|-----------------------------|
| CH0704A | Test Tube Rack - 6 Tubes    |
| CH0704B | 3 Test Tube Rack - 12 Tubes |



## Test Tube Racks with Drying Pegs

- Polished wood test tube racks with vertical pins for drying.
- Drilled holes are 25 mm
- 20 mm in diameter.

| Code    | Specs                                         |
|---------|-----------------------------------------------|
| CH0706A | 4 Test Tube Rack - 2 x 25 mm, 4 x 20 mm tubes |
| CH0706B | Test Tube Rack - 2 x 25 mm, 10 x 20 mm tubes  |



## Stands and Racks



1



2



3



4



5



6

### Test Tube Rack - 7 Place

- Polished wood test tube racks with dovetail construction can accommodate up to 7 x 13 mm tubes in drilled holes.

| Code    | Specs                      |
|---------|----------------------------|
| CH0703A | 1 Test Tube Rack - 7 Tubes |

### Test Tube Rack - 12 Place

- Polished wood test tube racks with dovetail construction can accommodate up to 12 x 18 mm tubes in drilled holes.

| Code    | Specs                          |
|---------|--------------------------------|
| CH0703B | 2 Test Tube Support - 12 Place |

### Test Tube Rack with Drying Pegs - 12 Place

- Polished wood test tube racks with 12 drying pegs can accommodate up to 12 x 25 mm tubes in drilled holes.

| Code    | Specs                                        |
|---------|----------------------------------------------|
| CH0703C | 3 Test Tube Rack with Drying Pegs - 12 Tubes |

### Test Tube Rack with Drying Pegs - 6 Place

- Polished wood test tube racks with 6 drying pegs can accommodate up to 6 x 19 mm tubes in drilled holes.

| Code    | Specs                                       |
|---------|---------------------------------------------|
| CH0705A | 4 Test Tube Rack with Drying Pegs - 6 Tubes |

### Test Tube Rack with Drying Pegs - 12 Place

- Polished wood test tube racks with 12 drying pegs can accommodate up to 12 x 25 mm tubes in drilled holes.

| Code    | Specs                                        |
|---------|----------------------------------------------|
| CH0705B | 5 Test Tube Rack with Drying Pegs - 12 Tubes |

### Test Tube Rack with Drying Pegs - 13 Place, 2 Tier

- Polished wood test tube racks with 6 drying pegs can accommodate up to 7 x 19 mm tubes.
- 6 x 25 mm tubes in drilled holes.

| Code    | Specs                                        |
|---------|----------------------------------------------|
| CH0705C | 6 Test Tube Rack with Drying Pegs - 13 Tubes |

## Test Tube Rack - 15 Place, Solid Block

- Solid wood block with 15 slotted holes can accommodate up to 15 x 13 mm tubes.

| Code      | Specs                     |
|-----------|---------------------------|
| CH0705D 1 | Test Tube Rack - 15 Tubes |



## Test Tube Rack with Drying Pegs - 6 Place

- Polished wood test tube rack with 6 drying pegs can accommodate up to 6 x 24 mm tubes in drilled holes.

| Code          | Specs                                     |
|---------------|-------------------------------------------|
| CH0002C-SPL 2 | Test Tube Rack with Drying Pegs - 6 Tubes |



## Test Tube Rack - 12 Place

- Polished wood test tube rack can accommodate up to 12 x 22 mm tubes in drilled holes.

| Code      | Specs                     |
|-----------|---------------------------|
| CH0004C 3 | Test Tube Rack - 12 Tubes |



## Test Tube Rack - 12 Place

- Polished wood test tube rack can accommodate up to 12 x 28 mm tubes in drilled holes.

| Code      | Specs                     |
|-----------|---------------------------|
| CH0005C 4 | Test Tube Rack - 12 Tubes |



## Wall Mounted Drying Rack - 90 Place

Wooden drying rack can accommodate up to 90 pieces, wall mounting brackets included.

- 6 rows with 12 x 3" pegs
- 3 rows with 6 x 5" pegs

| Code      | Specs                               |
|-----------|-------------------------------------|
| CH0001C 5 | Wall Mounted Drying Rack - 90 Place |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

# eisco

INDUSTRIAL

---

## EQUIPMENT

- Hot Plates
- Stirrers
- Water Bath



### Hot Plates with Magnetic Stirrer - Ceramic Top, Analog Control

Hot plates with chemical resistant ceramic coated aluminium plate and powerful heating elements, and spinning magnetic drives for fast and uniform heating of up to 5,000 ml of liquids. Includes 2 analog controls, metal rod and clamp.

- Plate Material: Ceramic coated aluminium
- Plate Dimensions: 190 x 190 mm
- Heater Power: 650 W
- Max. Plate Temperature: 400 degrees C
- Stirrer Speed: 100- 2000 rpm
- Max Stirring Capacity: 5000 ml
- Control: Analogue
- Accessories: Rod and vinyl coated 2 prong clamp
- Dimensions: 33 x 19.5 x 13 cm
- Weight: 4 kg

| Code                  | Specs                                   |
|-----------------------|-----------------------------------------|
| MGST1-V1              | Hot Plate, Input 110/120 V AC, 60 Hz    |
| MGST1-V2 <sup>1</sup> | Hot Plate, Input 220/240 V AC, 50/60 Hz |



### Hot Plates - Ceramic Top, Analog Control

Hot plates with chemically resistant ceramic coated aluminum plate and powerful heating elements for fast heating of up to 5,000 ml of liquids. Includes 1 analog control, metal rod and clamp.

- Plate Material: Ceramic coated aluminum
- Plate Dimensions: 190 x 190 mm
- Heater Power: 650 W
- Max. Plate Temperature: 400 degrees C
- Max Stirring Capacity: 5000 ml
- Control: Analogue
- Accessories: Rod and vinyl coated 2 prong clamp
- Dimensions: 33 x 19.5 x 13 cm
- Weight: 4 kg

| Code                  | Specs                                        |
|-----------------------|----------------------------------------------|
| HTPLT-V1              | Hot Plate Only, Input 110/120 V AC, 60 Hz    |
| HTPLT-V2 <sup>2</sup> | Hot Plate Only, Input 220/240 V AC, 50/60 Hz |





1

### Hot Plates with Magnetic Stirrer - Ceramic Top, Digital Display

Hot plates with chemical resistant ceramic coated aluminium plate, powerful heating elements and spinning magnetic drives for fast and uniform heating of up to 5,000 ml of liquids. Includes back lit LCD display, metal rod and clamp. Heating temperatures can be preset.

- Plate Material: Ceramic coated aluminium
- Plate Dimensions: 190 x 190 mm
- Heater Power: 650 W
- Max. Plate Temperature: 400 degrees C
- Stirrer Speed: 100- 2000 rpm
- Max Stirring Capacity: 5000 ml
- Control: Analogue
- Accessories: Rod and vinyl coated 2 prong clamp
- Dimensions: 33 x 19.5 x 13 cm
- Weight: 4 kg

| Code     | Specs                                           |
|----------|-------------------------------------------------|
| MGST2-V1 | Hot Plate Stirrer, Input 110/120 V AC, 60 Hz    |
| MGST2-V2 | Hot Plate Stirrer, Input 220/240 V AC, 50/60 Hz |



2

### Hot Plates with Magnetic Stirrer - Aluminium Top, Analog Control

Hot plates with aluminum plate, independently controlled heating elements and spinning magnetic drives allow uniform heating of up to 500 ml of liquids. Includes 2 analog controls, metal rod and clamp

- Plate Material: Aluminium
- Plate Dimensions: 190 x 190 mm
- Heater Power: 180 W
- Max. Plate Temperature: 280 degrees C
- Stirrer Speed: 0 - 1400 rpm
- Max Stirring Capacity: 500 ml
- Control: Analogue
- Accessories: Rod and vinyl coated 2 prong clamp
- Dimensions: 33 x 19.5 x 13 cm
- Weight: 4 kg

| Code        | Specs                                           |
|-------------|-------------------------------------------------|
| HTMST100-V1 | Hot Plate Stirrer, Input 110/120 V AC, 60 Hz    |
| HTMST100-V2 | Hot Plate Stirrer, Input 220/240 V AC, 50/60 Hz |

## Heating Mantles

- Electric heating mantles with spun metal casing.
- Painted with corrosion resistant coating.
- Incorporates a heating element suspended within an insulating ceramic fiber for safe heat transfer.
- Mantles have built-in power controller.
- Power and heating indication lights.
- Maximum temperature 450°C, 220 V AC.

| Code               | Specs                         |
|--------------------|-------------------------------|
| HMN25              | Heating Mantle 250ml - 150 W  |
| HMN50 <sup>1</sup> | Heating Mantle 500ml - 250 W  |
| HMN100             | Heating Mantle 1000ml - 350 W |
| HMN200             | Heating Mantle 2000ml - 450 W |



1

Ask us about **OEM** and  
Private Label Opportunities

## Heating Mantles - Digital Display

- Electric heating mantles with spun metal casing.
- Painted with corrosion resistant coating
- Incorporates a heating element suspended within an insulating ceramic fiber for safe heat transfer.
- Mantles have built-in power controller.
- Power and heating indication lights.
- Digital temperature display.
- Temperature sensor PT100, metal rod, and clamp.
- Maximum temperature 450°C.

| Code                 | Specs                                              |
|----------------------|----------------------------------------------------|
| HMDT250 <sup>2</sup> | Heating Mantle with Digital Display 250ml - 150 W  |
| HMDT500              | Heating Mantle with Digital Display 500ml - 250 W  |
| HMDT1000             | Heating Mantle with Digital Display 1000ml - 350 W |
| HMDT2000             | Heating Mantle with Digital Display 2000ml - 450 W |



2

GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS



1



2

### Unstirred Digital Water Bath

- Seamless stainless steel tank housed in a durable powder coated and chemical resistant outer case.
- 100 W power
- 100°C maximum temperature, 220 V AC.

| Code    | Specs                                    |
|---------|------------------------------------------|
| WBTH-4H | 1 Unstirred Digital Water Bath - 4 holes |
| WBTH-6H | 2 Unstirred Digital Water Bath - 6 holes |

### Water Still with Stainless Steel Stand

Water still provides high quality, pyrogen-free distilled water with low conductivity. The unit has built-in safety features including gate valve, flow switch, thermostat and fuse. The distilled output is cooled with a high efficiency coiled condenser. Temperature of distillate is 25 - 40°C (tested at ambient room temperature of 30° C).

- Output 4 Ltr. / Hour
- PH 5.5-6.0
- Cooling Water Required 1-2 Ltr /min.
- Minimum Pressure 3 p.s.i.
- Net Weight 6 kg

| Code    | Specs                                    |
|---------|------------------------------------------|
| CH0814  | 3 Water still with Stainless Steel Stand |
| CH0814A | 4 Spare Heater for above                 |



3



4



## LAB SAFETY

- Gloves
- Safety Goggles
- Lab Coats
- Aprons



## Gloves



### Disposable Gloves - Polyethylene

- Economical & light weight.
- Provide good protection while giving a bare handed touch.
- Supplied in Pack of 100.

| Code    | Specs                           |
|---------|---------------------------------|
| CH0612A | 1 Gloves Disposable - Polythene |

### Disposable Gloves - Vinyl

- Strong & thin
- Firm fitting gloves are designed for a skin-tight fit on either hand.
- Supplied in Pack of 100.

| Code    | Specs    |
|---------|----------|
| CH0612B | Small    |
| CH0612C | 2 Medium |
| CH0612D | Large    |

### Disposable Gloves - Latex

- Made from 100% highly purified rubber latex.
- Beaded cuff for greater protection.
- Non-sterile.
- Supplied in Pack of 100.

| Code    | Specs    |
|---------|----------|
| CH0612E | Small    |
| CH0612F | 3 Medium |
| CH0612G | Large    |

### Disposable Gloves - Nitrile

- Contains no latex, hence no allergic issues.
- Textured finish on fingers for better grip
- Ideal choice for chemical lab and use during dissection.
- Supplied in Pack of 100.

| Code    | Specs   |
|---------|---------|
| CH0612H | Small   |
| CH0612I | Medium  |
| CH0612J | 4 Large |

### Rubber Gloves

- General purpose, light weight.
- Anti-slip pattern, pure cotton flock lined.
- Supplied as pair.

| Code    | Size          |
|---------|---------------|
| CH0613A | 7" - 7.5"     |
| CH0613B | 8" - 8.5"     |
| CH0613C | 5 10" - 10.5" |

### Rubber Gloves

- Bi-color & 24 mil/0.60 mm.
- Flock lined & natural rubber gloves.
- Supplied as pair.

| Code                 | Size        |
|----------------------|-------------|
| CH0616A              | 7" - 7.5"   |
| CH0616B              | 8" - 8.5"   |
| CH0616C <sup>1</sup> | 9" - 9.5"   |
| CH0616D              | 10" - 10.5" |



### Neoprene Gloves

- Anti-slip pattern for good grip in wet and dry application.
- Good chemical resistance.
- Suitable for chemical industry, industrial cleaning & oil refining.
- Supplied as pair.

| Code                 | Size |
|----------------------|------|
| CH0617A              | 6.5" |
| CH0617B              | 7.5" |
| CH0617C              | 8.5" |
| CH0617D <sup>2</sup> | 10"  |



### Cut Resistant Gloves

- Dyneema® and synthetic fiber blend with PU coating.
- Suitable for glass and ceramic industry, sheet metal operation
- Injection moulding.
- Supplied as pair.

*These are cut resistant gloves - NOT CUT PROOF.*

| Code                 | Size   |
|----------------------|--------|
| CH0617E              | Small  |
| CH0617F <sup>3</sup> | Medium |
| CH0617G              | Large  |



### Rubber Gloves - Heavy Duty

- Extra thick for better protection
- Chlorinated for enhanced chemical resistance.
- Supplied as pair.

| Code                 | Size |
|----------------------|------|
| CH0618A <sup>4</sup> | 10"  |



### Gloves Leather

- Gauntlet style gloves designed to provide protection for the hand.
- Forearm against intermittent exposure to heat.

| Code                | Specs          |
|---------------------|----------------|
| CH0614 <sup>5</sup> | Gloves Leather |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS


## Hand Protectors & Safety Goggles



### Hand Protectors

- Heat & cold resistant silicon rubber protector.
- Allows easy handling of heated vessels like flask, bottles & beakers etc.
- Studded surface is used for firm grip, tested to 250°C.
- Pockets for fingers & thumb, ambidextrous.

| Code     | Specs             |
|----------|-------------------|
| CH0616HP | 1 Hand Protectors |



### Safety Goggles - Clear

- These safety glasses are perfect for fitting over most prescription lenses.
- Side vents prevent glasses from fogging & UV protected.
- Light weight polycarbonate, impact resistant & CE marked.
- Conforms to ANSI Z87.1 2010 side shields protection.

| Code    | Specs                    |
|---------|--------------------------|
| CH0619A | 2 Safety Goggles - Clear |



### Safety Goggles - Clear With Adjustable Temples

- These safety goggles offer better coverage around the sides, adjustable temples for universal fit.
- Hard coat lens with UV protection, light weight polycarbonate & impact resistant lens.
- CE marked. Conforms to ANSI --- side shields protection.

| Code  | Specs                                            |
|-------|--------------------------------------------------|
| G2612 | 3 Safety Goggles - Clear With Adjustable Temples |



### Safety Goggles - Clear, Vented

- These safety glasses are perfect for fitting over most prescription lenses.
- Vents keep the goggles from fogging, black elastic strap to fit all sizes.
- Designed to protect eyes from chemical splash and dust.
- CE marked, conforms to ANSI Z87.1, EN166 B.
- Side shield protection.

| Code     | Specs                            |
|----------|----------------------------------|
| CH0619DV | 4 Safety Goggles - Clear, Vented |



### Safety Goggles - Basic

- These safety glasses are perfect for fitting over most prescription lenses.
- Vents keep the goggles from fogging, green elastic strap to fit all sizes.
- Designed to protect eyes from chemical splash and dust.

| Code | Specs                    |
|------|--------------------------|
| G103 | 5 Safety Goggles - Basic |



### Safety Goggles - Basic, Indirect Vent

- These safety glasses are perfect for fitting over most prescription lenses.
- Vents keep the goggles from fogging, elastic strap to fit all sizes.
- Designed to protect eyes from chemical splash and dust.

| Code | Specs                                   |
|------|-----------------------------------------|
| G104 | 6 Safety Goggles - Basic, Indirect Vent |

## Lab Coats - White Polyester ①

- Modern styling, white polyester/cotton drill.
- Concealed press stud fastenings, backvent.
- Three large pockets.

| Code    | Size                                      |
|---------|-------------------------------------------|
| CH0622A | Lab Coats - White Polyester - Medium      |
| CH0622B | Lab Coats - White Polyester - Large       |
| CH0622C | Lab Coats - White Polyester - Small       |
| CH0622D | Lab Coats - White Polyester - Extra Large |



1

## Lab Coats - Navy Blue ②

- Modern styling, white polyester / cotton drill.
- Concealed press stud fastenings & backvent.
- Three large pockets.


| Code       | Size                                |
|------------|-------------------------------------|
| CH0622A-NV | Lab Coats - Navy Blue - Medium      |
| CH0622B-NV | Lab Coats - Navy Blue - Large       |
| CH0622C-NV | Lab Coats - Navy Blue - Small       |
| CH0622D-NV | Lab Coats - Navy Blue - Extra Large |



2

### Size Chart

| Code                  | Length (cm) | Chest (cm) | Arm Length (cm) |
|-----------------------|-------------|------------|-----------------|
| CH0622A<br>CH0622A-NV | 100         | 60         | 46              |
| CH0622B<br>CH0622B-NV | 105         | 64         | 50              |
| CH0622C<br>CH0622C-NV | 95          | 58         | 42              |
| CH0622D<br>CH0622D-NV | 110         | 67         | 54              |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Aprons



1



2




3



4



5



### Rubberised Aprons - Black 1

- Suitable for general lab work, comfortable and flexible.
- Waterproof and resistant to most chemicals.

| Code    | Specs                             |
|---------|-----------------------------------|
| RBAP-CH | Rubberised Apron - Black - Small  |
| RBAP-MD | Rubberised Apron - Black - Medium |
| RBAP-AD | Rubberised Apron - Black - Large  |

### Rubber Apron - Colored

- Suitable for general lab work, comfortable and flexible.
- Waterproof and resistant to most chemicals used by students in lab.

| Code         | Specs                                 |
|--------------|---------------------------------------|
| RBAP-CH-RD 2 | Rubber Apron - Colored - Small - Red  |
| RBAP-AD-BL 3 | Rubber Apron - Colored - Large - Blue |

### Vinyl Aprons 4

- These easy to wear transparent vinyl aprons are suitable for labs.
- Kitchens or shops.
- Protects clothes from acids and fluids.

| Code    | Specs             |
|---------|-------------------|
| VNAP-CH | Vinyl Apron Small |
| VNAP-AD | Vinyl Apron Large |

### Neoprene Aprons 5

- Neoprene aprons provides protection from a large range of chemicals such as acids, solvents, caustics and oils.
- Ideal for chemical plants, plating facilities & general lab work.

| Code     | Specs                |
|----------|----------------------|
| NUPAP-CH | Neoprene Apron Small |
| NUPAP-AD | Neoprene Apron Large |

### Size Chart

| Code       | Width (in) | Overall Length (in) |
|------------|------------|---------------------|
| RBAP-CH    | 24"        | 30"                 |
| RBAP-MD    | 27"        | 36"                 |
| RBAP-AD    | 27"        | 42"                 |
| VNAP-CH    | 27"        | 36"                 |
| VNAP-AD    | 36"        | 45"                 |
| NUPAP-CH   | 24"        | 30"                 |
| NUPAP-AD   | 27"        | 42"                 |
| RBAP-CH-RD | 24"        | 30"                 |
| RBAP-AD-BL | 27"        | 42"                 |

# eisco

INDUSTRIAL

---

## RUBBERWARE

- Bulbs
- Stoppers
- Rubber Tubing
- Natural Rubber Tubing
- Neoprene Tubing
- Nitrile Tubing





### Neoprene Tubing - For Burners

- Neoprene with reinforced ends, internal diameter of tubing 5 mm.
- Can accept inlet burner tube with dia. up to 8 mm.

| Code                 | Specs |
|----------------------|-------|
| CH0100D <sup>1</sup> | 45 cm |
| CH0100E              | 60 cm |
| CH0100F              | 90 cm |



Check out  
our new website  
[www.eiscoindustrial.com](http://www.eiscoindustrial.com)



### Pipette Filler - Rubber Bulb

- Rubber bulb filler for petting, suitable for noxious solutions.
- Microbial suspensions and other harmful and unpleasant solutions.
- Comprised of a bulb and three pinch valves controlling influx and efflux of pipette contents. Bulb dia. 54 mm.

| Code                 | Specs                    |
|----------------------|--------------------------|
| CH0580A <sup>2</sup> | Pipette Filler - Rubber  |
| CH0580D              | Pipette Filler - Silicon |

## Bulbs - Pipette

- Pear shaped, for use with pipettes etc.

| Rubber  |   | Silicon    | Capacity |
|---------|---|------------|----------|
| CH0230A | 1 | CH0230A-SI | 5 ml     |
| CH0230B |   | CH0230B-SI | 10 ml    |
| CH0230C |   | CH0230C-SI | 15 ml    |
| CH0230D | 2 | CH0230D-SI | 25 ml    |
| CH0230E |   | CH0230E-SI | 50 ml    |
| CH0230F |   | CH0230F-SI | 20 ml    |
| CH0230G |   | CH0230G-SI | 100 ml   |

## Bulbs - Teat

- For use with dropper pipettes.
- Capacity 1 ml approx.

| Code     |   | Specs           |
|----------|---|-----------------|
| CH0232A  |   | Bulbs - Rubber  |
| CH0232AA | 4 | Bulbs - Silicon |

## Bulbs - Rubber

- Heavy weight, for pipettes & medicine droppers.
- Dia. 5 - 6 mm.
- Capacity 2 ml approx.

| Code     |   | Specs           |
|----------|---|-----------------|
| CH0232B  |   | Bulbs - Rubber  |
| CH0232BB | 5 | Bulbs - Silicon |

## Bulbs - Latex

- Bulbs - Latex

| Code    |   | Capacity     |
|---------|---|--------------|
| CH0232C |   | 1 ml approx. |
| CH0232D | 6 | 2 ml approx. |



## Rubber Stoppers



1

### Rubber Stoppers 1

- Solid, superior quality rubber.
- Resistant to acid, alkali and ammonia.

| Code      | Bottom | Top     | Length | Pack Size |
|-----------|--------|---------|--------|-----------|
| CH0314A01 | 6 mm   | 9 mm    | 16 mm  | 10        |
| CH0314A02 | 8 mm   | 10.5 mm | 20 mm  | 10        |
| CH0314A   | 9 mm   | 11.5 mm | 20 mm  | 10        |
| CH0314B   | 11 mm  | 14 mm   | 24 mm  | 10        |
| CH0314C   | 13 mm  | 16 mm   | 24 mm  | 10        |
| CH0314D   | 15 mm  | 18 mm   | 24 mm  | 10        |
| CH0314E   | 17 mm  | 20 mm   | 26 mm  | 10        |
| CH0314F   | 18 mm  | 21 mm   | 26 mm  | 10        |
| CH0314G   | 19 mm  | 22 mm   | 28 mm  | 10        |
| CH0310H   | 21 mm  | 24 mm   | 28 mm  | 10        |
| CH0314I   | 23 mm  | 26 mm   | 28 mm  | 10        |
| CH0314J   | 25 mm  | 28 mm   | 28 mm  | 10        |
| CH0314K   | 27 mm  | 31 mm   | 32 mm  | 10        |
| CH0314L   | 29 mm  | 33 mm   | 32 mm  | 10        |
| CH0314M   | 31 mm  | 36 mm   | 35 mm  | 10        |
| CH0314N   | 33 mm  | 38 mm   | 38 mm  | 10        |
| CH0314N01 | 33 mm  | 41 mm   | 25 mm  | 10        |
| CH0314O   | 35 mm  | 45 mm   | 36 mm  | 10        |
| CH0314P   | 38 mm  | 42 mm   | 40 mm  | 10        |
| CH0314Q   | 40 mm  | 49 mm   | 40 mm  | 10        |
| CH0314R   | 44 mm  | 58 mm   | 45 mm  | 10        |
| CH0314S   | 48 mm  | 54 mm   | 33 mm  | 10        |
| CH0314T   | 51 mm  | 60 mm   | 57 mm  | 10        |
| CH0314U   | 54 mm  | 64 mm   | 25 mm  | 10        |
| CH0314X   | 75 mm  | 90 mm   | 39 mm  | 10        |



2

### Rubber Stoppers - Solid

- Assorted, 7 pieces of size 8, 9, 11, 13, 15 mm bottom dia. solid.
- 3 pcs. each of 17, 21, 23, 25, 29 mm bottom dia solid.
- Pack of 50.

| Code   | Specs                                |
|--------|--------------------------------------|
| CH0315 | Rubber Stoppers - Solid - Pack of 50 |

### Rubber Stoppers - Solid

- Assorted, 8 pieces of size 8, 9, 11, 13, 15, 17 mm bottom dia. 3 pcs. of size 21, 23 mm bottom dia.
- 2 pcs. each of 25, 29 mm dia. Pack of 58.

| Code   | Specs                              |
|--------|------------------------------------|
| CH0318 | Rubber Stoppers Solid - Pack of 58 |

## Rubber Stoppers - One hole and Two hole

- Resistant to acid, alkali and ammonia.
- Superior quality rubber.
- One hole and two hole stoppers available in following sizes.

| 1 Hole       | 2 Hole       | Bottom | Top     | Length | Pack Size |
|--------------|--------------|--------|---------|--------|-----------|
| CH0314A01/1H |              | 6 mm   | 9 mm    | 16 mm  | 10        |
| CH0314A02/1H |              | 8 mm   | 10.5 mm | 20 mm  | 10        |
| CH0314A/1H   | CH0314A/2H   | 9 mm   | 11.5 mm | 20 mm  | 10        |
| CH0314B/1H   | CH0314B/2H   | 11 mm  | 14 mm   | 24 mm  | 10        |
| CH0314C/1H   | CH0314C/2H   | 13 mm  | 16 mm   | 24 mm  | 10        |
| CH0314D/1H   | CH0314D/2H   | 15 mm  | 18 mm   | 24 mm  | 10        |
| CH0314E/1H   | CH0314E/2H   | 17 mm  | 20 mm   | 26 mm  | 10        |
| CH0314F/1H   | CH0314F/2H   | 18 mm  | 21 mm   | 26 mm  | 10        |
| CH0314G/1H   | CH0314G/2H   | 19 mm  | 22 mm   | 28 mm  | 10        |
| CH0310H/1H   | CH0310H/2H   | 21 mm  | 24 mm   | 28 mm  | 10        |
| CH0314I/1H   | CH0314I/2H   | 23 mm  | 26 mm   | 28 mm  | 10        |
| CH0314J/1H   | CH0314J/2H   | 25 mm  | 28 mm   | 28 mm  | 10        |
| CH0314K/1H   | CH0314K/2H   | 27 mm  | 31 mm   | 32 mm  | 10        |
| CH0314L/1H   | CH0314L/2H   | 29 mm  | 33 mm   | 32 mm  | 10        |
| CH0314M/1H   | CH0314M/2H   | 31 mm  | 36 mm   | 35 mm  | 10        |
| CH0314N/1H   | CH0314N/2H   | 33 mm  | 38 mm   | 38 mm  | 10        |
| CH0314N01/1H | CH0314N01/2H | 33 mm  | 41 mm   | 25 mm  | 10        |
| CH0314O/1H   | CH0314O/2H   | 35 mm  | 45 mm   | 36 mm  | 10        |
| CH0314P/1H   | CH0314P/2H   | 38 mm  | 42 mm   | 40 mm  | 10        |
| CH0314Q/1H   | CH0314Q/2H   | 40 mm  | 49 mm   | 40 mm  | 10        |
| CH0314R/1H   | CH0314R/2H   | 44 mm  | 58 mm   | 45 mm  | 10        |
| CH0314S/1H   | CH0314S/2H   | 48 mm  | 54 mm   | 33 mm  | 10        |
| CH0314T/1H   | CH0314T/2H   | 51 mm  | 60 mm   | 57 mm  | 10        |
| CH0314U/1H   | CH0314U/2H   | 54 mm  | 64 mm   | 25 mm  | 10        |
| CH0314X/1H   | CH0314X/2H   | 75 mm  | 90 mm   | 39 mm  | 10        |

## Rubber Stoppers - One Hole

- Assorted, containing 7 pieces each size 8, 9, 11, 13, 15 mm bottom dia. 3 pcs. each size 17, 21, 23, 25, 29 mm bottom dia.
- One hole. Pack of 50.

| Code   | Specs                      |
|--------|----------------------------|
| CH0316 | Rubber Stoppers - One Hole |

## Rubber Stopper - Two Hole

- Assorted, containing 10 pieces size 13, 15, 17, 18 mm bottom dia., 3 pcs. 21, 23 mm bottom dia, 2 pcs. 25, 29 mm dia.
- Two hole. Pack of 50.

| Code   | Specs                     |
|--------|---------------------------|
| CH0317 | Rubber Stopper - Two Hole |



1



2



3



4

GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Silicon Stoppers



### Silicone Stoppers - Translucent

1

| Code    | Bottom | Top     | Length |
|---------|--------|---------|--------|
| CH0319A | 6 mm   | 8 mm    | 16 mm  |
| CH0319B | 8 mm   | 10.5 mm | 20 mm  |
| CH0319C | 9 mm   | 11.5 mm | 20 mm  |
| CH0319D | 10 mm  | 12.5 mm | 20 mm  |
| CH0319E | 11 mm  | 14 mm   | 24 mm  |
| CH0319F | 13 mm  | 16 mm   | 24 mm  |
| CH0319G | 15 mm  | 18 mm   | 24 mm  |
| CH0319H | 17 mm  | 20 mm   | 26 mm  |
| CH0319I | 18 mm  | 21 mm   | 26 mm  |
| CH0319J | 19 mm  | 22 mm   | 28 mm  |
| CH0319K | 21 mm  | 24 mm   | 28 mm  |
| CH0319L | 23 mm  | 26 mm   | 28 mm  |
| CH0319M | 25 mm  | 28 mm   | 28 mm  |
| CH0319N | 27 mm  | 31 mm   | 32 mm  |
| CH0319O | 29 mm  | 31 mm   | 32 mm  |
| CH0319P | 31 mm  | 36 mm   | 35 mm  |
| CH0319Q | 33 mm  | 38 mm   | 38 mm  |
| CH0319R | 35 mm  | 45 mm   | 36 mm  |
| CH0319S | 38 mm  | 42 mm   | 40 mm  |
| CH0319T | 40 mm  | 49 mm   | 40 mm  |



Check out  
our new website  
[www.eiscoindustrial.com](http://www.eiscoindustrial.com)

## EISCO PREMIUM NEOPRENE STOPPERS

- Research grade for industrial and academic research labs
- Manufactured from long lasting neoprene
- Suitable for petroleum products, most organic acids and bases
- Oil resistant
- Suitable for high temperatures

### Neoprene Stoppers - Solid ①

- Manufactured from long lasting neoprene.
- Suitable for use with petroleum products, minerals, oils, fats as well as most inorganic acids and bases.
- Neoprene is oil resistant and suitable for temperatures up to 250°F/120°C.
- Solid stoppers available in the following sizes.

| Code        | Bottom | Top     | Length |
|-------------|--------|---------|--------|
| CH0319A-NUP | 6 mm   | 8 mm    | 16 mm  |
| CH0319B-NUP | 8 mm   | 10.5 mm | 20 mm  |
| CH0319C-NUP | 9 mm   | 11.5 mm | 20 mm  |
| CH0319D-NUP | 10 mm  | 12.5 mm | 20 mm  |
| CH0319E-NUP | 11 mm  | 14 mm   | 24 mm  |
| CH0319F-NUP | 13 mm  | 16 mm   | 24 mm  |
| CH0319G-NUP | 15 mm  | 18 mm   | 24 mm  |
| CH0319H-NUP | 17 mm  | 20 mm   | 26 mm  |
| CH0319I-NUP | 18 mm  | 21 mm   | 26 mm  |
| CH0319J-NUP | 19 mm  | 22 mm   | 28 mm  |
| CH0319K-NUP | 21 mm  | 24 mm   | 28 mm  |
| CH0319L-NUP | 23 mm  | 26 mm   | 28 mm  |
| CH0319M-NUP | 25 mm  | 28 mm   | 28 mm  |
| CH0319N-NUP | 27 mm  | 31 mm   | 32 mm  |
| CH0319O-NUP | 29 mm  | 31 mm   | 32 mm  |
| CH0319P-NUP | 31 mm  | 36 mm   | 35 mm  |
| CH0319Q-NUP | 33 mm  | 38 mm   | 38 mm  |
| CH0319R-NUP | 35 mm  | 45 mm   | 36 mm  |
| CH0319S-NUP | 38 mm  | 42 mm   | 40 mm  |
| CH0319T-NUP | 40 mm  | 49 mm   | 40 mm  |



## Neoprene Stoppers



### Neoprene Stoppers - One and Two Hole

- Manufactured from long lasting neoprene.
- Suitable for use with petroleum products, minerals, oils, fats as well as most inorganic acids and bases.
- Neoprene is resistant and suitable for temperatures up to 250°F/120°C.
- One hole and two hole stoppers available in the following sizes.

| 1 Hole         | 2 Hole         | Bottom | Top     | Length |
|----------------|----------------|--------|---------|--------|
| CH0319A-NUP/1H |                | 6 mm   | 8 mm    | 16 mm  |
| CH0319B-NUP/1H |                | 8 mm   | 10.5 mm | 20 mm  |
| CH0319C-NUP/1H |                | 9 mm   | 11.5 mm | 20 mm  |
| CH0319D-NUP/1H |                | 10 mm  | 12.5 mm | 20 mm  |
| CH0319E-NUP/1H |                | 11 mm  | 14 mm   | 24 mm  |
| CH0319F-NUP/1H | CH0319F-NUP/2H | 13 mm  | 16 mm   | 24 mm  |
| CH0319G-NUP/1H | CH0319G-NUP/2H | 15 mm  | 18 mm   | 24 mm  |
| CH0319H-NUP/1H | CH0319H-NUP/2H | 17 mm  | 20 mm   | 26 mm  |
| CH0319I-NUP/1H | CH0319I-NUP/2H | 18 mm  | 21 mm   | 26 mm  |
| CH0319J-NUP/1H | CH0319J-NUP/2H | 19 mm  | 22 mm   | 28 mm  |
| CH0319K-NUP/1H | CH0319K-NUP/2H | 21 mm  | 24 mm   | 28 mm  |
| CH0319L-NUP/1H | CH0319L-NUP/2H | 23 mm  | 26 mm   | 28 mm  |
| CH0319M-NUP/1H | CH0319M-NUP/2H | 25 mm  | 28 mm   | 28 mm  |
| CH0319N-NUP/1H | CH0319N-NUP/2H | 27 mm  | 31 mm   | 32 mm  |
| CH0319O-NUP/1H | CH0319O-NUP/2H | 29 mm  | 31 mm   | 32 mm  |
| CH0319P-NUP/1H | CH0319P-NUP/2H | 31 mm  | 36 mm   | 35 mm  |
| CH0319Q-NUP/1H | CH0319Q-NUP/2H | 33 mm  | 38 mm   | 38 mm  |
| CH0319R-NUP/1H | CH0319R-NUP/2H | 35 mm  | 45 mm   | 36 mm  |
| CH0319S-NUP/1H | CH0319S-NUP/2H | 38 mm  | 42 mm   | 40 mm  |
| CH0319T-NUP/1H | CH0319T-NUP/2H | 40 mm  | 49 mm   | 40 mm  |



Make sure to sign up  
for our **monthly** email

## Neoprene Stoppers - ASTM, Solid ①

- Manufactured according to sizes recommended by ASTM standards.
- Manufactured from long lasting neoprene.
- Suitable for use with petroleum products, minerals, oils, fats as well as most inorganic acids and bases.
- Neoprene is oil resistance and suitable for temperatures up to 250°F/120°C.
- Solid stoppers available in the following sizes.

| Code    | Bottom | Top     | Length |
|---------|--------|---------|--------|
| CH0321A | 8.2 mm | 12.7 mm | 25 mm  |
| CH0321B | 10 mm  | 15 mm   | 25 mm  |
| CH0321C | 13 mm  | 17 mm   | 25 mm  |
| CH0321D | 14 mm  | 19 mm   | 25 mm  |
| CH0321E | 16 mm  | 20 mm   | 25 mm  |
| CH0321F | 18 mm  | 24 mm   | 25 mm  |
| CH0321G | 20 mm  | 26 mm   | 25 mm  |
| CH0321H | 23 mm  | 27 mm   | 25 mm  |
| CH0321I | 24 mm  | 28 mm   | 25 mm  |
| CH0321J | 26 mm  | 32 mm   | 25 mm  |
| CH0321K | 27 mm  | 34 mm   | 25 mm  |
| CH0321L | 30 mm  | 37 mm   | 25 mm  |
| CH0321M | 31 mm  | 39 mm   | 25 mm  |
| CH0321N | 33 mm  | 41 mm   | 25 mm  |
| CH0321O | 36 mm  | 43 mm   | 25 mm  |
| CH0321P | 37 mm  | 45 mm   | 25 mm  |
| CH0321Q | 42 mm  | 50 mm   | 25 mm  |
| CH0321R | 45 mm  | 53 mm   | 25 mm  |
| CH0321S | 48 mm  | 56 mm   | 25 mm  |
| CH0321T | 50 mm  | 63 mm   | 25 mm  |
| CH0321U | 54 mm  | 64 mm   | 25 mm  |
| CH0321V | 58 mm  | 68 mm   | 25 mm  |
| CH0321W | 62 mm  | 75 mm   | 35 mm  |
| CH0321X | 75 mm  | 90 mm   | 39 mm  |
| CH0321Y | 83 mm  | 103 mm  | 39 mm  |
| CH0321Z | 90 mm  | 127 mm  | 50 mm  |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Neoprene Stoppers & Rubber Cones



1



2



3

### Neoprene Stoppers - ASTM, One and Two Hole

- Manufactured according to sizes recommended by ASTM standards.
- Manufactured from long lasting neoprene.
- Suitable for use with petroleum products, minerals, oils, fats as well as most inorganic acids and bases.
- Neoprene is oil resistant and suitable for temperatures up to 250°F/120°C.
- One hole and two hole stoppers available in the following sizes.

| 1 Hole     | 2 Hole     | Bottom | Top     | Length |
|------------|------------|--------|---------|--------|
| CH0321A/1H |            | 8.2 mm | 12.7 mm | 25 mm  |
| CH0321B/1H |            | 10 mm  | 15 mm   | 25 mm  |
| CH0321C/1H | CH0321C/2H | 13 mm  | 17 mm   | 25 mm  |
| CH0321D/1H | CH0321D/2H | 14 mm  | 19 mm   | 25 mm  |
| CH0321E/1H | CH0321E/2H | 16 mm  | 20 mm   | 25 mm  |
| CH0321F/1H | CH0321F/2H | 18 mm  | 24 mm   | 25 mm  |
| CH0321G/1H | CH0321G/2H | 20 mm  | 26 mm   | 25 mm  |
| CH0321H/1H | CH0321H/2H | 23 mm  | 27 mm   | 25 mm  |
| CH0321I/1H | CH0321I/2H | 24 mm  | 28 mm   | 25 mm  |
| CH0321J/1H | CH0321J/2H | 26 mm  | 32 mm   | 25 mm  |
| CH0321K/1H | CH0321K/2H | 27 mm  | 34 mm   | 25 mm  |
| CH0321L/1H | CH0321L/2H | 30 mm  | 37 mm   | 25 mm  |
| CH0321M/1H | CH0321M/2H | 31 mm  | 39 mm   | 25 mm  |
| CH0321N/1H | CH0321N/2H | 33 mm  | 41 mm   | 25 mm  |
| CH0321O/1H | CH0321O/2H | 36 mm  | 43 mm   | 25 mm  |
| CH0321P/1H | CH0321P/2H | 37 mm  | 45 mm   | 25 mm  |
| CH0321Q/1H | CH0321Q/2H | 42 mm  | 50 mm   | 25 mm  |
| CH0321R/1H | CH0321R/2H | 45 mm  | 53 mm   | 25 mm  |
| CH0321S/1H | CH0321S/2H | 48 mm  | 56 mm   | 25 mm  |
| CH0321T/1H | CH0321T/2H | 50 mm  | 63 mm   | 25 mm  |
| CH0321U/1H | CH0321U/2H | 54 mm  | 64 mm   | 25 mm  |
| CH0321V/1H | CH0321V/2H | 58 mm  | 68 mm   | 25 mm  |
| CH0321W/1H | CH0321W/2H | 62 mm  | 75 mm   | 35 mm  |
| CH0321X/1H | CH0321X/2H | 75 mm  | 90 mm   | 39 mm  |
| CH0321Y/1H | CH0321Y/2H | 83 mm  | 103 mm  | 39 mm  |
| CH0321Z/1H | CH0321Z/2H | 90 mm  | 127 mm  | 50 mm  |

### Rubber Cones

- Set of 7 rubber cone for use with Buchner funnels of various sizes

| Code   | Capacity               |
|--------|------------------------|
| CH0421 | Rubber cones, Set of 7 |

## EISCO PREMIUM TUBING

- Research grade for industrial and academic research
- Premium quality polypropylene racks
- Made from high quality natural rubber
- Resistant to acid and alkalis
- High pressure tubing
- Neoprene, nitrile, rubber, extra soft tubings

### Rubber Tubing 1

- Soft quality, resistant to acid and alkalis.
- Pack of 10 meters.
- Shore Hardness: 45A
- Tolerance:  $\pm 0.2\text{mm}$
- Temp.:  $-20^{\circ}\text{C}$  to  $+70^{\circ}\text{C}$



| Code      | Bore   | Wall Thickness |
|-----------|--------|----------------|
| CH0761A01 | 3 mm   | 1.0 mm         |
| CH0761A   | 4 mm   | 1.0 mm         |
| CH0761B   | 5 mm   | 1.5 mm         |
| CH0761C   | 6 mm   | 1.5 mm         |
| CH0761C01 | 6.5 mm | 1.5 mm         |
| CH0761D   | 7 mm   | 1.5 mm         |
| CH0761E   | 8 mm   | 2.0 mm         |
| CH0761F   | 9 mm   | 2.0 mm         |
| CH0761G   | 10 mm  | 2.0 mm         |
| CH0761H   | 12 mm  | 2.5 mm         |
| CH0761I   | 15 mm  | 3.0 mm         |
| CH0761J   | 12 mm  | 3.5 mm         |
| CH0761K   | 15 mm  | 3.5 mm         |
| CH0761L   | 16 mm  | 3.0 mm         |
| CH0761M   | 19 mm  | 3.0 mm         |
| CH0761N   | 20 mm  | 3.5 mm         |
| CH0761O   | 25 mm  | 3.5 mm         |



1

## Rubber Tubing



### Neoprene Rubber Tubing - Synthetic

- Synthetic neoprene rubber tubing.
- Pack of 10 meters.

| Code    | Bore  | Wall Thickness |
|---------|-------|----------------|
| NEOCH01 | 4 mm  | 1.0 mm         |
| NEOCH02 | 5 mm  | 1.5 mm         |
| NEOCH03 | 6 mm  | 1.5 mm         |
| NEOCH04 | 7 mm  | 1.5 mm         |
| NEOCH05 | 8 mm  | 2.0 mm         |
| NEOCH06 | 9 mm  | 2.0 mm         |
| NEOCH07 | 10 mm | 2.0 mm         |
| NEOCH08 | 12 mm | 2.0 mm         |

### Rubber Tubing - Nitrile

- Made of nitrile rubber.
- Pack of 5 meters.

| Code   | Bore  | Wall Thickness |
|--------|-------|----------------|
| NITB01 | 4 mm  | 1.0 mm         |
| NITB02 | 5 mm  | 1.5 mm         |
| NITB03 | 6 mm  | 1.5 mm         |
| NITB04 | 7 mm  | 1.5 mm         |
| NITB05 | 8 mm  | 2.0 mm         |
| NITB06 | 9 mm  | 2.0 mm         |
| NITB07 | 10 mm | 2.0 mm         |
| NITB08 | 12 mm | 2.0 mm         |

### Rubber Tubing - High Pressure

- High pressure, extra soft quality & Per meter.
- Shore Hardness: 50-55A
- Tolerance:  $\pm 0.2$ mm
- Temp.: -20°C to +70°C

| Code      | Bore    | Wall Thickness |
|-----------|---------|----------------|
| CH0762A   | 3 mm    | 4.5 mm         |
| CH0762B   | 4.5 mm  | 4.5 mm         |
| CH0762C   | 6 mm    | 4.5 mm         |
| CH0762D   | 8 mm    | 4.5 mm         |
| CH0762E   | 9 mm    | 4.5 mm         |
| CH0762E01 | 10 mm   | 6.0 mm         |
| CH0762F   | 12.5 mm | 6.5 mm         |

## Silicone Tubing <sup>1</sup>

- Seamless, non-toxic, translucent, very flexible, suitable for connecting to glass.
- Silicone tubing is suitable for use over a temperature range of -50°C to 250°C.
- Meets USP VI standards and FDA 21 CFR 177.2600 compliant.
- Shore Hardness: 45A.
- Sold by the meter.

| Code      | Bore  | Wall Thickness |
|-----------|-------|----------------|
| CH0764A   | 2 mm  | 1.0 mm         |
| CH0764B   | 3 mm  | 1.0 mm         |
| CH0764B01 | 3 mm  | 1.5 mm         |
| CH0764D   | 4 mm  | 1.0 mm         |
| CH0764E   | 4 mm  | 1.5 mm         |
| CH0764E   | 4 mm  | 2.0 mm         |
| CH0764F   | 5 mm  | 1.0 mm         |
| CH0764G   | 5 mm  | 1.5 mm         |
| CH0764I01 | 6 mm  | 1.0 mm         |
| CH0764I   | 6 mm  | 1.5 mm         |
| CH0764J   | 6 mm  | 2.0 mm         |
| CH0764J01 | 6 mm  | 3.0 mm         |
| CH0764M01 | 7 mm  | 1.5 mm         |
| CH0764M   | 7 mm  | 2.0 mm         |
| CH0764O01 | 8 mm  | 1.5 mm         |
| CH0764O   | 8 mm  | 2.0 mm         |
| CH0764P   | 8 mm  | 3.0 mm         |
| CH0764Q   | 9 mm  | 2.0 mm         |
| CH0764R   | 10 mm | 2.0 mm         |
| CH0764S   | 10 mm | 3.0 mm         |
| CH0764T01 | 12 mm | 2.0 mm         |
| CH0764T   | 12 mm | 3.5 mm         |
| CH0764V   | 14 mm | 3.0 mm         |
| CH0764W   | 15 mm | 4.0 mm         |
| CH0764X   | 16 mm | 3.0 mm         |
| CH0764X01 | 25 mm | 4.0 mm         |



## Hand Made Vacuum Tubing <sup>2</sup>

- High grade, handmade, extra soft rubber tubing.
- Extra heavy wall for vacuum work.
- Does not collapse under atmospheric conditions.
- Sold by the meter.

|         | Bore | Wall Thickness |
|---------|------|----------------|
| CH0763A | 6 mm | 7 mm           |
| CH0763B | 8 mm | 7 mm           |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS



## CERAMICS

- Evaporating Dishes
- Crucibles
- Funnels
- Mortar & Pestles
- Spotting Plates
- Caseroles



### Evaporating Dishes 1

- Porcelain, flat form, with spout.
- Glazed inside and outside except for rim.

| Code    | Capacity | O.D.   | Height |
|---------|----------|--------|--------|
| CH0070A | 35 ml    | 65 mm  | 16 mm  |
| CH0070B | 50 ml    | 75 mm  | 20 mm  |
| CH0070C | 100 ml   | 100 mm | 23 mm  |
| CH0070D | 200 ml   | 120 mm | 34 mm  |
| CH0070E | 300 ml   | 150 mm | 29 mm  |



1

### Evaporating Dishes 2

- Porcelain, round form, with spout & superior quality.
- Fully glazed except outside base.

| Code    | Capacity | O.D.   | Height |
|---------|----------|--------|--------|
| CH0071A | 35 ml    | 60 mm  | 24 mm  |
| CH0071B | 70 ml    | 75 mm  | 30 mm  |
| CH0071C | 80 ml    | 80 mm  | 30 mm  |
| CH0071D | 120 ml   | 94 mm  | 42 mm  |
| CH0071E | 150 ml   | 100 mm | 42 mm  |
| CH0071F | 250 ml   | 115 mm | 45 mm  |
| CH0071G | 385 ml   | 145 mm | 48 mm  |
| CH0071H | 525 ml   | 162 mm | 51 mm  |



2

### Evaporating Dishes 3

- Porcelain, with spout, deep form.
- Inside outside glazed.

| Code     | Capacity | O.D.   |
|----------|----------|--------|
| CH0072A  | 8 ml     | 25 mm  |
| CH00712B | 20 ml    | 50 mm  |
| CH0072C  | 40 ml    | 60 mm  |
| CH0072D  | 100 ml   | 75 mm  |
| CH0072E  | 175 ml   | 100 mm |
| CH0072F  | 275 ml   | 110 mm |
| CH0072G  | 450 ml   | 125 mm |
| CH0072H  | 850 ml   | 150 mm |
| CH0072I  | 1300 ml  | 175 mm |



3

### Evaporating Dishes- Silica 4

- Round bottom, with spout & translucent.

| Code    | Capacity | Dia.   | Height |
|---------|----------|--------|--------|
| CH0078A | 20 ml    | 55 mm  | 23 mm  |
| CH0078B | 45 ml    | 75 mm  | 27 mm  |
| CH0078C | 70 ml    | 87 mm  | 33 mm  |
| CH0078D | 85 ml    | 94 mm  | 36 mm  |
| CH0078E | 100 ml   | 100 mm | 45 mm  |
| CH0078F | 200 ml   | 113 mm | 47 mm  |



4

GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Crucibles | Evaporating Dishes

**Crucibles - Porcelain, Tall** 1

- Porcelain & tall
- Glazed inside and outside except base & with lid.

| Code    | Capacity | Height |
|---------|----------|--------|
| CH0074A | 10 ml    | 26 mm  |
| CH0074B | 15 ml    | 29 mm  |
| CH0074D | 30 ml    | 37 mm  |
| CH0074E | 40 ml    | 40 mm  |
| CH0074F | 50 ml    | 43 mm  |
| CH0074G | 100 ml   | 54 mm  |
| CH0074H | 250 ml   | 72 mm  |

**Crucibles - Porcelain, Squat** 2

- Porcelain & squat
- Glazed inside and outside except base.

| Code    | Capacity | Height |
|---------|----------|--------|
| CH0076A | 8 ml     | 20 mm  |
| CH0076B | 17 ml    | 26 mm  |
| CH0076C | 30 ml    | 30 mm  |
| CH0076D | 50 ml    | 37 mm  |
| CH0076E | 100 ml   | 46 mm  |
| CH0076F | 150 ml   | 52 mm  |
| CH0076G | 250 ml   | 60 mm  |

**Crucible - Silica** 3

- Low, glazed inside and outside.

| Code    | Capacity | Dia.  | Height |
|---------|----------|-------|--------|
| CH0080A | 15 ml    | 41 mm | 25 mm  |
| CH0080B | 25 ml    | 47 mm | 28 mm  |
| CH0080C | 50 ml    | 57 mm | 37 mm  |
| CH0080D | 80 ml    | 67 mm | 45 mm  |
| CH0080E | 100 ml   | 71 mm | 50 mm  |
| CH0080F | 150 ml   | 82 mm | 51 mm  |

### Buchner Funnels - Porcelain 1

- Buchner funnel with straight sides above the perforated plate.
- Thick stem.

| Code    | Dia.    |
|---------|---------|
| CH0462A | 3.5 cm  |
| CH0462B | 5 cm    |
| CH0462C | 6 cm    |
| CH0462D | 7.5 cm  |
| CH0462E | 10 cm   |
| CH0462F | 12.5 cm |
| CH0462G | 15 cm   |
| CH0462H | 20 cm   |
| CH0462I | 25 cm   |
| CH0462J | 30 cm   |



### Hirsch Funnels

- Glazed with fixed perforated plate.

| Code      | Dia.   |
|-----------|--------|
| CH0465A   | 3 cm   |
| CH0465B   | 4 cm   |
| CH0465C 2 | 5 cm   |
| CH0465D   | 7.5 cm |
| CH0465E   | 10 cm  |



### Gooch Crucibles

- Perforated bottom.
- Glazed inside and outside with the exception of outside bottom surface and rim.

| Code      | Capacity |
|-----------|----------|
| CH0509D 3 | 15 ml    |
| CH0509E   | 30 ml    |
| CH0509F   | 50 ml    |
| CH0509G   | 140 ml   |



GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

## Mortar &amp; Pestle | Spotting Plates | Caseroles

**Mortar and Pestle**

- Porcelain
- Grinding surface unglazed.

| Code                 | Outer Dia. | Capacity |
|----------------------|------------|----------|
| CH0540A <sup>1</sup> | 80 mm      | 60 ml    |
| CH0540B              | 105 mm     | 150 ml   |
| CH0540C              | 125 mm     | 275 ml   |
| CH0540D              | 155 mm     | 550 ml   |
| CH0540E              | 180 mm     | 900 ml   |
| CH0540F              | 200 mm     | 1500 ml  |
| CH0540G              | 230 mm     | 2000 ml  |

**Mortar and Pestle - Heavy pattern <sup>2</sup>**

- Made of all porcelain.

| Code    | Outer Dia. | Capacity |
|---------|------------|----------|
| CH0541A | 70 mm      | 50 ml    |
| CH0541B | 80 mm      | 70 ml    |
| CH0541C | 90 mm      | 150 ml   |
| CH0541D | 115 mm     | 275 ml   |
| CH0541E | 130 mm     | 400 ml   |
| CH0541F | 165 mm     | 750 ml   |
| CH0541G | 210 mm     | 1900 ml  |

**Spotting Plates**

- Made of porcelain
- Glazed on top and sides.

| Code                 | Specs                 |
|----------------------|-----------------------|
| CH0542A <sup>3</sup> | 6 cavity - Porcelain  |
| CH0542B <sup>4</sup> | 12 cavity - Porcelain |

**Caseroles <sup>5</sup>**

- Caseroles with lips
- Pouring spouts
- Porcelain handles.
- Fully glazed excluding rim.

| Code    | Capacity |
|---------|----------|
| CH0509A | 65 ml    |
| CH0509B | 200 ml   |
| CH0509C | 500 ml   |

# GLASSWARE

## A

### Adapters

|                                                            |    |
|------------------------------------------------------------|----|
| Bend Receiver Adapters with Vent.....                      | 61 |
| Cone/Screw Thread Adapters - Straight .....                | 61 |
| Expansion Adapters .....                                   | 59 |
| Multiple Adapters - Three Neck .....                       | 60 |
| Multiple Adapters - Two Neck, Angled .....                 | 60 |
| Multiple Adapters - Two Neck, Parallel .....               | 60 |
| Receiver Adapters - Bend with Vacuum .....                 | 62 |
| Receiver Adapters - Screw Thread .....                     | 71 |
| Receiver Adapters with Vacuum Connection - Screw Thread .. | 71 |
| Receiver Delivery Adapters .....                           | 61 |
| Reduction Adapters .....                                   | 59 |
| Splash Head Adapters - Pear Shape, Sloping.....            | 63 |
| Splash Head Adapters - Pear Shape, Vetricle .....          | 63 |
| Splash Heads Adapter - Pear Shape, Straight.....           | 62 |

## B

### Beakers

|                            |   |
|----------------------------|---|
| ASTM, Low Form .....       | 8 |
| ASTM, Tall Form.....       | 8 |
| Low Form with Spout .....  | 7 |
| Phillips Pattern .....     | 8 |
| Tall Form with Spout ..... | 7 |

### Bottles

|                                                             |    |
|-------------------------------------------------------------|----|
| BOD Bottles, Glass Stopper .....                            | 12 |
| Aspirator with Outlet for Stopcock .....                    | 14 |
| Aspirator with Outlet for Tubing.....                       | 14 |
| Gas Washing .....                                           | 13 |
| Gas Wash Woulff, Two Necks .....                            | 13 |
| Pipette Dropping .....                                      | 12 |
| Pipette Dropping, Amber .....                               | 12 |
| Weighing, Low Form .....                                    | 13 |
| Weighing, Tall Form .....                                   | 13 |
| Reagent Bottles - Graduated, Screw Cap .....                | 11 |
| Reagent Bottles - Narrow Mouth, Hexagonal Hollow Stopper .. | 10 |
| Reagent Bottles - Narrow Mouth, Polypropylene Stopper ..... | 10 |
| Reagent Bottles - Narrow Mouth, PP Stopper, Amber .....     | 10 |
| Reagent Bottles - Screw Cap, Amber .....                    | 11 |
| Reagent Bottles - Wide Mouth, Hexagonal Hollow Stopper .... | 11 |

### Burettes

|                                                                        |    |
|------------------------------------------------------------------------|----|
| Glass, Class-B ASTM with PTFE Key Stopcock .....                       | 36 |
| Glass with Glass Key Stopcock.....                                     | 33 |
| Glass Key Stopcock with Individual Works Certificate .....             | 82 |
| Glass with Needle Valve Stopcock.....                                  | 34 |
| Glass with Pinch Clip.....                                             | 33 |
| Glass with PTFE Key Stopcock .....                                     | 34 |
| Glass with PTFE Key Stopcock with Individual Works Certificate .....   | 82 |
| Modular, Glass with PTFE Key Stopcock .....                            | 35 |
| Schellbach, Glass with Glass Key Stopcock .....                        | 35 |
| Schellbach, Glass Key Stopcock with Individual Works Certificate ..... | 82 |

|                                                                                  |    |
|----------------------------------------------------------------------------------|----|
| Schellbach, Glass with PTFE Key Stopcock.....                                    | 35 |
| Schellbach, Glass with PTFE Key Stopcock with Individual Works Certificate ..... | 83 |
| Automatic - Glass with PTFE Key Stopcock .....                                   | 37 |
| Automatic - Glass with PTFE Key Stopcock with Individual Works Certificate ..... | 83 |
| Automatic - Glass with Intermediate PTFE Key Stopcock.....                       | 37 |
| Automatic - Schellbach, Glass with PTFE Key Stopcock.....                        | 38 |
| Micro Burettes .....                                                             | 36 |

## C

### Condensers

|                                            |    |
|--------------------------------------------|----|
| Allihn Bulb Condensers .....               | 57 |
| Allihn for Soxhlet Condenser .....         | 66 |
| Bulb Condensers - Screw Thread .....       | 70 |
| Dimroth Condenser .....                    | 58 |
| Graham Condenser - Coiled Distillate.....  | 58 |
| Liebig Condensers.....                     | 57 |
| Liebig Condensers - Jointed .....          | 57 |
| Liebig Condensers - Screw Thread .....     | 70 |
| Conical Measures .....                     | 44 |
| Crow Receivers.....                        | 44 |
| Crucibles - Gooch Type, Sintered Disc..... | 47 |
| Crystallizing Dishes without Spout .....   | 9  |
| Crystallizing Dishes with Spout.....       | 9  |

## D

|                      |    |
|----------------------|----|
| Dean and Stark ..... | 73 |
|----------------------|----|

## E

|                                                   |    |
|---------------------------------------------------|----|
| Evaporating Basins - Flat Bottom with Spout ..... | 9  |
| Extractors for Soxhlet Apparatus.....             | 66 |

## F

|                                   |    |
|-----------------------------------|----|
| Filtration Assembly - 47 mm ..... | 48 |
|-----------------------------------|----|

### Flasks

|                                                                               |    |
|-------------------------------------------------------------------------------|----|
| Boiling, Flat Bottom, Narrow Neck with Rim.....                               | 18 |
| Boiling, Flat Bottom, Wide Neck with Rim.....                                 | 19 |
| Boiling, Flat Bottom with Screw Thread Neck .....                             | 68 |
| Boiling, Flat Bottom with Screw Thread Neck .....                             | 68 |
| Boiling, Pear Shape, Two Neck with Joint .....                                | 21 |
| Boiling, Pear Shape with Joint.....                                           | 21 |
| Boiling, Round Bottom, Narrow Neck with Rim.....                              | 20 |
| Boiling, Round Bottom, Wide Neck, Beaded Rim.....                             | 21 |
| Distillation Flasks - Round Bottom, Oblique Neck with Screw Thread .....      | 69 |
| Distillation Flasks - Round Bottom, Tow Oblique Necks with Screw Thread ..... | 69 |
| Conical, Erlenmeyer .....                                                     | 17 |
| Conical Erlenmeyer, Hollow Glass Stopper .....                                | 16 |
| Conical, Erlenmeyer, Narrow Neck.....                                         | 15 |
| Conical, Erlenmeyer, Narrow Neck with Interchangeable Joint                   | 16 |
| Conical, Erlenmeyer, Wide Neck.....                                           | 15 |
| Conical, Erlenmeyer with Screw Thread Neck .....                              | 69 |
| Distilling, Pear Shape with Joint .....                                       | 22 |

GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

|                                                                                              |    |
|----------------------------------------------------------------------------------------------|----|
| Distilling, Round Bottom, 3 Neck .....                                                       | 22 |
| Distilling, Round Bottom, 3 Neck .....                                                       | 23 |
| Distilling, Round Bottom with Joint .....                                                    | 22 |
| Distilling, Round Bottom with Joint .....                                                    | 23 |
| Distilling with Side Arm .....                                                               | 24 |
| Evaporating .....                                                                            | 24 |
| Filtering .....                                                                              | 17 |
| Filtering, Buchner .....                                                                     | 18 |
| Florence - Boiling, Flat Bottom, Short Neck with Joint .....                                 | 19 |
| Florence - Boiling, Round Bottom with Joint .....                                            | 20 |
| <b>Flasks (cont.)</b> .....                                                                  |    |
| Iodine .....                                                                                 | 17 |
| Kjeldahl .....                                                                               | 24 |
| Kjeldahl, Plain .....                                                                        | 25 |
| Volumetric, Class-A Amber with Stopper .....                                                 | 28 |
| Volumetric, Class-A Amber with Stopper<br>with Individuals Works Certificate .....           | 86 |
| Volumetric Class-A .....                                                                     | 28 |
| Volumetric Class-A ASTM with Snap Cap .....                                                  | 32 |
| Volumetric Class-A ASTM with Solid Glass Stopper .....                                       | 30 |
| Volumetric Class-A ASTM with Solid Glass Stopper<br>with Individuals Works Certificate ..... | 87 |
| Volumetric Class-A USP with PP Stopper<br>with Individuals Works Certificate .....           | 87 |
| Volumetric Class-A ASTM with Stopper .....                                                   | 31 |
| Volumetric Class-A with Hollow Glass Stopper .....                                           | 29 |
| Volumetric Class-A with Stopper .....                                                        | 27 |
| Volumetric Class-A with Stopper<br>with Individual Works Certificate .....                   | 86 |
| Volumetric Class-B .....                                                                     | 28 |
| Volumetric Class-B ASTM with Glass Stopper .....                                             | 31 |
| Volumetric Class-B ASTM with Snap Cap .....                                                  | 32 |
| Volumetric Class-B ASTM with Solid Glass Stopper .....                                       | 30 |
| Volumetric Class-B with Hollow Glass Stopper .....                                           | 29 |
| Volumetric Class-B with Stopper .....                                                        | 27 |
| Fractional Distillation Apparatus .....                                                      | 80 |
| Fractional Distillation Apparatus - Secure Joint .....                                       | 72 |
| Fractional Vacuum Distillation Apparatus - Secure Joint .....                                | 72 |
| <b>Funnels</b> .....                                                                         |    |
| Buchner Funnels - Glass with Sintered Disc .....                                             | 47 |
| Buchner Funnels - Glass with Sintered Disk .....                                             | 48 |
| Dropping Funnel- Needle Valve Stopcock, Plastic Stopper .....                                | 50 |
| Dropping Funnel- Needle Valve Stopcock, Stem with Cone,<br>Plastic Stopper .....             | 51 |
| Dropping Funnel - Open with Glass Key Stopcock .....                                         | 49 |
| Dropping Funnel - Pressure Equalizing with Glass Key Stopcock<br>and Plastic Stopper .....   | 52 |
| Dropping Funnel - Pressure Equalizing with PTFE Key Stopcock and<br>Plastic Stopper .....    | 52 |
| Dropping Funnel - PTFE Key Stopcock, Plastic Stopper .....                                   | 50 |
| Dropping Funnel - PTFE Key Stopcock, Stem with Cone,<br>Plastic Stopper .....                | 51 |
| Dropping Funnels - Glass Key Stopcock, Plastic Stopper .....                                 | 49 |
| Dropping Funnels - Open with PTFE Key Stopcock .....                                         | 49 |
| Dropping Funnels - Screw Thread .....                                                        | 71 |
| Dropping Funnel - Glass Key Stopcock, Stem with Cone,<br>Plastic Stopper .....               | 50 |
| Filter Funnels .....                                                                         | 25 |

|                                                                                                   |    |
|---------------------------------------------------------------------------------------------------|----|
| Filter Funnels - Heavy .....                                                                      | 25 |
| Gilson Separating Funnel - Gas Key Stopcock, Plastic Stopper .....                                | 56 |
| Gilson Separating Funnel - PTFE Key Stopcock, Plastic Stopper .....                               | 56 |
| Powder Funnel .....                                                                               | 26 |
| Powder Funnels - Jointed .....                                                                    | 26 |
| Pressure Equalizing Separating Funnel - Pear Shaped, Glass Key<br>Stopcock, Plastic Stopper ..... | 55 |
| Pressure Equalizing Separating Funnel - Pear Shaped, PTFE Key<br>Stopcock, Plastic Stopper .....  | 55 |
| Separating Funnel - Pear Shaped, Glass Key Stopcock, Plastic<br>Stopper .....                     | 53 |
| Separating Funnel - Pear Shaped, Glass Key Stopcock, Plastic<br>Stopper, Stem with Cone .....     | 55 |
| Separating Funnel - Pear Shaped, PTFE Key Stopcock, Plastic<br>Stopper .....                      | 54 |
| Separating Funnel - Pear Shaped, PTFE Key Stopcock, Plastic<br>Stopper, Stem with Cone .....      | 54 |
| Separating Funnel - Pear Shaped, Rotaflow Stopcock, Plastic<br>Stopper .....                      | 54 |
| Squibb Separating Funnel - Glass Key Stopcock, Plastic Stopper<br>.....                           | 52 |
| Squibb Separating Funnel - PTFE Key Stopcock, Plastic Stopper<br>.....                            | 53 |
| Squibb Separating Funnel - Rotaflow Stopcock, Plastic Stopper<br>.....                            | 53 |

**G**

|                           |    |
|---------------------------|----|
| Gas Jars .....            | 26 |
| Gas Preparation Kit ..... | 65 |

**Graduated Cylinders - Glass**

|                                                                                       |    |
|---------------------------------------------------------------------------------------|----|
| Class-A ASTM with Hexagonal Base and Guard .....                                      | 43 |
| Class-A ASTM with Hexagonal Base and Guard<br>with Individual Works Certificate ..... | 85 |
| Class-A with Hexagonal Base .....                                                     | 40 |
| Class-A with Hexagonal Base with Individual<br>Works Certificate .....                | 84 |
| Class-A with Hexagonal Base USP with Individual<br>Works Certificate .....            | 84 |
| Class-A with Hexagonal Base and Stoppers .....                                        | 42 |
| Class-A with Hexagonal Base and Stoppers<br>with Individual Works Certificate .....   | 85 |
| Class-A with Round Base .....                                                         | 39 |
| Class-A with Round Base with Individual Works Certificate .....                       | 83 |
| Class-A with Round Base and Stoppers .....                                            | 41 |
| Class-B ASTM with Hexagonal Base and Guard .....                                      | 43 |
| Class-B with Hexagonal Base .....                                                     | 40 |
| Class-B with Hexagonal Base and Guard .....                                           | 41 |
| Class-B with Hexagonal Base and Stopper .....                                         | 42 |
| Class-B with Round Base .....                                                         | 39 |
| Class-B with Round Base and Stoppers .....                                            | 41 |

**H**

|                                                   |    |
|---------------------------------------------------|----|
| Hydro Distillation Apparatus .....                | 80 |
| Hydro Distillation Apparatus - Secure Joint ..... | 73 |

**K**

|                                       |    |
|---------------------------------------|----|
| Kipp's Gas Generator .....            | 65 |
| Kjeldahl Distillation Apparatus ..... | 79 |

## M

|                                      |    |
|--------------------------------------|----|
| Markham Distillation Apparatus ..... | 80 |
| Micro Distillation Apparatus .....   | 78 |

## N

|                         |    |
|-------------------------|----|
| Nessler Cylinders ..... | 44 |
|-------------------------|----|

## P

|                    |    |
|--------------------|----|
| Petri Dishes ..... | 14 |
|--------------------|----|

### Pipettes

|                                                                                   |    |
|-----------------------------------------------------------------------------------|----|
| Glass, Class-A .....                                                              | 45 |
| Glass, Class-A with Individuals Works Certificate .....                           | 88 |
| Pipettes, Glass, Class-B .....                                                    | 45 |
| Serological Pipettes - Glass, Class-A .....                                       | 46 |
| Serological Pipettes - Glass, Class-A<br>with Individuals Works Certificate ..... | 88 |
| Seological Pipettes - Glass, Class-B .....                                        | 46 |
| Plastic Screw Caps .....                                                          | 73 |
| Pycnometers .....                                                                 | 12 |

## R

|                      |    |
|----------------------|----|
| Rubber O Rings ..... | 73 |
|----------------------|----|

## S

|                                                                   |    |
|-------------------------------------------------------------------|----|
| Set 10 BU Organic Chemistry Kit .....                             | 74 |
| Set 27 BU/M Organic Chemistry Kit .....                           | 75 |
| Set 27 BU Organic Chemistry Kit - Secure Joint .....              | 76 |
| Set 27 X/M - Organic Chemistry Kit .....                          | 75 |
| Set 29 BU Organic Chemistry Kit .....                             | 74 |
| Set 34 BU Organic Chemistry Kit .....                             | 76 |
| Set 34 BU Organic Chemistry Kit in Grathells® Tray with Lid ..... | 77 |
| Set 46 BU Organic Chemistry Kit .....                             | 77 |
| Simple Distillation Apparatus .....                               | 79 |
| Simple Distillation Apparatus - Secure Joint .....                | 72 |
| Simple Organic Chemistry Set - in Grathells® Tray with Lid .....  | 78 |
| Simple Organic Chemistry Set - Micro .....                        | 78 |
| Simple Vacuum Distillation Apparatus - Secure Joint .....         | 72 |
| Single Cones .....                                                | 60 |
| Single Sockets .....                                              | 59 |
| Soxhlet Extraction Apparatus .....                                | 66 |
| Standard Distillation Apparatus .....                             | 79 |
| Still Heads - Plain with Thermometer Socket .....                 | 62 |
| Still Heads - Screw Thread .....                                  | 71 |
| Stopcock with standard ground joints .....                        | 14 |

### Stoppers

|                                                  |    |
|--------------------------------------------------|----|
| Flat Head .....                                  | 64 |
| Stoppers - Hexagonal .....                       | 64 |
| Stoppers - Hollow, Hexagonal with Flat End ..... | 63 |
| Stoppers - Penny Head .....                      | 64 |

## T

### Test Tubes

|                                   |    |
|-----------------------------------|----|
| Flat Bottom with Screw Cap .....  | 67 |
| Light Wall with Rim .....         | 67 |
| Rimless .....                     | 67 |
| Round Bottom with Screw Cap ..... | 67 |

## V

|                                     |    |
|-------------------------------------|----|
| Vigrex Columns - Screw Thread ..... | 70 |
|-------------------------------------|----|

## METALWARE

## B

### Base Retort Stands

|                                                                                 |     |
|---------------------------------------------------------------------------------|-----|
| 'A' Shaped Retort Stand Sets .....                                              | 117 |
| 'A' Shape Retort Stand Bases - Cast Iron with Leveling Screws .....             | 112 |
| 'A' Shape Retort Stand Bases - Cast Iron with Leveling Screws,<br>Nesting ..... | 112 |
| Deluxe Retort Stand Bases - Sheet Metal .....                                   | 112 |
| 'H' Shape Retort Stand Base - Metal .....                                       | 113 |
| Premium Retort Stand Bases - Cast Iron .....                                    | 112 |
| Retort Stand Bases - Triangular .....                                           | 112 |
| Retort Stand Sets - Cast Iron .....                                             | 117 |

### Bossheads

|                                           |     |
|-------------------------------------------|-----|
| Double Bossheads - Large .....            | 114 |
| Double Bosshead - Swivel .....            | 114 |
| Premium Bosshead - Die Cast .....         | 95  |
| Premium Bosshead - Heavy Duty .....       | 95  |
| Premium Bosshead - Metal .....            | 95  |
| Premium Bosshead - Square Head .....      | 95  |
| Premium Bosshead - Universal Swivel ..... | 95  |

|                      |        |
|----------------------|--------|
| Bunsen Burners ..... | 96-101 |
|----------------------|--------|

## C

### Clamps

|                                                                                       |     |
|---------------------------------------------------------------------------------------|-----|
| Burette Clamps .....                                                                  | 109 |
| Burette Clamps - Die Casted - Economy .....                                           | 109 |
| Burette Clamp with Bosshead - PVC Coated Round Jaws .....                             | 115 |
| Multipurpose Clamp Holder .....                                                       | 114 |
| Premium Clamp - Four Prong, Cork Lined .....                                          | 91  |
| Premium Clamp - Four Prong, Cork Lined, with Bosshead .....                           | 91  |
| Premium Clamp - Four Prong, Vinyl Coated .....                                        | 91  |
| Premium Clamp - Four Prong, Vinyl Coated, with Bosshead .....                         | 91  |
| Premium Clamp - Three Prong, Cork Lined .....                                         | 94  |
| Premium Clamp - Three Prong, Cork Lined, with Bosshead .....                          | 94  |
| Premium Clamp - Three Prong, Vinyl Coated .....                                       | 90  |
| Premium Clamp - Three Prong, Vinyl Coated .....                                       | 93  |
| Premium Clamp - Three Prong, Vinyl Coated, Dual Adjustment .....                      | 93  |
| Premium Clamp - Three Prong, Vinyl Coated, with Bosshead ..                           | 90  |
| Premium Clamp - Three Prong, Vinyl Coated, with Bosshead ..                           | 93  |
| Premium Clamp - Three Prong, Vinyl Coated, with Bosshead and<br>Dual Adjustment ..... | 93  |
| Premium Clamp - Two Prong, Cork Lined .....                                           | 94  |

GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

|                                                                                 |     |
|---------------------------------------------------------------------------------|-----|
| Premium Clamp - Two Prong, Cork Lined, with Bosshead .....                      | 94  |
| Premium Clamp - Two Prong, Vinyl Coated.....                                    | 92  |
| Premium Clamp - Two Prong, Vinyl Coated, Dual Adjustment .                      | 92  |
| Premium Clamp - Two Prong, Vinyl Coated, with Bosshead ....                     | 92  |
| Premium Clamp - Two Prong, Vinyl Coated, with Bosshead and Dual Adjustment..... | 92  |
| Premium Test Tube Clamp - Stainless Steel.....                                  | 124 |
| Premium Test Tube Clamp - Stainless Steel.....                                  | 124 |
| Suspension Clamps .....                                                         | 115 |
| Table Clamp.....                                                                | 116 |
| Thermometer Clamp with Bosshead .....                                           | 115 |
| Universal Clamp Holder - Jumbo .....                                            | 114 |
| Universal Clamp with Bosshead - PVC Coated Jaws .....                           | 115 |

**Clips**

|                                         |     |
|-----------------------------------------|-----|
| Hoffman's Clips.....                    | 102 |
| Hoffman's Clips - Heavy Duty .....      | 102 |
| Mohr's Clips .....                      | 102 |
| Burette Support Stand - Cast Iron ..... | 113 |
| Chattaway Spatula .....                 | 105 |
| Chemistry Lab Starter Set.....          | 120 |
| Clay Pipe Triangle .....                | 103 |

**L**

|                                         |     |
|-----------------------------------------|-----|
| Lab Essential Kit - Superior .....      | 120 |
| Laboratory Jacks - Stainless Steel..... | 108 |
| Laboratory Jacks - Standard.....        | 109 |
| Laboratory Tool Set - Basic.....        | 121 |
| Lab-TS .....                            | 121 |
| Lattice Clamps.....                     | 114 |

**M**

|                                |     |
|--------------------------------|-----|
| Metalware Set.....             | 119 |
| Metalware Set - Deluxe .....   | 118 |
| Metalware Set - Premium .....  | 118 |
| Metalware Set - Superior ..... | 119 |
| Measuring Scoops.....          | 26  |

**P**

|                                                         |     |
|---------------------------------------------------------|-----|
| Premium Test Tube Rack - Stainless Steel, 24 tubes..... | 124 |
|---------------------------------------------------------|-----|

**R****Rods**

|                           |     |
|---------------------------|-----|
| Aluminium Rods .....      | 113 |
| Stainless Steel Rods..... | 113 |
| Steel Rods .....          | 113 |

**Retort Stand Sets**

|                                    |     |
|------------------------------------|-----|
| Retort Stand Sets - Superior ..... | 117 |
| Retort Stand Sets - Value .....    | 117 |

**S****Spatulas**

|                                            |     |
|--------------------------------------------|-----|
| Micro Spatula Scoop - Stainless Steel..... | 104 |
| Micro Spatula Spoon - Stainless Steel..... | 104 |
| Palette Knife Spatula.....                 | 105 |

|                                                              |     |
|--------------------------------------------------------------|-----|
| Premium Double Sided Micro Spatula-Curved, Stainless Steel   | 107 |
| Premium Double Sided Micro Spatula Scoops - Stainless Steel  | 107 |
| Premium Double Sided Micro Spatulas - Stainless Steel .....  | 106 |
| Premium Double Sided Spatula-Wide Handle, Stainless Steel    | 106 |
| Semi Micro Arrowhead Spatula - Stainless Steel Dual End..... | 105 |
| Premium Nuffield Spatula - Wide Handle, Stainless Steel .... | 106 |
| Spatula & Scoops - Teflon Coated.....                        | 105 |
| Spatula Spoon - Stainless Steel, Wide Handle.....            | 104 |
| Spatula - Stainless Steel, Wide Handle, Dual End .....       | 104 |
| Spirit Lamp - Metal.....                                     | 103 |

**Spoons**

|                                                                                |     |
|--------------------------------------------------------------------------------|-----|
| Premium Combustion Spoon - Stainless Steel .....                               | 107 |
| Premium Double Sided Spoons-Wide Handle, Stainless Steel                       | 106 |
| Premium Fashion Spoon - Round Bottom, Stainless Steel with Wooden Handle ..... | 107 |
| Premium Micro Spoon Spatula - Stainless Steel.....                             | 107 |
| Premium Sodium Scoop - Brass with Wooden Handle .....                          | 107 |
| Scoop - Stainless Steel .....                                                  | 104 |
| Trowel Scoop - Stainless Steel, Dual End .....                                 | 105 |

**Stands**

|                                                     |     |
|-----------------------------------------------------|-----|
| Stands for Alcohol Burner .....                     | 110 |
| Stands for Alcohol Burner - Stainless Steel.....    | 110 |
| Tripod Stand - Circular.....                        | 110 |
| Tripod Stand - Circular with Removable Legs .....   | 111 |
| Tripod Stands - Circular .....                      | 110 |
| Tripod Stand - Triangular.....                      | 111 |
| Tripod Stand - Triangular with Removable Legs ..... | 111 |

**Support Rings**

|                                                         |     |
|---------------------------------------------------------|-----|
| Closed Support Rings .....                              | 116 |
| Closed Support Rings with Bossheads .....               | 116 |
| Closed Support Rings with Bossheads - Zinc Plated ..... | 116 |

**T**

|                                       |     |
|---------------------------------------|-----|
| Test Tube Holder - Brass.....         | 125 |
| Test Tube Holder - Vinyl coated ..... | 125 |

**Test Tube Racks - Epoxy Coated Steel Wire**

|                |     |
|----------------|-----|
| 20 Tubes ..... | 126 |
| 24 Tubes ..... | 126 |
| 36 Tubes ..... | 127 |
| 40 Tubes ..... | 126 |
| 40 Tubes ..... | 126 |
| 72 Tubes ..... | 127 |
| 80 Tubes ..... | 127 |
| 80 Tubes ..... | 127 |

**Tongs**

|                                                          |     |
|----------------------------------------------------------|-----|
| Beaker Tongs - Vinyl Coated, 3 Prong .....               | 125 |
| Crucible Tongs - General Use .....                       | 125 |
| Crucible Tongs - Straight, Stainless Steel .....         | 123 |
| Flask Tongs - Cork Lined, Heavy Duty .....               | 125 |
| Flask Tongs - Cork Lined, Stainless Steel.....           | 125 |
| Premium Beaker Tongs - Stainless Steel .....             | 123 |
| Premium Crucible Tongs - Stainless Steel, with Bow ..... | 123 |
| Premium Mercury Tongs - Stainless Steel .....            | 123 |
| Premium Flask Tongs - Stainless Steel .....              | 124 |
| Premium Flask Tongs - Stainless Steel, Angular.....      | 124 |

**U**

Ultimate Lab Starter Kit ..... 122

## W

### Wire Gauze

|                                           |     |
|-------------------------------------------|-----|
| Iron .....                                | 103 |
| Iron with Ceramic Center .....            | 103 |
| Stainless Steel.....                      | 103 |
| Stainless Steel with Ceramic Center ..... | 103 |

## PLASTICWARE

## B

### Beakers

|                                                 |     |
|-------------------------------------------------|-----|
| Polypropylene, Low Form .....                   | 130 |
| Beaker - Polypropylene, Low Form, Printed ..... | 130 |
| Beaker - TPX, Low Form, Printed .....           | 130 |
| Buchner Funnels - Polypropylene.....            | 143 |
| Burette Clamps - Polypropylene .....            | 133 |
| Burettes - Acrylic.....                         | 133 |

## C

|                                            |     |
|--------------------------------------------|-----|
| Centrifuge Tube Storage Box with Lid ..... | 134 |
| Conical Flasks - Polypropylene .....       | 138 |

## C (cont.)

|                                                              |     |
|--------------------------------------------------------------|-----|
| Conical Measures - Polypropylene, Hexagonal Base.....        | 138 |
| Conical Centrifuge Tubes - Polypropylene .....               | 134 |
| Conical Centrifuge Tubes - Polypropylene, self Standing..... | 134 |

## D

|                                                     |     |
|-----------------------------------------------------|-----|
| Desicators - Plain .....                            | 140 |
| Desicators - Vacuum .....                           | 140 |
| Desicator - Vacuum, Transparent.....                | 140 |
| Draining Basket .....                               | 149 |
| Dropping Bottles - LDPE, Euro Design .....          | 132 |
| Dual Sided Microcentrifuge Tube Rack with Lid ..... | 135 |

## F

### Flasks

|                                          |     |
|------------------------------------------|-----|
| Filtering Flasks - Polypropylene .....   | 138 |
| Volumetric Flasks - Polypropylene .....  | 142 |
| Floating Microcentrifuge Tube Rack ..... | 135 |

### Funnels

|                                     |     |
|-------------------------------------|-----|
| Funnel Holders - Polypropylene..... | 143 |
| Filter Funnels - Polyethylene ..... | 142 |
| Funnels - HDPE.....                 | 142 |

## G

|                                     |     |
|-------------------------------------|-----|
| Graduated Cylinder Protectors ..... | 137 |
|-------------------------------------|-----|

|                                                                    |     |
|--------------------------------------------------------------------|-----|
| Graduated Cylinders - Polypropylene, Class-B, Hexagonal Base ..... | 136 |
| Graduated Cylinders - TPX, Class-A, Hexagonal Base.....            | 137 |
| Graduated Cylinders - TPX, Class-B, Hexagonal Base .....           | 136 |

## J

|                   |     |
|-------------------|-----|
| Joint Clips ..... | 151 |
|-------------------|-----|

## K

|                                             |     |
|---------------------------------------------|-----|
| Kipp's Gas Generators - Polypropylene ..... | 144 |
|---------------------------------------------|-----|

## M

|                                                             |     |
|-------------------------------------------------------------|-----|
| Magnetic Stir Bars .....                                    | 150 |
| Magnet Retriever .....                                      | 151 |
| Measuring Jug - Polypropylene, Printed Graduations .....    | 139 |
| Measuring Jugs - Polypropylene.....                         | 139 |
| Measuring Jugs - Polypropylene, EURO Design .....           | 139 |
| Measuring Jugs - TPX, Printed Graduations .....             | 139 |
| Microcentrifuge Tubes .....                                 | 134 |
| Microcentrifuge Tube Rack - Polypropylene, Folding .....    | 149 |
| Microcentrifuge Tube Racks - Polycarbonate, Stackable ..... | 149 |
| Microcentrifuge Tube Storage Box with Lid .....             | 135 |

## P

|                                           |     |
|-------------------------------------------|-----|
| PCR Tube Rack with Lid .....              | 135 |
| PCR Tubes.....                            | 135 |
| Petri Dishes - Disposable, Sterile .....  | 141 |
| Petri Dishes - Reusable, Non Sterile..... | 141 |
| Petri Dish Racks .....                    | 141 |
| Pipette Stand - Rotary .....              | 145 |
| Pipette Stand - Vertical .....            | 145 |
| Pneumatic Troughs - Polypropylene .....   | 144 |

## R

|                                     |     |
|-------------------------------------|-----|
| Reagent Botles - Polypropylene..... | 131 |
| Rotatable Multi Tube Rack.....      | 149 |

## S

|                                                |     |
|------------------------------------------------|-----|
| Separating Funnel - Polypropylene.....         | 143 |
| Separatory Funnel Holder - Polypropylene ..... | 143 |
| Stoppers - Polypropylene .....                 | 151 |

## T

|                                                 |     |
|-------------------------------------------------|-----|
| Test Tube Baskets .....                         | 149 |
| Test Tube Peg Racks .....                       | 148 |
| Test Tube Rack - Polypropylene .....            | 146 |
| Test Tube Rack - Round, Polypropylene .....     | 146 |
| Test Tube Racks.....                            | 147 |
| Test Tube Racks - 3-Tier, White.....            | 147 |
| Test Tube Racks - Polypropylene, Blue.....      | 146 |
| Test Tube Racks - Polypropylene, Moulded .....  | 148 |
| Test Tube Racks - Polypropylene, Stackable..... | 148 |

GLASSWARE

METALWARE

PLASTICWARE

WOODWARE

EQUIPMENT

LAB SAFETY

RUBBERWARE

CERAMICS

|                                                      |     |
|------------------------------------------------------|-----|
| Test Tube Racks - Slant Racks.....                   | 148 |
| Tubing Connectors - Polypropylene, Cross .....       | 150 |
| Tubing Connectors - Polypropylene, L Shaped .....    | 150 |
| Tubing Connectors - Polypropylene, T & Y Shaped..... | 150 |

## U

|                                 |     |
|---------------------------------|-----|
| Universal Multi-Tube Rack ..... | 148 |
|---------------------------------|-----|

## W

### Wash bottles

|                                        |     |
|----------------------------------------|-----|
| Premium, Labeled, Wide Neck .....      | 132 |
| Premium, Narrow Neck .....             | 131 |
| Premium, Narrow Neck, Colored Top..... | 131 |
| Premium, Wide Neck, Colored Top.....   | 132 |
| Weighing Scoops - Polypropylene .....  | 144 |

## WOODWARE

## B

|                              |     |
|------------------------------|-----|
| Bottles Stands.....          | 153 |
| Burette Rack - Vertical..... | 153 |
| Burette Stands .....         | 153 |

## C

|                                           |     |
|-------------------------------------------|-----|
| Combination Funnel and Burette Stand..... | 155 |
|-------------------------------------------|-----|

## F

|                         |     |
|-------------------------|-----|
| Funnel Stand Clamp..... | 154 |
| Funnel Stands .....     | 154 |

## P

|                                 |     |
|---------------------------------|-----|
| Pipette Rack - Horizontal.....  | 153 |
| Pipette Stands - Vertical ..... | 154 |

## S

|                              |     |
|------------------------------|-----|
| Separating Funnel Stand..... | 154 |
|------------------------------|-----|

## T

|                                                 |     |
|-------------------------------------------------|-----|
| Test Tube Rack - 7 Place.....                   | 156 |
| Test Tube Rack - 12 Place.....                  | 156 |
| Test Tube Rack - 12 Place.....                  | 157 |
| Test Tube Rack - 12 Place.....                  | 157 |
| Test Tube Rack - 15 Place, Solid Block .....    | 157 |
| Test Tube Racks.....                            | 155 |
| Test Tube Racks - 4 Piece .....                 | 155 |
| Test Tube Racks with Drying Pegs.....           | 155 |
| Test Tube Rack with Drying Pegs - 6 Place.....  | 156 |
| Test Tube Rack with Drying Pegs - 6 Place.....  | 157 |
| Test Tube Rack with Drying Pegs - 12 Place..... | 156 |

|                                                         |     |
|---------------------------------------------------------|-----|
| Test Tube Rack with Drying Pegs - 12 Place.....         | 156 |
| Test Tube Rack with Drying Pegs - 13 Place, 2 Tier..... | 156 |

## W

|                                           |     |
|-------------------------------------------|-----|
| Wall Mounted Drying Rack - 90 Place ..... | 157 |
|-------------------------------------------|-----|

## EQUIPMENT

## H

|                                         |     |
|-----------------------------------------|-----|
| Heating Mantles .....                   | 161 |
| Heating Mantles - Digital Display ..... | 161 |

### Hot Plates

|                                                            |     |
|------------------------------------------------------------|-----|
| Ceramic Top, Analog Control .....                          | 159 |
| with Magnetic Stirrer - Aluminum Top, Analog Control ..... | 160 |
| with Magnetic Stirrer - Ceramic Top, Analog Control .....  | 159 |
| with Magnetic Stirrer - Ceramic Top, Digital Display ..... | 160 |

## W

|                                              |     |
|----------------------------------------------|-----|
| Unstirred Digital Water Bath .....           | 162 |
| Water Still with Stainless Steel Stand ..... | 162 |

## LAB SAFETY

## A

### Aprons

|                                 |     |
|---------------------------------|-----|
| Neoprene Aprons .....           | 168 |
| Rubber Apron - Colored.....     | 168 |
| Rubberised Aprons - Black ..... | 168 |
| Vinyl Aprons.....               | 168 |

## G

### Gloves

|                                        |     |
|----------------------------------------|-----|
| Cut Resistant Gloves .....             | 165 |
| Disposable Gloves - Latex .....        | 164 |
| Disposable Gloves - Nitrile.....       | 164 |
| Disposable Gloves - Polyethylene ..... | 164 |
| Disposable Gloves - Vinyl.....         | 164 |
| Gloves Leather.....                    | 165 |
| Neoprene Gloves .....                  | 165 |
| Rubber Gloves.....                     | 164 |
| Rubber Gloves.....                     | 165 |
| Rubber Gloves - Heavy Duty .....       | 165 |

## H

Hand Protectors ..... 166

## L

Lab Coats - Navy Blue..... 167

Lab Coats - White Polyester ..... 167

### Safety Goggles

Basic ..... 166

Basic, Indirect Vent..... 166

Clear ..... 166

Clear, Vented ..... 166

Clear With Adjustable Temples..... 166

## RUBBERWARE

## B

### Bulbs

Latex ..... 171

Pipette ..... 171

Rubber..... 171

Teat..... 171

## N

### Neoprene Stoppers

ASTM, One and Two Hole ..... 178

ASTM, Soild ..... 177

One and Two Hole..... 176

### Neoprene Stoppers (cont.)

Soild..... 175

## P

Pipette Filler - Rubber Bulb ..... 170

## R

Rubber Cones..... 178

### Rubber Stoppers

Rubber Stoppers ..... 172

Rubber Stoppers - One Hole..... 173

Rubber Stoppers - Solid ..... 172

Rubber Stoppers - Solid ..... 172

Rubber Stoppers - Solid ..... 173

Rubber Stopper - Two Hole ..... 173

Silicone Stoppers - Translucent..... 174

## T

### Tubing

Hand Made Vacuum Tubing..... 181

Neoprene Rubber Tubing - Synthetic ..... 180

Neoprene Tubing - For Burners ..... 170

Rubber Tubing ..... 179

Rubber Tubing - For Burner..... 170

Rubber Tubing - High Pressure ..... 180

Rubber Tubing - Nitrile..... 180

Silicone Tubing ..... 181

## CERAMICS

## B

Buchner Funnels - Porcelain ..... 185

## C

Caseroles..... 186

Crucible - Silica ..... 184

Crucibles - Porcelain, Squat..... 184

Crucibles - Porcelain, Tall..... 184

## E

Evaporating Dishes ..... 183

Evaporating Dishes ..... 183

Evaporating Dishes ..... 183

Evaporating Dishes- Silica..... 183

## G

Gooch Crucibles..... 185

## H

Hirsch Funnels ..... 185

## M

Mortar and Pestle ..... 186

Mortar and Pestle - Heavy pattern ..... 186

## S

Spotting Plates..... 186

|                  |     |
|------------------|-----|
| BUR-RCK 12.....  | 153 |
| CH0001C.....     | 157 |
| CH0002C-SPL..... | 157 |
| CH0004C.....     | 157 |
| CH0005C.....     | 157 |
| CH0070A.....     | 183 |
| CH0070B.....     | 183 |
| CH0070C.....     | 183 |
| CH0070D.....     | 183 |
| CH0070E.....     | 183 |
| CH0071A.....     | 183 |
| CH0071B.....     | 183 |
| CH0071C.....     | 183 |
| CH0071D.....     | 183 |
| CH0071E.....     | 183 |
| CH0071F.....     | 183 |
| CH0071G.....     | 183 |
| CH0071H.....     | 183 |
| CH0072A.....     | 183 |
| CH0072C.....     | 183 |
| CH0072D.....     | 183 |
| CH0072E.....     | 183 |
| CH0072F.....     | 183 |
| CH0072G.....     | 183 |
| CH0072H.....     | 183 |
| CH0072I.....     | 183 |
| CH0074A.....     | 184 |
| CH0074B.....     | 184 |
| CH0074D.....     | 184 |
| CH0074E.....     | 184 |
| CH0074F.....     | 184 |
| CH0074G.....     | 184 |
| CH0074H.....     | 184 |
| CH0076A.....     | 184 |
| CH0076B.....     | 184 |
| CH0076C.....     | 184 |
| CH0076D.....     | 184 |
| CH0076E.....     | 184 |
| CH0076F.....     | 184 |
| CH0076G.....     | 184 |
| CH0078A.....     | 183 |
| CH0078B.....     | 183 |
| CH0078C.....     | 183 |
| CH0078D.....     | 183 |
| CH0078E.....     | 183 |
| CH0078F.....     | 183 |
| CH0080A.....     | 184 |
| CH0080B.....     | 184 |
| CH0080C.....     | 184 |
| CH0080D.....     | 184 |
| CH0080E.....     | 184 |
| CH0080F.....     | 184 |
| CH0093QP-IP..... | 101 |
| CH0093QP-NG..... | 101 |
| CH0100D.....     | 170 |
| CH0100E.....     | 170 |
| CH0100F.....     | 170 |
| CH0100-LPG.....  | 170 |
| CH0137.....      | 130 |
| CH0137A.....     | 130 |
| CH0137B.....     | 130 |
| CH0137C.....     | 130 |
| CH0137D.....     | 130 |
| CH0137E.....     | 130 |
| CH0137F.....     | 130 |
| CH0137G.....     | 130 |
| CH0137H.....     | 130 |
| CH0138B/PR.....  | 130 |
| CH0138C/PR.....  | 130 |
| CH0138D/PR.....  | 130 |
| CH0138E/PR.....  | 130 |
| CH0138F/PR.....  | 130 |
| CH0167.....      | 153 |
| CH0167B.....     | 153 |
| CH0170A-N.....   | 131 |
| CH0170B-N.....   | 131 |
| CH0170C-N.....   | 131 |
| CH0170D-N.....   | 131 |
| CH0170D-W.....   | 131 |
| CH0171A.....     | 131 |
| CH0171B.....     | 131 |
| CH0171C.....     | 131 |
| CH0171D.....     | 131 |
| CH0172A-N.....   | 131 |
| CH0172A-W.....   | 131 |
| CH0172B-N.....   | 131 |
| CH0172B-W.....   | 131 |

|                 |     |
|-----------------|-----|
| CH0172C-N.....  | 131 |
| CH0172C-W.....  | 131 |
| CH0172D-N.....  | 131 |
| CH0172D-W.....  | 131 |
| CH0172E-N.....  | 131 |
| CH0172E-W.....  | 131 |
| CH0173A.....    | 132 |
| CH0173B.....    | 132 |
| CH0173C.....    | 132 |
| CH0179A.....    | 132 |
| CH0179B.....    | 132 |
| CH0179C.....    | 132 |
| CH0179D.....    | 132 |
| CH0179E.....    | 132 |
| CH0179F.....    | 132 |
| CH0179G.....    | 132 |
| CH0179H.....    | 132 |
| CH0179I.....    | 132 |
| CH0179J.....    | 132 |
| CH0179K.....    | 132 |
| CH0179L.....    | 132 |
| CH0181A.....    | 131 |
| CH0181B.....    | 131 |
| CH0181C.....    | 131 |
| CH0181D.....    | 131 |
| CH0182A.....    | 132 |
| CH0182B.....    | 132 |
| CH0182C.....    | 132 |
| CH0182D.....    | 132 |
| CH0182E.....    | 132 |
| CH0230A.....    | 171 |
| CH0230A-SI..... | 171 |
| CH0230B.....    | 171 |
| CH0230B-SI..... | 171 |
| CH0230C.....    | 171 |
| CH0230C-SI..... | 171 |
| CH0230D.....    | 171 |
| CH0230D-SI..... | 171 |
| CH0230E.....    | 171 |
| CH0230E-SI..... | 171 |
| CH0230F.....    | 171 |
| CH0230F-SI..... | 171 |
| CH0230G.....    | 171 |
| CH0230G-SI..... | 171 |
| CH0232A.....    | 171 |
| CH0232AA.....   | 171 |
| CH0232B.....    | 171 |
| CH0232BB.....   | 171 |
| CH0232C.....    | 171 |
| CH0232D.....    | 171 |
| CH0233A.....    | 133 |
| CH0233B.....    | 133 |
| CH0233C.....    | 133 |
| CH0256A.....    | 109 |
| CH0256B.....    | 109 |
| CH0258A.....    | 109 |
| CH0258B.....    | 109 |
| CH0262A.....    | 133 |
| CH0262B.....    | 133 |
| CH0264A.....    | 153 |
| CH0264B.....    | 153 |
| CH0270A.....    | 134 |
| CH0270B.....    | 134 |
| CH0270C.....    | 134 |
| CH0270D.....    | 134 |
| CH0271A.....    | 134 |
| CH0271B.....    | 134 |
| CH0272A.....    | 134 |
| CH0272B.....    | 134 |
| CH0284A.....    | 134 |
| CH0284B.....    | 134 |
| CH0284BB.....   | 134 |
| CH0285.....     | 135 |
| CH0286.....     | 135 |
| CH0287.....     | 135 |
| CH0288A.....    | 135 |
| CH0288B.....    | 135 |
| CH0289A.....    | 135 |
| CH0289B.....    | 135 |
| CH0290A.....    | 102 |
| CH0290B.....    | 102 |
| CH0298B.....    | 102 |
| CH0298C.....    | 102 |
| CH0298D.....    | 102 |
| CH0299A.....    | 102 |
| CH0299B.....    | 102 |
| CH0310H.....    | 172 |

|                     |     |
|---------------------|-----|
| CH0310H/1H.....     | 173 |
| CH0310H/2H.....     | 173 |
| CH0314A.....        | 172 |
| CH0314A01.....      | 172 |
| CH0314A01/1H.....   | 173 |
| CH0314A/1H.....     | 173 |
| CH0314A02.....      | 172 |
| CH0314A02/1H.....   | 173 |
| CH0314A/2H.....     | 173 |
| CH0314B.....        | 172 |
| CH0314B/1H.....     | 173 |
| CH0314B/2H.....     | 173 |
| CH0314C.....        | 172 |
| CH0314C/1H.....     | 173 |
| CH0314C/2H.....     | 173 |
| CH0314D.....        | 172 |
| CH0314D/1H.....     | 173 |
| CH0314D/2H.....     | 173 |
| CH0314E.....        | 172 |
| CH0314E/1H.....     | 173 |
| CH0314E/2H.....     | 173 |
| CH0314F.....        | 172 |
| CH0314F/1H.....     | 173 |
| CH0314F/2H.....     | 173 |
| CH0314G.....        | 172 |
| CH0314G/1H.....     | 173 |
| CH0314G/2H.....     | 173 |
| CH0314I.....        | 172 |
| CH0314I/1H.....     | 173 |
| CH0314I/2H.....     | 173 |
| CH0314J.....        | 172 |
| CH0314J/1H.....     | 173 |
| CH0314J/2H.....     | 173 |
| CH0314K.....        | 172 |
| CH0314K/1H.....     | 173 |
| CH0314K/2H.....     | 173 |
| CH0314L.....        | 172 |
| CH0314L/1H.....     | 173 |
| CH0314L/2H.....     | 173 |
| CH0314M.....        | 172 |
| CH0314M/1H.....     | 173 |
| CH0314M/2H.....     | 173 |
| CH0314N.....        | 172 |
| CH0314N01.....      | 172 |
| CH0314N01/1H.....   | 173 |
| CH0314N01/2H.....   | 173 |
| CH0314N/1H.....     | 173 |
| CH0314N/2H.....     | 173 |
| CH0314O.....        | 172 |
| CH0314O/1H.....     | 173 |
| CH0314O/2H.....     | 173 |
| CH0314P.....        | 172 |
| CH0314P/1H.....     | 173 |
| CH0314P/2H.....     | 173 |
| CH0314Q.....        | 172 |
| CH0314Q/1H.....     | 173 |
| CH0314Q/2H.....     | 173 |
| CH0314R.....        | 172 |
| CH0314R/1H.....     | 173 |
| CH0314R/2H.....     | 173 |
| CH0314S.....        | 172 |
| CH0314S/1H.....     | 173 |
| CH0314S/2H.....     | 173 |
| CH0314T.....        | 172 |
| CH0314T/1H.....     | 173 |
| CH0314T/2H.....     | 173 |
| CH0314U.....        | 172 |
| CH0314U/1H.....     | 173 |
| CH0314U/2H.....     | 173 |
| CH0314X.....        | 172 |
| CH0314X/1H.....     | 173 |
| CH0314X/2H.....     | 173 |
| CH0315.....         | 172 |
| CH0316.....         | 173 |
| CH0317.....         | 173 |
| CH0318.....         | 172 |
| CH0319A.....        | 174 |
| CH0319A-NUP.....    | 175 |
| CH0319A-NUP/1H..... | 176 |
| CH0319B.....        | 174 |
| CH0319B-NUP.....    | 175 |
| CH0319B-NUP/1H..... | 176 |
| CH0319C.....        | 174 |
| CH0319C-NUP.....    | 175 |
| CH0319C-NUP/1H..... | 176 |
| CH0319D.....        | 174 |
| CH0319D-NUP.....    | 175 |

|                     |     |
|---------------------|-----|
| CH0319D-NUP/1H..... | 176 |
| CH0319E.....        | 174 |
| CH0319E-NUP.....    | 175 |
| CH0319E-NUP/1H..... | 176 |
| CH0319F.....        | 174 |
| CH0319F-NUP.....    | 175 |
| CH0319F-NUP/1H..... | 176 |
| CH0319F-NUP/2H..... | 176 |
| CH0319G.....        | 174 |
| CH0319G-NUP.....    | 175 |
| CH0319G-NUP/1H..... | 176 |
| CH0319G-NUP/2H..... | 176 |
| CH0319H.....        | 174 |
| CH0319H-NUP.....    | 175 |
| CH0319H-NUP/1H..... | 176 |
| CH0319H-NUP/2H..... | 176 |
| CH0319I.....        | 174 |
| CH0319I-NUP.....    | 175 |
| CH0319I-NUP/1H..... | 176 |
| CH0319I-NUP/2H..... | 176 |
| CH0319J.....        | 174 |
| CH0319J-NUP.....    | 175 |
| CH0319J-NUP/1H..... | 176 |
| CH0319J-NUP/2H..... | 176 |
| CH0319K.....        | 174 |
| CH0319K-NUP.....    | 175 |
| CH0319K-NUP/1H..... | 176 |
| CH0319K-NUP/2H..... | 176 |
| CH0319L.....        | 174 |
| CH0319L-NUP.....    | 175 |
| CH0319L-NUP/1H..... | 176 |
| CH0319L-NUP/2H..... | 176 |
| CH0319M.....        | 174 |
| CH0319M-NUP.....    | 175 |
| CH0319M-NUP/1H..... | 176 |
| CH0319M-NUP/2H..... | 176 |
| CH0319N.....        | 174 |
| CH0319N-NUP.....    | 175 |
| CH0319N-NUP/1H..... | 176 |
| CH0319N-NUP/2H..... | 176 |
| CH0319O.....        | 174 |
| CH0319O-NUP.....    | 175 |
| CH0319O-NUP/1H..... | 176 |
| CH0319O-NUP/2H..... | 176 |
| CH0319P.....        | 174 |
| CH0319P-NUP.....    | 175 |
| CH0319P-NUP/1H..... | 176 |
| CH0319P-NUP/2H..... | 176 |
| CH0319Q.....        | 174 |
| CH0319Q-NUP.....    | 175 |
| CH0319Q-NUP/1H..... | 176 |
| CH0319Q-NUP/2H..... | 176 |
| CH0319R.....        | 174 |
| CH0319R-NUP.....    | 175 |
| CH0319R-NUP/1H..... | 176 |
| CH0319R-NUP/2H..... | 176 |
| CH0319S.....        | 174 |
| CH0319S-NUP.....    | 175 |
| CH0319S-NUP/1H..... | 176 |
| CH0319S-NUP/2H..... | 176 |
| CH0319T.....        | 174 |
| CH0319T-NUP.....    | 175 |
| CH0319T-NUP/1H..... | 176 |
| CH0319T-NUP/2H..... | 176 |
| CH0321A.....        | 177 |
| CH0321A/1H.....     | 178 |
| CH0321B.....        | 177 |
| CH0321B/1H.....     | 178 |
| CH0321C.....        | 177 |
| CH0321C/1H.....     | 178 |
| CH0321C/2H.....     | 178 |
| CH0321D.....        | 177 |
| CH0321D/1H.....     | 178 |
| CH0321D/2H.....     | 178 |
| CH0321E.....        | 177 |
| CH0321E/1H.....     | 178 |
| CH0321E/2H.....     | 178 |
| CH0321F.....        | 177 |
| CH0321F/1H.....     | 178 |
| CH0321F/2H.....     | 178 |
| CH0321G.....        | 177 |
| CH0321G/1H.....     | 178 |
| CH0321G/2H.....     | 178 |
| CH0321H.....        | 177 |
| CH0321H/1H.....     | 178 |
| CH0321H/2H.....     | 178 |
| CH0321I.....        | 177 |

|                 |     |
|-----------------|-----|
| CH0321I/1H..... | 178 |
| CH0321I/2H..... | 178 |
| CH0321J.....    | 177 |
| CH0321J/1H..... | 178 |
| CH0321J/2H..... | 178 |
| CH0321K.....    | 177 |
| CH0321K/1H..... | 178 |
| CH0321K/2H..... | 178 |
| CH0321L.....    | 177 |
| CH0321L/1H..... | 178 |
| CH0321L/2H..... | 178 |
| CH0321M.....    | 177 |
| CH0321M/1H..... | 178 |
| CH0321M/2H..... | 178 |
| CH0321N.....    | 177 |
| CH0321N/1H..... | 178 |
| CH0321N/2H..... | 178 |
| CH0321O.....    | 177 |
| CH0321O/1H..... | 178 |
| CH0321O/2H..... | 178 |
| CH0321P.....    | 177 |
| CH0321P/1H..... | 178 |
| CH0321P/2H..... | 178 |
| CH0321Q.....    | 177 |
| CH0321Q/1H..... | 178 |
| CH0321Q/2H..... | 178 |
| CH0321R.....    | 177 |
| CH0321R/1H..... | 178 |
| CH0321R/2H..... | 178 |
| CH0321S.....    | 177 |
| CH0321S/1H..... | 178 |
| CH0321S/2H..... | 178 |
| CH0321T.....    | 177 |
| CH0321T/1H..... | 178 |
| CH0321T/2H..... | 178 |
| CH0321U.....    | 177 |
| CH0321U/1H..... | 178 |
| CH0321U/2H..... | 178 |
| CH0321V.....    | 177 |
| CH0321V/1H..... | 178 |
| CH0321V/2H..... | 178 |
| CH0321W.....    | 177 |
| CH0321W/1H..... | 178 |
| CH0321W/2H..... | 178 |
| CH0321X.....    | 177 |
| CH0321X/1H..... | 178 |
| CH0321X/2H..... | 178 |
| CH0321Y.....    | 177 |
| CH0321Y/1H..... | 178 |
| CH0321Y/2H..... | 178 |
| CH0321Z.....    | 177 |
| CH0321Z/1H..... | 178 |
| CH0321Z/2H..... | 178 |
| CH0352A.....    | 137 |
| CH0352B.....    | 137 |
| CH0352C.....    | 137 |
| CH0352D.....    | 137 |
| CH0352E.....    | 137 |
| CH0352F.....    | 137 |
| CH0352G.....    | 137 |
| CH0354A/ES..... | 135 |
| CH0354B/ES..... | 135 |
| CH0354C/ES..... | 135 |
| CH0354D/ES..... | 135 |
| CH0354E/ES..... | 135 |
| CH0354F/ES..... | 135 |
| CH0354G/ES..... | 135 |
| CH0354I.....    | 135 |
| CH0354J.....    | 135 |
| CH0354K.....    | 135 |
| CH0354L.....    | 135 |
| CH0354M.....    | 135 |
| CH0354N.....    | 135 |
| CH0354O.....    | 135 |
| CH0354P.....    | 137 |
| CH0354Q.....    | 137 |
| CH0354R.....    | 137 |
| CH0354S.....    | 137 |
| CH0354T.....    | 137 |
| CH0354U.....    | 137 |
| CH0354V.....    | 137 |
| CH0355A.....    | 138 |
| CH0355B.....    | 138 |
| CH0355C.....    | 138 |
| CH0355D.....    | 138 |
| CH0355E.....    | 138 |
| CH0356A.....    | 139 |



|              |     |                 |     |                  |     |               |     |                  |     |
|--------------|-----|-----------------|-----|------------------|-----|---------------|-----|------------------|-----|
| CH0356B..... | 139 | CH0467B.....    | 143 | CH0622A.....     | 167 | CH0656A.....  | 113 | CH0710A.....     | 148 |
| CH0356C..... | 139 | CH0468.....     | 143 | CH0622A-NV.....  | 167 | CH0657A.....  | 113 | CH0710B.....     | 148 |
| CH0356D..... | 139 | CH0480A.....    | 143 | CH0622A-NV.....  | 167 | CH0657B.....  | 113 | CH0710C.....     | 148 |
| CH0356E..... | 139 | CH0480B.....    | 143 | CH0622B.....     | 167 | CH0657C.....  | 113 | CH0710D.....     | 148 |
| CH0356F..... | 139 | CH0480C.....    | 143 | CH0622B-NV.....  | 167 | CH0657D.....  | 113 | CH0710E.....     | 148 |
| CH0356G..... | 139 | CH0505A.....    | 144 | CH0622B-NV.....  | 167 | CH0658A.....  | 113 | CH0710F.....     | 148 |
| CH0356H..... | 139 | CH0505B.....    | 144 | CH0622C.....     | 167 | CH0658B.....  | 113 | CH0710G.....     | 148 |
| CH0356I..... | 139 | CH0505C.....    | 144 | CH0622C-NV.....  | 167 | CH0658C.....  | 113 | CH0711A.....     | 148 |
| CH0357A..... | 139 | CH0509A.....    | 186 | CH0622C.....     | 167 | CH0658D.....  | 113 | CH0711B1.....    | 148 |
| CH0357B..... | 139 | CH0509B.....    | 186 | CH0622C-NV.....  | 167 | CH0659A.....  | 113 | CH0711B2.....    | 148 |
| CH0357C..... | 139 | CH0509C.....    | 186 | CH0622C-NV.....  | 167 | CH0659B.....  | 113 | CH0711C.....     | 149 |
| CH0357D..... | 139 | CH0509D.....    | 185 | CH0622D.....     | 167 | CH0659C.....  | 113 | CH0711D.....     | 149 |
| CH0357E..... | 139 | CH0509E.....    | 185 | CH0622D-NV.....  | 167 | CH0659D.....  | 113 | CH0711E.....     | 149 |
| CH0357I..... | 139 | CH0509F.....    | 185 | CH0622D-NV.....  | 167 | CH0663D.....  | 114 | CH0712A.....     | 148 |
| CH0357J..... | 139 | CH0509G.....    | 185 | CH0622D-NV.....  | 167 | CH0664C.....  | 114 | CH0712B.....     | 148 |
| CH0357K..... | 139 | CH0518A.....    | 144 | CH0632A.....     | 104 | CH0664JM..... | 114 | CH00712B.....    | 183 |
| CH0357L..... | 139 | CH0518B.....    | 144 | CH0632A-TFC..... | 105 | CH0665A.....  | 114 | CH0712C.....     | 149 |
| CH0357M..... | 139 | CH0518C.....    | 144 | CH0632B.....     | 104 | CH0665B.....  | 114 | CH0713G-SUP..... | 124 |
| CH0364A..... | 140 | CH0521-SUP..... | 107 | CH0632C.....     | 104 | CH0666.....   | 114 | CH0727.....      | 125 |
| CH0364B..... | 140 | CH0525A.....    | 103 | CH0632F-SUP..... | 106 | CH0667.....   | 115 | CH0728A.....     | 125 |
| CH0364C..... | 140 | CH0525B.....    | 103 | CH0634A-SUP..... | 106 | CH0667A.....  | 115 | CH0728B.....     | 125 |
| CH0364D..... | 140 | CH0526A.....    | 103 | CH0634B-SUP..... | 106 | CH0668A.....  | 116 | CH0729A.....     | 125 |
| CH0366A..... | 140 | CH0526B.....    | 103 | CH0634C.....     | 104 | CH0668B.....  | 116 | CH0730C-SUP..... | 124 |
| CH0366B..... | 140 | CH0527A.....    | 103 | CH0634-TFC.....  | 105 | CH0668C.....  | 116 | CH0730D-SUP..... | 124 |
| CH0366C..... | 140 | CH0527B.....    | 103 | CH0635A.....     | 104 | CH0668D.....  | 116 | CH0732.....      | 125 |
| CH0366D..... | 140 | CH0527C.....    | 103 | CH0635A-TFC..... | 105 | CH0672A.....  | 116 | CH0733.....      | 125 |
| CH0367A..... | 140 | CH0527D.....    | 103 | CH0635B.....     | 104 | CH0672B.....  | 116 | CH0734E-SUP..... | 123 |
| CH0367B..... | 140 | CH0530A.....    | 103 | CH0635C.....     | 104 | CH0672C.....  | 116 | CH0734F-SUP..... | 123 |
| CH0367C..... | 140 | CH0530B.....    | 103 | CH0635D.....     | 104 | CH0674A.....  | 116 | CH0735E-SUP..... | 123 |
| CH0367D..... | 140 | CH0530C.....    | 103 | CH0635E.....     | 105 | CH0674B.....  | 116 | CH0735F-SUP..... | 123 |
| CH0372A..... | 141 | CH0535B.....    | 103 | CH0635E-TFC..... | 105 | CH0674C.....  | 116 | CH0736B-SUP..... | 123 |
| CH0372B..... | 141 | CH0535C.....    | 103 | CH0635H-SUP..... | 106 | CH0674D.....  | 116 | CH0736E-SUP..... | 107 |
| CH0372C..... | 141 | CH0540A.....    | 186 | CH0635I-SUP..... | 106 | CH0674E.....  | 116 | CH0737B-SUP..... | 107 |
| CH0372D..... | 141 | CH0540B.....    | 186 | CH0635J-SUP..... | 106 | CH0688A.....  | 115 | CH0738C-SUP..... | 124 |
| CH0372E..... | 141 | CH0540C.....    | 186 | CH0635K-SUP..... | 106 | CH0688B.....  | 115 | CH0738D-SUP..... | 123 |
| CH0374A..... | 141 | CH0540D.....    | 186 | CH0635O-SUP..... | 107 | CH0688F.....  | 116 | CH0740.....      | 125 |
| CH0374B..... | 141 | CH0540E.....    | 186 | CH0635P-SUP..... | 107 | CH0689A.....  | 154 | CH0761A.....     | 179 |
| CH0374C..... | 141 | CH0540F.....    | 186 | CH0635Q-SUP..... | 107 | CH0689B.....  | 154 | CH0761A01.....   | 179 |
| CH0376A..... | 141 | CH0540G.....    | 186 | CH0635R-SUP..... | 107 | CH0689C.....  | 154 | CH0761B.....     | 179 |
| CH0376B..... | 141 | CH0541A.....    | 186 | CH0635S-SUP..... | 107 | CH0689D.....  | 154 | CH0761C.....     | 179 |
| CH0421.....  | 178 | CH0541B.....    | 186 | CH0635-TFC.....  | 105 | CH0689E.....  | 154 | CH0761C01.....   | 179 |
| CH0422A..... | 138 | CH0541C.....    | 186 | CH0635T-SUP..... | 107 | CH0690.....   | 155 | CH0761D.....     | 179 |
| CH0422B..... | 138 | CH0541D.....    | 186 | CH0636A.....     | 105 | CH0694A.....  | 146 | CH0761E.....     | 179 |
| CH0431F..... | 138 | CH0541E.....    | 186 | CH0636A-TFC..... | 105 | CH0694B.....  | 146 | CH0761F.....     | 179 |
| CH0431G..... | 138 | CH0541F.....    | 186 | CH0636B.....     | 105 | CH0696.....   | 146 | CH0761G.....     | 179 |
| CH0431H..... | 138 | CH0541G.....    | 186 | CH0636B-TFC..... | 105 | CH0697.....   | 146 | CH0761H.....     | 179 |
| CH0452A..... | 142 | CH0542A.....    | 186 | CH0636C.....     | 105 | CH0702.....   | 155 | CH0761I.....     | 179 |
| CH0452B..... | 142 | CH0542B.....    | 186 | CH0637-SUP.....  | 106 | CH0702B.....  | 155 | CH0761J.....     | 179 |
| CH0452C..... | 142 | CH0580A.....    | 170 | CH0637-TFC.....  | 105 | CH0702C.....  | 155 | CH0761K.....     | 179 |
| CH0452D..... | 142 | CH0580D.....    | 170 | CH0639.....      | 105 | CH0703A.....  | 156 | CH0761L.....     | 179 |
| CH0452E..... | 142 | CH0589.....     | 153 | CH0640.....      | 105 | CH0703B.....  | 156 | CH0761M.....     | 179 |
| CH0452F..... | 142 | CH0590A.....    | 154 | CH0640B.....     | 105 | CH0703C.....  | 156 | CH0761N.....     | 179 |
| CH0454A..... | 142 | CH0590B.....    | 154 | CH0640C.....     | 105 | CH0704A.....  | 155 | CH0761O.....     | 179 |
| CH0454B..... | 142 | CH0594.....     | 141 | CH0642A.....     | 109 | CH0704B.....  | 155 | CH0762A.....     | 180 |
| CH0454C..... | 142 | CH0595.....     | 141 | CH0642B.....     | 109 | CH0705A.....  | 156 | CH0762B.....     | 180 |
| CH0454D..... | 142 | CH0612A.....    | 164 | CH0643.....      | 110 | CH0705B.....  | 156 | CH0762C.....     | 180 |
| CH0454E..... | 142 | CH0612B.....    | 164 | CH0643A.....     | 110 | CH0705C.....  | 156 | CH0762D.....     | 180 |
| CH0455A..... | 142 | CH0612C.....    | 164 | CH0643B.....     | 110 | CH0705D.....  | 157 | CH0762E.....     | 180 |
| CH0455B..... | 142 | CH0612D.....    | 164 | CH0643C.....     | 110 | CH0706A.....  | 157 | CH0762E01.....   | 180 |
| CH0461A..... | 144 | CH0612E.....    | 164 | CH0643D.....     | 110 | CH0706B.....  | 157 | CH0762F.....     | 180 |
| CH0461B..... | 144 | CH0612F.....    | 164 | CH0644A.....     | 110 | CH0707A.....  | 149 | CH0763A.....     | 181 |
| CH0461C..... | 144 | CH0612G.....    | 164 | CH0644B.....     | 110 | CH0707B.....  | 149 | CH0763B.....     | 181 |
| CH0461D..... | 144 | CH0612H.....    | 164 | CH0644-SPL.....  | 111 | CH0707C.....  | 149 | CH0764A.....     | 181 |
| CH0461E..... | 144 | CH0612I.....    | 164 | CH0645A-SPL..... | 111 | CH0708A.....  | 147 | CH0764B.....     | 181 |
| CH0461F..... | 144 | CH0612J.....    | 164 | CH0646A.....     | 111 | CH0708B.....  | 147 | CH0764B01.....   | 181 |
| CH0461G..... | 144 | CH0613A.....    | 164 | CH0646B.....     | 111 | CH0708C.....  | 147 | CH0764D.....     | 181 |
| CH0461H..... | 144 | CH0613B.....    | 164 | CH0648A.....     | 112 | CH0708D.....  | 147 | CH0764E.....     | 181 |
| CH0461I..... | 144 | CH0613C.....    | 164 | CH0648B.....     | 112 | CH0708E.....  | 147 | CH0764F.....     | 181 |
| CH0462A..... | 185 | CH0614.....     | 165 | CH0648C.....     | 112 | CH0708F.....  | 147 | CH0764G.....     | 181 |
| CH0462B..... | 185 | CH0616A.....    | 165 | CH0649A.....     | 112 | CH0708G.....  | 147 | CH0764I.....     | 181 |
| CH0462C..... | 185 | CH0616B.....    | 165 | CH0649B.....     | 112 | CH0708H.....  | 147 | CH0764I01.....   | 181 |
| CH0462D..... | 185 | CH0616C.....    | 165 | CH0649C.....     | 112 | CH0708I.....  | 147 | CH0764J.....     | 181 |
| CH0462E..... | 185 | CH0616D.....    | 165 | CH0653ARD1.....  | 117 | CH0708J.....  | 147 | CH0764J01.....   | 181 |
| CH0462F..... | 185 | CH0616HP.....   | 166 | CH0653BRD2.....  | 117 | CH0708K.....  | 147 | CH0764M.....     | 181 |
| CH0462G..... | 185 | CH0617A.....    | 165 | CH0653D1RD3..... | 117 | CH0708L.....  | 147 | CH0764M01.....   | 181 |
| CH0462H..... | 185 | CH0617B.....    | 165 | CH0653E1RD4..... | 117 | CH0708M.....  | 147 | CH0764O.....     | 181 |
| CH0462I..... | 185 | CH0617C.....    | 165 | CH0654F.....     | 112 | CH0709A.....  | 147 | CH0764O01.....   | 181 |
| CH0462J..... | 185 | CH0617D.....    | 165 | CH0654G.....     | 112 | CH0709B.....  | 147 | CH0764P.....     | 181 |
| CH0464A..... | 143 | CH0617E.....    | 165 | CH0654H.....     | 112 | CH0709C.....  | 147 | CH0764Q.....     | 181 |
| CH0464B..... | 143 | CH0617F.....    | 165 | CH0654I.....     | 112 | CH0709D.....  | 147 | CH0764R.....     | 181 |
| CH0464C..... | 143 | CH0617G.....    | 165 | CH0655A.....     | 112 | CH0709E.....  | 147 | CH0764S.....     | 181 |
| CH0465A..... | 185 | CH0618.....     | 165 | CH0655B.....     | 112 | CH0709F.....  | 147 | CH0764T.....     | 181 |
| CH0465B..... | 185 | CH0618A.....    | 165 | CH0655C.....     | 112 | CH0709G.....  | 147 | CH0764T01.....   | 181 |
| CH0465C..... | 185 | CH0618C.....    | 165 | CH0655E.....     | 112 | CH0709H.....  | 147 | CH0764V.....     | 181 |
| CH0465D..... | 185 | CH0619A.....    | 166 | CH0655F.....     | 112 | CH0709I.....  | 147 | CH0764W.....     | 181 |
| CH0465E..... | 185 | CH0619DV.....   | 166 | CH0655G.....     | 112 | CH0709J.....  | 147 | CH0764X.....     | 181 |
| CH0467.....  | 143 | CH0622A.....    | 167 | CH0656.....      | 113 | CH0709K.....  | 147 |                  |     |



# Index - Code | Alphabetical

GLASSWARE  
METALWARE  
PLASTICWARE  
WOODWARE  
EQUIPMENT  
LAB SAFETY  
RUBBERWARE  
CERAMICS

|             |     |
|-------------|-----|
| CH0764X01   | 181 |
| CH0765A     | 150 |
| CH0765B     | 150 |
| CH0765C     | 150 |
| CH0765D     | 150 |
| CH0765E     | 150 |
| CH0765EC    | 150 |
| CH0765F     | 150 |
| CH0765FC    | 150 |
| CH0765GC    | 150 |
| CH0765H     | 150 |
| CH0765I     | 150 |
| CH0765J     | 150 |
| CH0803A01   | 150 |
| CH0803B     | 150 |
| CH0803C     | 150 |
| CH0803D01   | 150 |
| CH0803E01   | 150 |
| CH0803F     | 150 |
| CH0803I     | 151 |
| CH0814      | 162 |
| CH0814A     | 162 |
| CH0869A     | 151 |
| CH0869B     | 151 |
| CH0869C     | 151 |
| CH0869D     | 151 |
| CH0869E     | 151 |
| CH0869F     | 151 |
| CH0873A     | 151 |
| CH0873B     | 151 |
| CH0873C     | 151 |
| CH0873D     | 151 |
| CH0873E     | 151 |
| CH0873F     | 151 |
| CH0971LP    | 98  |
| CH0971NG    | 98  |
| CH0972LP    | 98  |
| CH0972NG    | 98  |
| CH0973LP    | 98  |
| CH0973NG    | 98  |
| CH0974LP    | 98  |
| CH0974NG    | 98  |
| CH0991LP    | 99  |
| CH0991NG    | 99  |
| CH0992LP    | 99  |
| CH0992NG    | 99  |
| CH0993LP    | 99  |
| CH0993NG    | 99  |
| CH0994LP    | 100 |
| CH0994NG    | 100 |
| CH0995LP    | 100 |
| CH0995NG    | 100 |
| CH0996LP    | 100 |
| CH0996NG    | 100 |
| CH07336-SUP | 124 |
| CH182015A   | 126 |
| CH182015B   | 126 |
| CH182015C   | 126 |
| CH182015D   | 126 |
| CH182015E   | 127 |
| CH182015F   | 127 |
| CH182015G   | 127 |
| CH182015H   | 127 |
| CH1900A     | 82  |
| CH1900B     | 82  |
| CH1900C     | 82  |
| CH1900D     | 82  |
| CH1904A     | 82  |
| CH1904B     | 82  |
| CH1904C     | 82  |
| CH1904D     | 82  |
| CH1908A     | 82  |
| CH1908B     | 82  |
| CH1908C     | 82  |
| CH1908D     | 82  |
| CH1912A     | 83  |
| CH1912B     | 83  |
| CH1912C     | 83  |
| CH1912D     | 83  |
| CH1916A     | 83  |
| CH1916B     | 83  |
| CH1916C     | 83  |
| CH1916D     | 83  |
| CH1922A     | 83  |
| CH1922B     | 83  |
| CH1922C     | 83  |
| CH1922D     | 83  |


|         |    |
|---------|----|
| CH1922E | 83 |
| CH1922F | 83 |
| CH1922G | 83 |
| CH1922H | 83 |
| CH1922I | 83 |
| CH1926A | 84 |
| CH1926B | 84 |
| CH1926C | 84 |
| CH1926D | 84 |
| CH1926E | 84 |
| CH1926F | 84 |
| CH1926G | 84 |
| CH1926H | 84 |
| CH1926I | 84 |
| CH1930A | 84 |
| CH1930B | 84 |
| CH1930C | 84 |
| CH1930D | 84 |
| CH1930E | 84 |
| CH1930F | 84 |
| CH1930G | 84 |
| CH1930H | 84 |
| CH1930I | 84 |
| CH1932A | 85 |
| CH1932B | 85 |
| CH1932C | 85 |
| CH1932D | 85 |
| CH1932E | 85 |
| CH1932F | 85 |
| CH1932G | 85 |
| CH1932H | 85 |
| CH1932I | 85 |
| CH1936A | 85 |
| CH1936B | 85 |
| CH1936C | 85 |
| CH1936D | 85 |
| CH1936E | 85 |
| CH1936F | 85 |
| CH1936G | 85 |
| CH1936H | 85 |
| CH1936I | 85 |
| CH1952A | 86 |
| CH1952B | 86 |
| CH1952C | 86 |
| CH1952D | 86 |
| CH1952E | 86 |
| CH1952F | 86 |
| CH1952G | 86 |
| CH1952H | 86 |
| CH1952I | 86 |
| CH1952J | 86 |
| CH1952K | 86 |
| CH1960A | 86 |
| CH1960B | 86 |
| CH1960C | 86 |
| CH1960D | 86 |
| CH1960E | 86 |
| CH1960F | 86 |
| CH1960G | 86 |
| CH1960H | 86 |
| CH1960I | 86 |
| CH1960J | 86 |
| CH1960K | 86 |
| CH1964A | 87 |
| CH1964B | 87 |
| CH1964C | 87 |
| CH1964D | 87 |
| CH1964E | 87 |
| CH1964F | 87 |
| CH1964G | 87 |
| CH1964H | 87 |
| CH1964I | 87 |
| CH1964J | 87 |
| CH1964K | 87 |
| CH1972A | 87 |
| CH1972B | 87 |
| CH1972C | 87 |
| CH1972D | 87 |
| CH1972E | 87 |
| CH1972F | 87 |
| CH1972G | 87 |
| CH1972H | 87 |
| CH1972I | 87 |
| CH1972J | 87 |
| CH1990A | 88 |
| CH1990B | 88 |
| CH1990C | 88 |

|             |     |
|-------------|-----|
| CH1990D     | 88  |
| CH1990E     | 88  |
| CH1990F     | 88  |
| CH1990G     | 88  |
| CH1990H     | 88  |
| CH1990I     | 88  |
| CH1990J     | 88  |
| CH1990K     | 88  |
| CH1996A     | 88  |
| CH1996B     | 88  |
| CH1996C     | 88  |
| CH1996D     | 88  |
| ECH0139B/PR | 130 |
| ECH0139C/PR | 130 |
| ECH0139D/PR | 130 |
| ECH0139E/PR | 130 |
| ECH0139F/PR | 130 |
| G103        | 166 |
| G104        | 166 |
| G2612       | 166 |
| HMDT250     | 161 |
| HMDT500     | 161 |
| HMDT1000    | 161 |
| HMDT2000    | 161 |
| HMN25       | 161 |
| HMN50       | 161 |
| HMN100      | 161 |
| HMN200      | 161 |
| HTMST100-V1 | 160 |
| HTMST100-V2 | 160 |
| HTPLT-V1    | 159 |
| HTPLT-V2    | 159 |
| MGST1-V1    | 159 |
| MGST1-V2    | 159 |
| MGST2-V1    | 160 |
| MGST2-V2    | 160 |
| NEOCH01     | 180 |
| NEOCH02     | 180 |
| NEOCH03     | 180 |
| NEOCH04     | 180 |
| NEOCH05     | 180 |
| NEOCH06     | 180 |
| NEOCH07     | 180 |
| NEOCH08     | 180 |
| NITB01      | 180 |
| NITB02      | 180 |
| NITB03      | 180 |
| NITB04      | 180 |
| NITB05      | 180 |
| NITB06      | 180 |
| NITB07      | 180 |
| NITB08      | 180 |
| NUPAP-AD    | 168 |
| NUPAP-AD    | 168 |
| NUPAP-CH    | 168 |
| NUPAP-CH    | 168 |
| PRBH01      | 95  |
| PRBH02      | 95  |
| PRBH03      | 95  |
| PRBH04      | 95  |
| PRBH05      | 95  |
| PRCMP01     | 90  |
| PRCMP01BH   | 90  |
| PRCMP02     | 91  |
| PRCMP02BH   | 91  |
| PRCMP03     | 91  |
| PRCMP03BH   | 91  |
| PRCMP04     | 92  |
| PRCMP04BH   | 92  |
| PRCMP05     | 92  |
| PRCMP05BH   | 92  |
| PRCMP06     | 93  |
| PRCMP06BH   | 93  |
| PRCMP07     | 93  |
| PRCMP07BH   | 93  |
| PRCMP08     | 94  |
| PRCMP08BH   | 94  |
| PRCMP09     | 94  |
| PRCMP09BH   | 94  |
| RBAP-AD     | 168 |
| RBAP-AD     | 168 |
| RBAP-AD-BL  | 168 |
| RBAP-AD-BL  | 168 |
| RBAP-CH     | 168 |
| RBAP-CH     | 168 |
| RBAP-CH-RD  | 168 |
| RBAP-CH-RD  | 168 |

|            |     |
|------------|-----|
| RBAP-MD    | 168 |
| RBAP-MD    | 168 |
| SSJCT-10   | 108 |
| SSJCT-20   | 108 |
| SSJCT-30   | 108 |
| SSJCT-40   | 108 |
| SSJCT-50   | 108 |
| SSJCT-60   | 108 |
| STD-RD     | 117 |
| STD-RD1    | 117 |
| STDRDAB01  | 117 |
| STDRDAB02  | 117 |
| STDRDAB03  | 117 |
| STD-RDSP02 | 117 |
| STD-RDSP03 | 117 |
| STD-RDVL02 | 117 |
| STD-RDVL03 | 117 |
| VNAP-AD    | 168 |
| VNAP-AD    | 168 |
| VNAP-CH    | 168 |
| VNAP-CH    | 168 |
| WBTH-4H    | 162 |
| WBTH-6H    | 162 |



NOW SERVING OVER 110 COUNTRIES



[www.eiscoindustrial.com](http://www.eiscoindustrial.com)

# eisco

I N D U S T R I A L




Glassware



Metalware



Plasticware



Woodware



Equipment



Rubberware



Ceramics



EISCO SCIENTIFIC LLC

788 Old Dutch Rd

Victor, NY 14564

Phone: 585.563.7963

Web: [www.eiscoindustrial.com](http://www.eiscoindustrial.com)

Email: [sales@eiscoindustrial.com](mailto:sales@eiscoindustrial.com)

eisco

I N D U S T R I A L