

# CairMax Duo

A versatile foam air support surface with the ability to alternate between static or dynamic support, in response to the needs of the client. With a focus on support and versatility, the mattress functions as a static pressure relieving surface without the control unit, however when clinically required, can 'step up' to provide dynamic support with the control unit connected.

#### Versatility

Mattress and control unit can be purchased as a system or separately – place a CairMax Duo mattress on every bed frame to provide all clients with superior static support then add a control unit for increased therapy as needed.

### Simplicity

This 'plug and play' system is very simple to use and requires minimum training to operate, making it ideal for community or aged care use.

## Safety

No more concerns over correct pressure level settings with our AutoCair feature – automatic weight sensing technology that automatically adjusts the pressure settings to keep the client within a safe pressure range for effective therapy.

## Comfort

A durable foam core combined with alternating air cells embedded with castellated foam, and a memory foam topper, delivers comfort and support for clients up to 230 kilograms.

CairMax Duo	Standard	2000 x 880 x 180 mm	FAMCM-R01
System	King Single	2000 x 1050 x 180 mm	FAMCM-R01K
CairMax Duo,	Standard	2000 x 880 x 180 mm	FAMCM-M01
Mattress Only	King Single	2000 x 1050 x 180 mm	FAMCM-M01K
CairMax Duo, Control Unit Only			APMBL-CU01



