



Easy to Use



Clear Window



Infection Control

The EyePro® is constructed from a clear polyurethane film with a non woven fabric border. It is intended to aid in keeping the eye closed during medical procedures.

GENTLE ADHESIVE BORDER


Allows for easy handling and excellent conformity to the eye socket

TRANSPARENT WINDOW

To easily monitor eyelid closure

STERILE WRAPPED

Decreases the risk of cross contamination



Suitable for:

Patients undergoing anaesthesia

Patients being monitored in an unconscious state

Reorder	Pcs	Size	Pcs/Unit	Units/Ctn
EP010	2pcs	9.3cm x 3.7cm	50	10

