



Highly Absorbent



High Fluid Retention



Sideways Swelling

Impervia Neosorb Superabsorbent Dressings quickly absorb and binds large amounts of wound exudate into its 'heart' of superabsorbent polymers. The outer 'envelope' is ultrasonically bonded (adhesive-free construction) increasing conformability and minimising risk of possible skin reactions.

HIGH FLUID RETENTION


Can hold more than 30 times its dry weight in fluid*

DRY CONTACT SURFACE

Free draining wound contact surface encouraging exudate to be bound in the SAP absorbent core

SIDEWAYS SWELLING

Fluid uptake wicks outward, not upward, providing minimal change in the thickness of the dressing



Suitable for:

Heavily exuding wounds

Exudate management under compression therapy

Protection of clothing and linen

Reorder	Size	Pcs/Unit	Units/Ctn
NWC040	10cm x 10cm	20	20
NWC046	10cm x 20cm	20	20
NWC041	15cm x 15cm	20	20
NWC042	15cm x 20cm	20	20
NWC043	20cm x 20cm	20	20
NWC044	20cm x 30cm	20	10
NWC045	20cm x 45cm	5	20
NWC047	30cm x 45cm	5	10



*In-house testing