

Volumat Line of Administration Sets

for use with Agilia® Infusion System

Standard Sets (No secondary access)

Order Number	Product Code	Description	Priming Volume (mL)	Nominal Length (inch)	Drops (mL)	Tubing ID (inch)	Units / Case	Available
M46441360	VL ST10-O	Standard Set	25	112	20	0.118	30	Yes
M46441375	VL ST10-XO	Standard Set, Non-siliconized Spike	25	112	20	0.118	30	Yes
M46441960	VL ST02-O	Standard Set, Y-site	25	112	20	0.118	30	Yes
M46441985	VL ST02-XO	Standard Set, Y-site, Non-siliconized Spike	25	112	20	0.118	30	Yes

Primary Sets

Order Number	Product Code	Description	Priming Volume (mL)	Nominal Length (inch)	Drops (mL)	Tubing ID (inch)	Units / Case	Available
M46445660	VL PR42-11	Primary Set, 2 Y-sites	25	112	20	0.118	30	Yes
M46445690	VL PR42-12	Primary Set, 3 Y-sites	25	112	20	0.118	30	Yes
M46445780	VL PR44-12	Primary Set, 3 Y-sites	25	112	20	0.118	30	Yes

Filter Sets

Order Number	Product Code	Description	Priming Volume (mL)	Nominal Length (inch)	Drops (mL)	Tubing ID (inch)	Units / Case	Available
M46442460	VL PNO2-1	Standard Set, Y-site, 1.2-micron Filter, not made with PVC	25	114	20	0.118	30	Yes
M46444160	VL ON72-1	Primary Set, 2 Y-sites, 0.2-micron Filter, not made with PVC	26	114	20	0.118	30	Yes
M46444175	VL ON72-X1	Primary Set, Y-site, 0.2-micron Filter, Non-siliconized spike, not made with PVC	26	114	20	0.118	30	Yes
M46445860	VL PR72-11	Primary Set, 2 Y-sites, 0.2-micron Filter	27	112	20	0.118	30	Yes

Burette Sets

Order Number	Product Code	Description	Priming Volume (mL)	Nominal Length (inch)	Drops (mL)	Tubing ID (inch)	Units / Case	Available
M46445260	VL PA92-1	Burette Set, 20 drop factor	25	118	20	0.118	10	Yes
M46445270	VL PA92-2	Burette Set, 60 drop factor	20	118	60	0.118	10	Yes

Specialty Sets

Order Number	Product Code	Description	Priming Volume (mL)	Nominal Length (inch)	Drops (mL)	Tubing ID (inch)	Units / Case	Available
M46443450	VL SP62-1	Standard Set, Y-site, not made with PVC	25	112	20	0.118	30	Yes
M46445475	VL SP40-11	Primary Set, 2 Y-sites, not made with PVC, Reduced length	20	81	20	0.118	30	Yes
M46445670	VL PR42-21	Primary Set, 3 Y-sites	25	112	20	0.118	30	Special Order

Transfusion Sets

Order Number	Product Code	Description	Priming Volume (mL)	Nominal Length (inch)	Drops (mL)	Tubing ID (inch)	Units / Case	Available
M46442860	VL TROO-O	Transfusion-Single Spike	25	112	20	0.118	40	Special Order
M46443160	VL SP22-O	Transfusion set, Dual Spike	30	114	20	0.118	40	Yes

Secondary Sets

Order Number	Product Code	Description	Priming Volume (mL)	Nominal Length (inch)	Drops (mL)	Tubing ID (inch)	Units / Case	Available
M77460051	VL SLOO-XO	Secondary Set, Non-siliconized Spike	12	35	20	0.118	40	Yes
M77460052	VL SLOO-XP	Secondary Set, Non-siliconized Spike, not made with PVC	10.5	27	20	0.118	40	Yes

Refer to the Instructions for use for a complete list of warnings and precautions associated with the use of this product.

For more information contact your sales representative or call customer service at 888-386-1300 or visit www.agiliasystem.com

Copyright © 2021 Fresenius Kabi AG. All rights reserved. All trademarks shown are the property of their respective owners.
0102-VL-01-08/21

Fresenius Kabi USA
Three Corporate Drive
Lake Zurich, Illinois 60047
855-FK-INFUSION
855-354-6387
www.fresenius-kabi.us/agilia