



Medical Equipment  
Furniture & Uniforms

**Product Guide 2021**

[www.medeleg.com.au](http://www.medeleg.com.au)



**Brands** we are proud to **supply**



**Medeleq Pty Ltd was founded in 1980, by Michael and Denise Cutajar, to supply Medical Electronic Equipment to Hospitals and Allied Professionals in Australia.**

**Our beginnings were at Coolangatta, Qld and within 18 months were also located in Brisbane and Sydney. In 1987 the Queensland offices were consolidated into Medeleq's own premises at Nerang, and now in our larger premises and showroom at Arundel.**

We hold Quality Assurance Certification to ISO 9001: 2015 and in 1998 Business Queensland placed our company in the top 450 Privately Owned Companies in the state. Of paramount importance to Medeleq is the relationship between our clients and ourselves. We also place major importance on the ongoing training of clients in any product they purchase from us. We are constantly updating our own knowledge to enable us to support our business and products.

Medeleq has adopted a policy of supplying only those products which clearly demonstrate their quality and which are suitable for their intended purpose. We place particular emphasis on expertise, experience, capability, reliability and quality.

We believe that a high standard of quality in all our services is essential to the achievement of company objectives and the continued growth and success of the company. The staff at Medeleq have an extensive knowledge of the Australian Medical Electronic market. Our backgrounds are Medical, Nursing, Scientific, Biomechanics (Sports science), Technical and Business Management.



Individual Product **Brochures** available on request

## Contents

## Page

ADE Medical Scales	4
Baby Scales	6
Column Scales	8
Bariatric Scales	10
Chair & Wheel Chair Scales	11
Vital Signs Monitoring ( <i>Creative Medical</i> )	12
Accessories	14
Mobile Trolleys and Wall Mounts	15
Pulse Oximeters	16
Veterinary Products	18
Medical Assessment Products	20
Medical Suction	22
Portable Suction	24
Vac Sax	26
Silentia Privacy Screens	28
Francehospital Medical Trolleys	34
Francehospital Accessories	37
Stainless Steel Trolleys & Equipment	38
Wire shelving products (Chrome & Stainless)	40
Wire Baskets (Chrome & Stainless)	41
Patslide	42
Rehabilitation Stairs	43
Cherokee Medical Uniforms	44

# Scales & Measurement

These devices meet the highest quality standards, they are precise, reliable and practical in use. We combine the technical know-how of our ADE designers and engineers with the experience of doctors and nurses in their daily practice. Certified by the ISO 9001 standard.

**Quality you can trust.**

## ADE M320600

Premium Hospital Scale, with advanced features. This is a strong and sturdy scale that has an oversized platform.

- Calibratable
- Mother / Child function
- Hold / Auto-Off, Tare, BMI
- Battery Operation
- Capacity 250kg / 50g < 50 kg > 100g
- Dim: 400 x 350 x 44mm



Scales & Measurement

# Floor Scales



# ADE Scales



## Mechanical Floor Scale

### ADE M308800

- Reliable automatic dial
- Large weighing platform covered with non slip mat
- Clear and easy to read dial

Capacity 160kg/500g  
Dimensions 325 x 480 x 108mm



## Hospital Grade Floor Scale

### ADE M322600

- Strong and sturdy construction with oversize platform
- Turning display
- Mother/Child function
- Mains or battery operated

Capacity 250kg/100g Dimensions 400 x 350 x 44mm



## Electronic floor scale

### ADE BE1509

- Solid safety glass housing and silicon mat
- Auto-hold function
- Large LCD step-on operation
- Easy to clean

Capacity 180kg/100g  
Dimensions 302 x 302 x 18 mm



## Ultimate Floor Scale

### ADE M320600

- All the features of the Hospital Grade Scale
- Plus BMI, hold/auto, tare and calibratable

Capacity 250kg/50g <50kg> 100g Dimensions 400 x 350 x 44mm

# Baby Measurement



## Baby Toddler Scale

### MEDBAK

- Removable tray for toddler to stand
- Direct reading on backlit display
- Auto-off, battery operated

Capacity 20kg/5g



## Electronic Baby Hanging Scale

### M111600

- Electronic hanging type baby scale
- Large LCD, 15 mm
- Automatic zero setting and switch off
- Tare function
- Small and practical device made out of steel with integrated length measuring tape, hand - grip and hook for weighing trouser
- Battery operated

Capacity 20kg/20g



# Baby Scales



# ADE Scales

## Advanced Hospital Baby Scale

### M118600

- Contoured tray provides secure space for all infants
- Fine weight graduation range and auto-hold function offer precise and fast weighing results
- Easy to read backlit LCD display

Capacity 20kg/5g <10kg>10g



## Electronic Baby Scale with Integrated Length Measure

### M118600-01

- All the features of M118600
- The fully integrated digital length measuring device will take precise readings of the infant up to 80cm by 2mm

Capacity 20kg/5g <10kg>10g

## Baby Measure Accessories



## Baby Measure Mat

### ADE MZ10027-1

- Easy to clean, skin friendly, roll up, compact

Range 100-1000mm/5mm



## Baby Measure Rod

### ADE MZ10028

- Aluminium with foldable measuring tongues
- Support feet for table measurement

Range 100-1000mm/2mm



## Sliding Baby Scale

### ADE MZ10040

- Perfectly suited for mobile & stationary use
- Running foot stop guided by a central rail measures accurately to 1mm

Range 100-1000mm/1mm



## Circumference Measure Tape

### ADE MZ10021

- Auto rewind printed both sides, measures both circumference and length

Range 50-1500mm/1mm



## Disposable Measuring Tapes

### MED5521BMT

- Disposable tapes for circumference

Range to 59cm/1mm



## Baby Scale Carry Case

### ADE MZ10031

- Padded carry bag, suits M112600, MEDBAK, MEDBAY



## Electronic Baby Scale

### M112600

- Light-weight weighing scale ideal for GPs and home visiting nurses
- Clear 30 mm display and three key operation make it easy to use
- Operated by standard AA batteries and mains

Capacity 20kg/5g



## Mechanical Baby Scale

### M108800

- Quick to use baby and toddler weighing scale
- Safe and comfortable platform
- Clear weight indication by needle readout
- Manual zero setting
- Robust and simple to handle perfect for general use

Capacity 20kg/50g <10kg>100g

# Column Scales



## MED-SKL06B Ultrasonic Height & Weight with BMI

- Self-checking function
- Text-to-speech function
- 5" Touch Screen
- Ultrasonic Height measure
- Completely automatic height, weight measurement and body shape, display of countdown function.
- External Printer

Capacity 300kg/50g



## Electronic Medical Column Scale

FM650PR

- 310 x 310mm Platform
- 25 User memory
- Tare, Hold, Auto Switch Off
- Mains charger / Internal battery
- 4 Castors

Capacity 250kg/100g



## TCS Series Column Scale

MED HH

- Large Platform 440 x 460mm
- Adjustable Feet
- BMI Function
- Tare, Hold, Auto Switch Off
- Integrated Height Measure

Capacity 300kg/50g



# Column Scales & Height Measuring Stadiometers

## Round Dial Column Scale

ADE M306800

- 270 x 370mm Platform
- Integrated Height Measure

Capacity 150kg/500g



## Electronic Column Scale

ADE M320600-01

- Extra wide, flat and slip resistant platform
- Optional Height Measure
- Auto Zero & Switch Off
  - Tare, Hold, BMI, mother & child function
  - Turning display, low battery indication

Capacity 250kg/100g



## Telescopic Measure Rod

MZ10023-1

Wall or column scale mounted.

Range 60-210cm  
Graduation 1mm



## Electronic Measure Rod

MZ10032

Wall mounted.

Range 110-200cm  
Graduation 1mm



## Free Standing Stadiometer

MZ10038

Metal base for free standing.

Range 60-210cm  
Graduation 1mm



## Roll Up Height Measure

MZ10017

Wall mounted measure tape.

Range 0-220cm  
Graduation 1mm



# Bariatric Scales

## Bariatric Wide Platform Scale

**ADE M319660**

- Big platform, remote cable display
- BMI, hold/auto, tare and calibratable
- Mains or battery operated

Capacity 300kg/100g  
Dimensions 520 x 520 x 75mm



## Single Sided Handrail Scale

**ADE M319660-02**

- Large Platform 520 x 520mm
- BMI Function
- Tare, Hold Auto Switch Off
- Mains charger / Internal battery
- 2 Rear Castors

Capacity 300kg/100g

## Electronic Handrail Scale

**ADE M319660-01**

- Large Platform 520 x 520mm
- BMI Function
- Tare, Hold Auto Switch Off
- Mains charger / Internal battery
- 2 Rear Castors

Capacity 300kg/100g



# Chair & Wheelchair Scales

## HL Series Wide Chair Scale

**MED HL**

- Wide 530mm seat
- Lift up arm & footrests. Individual swivelling and braking wheels

Capacity 250kg/50g



## Wheelchair Scale

**MEDTCS-RWB**

- 1m x 1m platform with integrated 20cm ramps
- Fold up handrail
- Handles & wheels for moving

Capacity 500kg/200g Dimensions 1.2 X 1.4m



## Additional Scales

### Electronic Table Scale

**HW945-3**

- Stainless weighing plate 260 x 215mm with counting function

Capacity 3kg/0.1g



### Bed Weighing Scale

**M601620**

- Two independent weighbridges with built-in ramp

Capacity 500kg/100g



### Electronic Table Scale

**HW945-30**

- Stainless weighing plate 260 x 215mm with counting function

Capacity 30kg/1g

# Vital Signs Monitors



## Vital Signs Monitor PC900 PLUS

- 7" high resolution TFT with touch screen
- Up to 12,000 NIBP and 2,000 groups SpO<sub>2</sub> storage
- Trend Graphs • Auto, Manual, Stat
- NIBP measuring modes
- Built in rechargeable LI battery – 6 hour monitoring
- USB connection for data output and software upgrade
- SpO<sub>2</sub> NIBP standard

*Optional Infrared Ear Thermometer, ECG, Temp & Thermal Printer*

## Vital Signs Monitor PC900Pro SN

- 3.5" High resolution LCD display
- Tabletop design with user-friendly operation
- Up to 12,000 NIBP and 2,000 groups SpO<sub>2</sub> storage
- Trend Graphs
- Auto, Manual, Stat
- NIBP measuring modes
- Built in rechargeable LI battery – 6 hour monitoring
- USB connection for data output & software upgrade
- SpO<sub>2</sub> NIBP standard

*Optional ECG, Temp & Thermal Printer*



## Vital Signs Monitor PC900Pro SNET

As PC900Pro SN with ECG & Temp.



## Patient Monitor PC3000

- Lightweight and easy to operate
- 7" high resolution colour TFT display seven channels of wave forms – 5 channel of real time and two hours of trend screen
- Auto analysis of 20 types arrhythmia
- Manual measurement of S-T segment while ECG waveform frozen or automatic measurement while monitoring
- Large parameter storage capability

*Optional - Thermal Printer which is remote to main unit using USB connection Cable*

# Multi Parameter Monitors



## Bedside Monitor K10

- 10.4" high resolution colour TFT touch screen
- ECG/RESP/SpO<sub>2</sub>/NIBP/S-T/Temp/PR
- NIBP customized 5-phase multi-cycle setting
- 8-waveform simultaneous display for different views history data records reviewed in waveform, lists and trend graph
- Data export and software upgrade via USB Parameter box for extended functions

*Optional: 2-IBP, EtCO<sub>2</sub>, CMS, C/O*

## Bedside/ICU Monitor K12

- 12.1" high resolution colour TFT touch screen
- ECG/RESP/SpO<sub>2</sub>/NIBP/S-T/Temp/PR
- NIBP customized 5-phase multi-cycle setting 8-waveform simultaneous display for different views
- History data records reviewed in waveform, lists and trend graph
- Data export and software upgrade via USB Parameter box for extended functions

*Optional: 2-IBP, EtCO<sub>2</sub>, CMS, 12-lead ECG*



## The K Series features:



USB connection for data output and software upgrade.



## Bedside/ICU Monitor K15

- 15" high resolution colour TFT touch screen
- ECG/RESP/SpO<sub>2</sub>/NIBP/S-T/Temp/PR
- NIBP customized 5-phase multi-cycle setting 5 types of calculations: **Oxygenation, Renal Function, Ventilation, Medication and Hemodynamic**
- 8-waveform simultaneous display
- Multiple customized display
- History data records reviewed in waveform, lists and trend graph
- Data export and software upgrade via USB Parameter box for extended functions

*Optional: 2-IBP, EtCO<sub>2</sub>, CMS, 12-lead ECG*



### NIBP Cuffs

- 15020009 - Baby 6-11cm
- 15020010 - Infant 10-19cm
- 15020011 - Paed 18-26cm
- 15020013 - Adult 25-35cm
- 15020014 - Large Adult 32-47cm

- Creative NIBP Cuffs are durable and reusable
- Fits most brand monitors on today's market
- All connecting tubing has Creative Quick Fit Metal connectors making them interchangeable with each other



**Neonate Wrap Sensor**

- Reusable L-type supplied with 2 x disposable wraps
- USB connection

### ADCUFF & Bladder Bariatric

- 845-12BXBD - Bariatric 44-66cm
- Bariatric 1 tube, L/F.
- Burgundy colour. 44-66 cms
- With connector



**Disposable Wraps**  
3101-000026

- Wrap for Toe, L = 205mm
- Packet of 5




**Paediatric Sensor Clip**  
15040055

- Type for SP20 HH Oximeter



**Paediatric Soft Sensor**  
15040051

- Type for SP20 HH Oximeter



**Adult Sensor Clip**  
15040022

- Type for SP20 HH Oximeter



**Disposable Wraps**  
3101-000027


- Auxiliary Wrap, L = 200mm
- Packet of 5

Coriton Monitor Mounting:  
High quality Medical Mounting Solutions.  
**Good Health Support.**



### Patient Monitor Wall Mount with Basket

- WM001-1-BS004**
- Load Capacity 30kg
  - Standard channel length: 33cm
  - Vertical position lock with knob lock



### Mobile Roll Stand for Monitors

- MDQCMC**
- Heavy duty 5 star base with lockable wheels
  - Accessory basket and push handle standard
  - Fits all models



### Equipment Trolley with Drawer

- TR800**
- Heavy duty lockable wheels

### Laptop Trolley Adjustable Height

- TR800-100-5**
- Foot pedal with gas spring assisted control for variable height adjustment
  - 4" diameter casters – 4 locking
  - Integrated adjustable cable management





# Pulse Oximeters



# Pulse Oximeters



## SP-20 Handheld Pulse Oximeter

Applicable for Adult, Paediatric and Neonate use.

- Measuring SpO<sub>2</sub>, Pulse Rate
- 3.5 inch TFT display of P1%, plethysmograph & pulse bar
- Adult, pediatric and neonate measurement
- Lithium battery up to 18 hours
- 500 hours SpO<sub>2</sub> data storage
- PC software for data upload

NSW SOA PTM1804



Accurate, reliable and sensitive monitoring against low perfusion, motional situation and weak pulse signal

## Handheld Oximeter PC66B2

- Hi-Res Colour LCD Display
- SpO<sub>2</sub>, PR, PI, Plethysmogram & Pulse Bar
- Patented **MoveOxy™ Technology** for removal of artefact and anti-motion performance
- 3 x AAA Batteries
- Adjustable alarm limit, pitch tone and key tone level

Supplied with Adult and Paediatric clip type sensors



## Handheld Capnograph & Oximeter PC900B

- 3.5" Colour display with audio and visual alarm
- Monitoring SpO<sub>2</sub>, PR, EtCO<sub>2</sub>/CO<sub>2</sub> and RR
- Temp. Pressure Compensation and Zero Adjustment System Balance Gas compensation minimises influence against other gases
- BT/PS compensation available
- Applicable to Intubate. Non intubated patients
- Water filter specially designed for high humidity
- Patented **MoveOxy™ Technology**
- 24 hour graphic for trends

Compact Design Specially For Emergency, Transport & Outdoor Rescue



## Finger Sensor PC60E

- Drop resistant
- 4 Direction display
- OLED colour display of SpO<sub>2</sub>, PR, PI waveform and Pulse bar
- 30+ hours operation x 2 AAA Batteries
- Audio Visible Alarm
- Battery Indication and Low Battery Alarm



## Finger Tip Oximeter PC60B1

- Splash proof & drop resistant
- OLED colour display of SpO<sub>2</sub>, PR, PI waveform and Pulse bar
- 30+ hrs operation x 2 AAA Batteries
- Audio Visible Alarm
- Battery Indication & Low Battery Alarm



## Finger Sensor PC60E+

- with 15040014 Paediatric Sensor



## Finger Pulse Oximeter PC60F

- 4 direction display
- **MoveOxy™ Technology**,
- Waterproof and Drop Resistant.
- Wireless Connection

## Smart Watch CMSW

- A1-ECG Analysis, BP, HR & Alarms
- Sleep Monitoring
- Waterproof & Multisport modes
- Personal reminder.
- Battery life 7 days



## Smart Band CMSB

- A1-ECG Analysis, BP, HR & Alarms
- Sleep Monitoring
- Waterproof & Multisport modes
- Personal reminder.
- Battery life 7 days



## Paediatric Finger Oximeter PC60D2-P

- Special design for paediatric
- OLED display of SpO<sub>2</sub>, PR, PI, Pulse Bar & waveform
- 4 direction display
- Audio & visible alarm
- Automatic power on/off
- Low battery indication



## Creative Wrist Oximeter PC-68B

- Display of SPO<sub>2</sub>, RR, Pulse Bar, PI & waveform.
- Wearable design for continuous monitoring.
- Built-in Wireless module
- Data storage
- PC analysis software
- Audio & variable alarm



## Cuff with Finger Sensor PC-102

- Creative Vital Signs Spot Check Monitor
- Measures Blood Pressure & SPO<sub>2</sub>
- Rechargeable Battery.
- Built-in Wireless module for real time transfer.
- PC analysis software



# Veterinary Equipment



## Veterinary Scale

### ADE EHR4

- Hygienic Stainless Steel Platform 900 x 550mm
- Automatic zero setting, tare, unit switch kg/lbs, low battery indication, auto switch off, hold function

Capacity 150kg/50g



## Handheld Oximeter

### SP20 VET

- 3.5" auto rotation display of SpO<sub>2</sub>, PR, Pl% & Plethysmogram
- Patented **MoveOxy™ Technology** for removal of artefact and anti-motion performance
- Various probes
- Adjustable alarm limit, pitch, key tone level and brightness
- Up to 1000 hours data storage



## Electronic Small Animal Scale

### M112600

- Light-weight
- Clear 30 mm display and easy to use three key operation
- Operated by standard AA batteries and mains

Capacity 20kg/5g



## Vital Signs Monitor

### PC900PRO SN VET

- 3.5" High resolution LCD display
- User Friendly Operation
- Up to 12,000 NIBP and 2,000 groups SpO<sub>2</sub> storage
- Auto, Manual, Stat, Blood Pressure Measuring
- Built in rechargeable Li-ion battery
- USB connection for data output
- SpO<sub>2</sub> NIBP standard

Optional ECG, Temp & Thermal Printer



## Proscope™ 676

- Chrome-plated zinc combination chestpiece with 20mm diaphragm and bell
- Slip on gray non-chill bell ring for patient comfort
- Adjustable aluminum binaurals with PVC eartips
- Flexible 21" PVC tubing is made without phthalates
- Weighs 75g
- Overall length 78.5cm



## Proscope™ 675 Dual-Head Stethoscope

- Lightweight aluminum combination chestpiece with 1.4" diaphragm and 1" non-chill bell & rim
- Adjustable aluminum binaurals with PVC eartips
- Flexible 21" PVC tubing is made without phthalates
- Weighs 75g
- Overall length 80cm



## Multi Parameter Monitor

### UP7000 VET

- 12.1" High resolution colour touch screen
- Customised display of parameters
- 9 wave-form display simulant
- 5 kind calculation - Oxygenation, Renal Function, Ventilation, Medication and Hemodynamic



## Handheld Capnograph & Oximeter

### PC900B VET

- 3.5" Display of EtCO<sub>2</sub>, SpO<sub>2</sub>, PR, RR & Plethysmogram
- Various Probe Options
- Adjustable alarm limit, pitch, key tone level & brightness
- Up to 1000 hours data storage



# Medical Assessment Products



## ADE Pocket Oscope & Ophthalmoscope

- Fiber optics for cool, obstruction-free light transmission
- Choose from 2.5V xenon or LED lamp for true tissue color
- Wide-angle swivel 2.5X viewing lens permits instrumentation
- Insufflator port (fitting and insufflator sold separately) for pneumatic otoscopy
- AA battery handle (batteries included) with on/off power switch integrated into pocket clip



## Adscope Lite 619 Ultra-lite Clinician Stethoscope

- Combination diaphragm/bell adult chestpiece precisely machined from aluminum alloy
- Ultra-lightweight, only 105g
- Extra-large 1 3/8" diameter bell and highly sensitive 1 3/4" diaphragm
- Enhanced acoustic response; use light pressure for low frequencies, firm pressure for high frequencies
- Clinician headset boasts reinforcing yoke molded into the flexible 22" PVC tubing
- Made without natural rubber latex or phthalates
- Overall length 30"
- 18 tubing colors and 4 satin metal finishes
- Adsoft™ Plus Snap-on silicone eartips
- Non-chill bell and diaphragm retaining rim for patient comfort
- Includes both large and small Adsoft™ Plus eartips and scope ID tag



## ADC 760 Aneroid Sphygmomanometers

- Contemporary black enamel 300mmHg manometer provides day-in, day-out dependability exceeds industry standards
- Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system, prevents mis-cuffing
- Chrome-plated brass air release valve
- Leatherette carrying case with durable nylon zipper



## ADC Pro's Combo 768/641 Kit Adult

- Prophyg 768 Blood Pressure Sphygmomanometer
- Adscope 641 Sprague Stethoscope with Nylon Carrying Case

## Forehead Thermometer

- No touch
- Power: 2 x AAA battery
- Weight: about 100g (including batteries)
- Sizes: 168.5\*36\*48mm
- Measuring range: 32-42.9 degree
- Accuracy: +/- 0.2 degree



## The Metalite II™ 353

- Aluminum housing for lightweight comfort
- LED lamp for bright white illumination
- Locking push-toggle activated
- Pocket clip for convenience
- Measures 5 1/4" long by 1/2" diameter
- Two AAA batteries included



## ADC Proscope 665 Adult Stethoscope

- Single patient use scope
- ABS construction throughout is super lightweight
- May be ETO sterilized and safely incinerated
- MRI compatible and ultra-sensitive diaphragm for outstanding acoustic performance
- Individually packaged in a sealed polybag
- Superlightweight and durable
- Overall length 32"



## Adscope 603 Clinician Stethoscope

- Combination diaphragm/bell adult chestpiece from surgical stainless steel for outstanding performance and rugged durability
- Light pressure for low, bell-type frequency, firm pressure for high, diaphragm-type frequency
- Extra large bell (1 3/8" diameter) for unsurpassed low-frequency response
- Adsoft™ Plus silicone eartips with new snap-on connection for the ultimate in wearing comfort and acoustic seal



## Our range of Suction Units are:

- Smart designed
- High Quality
- Easy to manoeuvre
- Light weight
- Low noise
- Powerful suction

### DF-660

The Surgical Suction Unit has been designed for professional aspiration of body fluids, tissues or bones of the patient during or after surgery with the latest innovation in high-vacuum suction.

- 240v Mains Operation
- Oil-Less Piston Vacuum Pump
- **120LPM Flow Rate**
- 650mm Hg Max Pressure



### DF-770

The Portable Suction Unit provides reliable high vacuum for a variety of applications in hospitals, clinics and medical practices. A truly innovative and uniquely designed medical suction unit, offering new standards of hygiene, ease of use & low noise levels.

- 240v Mains Operation
- Oil-Less Vacuum Pump
- 40 liter / Minute Flow Rate
- 550mm Hg Max Pressure



Also available in the **Doctor's Friend range.**

Medeleq has been at the forefront of Medical Suction for over 35 years. That is why we depend on one of the worlds leading manufacturers of portable suction products. Quality, reliability and one of the most comprehensive ranges available in Australia.

# Portable Suction

## Sam E.P.S. Emergency Portable Suction Unit SAMEPS/01060601

- Hi Performance European manufactured Suction Device for Emergency use
- 12-28V DC, 240vMains or Internal Battery Operation
- 32LPM Flow Rate with 600mm mmHg Maximum pressure
- Lightweight 4.7kg including collection bottle
- With or without Carry bag
- 125min Battery on full power, or up to 11 hours on low settings
- Digital Display and Onboard lighting
- Data Logging and service light for Hospital Biomedical records
- Anti Bacterial coating available for extra infection Control requirements



**MGE**  
MEDICAL SUPPLIES

## Medical Professional Suction DF-506K

- 240v Mains Operation
- Single Piston Pump Design
- 40LPM Flow Rate
- 720mm Hg Max Pressure



**DF-506M**  
With the addition of:  
Integrated handle  
and wheels



## General Purpose Medical Suction Pump DF-760

- 240v Mains Operation
- Twin Diaphragm Pump Design
- 35LPM Flow Rate
- 650~700mm Hg Max Pressure



**DF-760B**  
With the addition of:  
Mounted on  
Trolley for  
easy handling

## Home Use or Light Suction DF-750

- 240v Mains Operation
- Twin Diaphragm Pump Design
- >20LPM Flow Rate
- 650~700mm Hg Max Pressure



## Portable Suction Unit JX820-D2

- Mains or Battery
- Long Life Piston Pump Design
- >20LPM Flow Rate
- 600mm Hg Max Pressure

# Suction Liner System



# BactiClear<sup>®</sup>

ANTIMICROBIAL

**Complete range of Disposable Suction Liners for all Hospital areas.**

Medeleq is the Exclusive Australian importer of VacSax UK Suction products.

VacSax now has the largest and most comprehensive product range available in the disposable suction market.



**SIZES AVAILABLE:**  
1 Litre, 2 Litre  
& 3 Litre  
Cascade or Shutoff



**VacSax BactiClear<sup>®</sup> Antimicrobial Disposable Suction Liner System.** An innovation unique to us, proven to reduce the growth of harmful bacteria, including MRSA and E.coli.

- ✓ **FREE STANDING LINERS**
- ✓ **GRADUATED**
- ✓ **HEAVY DUTY**
- ✓ **SIMPLE TO USE**
- ✓ **MINIMAL ENVIRONMENTAL IMPACT**
- ✓ **PVC FREE**



Pipeline Protector 9510-100



Fumevax in-line Smoke Filter 9510-080



# The Original Privacy Screen

3 easy and important facts to consider when comparing privacy screen systems

The Easy Concept consists of our three trademarks: EasyClean™ EasyReturn™ and EasyClick™.

- These key features make our wide range of products easily installed, cleaned and handled - contributing to optimal hygiene and easily-achieved patient privacy.
- The products within the Easy Concept are offered in a wide range of heights, lengths, combinations and colours.
- Our screens come in both fixed and mobile, as well as folding and frosted to meet the needs of an ever-changing environment.


## 1. EasyClean

Optimises infection prevention



## 2. EasyReturn

Quick and easy to fold back



## 3. EasyClick

Simple to connect and combine our screens



## Easy, reliable and flexible

Our Silentia Folding Screen, invented in 1990, is the essential piece of the Silentia Screen System. It's accordion design and hard surface make it ideal in any care setting supporting patient privacy and most importantly, instant cleaning.



### Fixed Folding

- Comes in a wide range of lengths, heights and colors
- Our standard screen is ordered with 25cm panels starting at 0.75m and a maximum length at 3.75m. Standard heights are 1.45m, 1.65m and 1.85m.



### Mobile Folding Screens

- Can be utilised as a freestanding privacy screen. ideal during bed side procedures or any situation where mobility is a requirement.
- The patented wheelbase design takes up a minimal amount of space without interfering with stability.



### Mobile Bed-End


- Extremely versatile trifold mobile privacy screen.
- Our mobile screens can be used separately or in combination with wall-mounted screens.
- Can be extended into two length sizes: 1.5m and 2m.




White



Forest



Ocean



Tropical




Savannah



Birch



Lime



Mandarin



Daylight

# Photo Panels

**Silentia's PhotoPanels help create not only privacy, but peace of mind and a comfortable feel.**

Studies show, that images and colors can have a positive influence on people by lowering their blood pressure, stimulating positive thoughts and relieving stress, pain and anxiety, ideal for elder care, hospitals and other healthcare settings.

**Relaxing, Remembering, Protecting**  
**PhotoPanels** has been developed as a professional helping tool to stimulate people in different situations by creating peace and calm, security and protection or memories of recognisable places. They offer at the same time optimal hygiene and privacy without feeling closed in.

**Find images & examples in our photo archive** at our website photo archive you can choose from a large variety of imagery. You can also use your own photo.

**Find your image at [silentia.eu](http://silentia.eu)**



# Story Panels & Colours



**Silentia's StoryPanels™** offer a new way to give children a fun and cheerful atmosphere by creating a positive distraction and engaging them with bright colours, interesting shapes and playful figures. Helping them to get lost in their imagination and to take their mind off their current circumstance.







**Shower Privacy Screens**

- Water Resistant
- Consists of 2 panels 750mm x 1450 mm
- Daylight Panels
- One side attaches to wall with 0-180 degrees opening

**DayLight Duo Screens**

Daylight Screens available in **Daylight Top** with **coloured panels**, or **frosted** sections which come in pairs with coloured panels. Also available in total **Daylight** panels.



- ClearPanels are a completely transparent, mobile screen that makes it easy to monitor & communicate with patients.
- ClearPanels meet the healthcare sector's strict hygiene requirements with smooth, easily accessible surfaces that withstand disinfectants used in any healthcare setting.
- A screen solution that is easy to clean, move and store when not in use.
- Section panels are made of 4mm thick polycarbonate
- Height either 1.55m. or 1.85m. Length: 1.5m



Transparent • Easy to clean • Easy to move • Easy to store

# francehospital® Medical Trolleys

## Innovative Systems for Medical and Ward Trolleys and Carts

*Innovative Designs feature:*

- ✓ Colour coded drawer handles
- ✓ Soft-closing drawers
- ✓ Configurable
- ✓ Pull-out shelf
- ✓ Various drawer sizes
- ✓ Drawer dividers
- ✓ Large wheels



### Perfectly Engineered Trolleys

To increase the quality of life in the wards, why not equip yourself with medication, therapy and emergency trolleys that are **comfortable, efficient, ergonomic and... silent?**



PRE9AAA-G

#### PRECISO Modular series

*Designed & manufactured in Europe to perfectly meet the needs of modern hospital hygiene, security & functionality requirements as well as the latest décor trends.*

- Aluminium & Composite plastic design for strength, light weight & long term durability
- 125mm Premium Tente Wheels
- Unique Covered Ball bearing drawer rails with soft-close function
- Lift-out Composite plastic drawers for ease of re-stocking or remote storage (very handy)
- Integrated push handles on 3 sides
- Upper & lower, molded bumper system for greater wall & trolley protection
- Keyed central locking is standard, breakaway lock or digital keypad options also available
- Fully configurable to any healthcare requirement



#### Integration of the finest materials into an exceptional Emergency Trolley platform

##### PRE9AER-G

*The Preciso Emergency Trolley configuration offers colour coded panels, defib tray, O2 cylinder holder, cardiac board, suction pump tray and IV Pole into an already feature packed platform. Ready to Go in Standard Configuration, but fully customisable to suit individual Hospital needs.*

#### PRECISO Anaesthetic PRE9AV-1

- 5 drawers 2 x 3", 2 x 6", 1 x 9"
- Overbridge with 3 accessory rails: 2 rows dispensers, 3rd with plastic tray and disinfectant dispenser
- IV Holder and kidney dish holder
- Central lock system
- LH - rail with glove box & sharps container
- RH - Maxi Wally waste bin
- Infinite Drawer Dividers
- ISO Rails



#### PRESTO 400 Medium Trolley PS-404-BL

- ABS worktop with three raised edges and anti-drip front rim
- Integrated push handle on three sides
- Super-lightweight composite aluminium/plastic panels
- Equipped with FH-drawers: transportable and mounted on square guides
- Overbridge available
- ABS base with perimeter bumper;
- High quality wheels: Ø 125 mm, double ball bearing, swiveling, non-marking, 2 with brake.

Size 662 x 594 x 1015 H



#### PRESTO 600 Large Trolley PS605-BL

- ABS worktop with three raised edges and anti-drip front rim
- Integrated push handle on three sides
- Super-lightweight composite aluminium/plastic panels
- Equipped with FH-drawers: transportable and mounted on square guides
- Overbridge available
- ABS base with perimeter bumper;
- High quality wheels: Ø 125 mm, double ball bearing, swiveling, non-marking, 2 with brake.

Size 855 x 594 x 1015 H



#### Perfectly Engineered Wheels

PRESTO medication and therapy trolleys and PERSOLIFE emergency trolleys feature *Nerviani Wheels*. These are high quality, carefully engineered and constructed wheels, that help to ensure a trolley that rolls smoothly through wards without rattling and clanking.

# Accessories



### Persolife Emergency PL60C5

- Orange lateral & back walls
  - 5 Drawers - 2 x 76mm, 2 x 156mm, 1 x 236mm
  - Lateral Pull-out shelf for extra working surface
  - LH: 2 tilt-out bins and storage compartment
  - RH: waste bin
  - Central Lock by cardiac board and seal system
  - Rear: Catheter holder, Oxygen tank compartment
  - Included: IV Pole & DEFIB holder
- 855mm x 594mm x 1050mm Capacity: 300kgs**



### Essential Clean Linen Trolley ES10GR-AK1

- Bed change (clean linen distribution, dirty linen and waste collection)
- Daily multiple activities thanks to its huge features and versatility
- Double ball-bearing wheels 100 mm
- Surface in turquoise ABS plastic with round corners
- Aluminum structure



### Persocar Trolley PT60

- 5 soft-closing drawers: 2 x 76mm, 2 x 156mm, 1 x 236mm
  - Full-width frontal handle, label holder & I D Tag
- Capacity: 300kg - 45kg/drawer**

### ERGODOC Patient File Trolley ED1

- Powder coated steel, with moulded plastic fronts
  - Large lockable purple work top (580x480mm)
  - Integrated push-handles on both sides
  - LH 2 compartments, 1 tilt-out bin, 1 hand disinfectant bottle holder
  - RH 3 x gloves boxes compartment
  - 125 mm wheels on double set of ball bearings, non-marking grey rubber, 2 with brakes
- Capacity about 35 patient files**



### Persocar Compact Trolley PT40L4

- 4 soft-closing drawers 3 x 155mm, 1 x 235 mm
  - Full-width frontal handle, label holder & I D Tag
- Capacity: 300kg - 45kg/drawer**



### Drawer Organiser Bins FH-V6

Available in 2 widths, 2 heights, 2 lengths. ABS Plastic, dividers separate.



### RED & YELLOW Plastic Breakable Security Seal SG-RED SG-YEL

Pack 100 pcs, Security number



### MEDSOPPY PCSOP21

Comb style dividers  
2 long 2 short



### Dividers for Organiser Boxes FH-DIV

Available in 2 widths & 2 heights



### SPG-T1

Wire support for gloves box, to hang on the accessories rail



### PCMPC

Support for laptop in coated steel with Plexiglas base



### Multiple Drawer Unit

**PCMBC54**  
5+4 - 600x140x370  
**PCMBC55**  
5+5 - 600x110x330



### Lable Holder PCMLH

Holds 18  
Includes bracket for mounting.



### IV Pole

Stainless steel, height adjustable, with 2 hooks  
**APF**  
PERSOCAR | PERSOCARCOLOR SHANGO  
**APF-CL**  
EVOLIGHT | COLORLIGHT



### Universal Support For Sharps Containers SUPUN

Wire Rilsan coated steel with fixing strap. Can be fixed with clamp (on request) For containers: from 0.6 to 7 Litres



### Maxi WALLY Waste Bin MWALI

Hanging waste collector with knee opening system. Internal removable bin, with flats for fixing the waste bag  
Capacity: 20 litres, 350 x 250 x 520mm



### Defibrillator Shelf DFS5

Easily securable for transport



### Patient Backboard (Cardiac Board) CB5



### Oxygen Cylinder Holder OC5



### Suspended Chart

**C1-FH**  
318 x 245mm (DIN A4)  
**C3-FH**  
465 x 390mm for X-rays / foolscap size



### CRUSHY Pill Crusher & Splitter PT-V

ABS plastic in violet  
Pack: 12 pcs in a display box



### Wire Basket CLT34 on SABU Lateral Accessories Rail

# Stainless Steel Trolleys

## Six Drawer Specialty Trolley SSTROLST006

Size 900 x 500 x 900mm  
 • 6 Drawer Specialty Trolley with 2 Shelves



## Seven Drawer Specialty Trolley SSTROLSCT005

Size 900 x 500 x 900  
 • 5 Horizontal Drawers  
 • 2 Large Half Drawers  
 • Push Handle LHS, Pull Up Shelf RHS



## 3-Tier Trolley SSTROLTM

SSTROLTM001: 490 x 490 x 900mm  
 SSTROLTM002: 600 x 490 x 900mm  
 SSTROLTM003: 750 x 490 x 900mm  
 SSTROLTM004: 800 x 490 x 900mm  
 SSTROLTM005: 900 x 490 x 900mm  
 SSTROLTM006: 1000 x 490 x 900mm



## Preparation Table SSTROLTG1500

SSTROLTG1500: 1500 x 600 x 900mm  
 SSTROLTG2400: 2400 x 600 x 900mm



## Flat Top Trolley SSTROLTE

SSTROLTE001: 490 x 490 x 900mm  
 SSTROLTE002: 600 x 490 x 900mm  
 SSTROLTE003: 750 x 490 x 900mm  
 SSTROLTE004: 800 x 490 x 900mm  
 SSTROLTE005: 900 x 490 x 900mm  
 SSTROLTE006: 1000 x 490 x 900mm



## Trolley with Rails SSTROLTA

SSTROLTA001: 490 x 490 x 900mm  
 SSTROLTA002: 600 x 490 x 900mm  
 SSTROLTA003: 750 x 490 x 900mm  
 SSTROLTA004: 800 x 490 x 900mm  
 SSTROLTA005: 900 x 490 x 900mm  
 SSTROLTA006: 1000 x 490 x 900mm



## Single Drawer Trolley SSTROLTB

SSTROLTB001: 490 x 490 x 900mm  
 SSTROLTB002: 600 x 490 x 900mm  
 SSTROLTB003: 750 x 490 x 900mm  
 SSTROLTB004: 800 x 490 x 900mm  
 SSTROLTB005: 900 x 490 x 900mm  
 SSTROLTB006: 1000 x 490 x 900mm

## IV Cannulation Trolley SSTROLWBT4

Size 590 x 440 x 1000mm  
 • Specialty Basket Trolley developed for IV Cannulation and other uses  
 • 2 x 100mm, 2 x 200mm Baskets  
 • Bracket for Sharps on RHS (container not included)



## Wire Basket Trolleys SSTROLWBT

• Several Sizes available including:  
 WBT1: 460 x 610 x 1140 WBT2: 1000 x 600 x 1350 WBT3: 1000 x 400 x 730  
 • Chrome Wire Frame and Baskets  
 • Stainless Steel Top



## Two Drawer Trolley SSTROLTC

SSTROLTC001: 490 x 490 x 900mm  
 SSTROLTC002: 600 x 490 x 900mm  
 SSTROLTC003: 750 x 490 x 900mm  
 SSTROLTC004: 800 x 490 x 900mm  
 SSTROLTC005: 900 x 490 x 900mm  
 SSTROLTC006: 1000 x 490 x 900mm



## Three Drawer Trolley SSTROLTN

SSTROLTN001: 490 x 490 x 900mm  
 SSTROLTN002: 600 x 490 x 900mm  
 SSTROLTN003: 750 x 490 x 900mm  
 SSTROLTN004: 800 x 490 x 900mm  
 SSTROLTN005: 900 x 490 x 900mm



## Four Drawer Trolley SSTROLTF

SSTROLTF001: 490 x 490 x 900mm  
 SSTROLTF002: 600 x 490 x 900mm  
 SSTROLTF003: 750 x 490 x 900mm



## Linen Trolley SSTROLGA

SSTROLGA001: Single Linen Skip  
 SSTROLGA002: Double Linen Skip  
 SSTROLGA003: Triple Linen Skip




## Linen Trolley with Lid SSTROLGA001

SSTROLGA001: Single Linen Skip with Lid  
 SSTROLGA002: Double Linen Skip with Lid  
 SSTROLGA003: Triple Linen Skip with Lid

## Stackable IV Pole SSTROLSIVP

• Stainless Steel IV Pole  
 • Weighted 5 Star Base  
 • 4 Hooks  
 • Height Adjustable 1300 - 2400 mm



## Five Drawer Trolley SSTROLTH

SSTROLTH001: 490 x 490 x 900mm  
 SSTROLTH002: 600 x 490 x 900mm  
 SSTROLTH003: 750 x 490 x 900mm



## Two Drawer Trolley Side by Side SSTROLTCW

SSTROLTCW004W: 490 x 490 x 900mm  
 SSTROLTCW005W: 900 x 490 x 900mm



## Four Drawer Trolley Side by Side SSTROLTD

SSTROLTD004: 490 x 490 x 900mm  
 SSTROLTD005: 900 x 490 x 900mm



## Mayo Table SSTROLTMAYO

• Stainless Steel Mayo Table  
 • With round tube Upright post  
 • Tray Size 600 x 400 mm



## Kick Bucket and Stand SSTROLKBS

• Stainless Steel Kick Bucket and Stand



## Step Stools SSTROLSTEP

# Wire Shelving & Baskets

**Wire shelving** is an ideal solution to the storage needs of Hospitals and Healthcare environments; trying to meet modern infection control guidelines.

Medeleq's wire shelving is manufactured to accommodate heavy loading and ease of access.

Product is available in stainless steel or chrome plated finish. We can supply basic shelving, right through to larger, high density compactus systems.

High Quality **Wire Storage Baskets** for many areas of Healthcare and Hospital.

Available with a Chrome Wire or Stainless Steel finish.

Various Sizes and configurations. Wall mountable or freestanding.

Ideal to compliment our Storage Shelving.




## Wire Compactus Systems

- Fixed Ends with rolling centre sections
- Various lengths
- Available in Chrome Wire or Stainless Steel
- Shelf lengths: 900mm, 1500mm, 1800mm
- Shelf widths: 450mm, 530mm, 610mm


## WIRE SHELVING SYSTEMS

- Mobile or fixed units of various lengths, heights, shelf sizes and number
- Available in Chrome Wire or Stainless Steel
- Shelf lengths: 900mm, 1200mm, 1500mm, 1800mm
- Shelf widths: 450mm, 530mm, 610mm



## LOUVRE WALL PANEL

- Roll/slide/hang along wall mounted track system allowing for removal & access to wall for cleaning
- Available in kits ranging from 2 to 5 Panels wide, with or without sliding space



**Extra Small**  
**FHWB-05-S**  
120mm D x 100mm  
W x 80mm H



**Small**  
**FHWB-10-S**  
240mm D x 130mm  
W x 135mm H



**Small**  
**FHWB-20-S**  
240mm D x 130mm  
W x 135mm H



**Small Hooked**  
**FHWB-25-S**  
120 mm D x 200mm W x 100mm H



**Medium**  
**FHWB-30-S**  
250mm D x 205mm W x 160mm H



**BASKET DIVIDER**  
FHWB-SPACS-S, FHWB-SPACL-S  
FHWB-40-S and FHWB-60-S



**Large**  
**FHWB-40-S**  
420mm D x 270mm W x 160mm H

**Small Tall**  
**FHWB-55-S**  
150mm D x 350mm W x 250mm H



**Extra Large**  
**FHWB-60-S**  
420mm D x 450mm W x 215mm H



**Long Wall Mountable**  
**HWB-45-S**  
Direct Wall Mount  
380mm D x 160mm W x 80mm H


## Accessories



Push Handle to suit Shelf Widths



Shelf Ledges  
130mm



Shelf Dividers  
150mm



Label Holders



Mobile or Fixed

## PATSLIDE®

*Patient Transfer made easy! in 3 safe steps...*




**All Purpose** 1525 L x 635 mm W  
**Theatre** 1525 L x 560 mm W  
**Slimline** 1400 L x 480 mm W  
**Duo-set (2pc)** 750 L x 635 mm W each

*Hanging hole for storage on wall mounted hook*

- Made from a specially formulated thermoplastic, the PATSLIDE static insulative and X-ray translucent
- Upper surface provides low resistance and the slide flexes to facilitate easy movement of the patient
- Easily cleaned by washing with a detergent, or any other aqueous based disinfectant
- Smooth routed edges to protect both patient and operator

*Handle cut-outs to aid board movement during patient transfer*

*Rubber stripping to increase grip and minimise board movement*



**Available in 3 configurations.**

Side Hand Rails are approx. 0.85m above the flat surface (fixed to outside of stairs)  
 Optional • Dual Height • Second Rail 650mm





**90 Degree – 3 x end steps, 2 side steps**  
 Overall length: 1.35m, width 1.05m  
 Main Platform: 0.75m x 0.75m  
 3 x end steps: 0.75m long, width 0.30m, height 0.13m  
 2 x side step: 0.75m long, width 0.30m, height 0.195m

**Straight Through – 3 x end steps, 2 x back steps**  
 Overall length: 1.65m, width 0.75m  
 Main Platform: 0.75m x 0.75m  
 3 x end steps: 0.75m long, width 0.30m, height 0.13m  
 2 x back step: 0.75m long, width 0.30m, height 0.195m



**Dead end – 3 end steps**  
 Overall length: 1.35m, width 0.75m  
 Main Platform: 0.75m x 0.75m  
 3 x end steps: 0.75m long, width 0.30m, height 0.13m



## One Brand Eight Fabrics

All designed with your comfort in mind.

# cherokee

INSPIRED COMFORT



Office and Showroom: 23/9-15 Sinclair Street, Arundel QLD 4214

Phone (07) 5529 0753  
or Free Call 1800 677 712

General enquiries [sales@medeleq.com.au](mailto:sales@medeleq.com.au)  
Shop online now [www.medeleg.com.au/Shop](http://www.medeleg.com.au/Shop)

Distribution all states

[www.medeleg.com.au](http://www.medeleg.com.au)