

Healthcare  
furniture,  
enhancing  
lives.

# Caremed Alrick.

Overbed Table Range



All our overbed tables have  A smooth surface,  
 Non-absorbent and free from gaps, so easy to clean


### OBTPL-NATOAK PREMIUM TILT TOP TABLE



#### Specifications

Tilt Area	385 x 560mm
Fixed Area	250 x 385mm
Tilt Angles	0°/18°/25°/32°/38°/43°
Frame	Powder Coated Steel
Height Range	785mm to 1050mm
Maximum Load	20 kg (Evenly Distributed)
Castors	4 x 50mm, 2 x Lockable

#### Colours



Natural Oak

#### Features



One Handed Operation



Tilt Top



Round Base



Locking Castors

### A6010 TIMBER TOP OVERBED TABLE




#### Specifications

Top Area	800 x 400mm
Frame	Powder Coated Steel
Height Range	785mm to 1050mm
C-Base Length	610mm (internal)
Maximum Load	20 kg (Evenly Distributed)
Castors	4 x 50mm, 2 x Lockable

#### Stain Colours



Cinnamon Beech Chestnut Cognac




Saddle Cocoa Natural Burnt Walnut

#### Features



Round Base



Locking Castors


### A6010 LAMINATE OVERBED TABLE



#### Specifications


Top Area	800 x 400mm
Frame	Powder Coated Steel
Height Range	785mm to 1050mm
C-Base Length	610mm (internal)
Maximum Load	20 kg (Evenly Distributed)
Castors	4 x 50mm, 2 x Lockable

#### Laminate Colours




Natural Oak (Standard)

#### Top Optional Colours



Oyster Grey Artisan Oak Soft Walnut



Jamaican Walnut European Walnut Shannon Oak South West Jarrah

#### Features



Round Base



Locking Castors

VIC  
2/77 Ricky Way  
Epping VIC 3076

NSW  
104 Loughnan Street  
Lake Cargelligo NSW 2672

UK/EU  
Port Stanley House  
Falkland Close Coventry CV4 8AU

T 1800 282 001  
sales@caremed.healthcare  
www.caremed.healthcare