



AAMI
 Level 1



Splash
 Resistant



Various
 Colours

Provides arm and body protection, ensuring coverage at all times.
 AAMI level 1 barrier protection.

STITCHED NECKLINE AND CUFF


Fold over stitched neckline and cuff for clean and protective finish

NECK AND WAIST TIES


Adjustable neckline and back ties to quickly apply and remove

SONICALLY BONDED


No perforation or chance of strike through



Blue



White



Yellow

Suitable for:

- Hospital and medical
- Emergency services
- Mortuary services
- Aged care

Reorder	Size	Feature	Colour	Pcs/Unit	Units/Ctn
GN001	Regular	AAMI Level 1	Blue	50	1
GN002	Regular	AAMI Level 1	White	50	1
GN003	Regular	AAMI Level 1	Yellow	50	1