



Anti-Static



Comfortable



Low Linting

Owear Theatre Caps are constructed from a non-woven, spun bond polypropylene fabric.

The ties are 18cm in length to allow for easy application and help to conform to the head for a comfortable fit.

ASSORTED COLOURS

Wide variety of colours to choose from

SINGLE USE

To aid in maintaining infection control



Suitable for:


Hospital and medical

Emergency rooms


Operating theatres

Delivery suites


ICU




Blue



Green



Pink



White



Dark Blue

Reorder	Size	Colour	Sterility	Pcs/Unit	Units/Ctn
TC005	One Size	Blue	Non-Sterile	250	8
TC006	One Size	Green	Non-Sterile	250	8
TC007	One Size	Pink	Non-Sterile	250	8
TC008	One Size	White	Non-Sterile	250	8
TC009	One Size	Dark Blue	Non-Sterile	250	8