

Level 3 Visor Mask

Anti-Fog - Anti-Glare - Surgical Ties - Indigo



PROTECTION STANDARDS



BFE (%): ≥ 98
 PFE (%): ≥ 98
 ΔP (mmH₂O/cm²): <6.0
 Fluid Resistance: 160mm Hg
 Flame Spread: Class 1

Variation in lot results exist.
 All results must meet or exceed
 ASTM F2100 Requirements.

Based on ASTM F2100 standards, this mask provides the highest level of protection during surgical procedures. The attached visor and advanced fluid-resistant and filtration capabilities of the mask provide all-in-one protection for procedures with a high risk of fluid splash or surgical debris, such as orthopedic surgery or trauma.

LONG SURGICAL TIES

Long ties allow for easy donning and doffing, with flexibility to adjust the tightness and fit of the mask over a surgical or bouffant cap

ANTI-GLARE VISOR

Designed with wrap-around wings, anti-glare film and anti-fog properties, providing all-in-one splash protection with an unobstructed view

WIDE NOSEPIECE WITH ANTI-FOG FOAM

Constructed with high-quality aluminum, anti-fog foam and a black anti-glare strip, providing enhanced visibility and comfort

ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the clinician's skin



Suitable for:

- Hospital and medical
- Emergency rooms
- ICU and trauma areas
- Areas with high risk of fluid splash

| Reorder | Size | Feature | Pcs/Unit | Units/Ctn |
|----------|----------|---------|----------|-----------|
| PG4-5073 | One Size | Level 3 | 25 | 4 |



| Reorder | Description | Colour | Anti-fog Feature | Fluid Resistance ¹ | BFE ² | Pcs/Box | Boxes/Case | Pcs/Case |
|--|---|-----------|------------------|-------------------------------|------------------|---------|------------|----------|
| LEVEL 2 BARRIER PROTECTION - FLUID RESISTANCE 120MM HG | | | | | | | | |
| PG4-1592 | ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Procedure Ear Loops | Blue | Foam | 120mm Hg | ≥ 98% | 25 | 4 | 100 |
| PG4-2592 | ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Surgical Ties | Blue | Foam | 120mm Hg | ≥ 98% | 25 | 4 | 100 |
| PM4-305 | ASTM Level 2 Mask / Surgical Ties | Blue | - | 120mm Hg | ≥ 98% | 50 | 6 | 300 |
| PM4-306 | ASTM Level 2 Mask / Procedure Ear Loops | Blue | - | 120mm Hg | ≥ 98% | 50 | 10 | 500 |
| PM4-311 | ASTM Level 2 Anti-Fog Strip Mask / Surgical Ties | Blue | Film | 120mm Hg | ≥ 98% | 50 | 6 | 300 |
| PM4-312 | ASTM Level 2 Anti-Fog Strip Mask / Procedure Ear Loops | Blue | Film | 120mm Hg | ≥ 98% | 50 | 10 | 500 |
| LEVEL 3 BARRIER PROTECTION - FLUID RESISTANCE 160 MM HG | | | | | | | | |
| PG4-1073 | ASTM Level 3 Anti-fog Foam Mask / Procedure Ear Loops | Indigo | Foam | 160mm Hg | ≥ 98% | 50 | 10 | 500 |
| PG4-1273 | ASTM Level 3 Mask / Procedure Ear Loops | Blue | - | 160mm Hg | ≥ 98% | 50 | 10 | 500 |
| PG4-1333 | ASTM Level 3 Soft Mask / Procedure Ear Loops / Irritant-Free | White | - | 160mm Hg | ≥ 98% | 50 | 10 | 500 |
| PG4-1563 | ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Procedure Ear Loops | Blue | Foam | 160mm Hg | ≥ 98% | 25 | 4 | 100 |
| PG4-2333 | ASTM Level 3 Soft Mask / Surgical Ties / Irritant-Free | White | - | 160mm Hg | ≥ 98% | 50 | 6 | 300 |
| PG4-3033 | ASTM Level 3 Irritant-Free Duckbill Mask / Surgical Ties | White | - | 160mm Hg | ≥ 98% | 50 | 6 | 300 |
| PG4-4053 | ASTM Level 3 Anti-Fog Tape Mask / Surgical Ties | Green | Adhesive Tape | 160mm Hg | ≥ 98% | 50 | 6 | 300 |
| PG4-4073 | ASTM Level 3 Anti-fog Foam Mask / Surgical Ties | Indigo | Foam | 160mm Hg | ≥ 98% | 50 | 6 | 300 |
| PG4-5073 | ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Surgical Ties | Indigo | Foam | 160mm Hg | ≥ 98% | 25 | 4 | 100 |
| PG4-6003 | ASTM Level 3 High-Filtration Mask / Surgical Ties | Grey | - | 160mm Hg | ≥ 98% | 50 | 6 | 300 |
| PAEDIATRIC LINE | | | | | | | | |
| PG4-1680 | Paediatric Mask / Dinosaur Pattern | Dinosaurs | - | - | - | 50 | 10 | 500 |

1. Resistance to Penetration by Synthetic Blood. ASTM 1862

2. Bacterial Filtration Efficiency. ASTM F2101