



Pacific Medical Catalogue

TABLE OF CONTENTS

Step Up Stools 2

STEP UP STOOLS

SINGLE STEP STOOL



Perfect for any examination couch this sturdy single step stool is made from 1" stainless steel tubing with reinforced rubber tips & a non slip ribbed rubber platform.

Specifications:

Stainless steel frame
Non-slip ribbed rubber platform
Reinforced rubber tips
Manufactured with 1" steel tubing
Dimensions: 455x245x200 (LxWxH)
Weight Capacity: 150KG

ORDER CODE: SST

[Click For More Information](#)

DOUBLE STEP STOOL

Perfect for any examination couch this sturdy double step stool is made from 1" stainless steel tubing with reinforced rubber tips & a non slip ribbed rubber platform.

Specifications:

Stainless steel frame
Non-slip ribbed rubber platform
Reinforced rubber tips
Manufactured with 1" steel tubing
Overall Dimensions: 400x550x470mm (LxWxH)
First Step: 250mm
Second Step: 220mm
Weight Capacity: 150KG

ORDER CODE: DST

[Click For More Information](#)



INDEX

D

Double Step Stool 2

S

Single Step Stool 2



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au