



Latex Free



Comfortable



Various Colours

Owear Round Caps are constructed from a non-woven, spun bond polypropylene fabric, which is anti-static, does not shed lint, and is latex free. The opening is elasticised and secured with a blanket stitch for a reliable and long lasting wear.

ASSORTED COLOURS

Wide variety of colours and patterns to choose from

ELASTICISED OPENING

Free of natural rubber latex and comfortable for all day wear

SINGLE USE

To maintain infection control




Suitable for:


Hospital and medical

Aged care


Patient




Blue




Green



Pink




White



Yellow



Christmas



Red



Dark Blue

Reorder	Size	Colour	Pcs/Unit	Units/Ctn
RC001	Regular	Blue	250	4
RC001L	Large	Blue	250	4
RC002	Regular	Green	250	4
RC003	Regular	Pink	250	4
RC004	Regular	White	250	4
RC005	Regular	Yellow	250	4
RC006X	Regular	Christmas	250	1
RC007	Regular	Red	250	4
RC008	Regular	Dark Blue	250	4