

ri-thermo® N

Simple, safe, reliable



Measuring range	0°C to 100.0°C (32.0°F to 212.0°F)
Measuring time	Approx. 1-2 seconds
Measuring accuracy	± 0.2°C, 32.0-42.2°C (± 0.4°F, 89.6-108.0°F)
Display	Liquid crystal display with 0.1°C (0.1°F) display units
Alarm	a. The device is switched on (ON) and ready to measure: 1 short beep. b. The measuring process has finished: 1 long beep. c. System error or fault: 3 short beeps
Memory	a. The last measured value appears when the device is switched on. b. 12 measurements can be retrie- ved in memory mode.
Operating temperature	5°C - 40°C (41.0°F - 104°F)
Storage/ transport temperature	-25°C - +55°C (-13°F - 131°F)
Automatic switch-off	About 1 minute after the last mea- surement. To prevent misuse, the device switches itself off automatically when the START button is pressed and held down for more than 15 seconds.
Battery	CR2032 BATTERY (X1) – at least 1000 measurements
Size	153 mm (L) x 31 mm (W) x 40 mm (H)
Weight	53g (with battery), 50g (without battery)
Standards	Meets ASTM E-1965 requirements

ri-thermo® N infrared thermometer

Is the new ergonomic infrared thermometer for uncomplicated use in medical practices, hospitals and in the home. Its innovative infrared technology not only allows **ri-thermo® N** to measure the body temperature of children and adults easily in a matter of seconds but also means that it can measure the temperature of surfaces, rooms and bathwater. The temperature values are displayed in a range from 0 to 100 °C or on the corresponding scale in ° Fahrenheit. Its large LCD display with backlighting makes the values easy to read. Even in the dark. Displays up to 12 stored temperatures. Temperature curves can also be plotted retrospectively or different measured values compared, so patients large and small can be monitored even more closely. The **ri-thermo® N** sounds a handy acoustic signal when it has finished measuring and if it detects a fever, has an energy-saving automatic switch-off function and comes with a hygienic storage container with probe covers. **ri-thermo® N** is of the usual high Riester quality, is clinically validated and is supplied with batteries.

- **ri-thermo® N** in a plastic box with 25 probe covers.
- Pack of 25 disposable probe covers.
- Pack of 100 disposable probe covers.
- Pack of 1000 disposable probe covers

Art. No.:	1805
Art. No.:	12700
Art. No.:	12701
Art. No.:	12702

Riester

Rudolf Riester GmbH | P.O. Box 35 | Bruckstr. 31 | D- 72417 Jungingen | Germany
Tel.: (+49) +7477-9270-0 | Fax: (+49) +7477-9270-70 | info@riester.de | www.riester.de

ECOMED
Serving The HealthCare Industry

Distributed by



 1800 67 5432

 1800 67 5433

 sales@ecomед.com.au

 www.ecomed.com.au