

seca 310

Mechanical circular dial scales for use far from any hospital

Tough and reliable, the seca 310 has been designed especially for work in areas where there is no access to hospitals and baby clinics. Weighing just over one kilo, with shatter-proof glass dial, corrosion-resistant metal casing and internal parts and robust mechanics, it is ideal for use anywhere in the world. The seca 310 has a fine graduation of either 25 or 50 g with a capacity of up to 15 or 25 kg. The suspension hooks are made of galvanized steel and are extremely durable and safe. The sling seat seca 410 is included in the basic equipment delivered with the circular dial scales seca 310.

- Capacity: 15 or 25 kg
- Graduation: 25 or 50 g
- With sling seat seca 410
- Optional: Carrying case seca 411

Accessories



seca 410

Sling seat for the circular dial scales seca 310

The ergonomically shaped sling is especially designed for weighing babies with the circular dial scales seca 310. The sling positions babies and toddlers comfortably and safely for weighing. The integrated padding can be adjusted individually to the size of the child. Made of durable nylon material, the sling seat is stable, hygienic and easy to clean.

Accessories

seca 411

Carrying case for seca 310 and seca 410

The handy carrying case seca 411 with strap provides sufficient room and protection for the circular dial scales seca 310 and the sling seat seca 410.



ECOMED

Serving the HealthCare Industry

T 1800 67 5432

F 1800 67 5433

E enquiries@ecomед.com.au

W ecomед.com.au