EMR-validated handrail scale with ID-Display[™] and optional height measurement



seca Scale-up Line

- + Display shows user and patient names or IDs for increased patient safety
- + Data transmission controlled by touchscreen
- + Safe weighing with a stable handrail
- + Large platform made of tested seca Bearclaw Glass® for absolute breakage resistance and perfect hygiene
- + EMR-validated: transmit measured data directly to any Electronic Medical Record (EMR) system



Scale

2011111116320

EMR-validated handrail scale with ID-Display[™] and optional height measurement



Robust, ergonomic handrail for patient stability

The elegantly-curved handrail is integrated with the cast-iron base and is part of the weight measurement, while providing added safety for patients who need it. This allows patients with restricted mobility to step onto the scale easier, the risk of falling is considerably reduced, and improves user workflow. The handrail is ergonomically designed to offer the ideal grip for patients, no matter which height.



Spacious and hygienic platform

The platform of the seca Scale-up Line consists of absolutely break-proof seca Bearclaw Glass®. The glass not only looks good, it also ages significantly slower than e.g. plastic and is particularly hygienic. It can be easily and quickly cleaned and has germ-reduction rates similar to stainless steel, with an anti-slip coating for safe standing. Transport is very easy due to the smooth-running, robust rubber wheels.



ID-Display[™] ensures easy and user-friendly data transmission to EMR systems

seca Scale-up Line is developed especially for perfect integration of measured data into the EMR. Weight and data such as time of measurement are stored centrally and accessible for every user. The ID-Display™ shows name or patient ID with every measurement, which are gathered by an additional barcode scanner. This makes the integration process transparent, safe and even more user-friendly. For more information please visit integration.seca.com

Order code ON SHM INT NN

+ Flat scale seca 654

KG-only available

+ Handrail seca 455

Capacity 800 lbs / 360 kg or 360 kg Graduation 0.1 lb / 50 g23.6 x 2.8 x 18.4 " Dimensions platform (W x H x D) 600 x 70 x 467 mm **Functions** Display of user and patient name or ID, TARE, Pre-TARE, mother+child function, HOLD, Auto-HOLD, Auto-BMI, CLEAR, Auto-CLEAR, lbs/kg switch-over, inch/cm switch-over, damping, standby mode, Wi-Fi, Ethernet Accessories Handrail seca 455, Ultrasonic height measurement seca 257, Barcode scanner mount seca 463, Panda seca 459, Large transport wheel set seca 468