

*Caring
for Every Body*



Bandages

1

Gauzes & Swabs

2

Tapes

3

Dressings

4

Sterilizing Consumables

5

Apparel

6

Miscellaneous

7

First Aid

8

Waste Bags

9

Advanced Wound Care Dressings

10

Compression Range

11

Surgical Consumables

12

Community Support

13

Disclaimer: All product information is correct at date of publication

Date Published: November 2019

All Sentry Medical Latex Free Products are "Not made with natural rubber latex"

Sentry Medical Pty Ltd
22 Peter Brock Drive Eastern Creek NSW 2766 | PO Box 6820 Wetherill Park NSW 1851
T: 1300 995 999 Int T: +612 8884 3995 | F: 1300 995 998 Int F: +612 8884 3998
E: sales@sentrymedical.com.au | W: www.sentrymedical.com.au



Management
System
EN ISO
13485:2012



www.tuv.com
ID 000041040

Bandages



CoWrap® Cohesive Bandage Tan

Non Sterile | Latex Free

CoWrap® Cohesive Bandages are constructed from a blend of non-woven fibres and methacrylate adhesive.

Ideal for:

- support
- securing primary dressings
- compression



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
COB001	2.5cm x 4.5m Stretched Tan	60	10
COB002	5cm x 4.5m Stretched Tan	36	10
COB003	7.5cm x 4.5m Stretched Tan	24	10
COB004	10cm x 4.5m Stretched Tan	18	10
COB005	15cm x 4.5m Stretched Tan	12	10

CoWrap® Cohesive Bandage Blue

Non Sterile | Latex Free

CoWrap® Cohesive Bandages are constructed from a blend of non-woven fibres and methacrylate adhesive. Blue colour allows easy detection as required within food or veterinary industries.

Ideal for:

- support
- securing primary dressings
- compression



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
COB051	2.5cm x 4.5m Stretched Blue	60	10
COB052	5cm x 4.5m Stretched Blue	36	10
COB053	7.5cm x 4.5m Stretched Blue	24	10
COB054	10cm x 4.5m Stretched Blue	18	10
COB055	15cm x 4.5m Stretched Blue	« 12	10

«Indent Item

Conforming Gauze Bandage

Non Sterile | No Clips

Sentry Medical Conforming Gauze Bandages are constructed from a blend of Polyamide yarn and Rayon.

These bandages are

- soft
- light weight
- porous, giving high comfort and air circulation



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
CFB001	2.5cm x 1.5m Unstretched	12	120
CFB002	5cm x 1.5m Unstretched	12	60
CFB003	7.5cm x 1.5m Unstretched	12	40
CFB004	10cm x 1.5m Unstretched	12	30
CFB005	15cm x 1.5m Unstretched	12	20

Bandafix® Nasal Sling

Non Sterile

Bandafix® Nasal Slings are a high quality elastic net with absorbent viscose wadding.

- quick and easy to apply
- high absorption capacity
- one size fits all

Ideal for

- post-operative nasal procedures



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NB001	Bandafix® Nasal Sling	25	24

Bandages

MEDICREPE®

Triangular Bandages

Non Sterile

Sentry Medical Triangular Bandages are available in cotton and non-woven material.

Cotton range:

- washable
- fast edges to prevent fraying
- good air circulation

Non-woven range:

- disposable
- cost effective



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
TB001	Cotton Calico 90cm x 90cm	200	1
TB002	Non-Woven 110cm x 110cm	400	1
TB002H	Cotton Calico 110cm x 110cm	200	1

MEDICREPE® All Cotton Crepe Bandage Light

Non Sterile | No Clips

The MEDICREPE® (All Cotton Crepe Bandage) is lightweight

Ideal for

- retention bandaging
- securing dressings



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ACL002	5cm x 1.6m Unstretched No Clips	12	60
ACL003	7.5cm x 1.6m Unstretched No Clips	12	40
ACL004	10cm x 1.6m Unstretched No Clips	12	30
ACL005	15cm x 1.6m Unstretched No Clips	12	20

MEDICREPE® All Cotton Crepe Bandage Medium

Non Sterile | No Clips

The MEDICREPE® (All Cotton Crepe Bandage) is medium weight

Ideal for

- retention bandaging
- securing dressings.
- conforms to BP standards



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ACC001	2.5cm x 1.6m Unstretched No Clips	12	120
ACC002	5cm x 1.6m Unstretched No Clips	12	60
ACC003	7.5cm x 1.6m Unstretched No Clips	12	40
ACC004	10cm x 1.6m Unstretched No Clips	12	30
ACC005	15cm x 1.6m Unstretched No Clips	12	20
ACC0051	15cm x 1.6m Unstretched No Clips (12) BULK	20	1

MEDICREPE® All Cotton Crepe Bandage Medium

Sterile | No Clips

The MEDICREPE® (All Cotton Crepe Bandage) is medium weight and provided sterile.

Ideal for

- retention bandaging
- securing dressings
- post-operative use
- conforms to BP standards



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ACS001	2.5cm x 1.6m Unstretched No Clips	150	1
ACS002	5cm x 1.6m Unstretched No Clips	120	1
ACS003	7.5cm x 1.6m Unstretched No Clips	80	1
ACS004	10cm x 1.6m Unstretched No Clips	70	1
ACS005	15cm x 1.6m Unstretched No Clips	50	1

Bandages

MEDICREPE®

MEDICREPE® Elastic Crepe Bandage Medium

Non Sterile | Clips | Overlocked

The **MEDICREPE®** (Elastic Crepe Bandage) Medium range is constructed from a blend of cotton and polyurethane fibres (spandex), giving excellent stretch and memory, with overlocked fast edges. Conforms to BP standards.

These bandages are

- washable
- cost effective

Ideal for:

- support bandaging



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
WB001	2.5cm x 1.5m Unstretched	12	120
WB002	5cm x 1.5m Unstretched	12	60
WB003	7.5cm x 1.5m Unstretched	12	40
WB004	10cm x 1.5m Unstretched	12	30
WB005	15cm x 1.5m Unstretched	12	20

MEDICREPE® Elastic Crepe Bandage Heavy

Non Sterile | Clips

The **MEDICREPE®** (Elastic Crepe Bandage) Heavy range is constructed from a blend of cotton and polyurethane fibres (spandex), giving excellent stretch and memory, with overlocked fast edges. Conforms to BP standards.

These bandages are

- washable
- cost effective

Ideal for:

- strong support bandaging



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
WB002H	5cm x 1.5m Unstretched	12	60
WB003H	7.5cm x 1.5m Unstretched	12	40
WB004H	10cm x 1.5m Unstretched	12	30
WB005H	15cm x 1.5m Unstretched	12	20
WB007H	5cm x 2.3m Unstretched	12	60
WB008H	7.5cm x 2.3m Unstretched	12	40
WB009H	10cm x 2.3m Unstretched	12	30
WB010H	15cm x 2.3m Unstretched	12	20

Gauzes & Swabs



Gauze Swabs

Sterile | 8 Ply

Sentry Medical Gauze Swabs are 100% cotton and are french folded to reduce fraying and linting. Available in

- plain white and green
- varying sizes
- different quantity sterile packs

Also with x-ray detectable thread for surgical procedures (see ThreadX® range, page 11-1 to 11-4).



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
GSG2425	10cm x 10cm x 5pcs Green	50	9
GSR2223	5cm x 5cm x 3pcs White	50	30
GSR2225	5cm x 5cm x 5pcs White	50	24
GSR2322	7.5cm x 7.5cm x 2pcs White	50	24
GSR2323	7.5cm x 7.5cm x 3pcs White	50	24
GSR2325	7.5cm x 7.5cm x 5pcs White	50	24
GSR2422	10cm x 10cm x 2pcs White	50	16
GSR2423	10cm x 10cm x 3pcs White	50	8
GSR2425	10cm x 10cm x 5pcs White	50	8

Gauze Swabs

Non Sterile | 8 Ply

Sentry Medical Gauze Swabs are 100% cotton and are french folded to reduce fraying and linting. Available in

- white and green
- varying sizes

Also with x-ray detectable thread for surgical procedures (see ThreadX® range, page 11-1 to 11-4).



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
GS001	7.5cm x 7.5cm 100pcs	1	40
GS002	10cm x 10cm 100pcs	1	40
GS003	10cm x 10cm Green 100pcs	1	40
GS004	5cm x 5cm 100pcs	1	160

Ribbon Gauze

Sterile

Sentry Medical Ribbon Gazes are constructed from 100% cotton gauze, making them soft, gentle, and absorbent. Also available with x-ray detectable thread (see ThreadX® range, page 11-1 to 11-4).

- Ideal for packing
- deep wounds
 - undermining and tunnels



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
GSV001	10cm x 3m	100	1

Gauzes & Swabs



Cotton Balls

Sterile | Non Sterile

Sentry Medical Cotton Balls are 100% pure cotton, that are both gentle and absorbent. Cotton Balls are available in

- small and large sizes
- sterile packets
- non sterile bulk packs



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
CB001	Small 100pcs	1	40
CB001L	Large 100pcs	1	20
CB002S	Large Sterile 5pcs	50	12

Cotton Tip Applicators

Sterile | Non Sterile

Sentry Medical Cotton Tip Applicators are available

- sterile and non sterile
- in different lengths and packet sizes
- with a choice of wooden and plastic stems



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
CTA001	7.5cm Double Ended Plastic 100pcs	150	1
CTA004	15cm Single End Wooden 100pcs	24	10
JTA001	18cm Jumbo Single End Plastic 20pcs	20	10

Non-Woven Balls

Non Sterile | Latex Free

Sentry Medical Non-Woven Balls are constructed from a polyester and rayon blend.

They have four times the absorbency of conventional cotton gauze and are low linting. Also available with x-ray detectable thread (see ThreadX® range, page 11-1 to 11-4).



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NWB001	Non-Woven Balls 100pcs	1	20

Gauzes & Swabs



Non-Woven Swabs

Sterile | 4 Ply

Sentry Medical Non-Woven Swabs are constructed from a polyester and rayon blend. They have four times the absorbency of conventional cotton gauze and are low linting.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NWS2322	7.5cm x 7.5cm x 2pcs	50	40
NWS2323	7.5cm x 7.5cm x 3pcs	50	10
NWS2425	10cm x 10cm x 5pcs	50	6

Non-Woven Swabs

Non Sterile | 4 Ply

Sentry Medical Non-Woven Swabs are constructed from a polyester and rayon blend. They have four times the absorbency of conventional cotton gauze and are low linting.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NWS001	5cm x 5cm 100pcs	1	160
NWS002	7.5cm x 7.5cm 100pcs	1	40
NWS003	10cm x 10cm 100pcs	1	40

Non-Woven Split Drainage Swabs

Sterile | 4 Ply

Sentry Medical Non-Woven Split Drainage Swabs are constructed from a polyester and rayon blend. They have four times the absorbency of conventional cotton gauze and are low linting.



Cut with a Y split for easy placement around drains, they require no cutting and are ready to use.

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NWD1321	7.5cm x 7.5cm x 2pcs	50	10
NWD1421	10cm x 10cm x 2pcs	50	10

Tapes



SofSecure® Silicone Fixation Tape

Non Sterile | Latex Free

Sentry Medical SofSecure® Silicone Fixation tape provides:

- gentle, secure and water resistant fixation for dressings and devices
- scissor-free use with “Easy Tear” perforations every 1cm
- 3 widths to cover a variety of applications



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ISB160	1.9cm x 1.5m	6	20
ISB163	2.5cm x 5m	6	20
ISB164	5cm x 1.5m	6	20

Paper Tape

Non Sterile | Latex Free | Individually Wrapped

A low allergy paper tape that provides general purpose medical securement. Available in a variety of widths, with and without dispenser.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
HYP001	1.25cm x 9.1m	24	30
HYP002	2.5cm x 9.1m	12	30
HYP003	5cm x 9.1m	6	80
HYP004	7.5cm x 9.1m	4	30
HYP100	2.5cm x 9.1m with Dispenser	12	30
HYP101	1.25cm x 9.1m with Dispenser	24	30

Zinc Oxide Tape

Non Sterile | Latex Free

Sentry Medical Zinc Oxide Tapes are constructed from plain cotton weave fabric with an adhesive containing Zinc Oxide.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ZAT002	1.25cm x 5m	24	30
ZAT001	2.5cm x 5m	12	30

Tapes



AsGUARD® Flex Non-Woven Roll

Non Sterile | Latex Free

AsGUARD® Flex Non-Woven Roll is self-adhesive, non-woven fabric fixation tape. Presented in a dispenser box.

Provides

- comfortable, conforming adhesion
- simple, secure fixation of dressings and access devices



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ISD020	2.5cm x 10m	1	30
ISD021	5cm x 10m	1	24
ISD022	7.5cm x 10m	1	24
ISD023	10cm x 10m	1	12
ISD024	15cm x 10m	1	12
ISD025	20cm x 10m	1	6

Dressings



AsGUARD® Clear Film Dressing and Roll

Transparent Film

AsGUARD® Clear Film Dressing and Roll is a transparent adhesive polyurethane film with high moisture vapour permeability which

- adheres to the skin
- forms a flexible, comfortable layer
- is a water resistant barrier
- allows skin to breath
- assists in the healing process



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ISB001	4cm x 5cm Sterile Dressing	50	10
ISB002	6cm x 7cm Sterile Dressing	^A 50	10
ISB003	10cm x 12cm Sterile Dressing	^A 25	10
ISB004	15cm x 20cm Sterile Dressing	^A 25	10
ISB010	5cm x 10m Non Sterile Roll	1	24
ISB011	10cm x 10m Non Sterile Roll	1	12

^A“Date of application” label included

AsGUARD® EXCEED

Sterile | High MVTR Film Dressing

AsGUARD® EXCEED High MVTR Film Dressing

- provides water and microbial protection to any type of IV access site
- may assist in minimising the risk of CLABSI when managing PICC and Central lines
- with stabilizing outer layer facilitates one handed application
- includes fixation strips to secure IV lines and ports
- insertion or dressing change date label ensuring maintenance of procedure documentation



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ISB202	6cm x 7cm	50	10
ISB203	10cm x 12cm	25	10

AsGUARD® Clear + Film Island Dressing

Sterile | Transparent | Non-Adherent Pad

AsGUARD® Clear + Film Island Dressings combine an absorbent pad with a transparent adhesive polyurethane film. These dressings are

- flexible
- water resistant
- gentle on the skin

Ideal for

- simple and effective management of minor wounds



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ISC001	4cm x 5cm	50	10
ISC002	6cm x 7cm	^A 50	10
ISC003	10cm x 12cm	^A 25	10
ISC004	15cm x 20cm	^A 25	10
ISC025	10cm x 25cm	25	10
ISC026	10cm x 30cm	25	10
ISC028	10cm x 20cm	25	10

^A“Date of application” label included

Dressings



AsGUARD® FLEX + Non-Woven Island Dressing

Sterile | Adhesive | Non-Adherent Pad

AsGUARD® FLEX + Non-Woven Island Dressings combine an absorbent pad with a soft and comfortable non-woven adhesive backing.

Ideal for

- simple and effective management of minor wounds



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ISD001	7.5cm x 5cm	50	10
ISD002	8cm x 10cm	50	10
ISD003	15cm x 10cm	20	10
ISD004	20cm x 10cm	20	10
ISD006	10cm x 25cm	25	10
ISD008	10cm x 30cm	25	10
ISD030	8.25cm x 6cm	50	10

AsGUARD® FLEX + Roll


Non Sterile | Latex Free

AsGUARD® FLEX + Roll is a “cut-to-required-length” roll of self-adhesive non-woven fabric fixation tape with low adherent wound pad, presented in a dispenser box.


The measurement guide on release paper makes sizing of dressing easy.

Ideal for

- simple dressing of minor wounds, where sterile dressing not required.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ISD011	5cm x 10m	1	24
ISD012	7.5cm x 10m	1	24
ISD013	10cm x 10m	1	12



AsGUARD® FLEX + Haemoplast

Sterile | Latex Free

AsGUARD® FLEX + Haemoplast is a non-woven adhesive fabric with absorbent pressure pad.

Ideal as

- a dressing for injection, venepuncture and IV access sites.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ISD100	AsGUARD® Flex + Haemoplast	100	40



Dressings



AsGUARD® I.V. Dressing

Sterile | Transparent

AsGUARD® I.V. Dressing is a transparent adhesive polyurethane film dressing. The 'U' shape cut out allows easy IV line access.

Ideal for

- dressing peripheral intravenous access sites



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
IVD002	6cm x 7cm	50	10



AsGUARD® SurgiStrip®

Sterile | Reinforced Skin Closure Strips

Sentry Medical AsGUARD® SurgiStrip® Skin Closures provides

- superior strength
- exceptional porosity, and
- high quality, medical grade adhesive, giving maximum support and protection to your wounds



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
C-01	3mm x 75mm, (5's)	50	20
C-02	6mm x 75mm (3's)	50	20
C-03	6mm x 38mm (6's)	50	20
C-04	6mm x 100mm (10's)	50	20
C-05	13mm x 100mm (6's)	50	20
C-06	25mm x 125mm (4's)	25	20

Dressings



IMPERVIA® Absorbent Dressing

Sterile | Impervious Backing

IMPERVIA® absorbent sterile dressings are designed to manage moderate to heavily exuding wounds. These dressings are a highly absorbent inner layer of cellulose blend with blue impervious backing enclosed by a soft non-woven fabric.

Ideal for use

- in direct contact with wound, or as a secondary dressing
- where protection of clothing and linen is desired



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NWC009	10cm x 10cm	50	16
NWC010	10cm x 20cm	25	8
NWC011	20cm x 20cm	25	8

Non-Woven Combine Dressing

Sterile

Sentry Medical Combine Dressings consist of a highly absorbent inner, enclosed in a soft and comfortable non-woven covering fabric with soft end seals.

Ideal for

- use as primary or secondary dressing
- medium to heavy exuding wounds



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NWC001	10cm x 10cm Sterile	50	16
NWC002	10cm x 20cm Sterile	50	8
NWC003	20cm x 20cm Sterile	50	5

Non-Woven Combine Roll Dressing

Non Sterile

Sentry Medical Combine roll dressings consist of a highly absorbent inner, enclosed in a soft and comfortable non-woven fabric. The roll can be cut to sizes that suit.

Ideal for

- medium to heavy exuding wounds
- padding or protection for major wounds
- primary or secondary dressing



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NWC100	10cm x 10m	10	1
NWC101	20cm x 10m	10	1

Dressings



Universal Multi Trauma Combine Dressing

Sterile

Sentry Medical Multi Trauma Dressings consist of a highly absorbent inner, enclosed in a soft and comfortable non-woven fabric with soft end seals.



Ideal for

- medium to heavy exuding wounds
- padding or protection for major wounds
- use as primary or secondary dressing

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
UD001	Large 20cm x 90cm	100	1
UD002	Small 20cm x 40cm	100	1

Wound Dressing

Sterile

Sentry Medical Wound Dressings are constructed of thick, highly absorbent cotton pad attached to a cotton gauze bandage.



Ideal for

- First Aid emergency response

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
WD001	Small No. 13	12	48
WD002	Medium No. 14	12	22
WD003	Large No. 15	12	18

Burns Dressing

Sterile

Sentry Medical Burns Dressings have on one side a thin, non-adherent layer: prevents adherence to burns, allows ease of removal.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
BD001	75cm x 220cm	50	1
BD002	75cm x 75cm	100	1

Apply to the burn site and wet thoroughly.
First aid situations:

- provides cool, moist and soothing relief
- protects the affected wound site
- ideal in First Aid situations

Burns Pad

Sterile

Sentry Medical Burns Pads are constructed with 100% cotton, with inner wadding covered by a layer of gauze. Burns Pads are a secondary dressing, used for



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
BD003	50cm x 76cm	20	1

Ideal for use as padding while awaiting medical review.

- protecting the wound
- comfort of the patient

Dressings



Mediam® 135 Low Adherent Dressing- Lite

Sterile

Mediam® 135 Low Adherent Dressing- Lite consists of a thin absorbent inner pad encased between two layers of highly embossed polyethylene film. This film minimises adherence of dressing to the wound.



Ideal for

- lightly exuding and small wounds
- First Aid applications

REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
NAL001	7.5cm x 5cm	100	8
NAL002	7.5cm x 7.5cm	100	8
NAL003	7.5cm x 10cm	100	6
NAL004	7.5cm x 20cm	100	4

Mediam® 210 Low Adherent Dressing

Sterile

Mediam® 210 Low Adherent Dressing consists of a heavier gauge, highly absorbent inner pad encased between two layers of highly embossed polyethylene film. This film minimises adherence of dressing to the wound.



Ideal for

- light to moderate exuding wounds

REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
NAD001	5cm x 5cm	50	10
NAD002	7.5cm x 7.5cm	50	10
NAD003	10cm x 10cm	50	8
NAD004	10cm x 20cm	50	8

Dressings



Non Woven Eye Pad

Sealed Edge

Sentry Medical Eye Pads are constructed from absorbent inner, and are covered with non-woven outer layer with sealed edges.

- prevents linting



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
EP001	6cm x 8cm Sterile	50	24
EP002	6cm x 8cm Non Sterile Bulk	2000	1

EyePro® Eyelid Occlusion Dressing

Sterile | Latex Free

The **EyePro®** is constructed from a non-woven fabric and film laminate.

Ideal for use on patients

- undergoing anaesthesia
- being managed in an unconscious state

The dressing is designed to

- position the transparent film over the closed eyelid margin
- provide total closure of the eye
- protect the eye from dryness and accidental trauma



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
EP010	EyePro™ Eyelid Occlusion Dressing (2pcs)	50	10



SENTURIAN® Dialysis Pack

Sterile

Dialysis Packs are intended to be used by trained personnel in a dialysis arena.

The components include:

- 1 x PE/Cellulose Sterile Field, Blue, (60 x 60cm)
- 1 x White Dressing Tray Type 3 (140mm x 130mm x 18mm)
- 2 x Dressing Towel (40cm x 40cm)
- 1 x Water Proof Drape (330mm x 400mm)
- 15 x Non Woven balls
- 10 x Gauze Swabs 10cm x 10cm (8 ply)
- 3 x Disposable Plastic Forceps: 2 x Blue, 1 x Yellow, 11cm

Sterile field fold method type 4

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
WDP030	SENTURIAN® Dialysis Pack	« 120	1

«Indent Item



Dressings



SENTURIAN® Basic Dressing Pack

Sterile

SENTURIAN® Basic Dressing Packs consists of

- absorbent inner layer sterile field with impervious back
- tip resistant tray
- 3 forceps - 2 blue, 1 yellow
- dressing towel
- cleansing components of choice

The dressing packs come in a variety of presentations with options of included components

- cotton or non-woven balls and swabs
- tray with or without swab drainage area
- tear or peel open outer packaging




REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
WDP004	Type 4 Basic Dressing Pack (Tear) 1 x PE/Non-woven Sterile Field, Blue, 50cm x 50cm 1 x White Dressing Tray plain , 140mm x 130mm x 18mm 1 x Dressing Towel, 40cm x 40cm 5 x Cotton Balls 2 x Gauze Swabs, 7.5cm x 7.5cm, 8 ply 3 x Disposable Plastic Forceps: 2 x Blue, 1 x Yellow, 11cm	160	1
WDP005	Type 5 Basic Dressing Pack (Tear) 1 x PE/Non-woven Sterile Field, Blue, 50cm x 50cm 1 x White Dressing Tray plain , 140mm x 130mm x 18mm 1 x Dressing Towel, 40cm x 40cm 5 x Non Woven Balls 2 x Non Woven Swabs, 7.5cm x 7.5cm, 4 Ply 3 x Disposable Plastic Forceps: 2 x Blue, 1 x Yellow, 11cm	160	1
WDP006	Type 6 Basic Dressing Pack (Tear) 1 x PE/Non-woven Sterile Field, Blue, 50cm x 50cm 1 x White Dressing Tray graded swab drainage area, 140mm x 130mm x 18mm 1 x Dressing Towel, 40cm x 40cm 5 x Non Woven Swabs, 7.5cm x 7.5cm, 4 Ply 3 x Disposable Plastic Forceps: 2 x Blue, 1 x Yellow, 11cm	160	1
WDP007	Type 7 Basic Dressing Pack (Peel) 1 x PE/Non-woven Sterile Field, Blue, 50cm x 50cm 1 x White Dressing Tray graded swab drainage area, 140mm x 130mm x 18mm 1 x Dressing Towel, 40cm x 40cm 6 x Gauze Swabs, 7.5cm x 7.5cm, 8 Ply 3 x Disposable Plastic Forceps: 2 x Blue, 1 x Yellow, 11cm	160	1
WDP009	Type 9 Basic Dressing Pack (Peel) 1 x PE/Non-woven Sterile Field, Blue, 50cm x 50cm 1 x White Dressing Tray graded swab drainage area 140mm x 130mm x 18mm 1 x Dressing Towel, 40cm x 40cm 6 x Non Woven Swabs ,7.5cm x 7.5cm 3 x Disposable Plastic Forceps: 2 x Blue, 1 x Yellow, 11cm	160	1
WDP010	Type 10 Basic Dressing Pack (Tear) 1 x PE/Non-woven Sterile Field, Blue, 50cm x 50cm 1 x White Dressing Tray graded swab drainage area 140mm x 130mm x 18mm 1 x Dressing Towel, 40cm x 40cm 6 x Non Woven Swabs, 7.5cm x 7.5cm, 4 Ply 3 x Disposable Plastic Forceps: 2 x Blue, 1 x Yellow, 11cm	160	1
WDP011	Type 11 Basic Dressing Pack (Tear) 1 x PE/Non-woven Sterile Field, Blue, 50cm x 50cm 1 x White Dressing Tray plain , 140mm x 130mm x 18mm 1 x Dressing Towel, 40cm x 40cm 5 x Gauze Swabs, 7.5cm x 7.5cm 3 x Disposable Plastic Forceps: 2 x Blue, 1 x Yellow, 11cm	160	1

Sterilizing Consumables



Sentry Medical imports and distributes the Bolsaplast film laminate range from Europe. CE marked, manufactured and labelled according to EN868 and ISO11607 whilst complying with AS4187:2014. It is produced under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilization process. Produced with 60GSM medical grade paper and PET/PP blend for the film, it offers extremely low fibre shedding and is resistant to film tear. It is available in heat sealable and self sealable pouches, flat and gusseted reels. It is also available with a Tyvek® backing for plasma sterilization.


Bolsaplast Flat Rolls



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
MRP1040	50mm x 200m	1	1
MRP1041	75mm x 200m	1	1
MRP1042	100mm x 200m	1	1
MRP1043	125mm x 200m	«	1
MRP1044	150mm x 200m	1	1
MRP1045	200mm x 200m	1	1
MRP1046	250mm x 200m	1	1
MRP1047	300mm x 200m	1	1
MRP1048	400mm x 200m	1	1

«Indent Item

Bolsaplast Gusset Rolls



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
MRF1060	100mm 40mm x 100m	1	1
MRF1061	150mm 40mm x 100m	1	1
MRF1062	200mm x 60mm x 100m	1	1
MRF1063	250mm x 60mm x 100m	«	1
MRF1064	300mm x 80mm x 100m	«	1
MRF1065	400mm x 80mm x 100m	«	1

«Indent Item

Bolsaplast Tyvek® Flat Rolls



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
MRT1170	50mm x 70m	«	1
MRT1171	75mm x 70m	«	1
MRT1172	100mm x 70m	1	1
MRT1174	150mm x 70m	1	1
MRT1175	200mm x 70m	1	1
MRT1176	250mm x 70m	«	1
MRT1177	300mm x 70m	«	1
MRT1178	400mm x 70m	«	1

«Indent Item

Sterilizing Consumables



Bolsaplast Heat Sealable Pouches



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
MBO1000	75mm x 210mm	2500	1
MBO1001	75mm x 280mm	2400	1
MBO1003	100mm x 200mm	2400	1
MBO1004	100mm x 240mm	2400	1
MBO1005	100mm x 280mm	2400	1
MBO1006	100mm x 420mm	2400	1
MBO1007	150mm x 210mm	2000	1
MBO1008	150mm x 280mm	2000	1
MBO1009	150mm x 320mm	2000	1
MBO1010	150mm x 360mm	2000	1
MBO1011	200mm x 320mm	1800	1
MBO1012	200mm x 360mm	1800	1
MBO1013	200mm x 420mm	1200	1
MBO1014	250mm x 360mm	1200	1
MBO1015	250mm x 420mm	1200	1
MBO1016	250mm x 500mm	1200	1
MBO1017	300mm x 420mm	600	1
MBO1018	300mm x 500mm	600	1
MBO1019	400mm x 580mm	600	1
MBO1020	100mm x 550mm	« 2000	1
MBO1021	150mm x 400mm	1200	1
MBO1022	100mm x 580mm	« 2000	1
MBO1023	200mm x 500mm	« 1200	1
MBO1024	300mm x 580mm	« 700	1

«Indent Item

Bolsaplast Tyvek® Heat Sealable Pouches



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
MBT1160	75mm x 210mm	« 1000	1
MBT1161	100mm x 280mm	« 1000	1
MBT1162	150mm x 320mm	« 1000	1
MBT1163	200mm x 420mm	« 350	1
MBT1164	250mm x 500mm	« 350	1
MBT1166	400mm x 580mm	« 350	1

«Indent Item

Bolsaplast Heat Sealable Gusset Pouches



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
MBF1202	200mm x 55mm x 400mm	1	1

«Indent Item

Sterilizing Consumables



bolsaplast

Bolsaplast Self Sealable Pouches



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
MBA1080	75mm x 200mm	200	6
MBA1081	90mm x 200mm	«	200
MBA1082	90mm x 230mm	200	6
MBA1084	100mm x 240mm	«	200
MBA1085	130mm x 250mm	200	6
MBA1086	130mm x 330mm	200	6
MBA1087	140mm x 280mm	«	200
MBA1088	200mm x 330mm	200	4
MBA1089	250mm x 400mm	200	4
MBA1090	300mm x 390mm	400	1
MBA1091	300mm x 450mm	«	400

«Indent Item

Dust Cover Bags

Sentry Medical Dust Cover Bags are constructed from varying grades of polythene, making them highly durable with a slightly rubberised finish. These bags can also be heat sealed using an impulse bar sealer. Dust Cover Bags provide an ideal storage environment for instruments and trays that have undergone the sterilization process, and need protection from dust whilst in storage.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
DC001	152mm x 254mm	100	20
DC002	203mm x 305mm	100	15
DC003	254mm x 380mm	100	15
DC004	305mm x 508mm	100	15
DC005	406mm x 508mm	100	5
DC006	450mm x 762mm	100	4
DC007	610mm x 762mm	100	4

ISP® Self Sealable Dust Cover Bags

The ISP® Self Sealable Dust Cover Bags are constructed from a polyester and polypropylene blend. These bags have an easy peel opening, adhesive sealing tab and an area for clear labelling. Self Sealable Dust Cover Bags provide an ideal storage environment for instruments and trays that have undergone the sterilization process, and need protection from dust whilst in storage.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
	150mm x 250mm	100	5
	200mm x 300mm	100	5
	250mm x 375mm	100	5
	300mm x 500mm	100	5
	400mm x 550mm	100	5
	400mm x 750mm	50	5
	600mm x 750mm	50	5

Instrument Protectors

The ISP® Medical Instrument Protector is used during the sterilization process to reduce the risk of sharp instrument piercing of the sterilization packaging material. We offer 2 sizes.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
3FSZB730802	Small 100 pcs	1	1
3FSZB730804	Medium 100 pcs	1	1

Sterilizing Consumables



Sentry Medical imports and distributes the INTERSTER® range of specialised CSSD products. These products meet the applicable ISO and CE requirements. This range includes challenge devices, seal checks, steam monitoring strips and other CSSD consumables. INTERSTER® use class 6 ink technology across their entire steam sterility assurance range.

ISP® Bowie & Dick Autoclave Test Pack

Type 2

The ISP® Bowie & Dick Autoclave Test Pack aligns with the requirement of AS4187:2014 as a steam penetration, air removal test.

We have 2 versions of Bowie & Dick Test Packs available. Both are 7kg mimic. One test pack has obtained a BSI Kitemark (independent certification to ISO11140:1, ISO11140:4)



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
3FSKS610804	Bowie & Dick Autoclave Test	20	1
3FSKS610860	Bowie & Dick Autoclave Test, Kitemark	20	1

Conforms to ISO11140.

ISP® Steam Emulating Indicator

Type 6 | 3.5min | 134 degree

The ISP® Steam Emulating Indicator is a test strip that is reflective of time, temperature + steam quality.

The Type 6 chemical indicator is non-toxic and lead-free, and the strip is laminated to prevent the ink coming in contact with the contents of the tray. One version we offer is now certified with BSI Kitemark (independent certification to ISO11140).



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
3FSKS645818	3.5min/134 °C Chemical Indicator	250	1
3FSKS645820	3.5min/134 °C Self Adhesive Indicator	250	1
3FSKS645834	3.5min/134 °C Kitemark Indicator	250	1

Conforms to ISO11140-1.

ISP® Plasma Indicator

Type 1

The ISP® Plasma Indicator complies with ISO standards and has lead free ink.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
3FSKS645808	Plasma Indicator	« 250	1

«Indent Item

ISP® Seal Check

The INTERSTER® Seal Check® is a standardized method in accordance with ISO 11607-2, and can be used for every sterilization pouch or reel in combination with any sealing machine with a rotating printing wheel. The Seal Check® displays a good seal and also shows an indication of possible sealing problems.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
3FSZB830802	180mm 250pcs	1	25

It is for use with Sterilization Pouches or Reels that are made of a medical grade paper and a laminate.

An incorrect temperature or pressure will be visible on the Seal Check®.

Sterilizing Consumables



Liquid Control Dye Penetration Test

Liquid Control Dye Penetration Tests are a new generation of ink tests for routine testing of your different sealing devices. Liquid Control Dye Penetration Tests are specially developed to check the integrity of the seals, produced with different heat sealers (for steam sterilisation paper and low temperature sterilisation using Tyvek® packaging)



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
3FSZB830816	Liquid Dye Test I Laminate/Medical Paper	1	30
3FSZB830818	Liquid Dye Test II Laminate/Tyvek®	1	30

- Test I suitable for film/medical grade paper packaging
- Test II suitable for film/Tyvek® packaging
- Easy, safe and quick test according to the ISO 11607-1 as written in Appendix B.

ISP® Sterilization Marker Pen

The ISP® Sterilization Environment Utility Marker is used for marking sterilization packaging material, with special autoclave low odour and non toxic ink.

The ink is suitable for steam, gamma, formaldehyde, plasma and ethylene oxide sterilization processes. The pen is designed for light touch writing to avoid damage to the packaging material.




REORDER	PENDING CODE	DESCRIPTION	PCS/UNIT	UNITS/CTN
3FSZB710824		Black 25pcs	1	15
33005030	3FSZB710826	Blue 25pcs	«	15
33005032	3FSZB710828	Red 25pcs	1	15

The solid writing point of the marker ensures a constant writing quality and long durability. The ink dries immediately and does not fade.

«Indent Item

ISP® Helix Test Type 2+6 Indicator Strip

The ISP® Helix Test System is the ideal process challenge device to perform the strongest challenge when sterilizing instruments. It can also be used to ensure maximum steam penetration where hollow instruments or instruments with small lumina are concerned. It is easy to interpret with a clear and even colour change. Complies with EN 867:5



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
3FSKS630802	100 Strips 3.5min Type 6 Ink for Helix	1	1

Non toxic Ink and Adhesive back

ISP® Steam Indicator Tape

The ISP® Steam Indicator tapes are intended for sealing packages to be steam sterilized.

They indicate exposure to saturated steam as a Type 1 indicator.

The non toxic Ink indicator has a light to dark color change in accordance with ISO11140-1 (Type 1)

The buff coloured tape is intended for use on Linen/Woven material.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
3FKLB410802	19mm x 50m (buff)	1	48
3FKLB410804	25mm x 50m (buff)	1	36
3FKLB410806	19mm x 50m (blue)	1	48
3FKLB410808	25mm x 50m (blue)	1	36

The blue coloured tape has a higher level of adhesion and is intended for use on Non Woven material.

Conforms to ISO11140-1.

Sterilizing Consumables





STERIFIT® Sterilization Wraps

STERIFIT® Sterilization Wraps have been developed with the latest SMMMS bonding technology, creating a more reliable and effective sterile barrier. They are suitable for steam, ethylene oxide, and plasma sterilization processes. Dual colouring allows for easy detection of strike through.

Benefits of STERIFIT® include:

- Constructed with Penta-Tech Layering
- Effective Sterile Barrier
- High Resistance to Tearing and Abrasion
- Highly Drapeable
- Available in Varying Strengths



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
STS001	Standard 30cm x 30cm	« 1000	1
STS002	Standard 43cm x 43cm	500	1
STS003	Standard 61cm x 61cm	500	1
STS004	Standard 71cm x 71cm	300	1
STS005	Standard 91.5cm x 91.5cm	150	1
STS006	Standard 100cm x 106cm	150	1
STS007	Standard 100cm x 140cm	100	1
STS008	Standard 121.5cm x 121.5cm	100	1
STS009	Standard 130cm x 152.5cm	« 100	1
STS010	Standard 137cm x 183cm	100	1
STS011	Standard 213.5cm x 213.5cm	« 50	1
STS012	Standard 100 x 150cm	« 50	1
STH001	Heavy Duty 71cm x 71cm	200	1
STH002	Heavy Duty 91.5cm x 91.5cm	150	1
STH003	Heavy Duty 106cm x 120cm	50	1
STH004	Heavy Duty 106cm x 140cm	« 50	1
STH005	Heavy Duty 121.5cm x 121.5cm	50	1
STH006	Heavy Duty 130cm x 152.5cm	50	1
STH007	Heavy Duty 137cm x 183cm	50	1
STU001	Ultra Heavy Duty 106cm x 120cm	50	1
STU002	Ultra Heavy Duty 106cm x 140cm	« 50	1
STU003	Ultra Heavy Duty 106cm x 160cm	« 50	1
STU004	Ultra Heavy Duty 121.5cm x 121.5cm	50	1
STU005	Ultra Heavy Duty 137cm x 183cm	40	1
STU006	Ultra Heavy Duty 137cm x 137cm	« 40	1

«Indent Item



Sterilizing Consumables



STERIFIT® 2 Dual Layer Sterilization Wraps

STERIFIT® 2 Dual layer wrap, for single wrap technique.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
STS102	STERIFIT® 2 Standard 60 x 60cm	250	1
STS103	STERIFIT® 2 Standard 76 x 76cm	150	1
STS104	STERIFIT® 2 Standard 91 x 91cm	75	1
STS106	STERIFIT® 2 Standard 121 x 121cm	50	1
STS109	STERIFIT® 2 Standard 101 x 101cm	75	1
STH101	STERIFIT® 2 Heavy Duty 91 x 91cm	75	1
STH102	STERIFIT® 2 Heavy Duty 101 x 101cm	75	1
STH103	STERIFIT® 2 Heavy Duty 121 x 121cm	25	1
STU102	STERIFIT® 2 Ultra Heavy Duty 121 x 121cm	25	1
STU103	STERIFIT® 2 Ultra Heavy Duty 137 x 183cm	25	1



Sterilization Transport Tray

Sentry Medical Sterilization Transportation Trays make transportation and storage of large or multiple sterile packs easier. The trays are used to prevent damage to the **STERIFIT®** wrapped packs before, during and after the sterilization process.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
STT001	560mm x 300mm x 35mm «	20	1

The trays can also be conveniently placed inside the autoclave during sterilization, making transportation of heavy or multiple packs easier. The trays are durable and are reusable.

«Indent Item

Apparel




Non-woven, spunbond polypropylene is used to make all OWEAR® range products. The weave of these products have an all-round consistent tensile strength, yet are soft enough to be worn all day long. Unlike many other types of apparel, OWEAR® allows greater airflow and breathability to reduce temperature build-up.

OWEAR® Berets & Bouffant Caps



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
RC001	Blue	250	4
RC001L	Blue 24inch	250	4
RC002	Green	250	4
RC003	Pink	250	4
RC004	White	250	4
RC005	Yellow	250	4
RC006X	Christmas Theme	250	1
RC007	Red	250	4
RC008	Dark Blue	250	4
RCD001	Round Cap Dispenser	1	5

OWEAR® Headscarves



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
SF001	Blue	250	2
SF002	Green	250	2
SF003	Pink	250	2

OWEAR® Hood

Ties



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
HD001	Blue Regular	250	4
HD002	Blue Small	250	4

OWEAR® Theatre Cap

Tie Back



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
TC005	Blue	250	8
TC006	Green	250	8
TC007	Pink	250	8
TC008	White	250	8
TC009	Dark Blue	250	8

Apparel



OWEAR® Isolation Gowns

Non Sterile | Latex Free

* Tested for Chemotherapy Use

^ Anterior & Sleeves Only

+ Back of Gown Only

When Donned with Gloves

All OWEAR® impervious and chemically tested Isolation Gowns comply with the ISO16603:2004, ISO16604:2004 and ISO6529:2006 standards, as well as meet AAMI Liquid barrier performance levels (excluding GN015P and GN016P). The gowns come in a variety of colours, and are designed for specific procedures to encompass comfort, functionality and protection. All gowns are latex free and have either cotton cuffs (Sof-Cuff) or elastic sonically bonded cuffs.

REORDER	DESCRIPTION	AAMI Level	PCS/ UNIT	UNITS/ CTN	IMPERVIOUS	SPLASH RESISTANT	CHEMICALLY TESTED*
GN001	Blue Splash Resistant	1	50	1		•	
GN001P	Blue Impervious Sof-Cuff	2	50	1	•#		•
GN002	White Splash Resistant	1	50	1		•	
GN003	Yellow Splash Resistant	1	50	1		•	
GN003P	Yellow Impervious Sof-Cuff	2	50	1	•#		•
GN004P	Purple Impervious Oncology	2	100	1	•		•
GN0041P	Purple Impervious Oncology Sof-Cuff	2	50	1	•#		•
GNX004	Mid Blue Short Sleeve Patient Modesty		100	1		•	
GN005P	Blue Impervious Sonically Bonded	2	50	1	•		•
GN006P	Yellow Impervious Sonically Bonded	2	50	1	•		•
GNX006	Mid Blue Short Sleeve Patient Modesty Bariatric		50	1		•	
GN008P	Blue Impervious Extra Long Sof-Cuff	2	50	1	•#		•
GN009P	White Polybreath Impervious Sonically Bonded	2	50	1	•		•
GN013P	Blue Hybrid Regular Sof-Cuff	2	50	1	•^#	•+	•^
GN015P	Blue Impervious Thumb Loop Regular		75	1	•		
GN016P	Blue Impervious Thumb Loop XL		75	1	•		
GN023P	Hybrid Gown Sof-Cuff Yellow Medium	2	50	1	•^#	•+	•^
GN024P	Hybrid Gown Sof-Cuff Yellow Large	2	50	1	•^#	•+	•^
GN025P	Hybrid Gown Sof-Cuff Yellow X-Large	2	50	1	•^#	•+	•^
GN028P	Hybrid Gown Sof-Cuff Blue Medium	2	50	1	•^#	•+	•^
GN029P	Hybrid Gown Sof-Cuff Blue Large	2	50	1	•^#	•+	•^
GN030P	Hybrid Gown Sof-Cuff Blue X-Large	2	50	1	•^#	•+	•^



Apparel



OWEAR® Aprons

4 Sizes | Individually Folded



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
AP001	610mm x 1220mm White	100	10
AP002	810mm x 1320mm White	100	6
AP003	810mm x 1420mm White	100	6
AP004	810mm x 1600mm White	100	6
AP006	810mm x 1420mm White	100	10
DISP-AP2	Acrylic Apron Dispenser	1	6

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
AP022	810mm x 1320mm Pink	100	6
AP032	810mm x 1320mm Blue	100	6

OWEAR® Briefs & Pelvic 'V'



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
BF001	Brief White	100	8
BF002	Brief White Large	75	8
PV001	Pelvic V White	100	8

OWEAR® Overshoe Protective Covers



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
OSN001	Non Slip	100	6
OSN002	Non Slip X-Large	100	6
OSR002	Regular	100	6
OSW001	Waterproof Plastic	100	10
OSW003	Waterproof Non Slip Overboot	100	6
OSW004	Waterproof Non Slip Below Knee	100	6

Apparel



Sentry Medical is the proud supplier of the PRIMAGARD® range of surgical masks in Australia. Developed in Canada and distributed throughout North America, all surgical masks are independently tested to comply with the requirements of AS4381:2015. The masks are available in a variety of protection levels; 16kPa (120mmHg) for procedures with moderate risk of fluid exposure and 21kPa (160mmHg) for procedures with high risk of fluid exposure.

The masks have longer ties for easier fitting, and are made from low allergenic components for comfort.

PRIMAGARD® Procedure 160 Mask

Level 3 | Latex Free | Anti-Fog | Earloop

- Anti-fog strip made of soft, open-celled foam, prevents moist air from rising up to fog glass.
- Fluid resistant to 21kPa (160mmHg), highest level of fluid protection based on AS4381:2015 standard.
- Suitable for procedures where there is a high risk of fluid splash, such as trauma, the ER or the ICU.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-1073	Anti-Fog, Earloop, Indigo, 4 Ply	50	10

- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 3 Barrier

PRIMAGARD® Procedure 160 Mask

Level 3 | Latex Free | Anti-Fog | Anti-Glare | Visor | Earloop

- Attached visor has anti-fog properties and wrap-around wings to protect the eyes from splash from the front and sides.
- Fluid resistant to 21kPa (160mmHg), highest level of fluid protection based on AS4381:2015 standard.
- Suitable for procedures where there is a high risk of fluid splash, such as trauma, the ER or the ICU.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-1563	Visor, Anti-Fog, Anti-Glare, Earloop, Indigo, 4 Ply	25	4

- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 3 Barrier.

PRIMAGARD® Surgical 160 Mask

Level 3 | Latex Free | Adhesive Nose Strip | Anti-Fog | Tieback

- Anti-fog strip made of skin-safe adhesive tape, prevents moist air from rising up to fog glasses and eye protection.
- Fluid resistant to 21kPa (160mmHg), highest level of fluid protection based on AS4381:2015 standard.
- Suitable for surgical procedures where there is a high risk of fluid splash or surgical debris, such as orthopedic surgery or trauma.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-4053	Adhesive Nose Strip, Anti-Fog, Tieback, Green, 4 Ply	50	6

- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 3 Barrier.

Apparel



PRIMAGARD® PRIMASOFT® Surgical 160 Mask

Level 3 | Latex Free | Softbill | Tieback

- PRIMASOFT® masks are silky, smooth and irritant-free, ideal for those with sensitivities to other mask.
- Softbill-styleQ keeps the mask away from the wearer's mouth, for increased comfort.
- Fluid resistant to 21kPa (160mmHg), highest level of fluid protection based on AS4381:2015 standard.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-3033	Softbill, Irritant Free, Tieback, White, 4 Ply	50	6

- Suitable for surgical procedures where there is a high risk of fluid splash or surgical debris, such as orthopedic surgery or trauma.
- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 3 Barrier.

PRIMAGARD® PRIMASAURUS™ Surgical 160 Mask

Level 3 | Latex Free | Tieback

- Fluid resistant to 21kPa (160mmHg), highest level of fluid protection based on AS4381:2015 standard.
- Suitable for surgical procedures where there is a high risk of fluid splash or surgical debris, such as orthopedic surgery or trauma.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-2083	Tieback, Dinosaur Print, 4 Ply	50	6

- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 3 Barrier.

PRIMAGARD® Surgical 160 Mask

Level 3 | Latex Free | Anti-Fog | Tieback

- Anti-fog strip made of soft, open-celled foam, prevents moist air from rising up to fog glasses and eye protection.
- Fluid resistant to 21kPa (160mmHg), highest level of fluid protection based on AS4381:2015 standard.
- Suitable for surgical procedures where there is a high risk of



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-4073	Anti-Fog, Tieback, Indigo, 4 Ply	50	6

- fluid splash or surgical debris, such as orthopedic surgery or trauma.
- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 3 Barrier.

PRIMAGARD® Surgical 160 Mask

Level 3 | Latex Free | Anti-Fog | Anti-Glare | Visor | Tieback

- Attached visor has anti-glare film and anti-fog properties, and wrap-around wings to protect the eyes from splash from the front and sides.
- Fluid resistant to 21kPa (160mmHg), highest level of fluid protection based on AS4381:2015 standard.
- Suitable for surgical procedures where there is a high risk of



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-5073	Anti-Fog, Anti-Glare, Visor, Tieback, Indigo, 4 Ply	25	4

- fluid splash or surgical debris, such as orthopedic surgery or trauma.
- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 3 Barrier.

Apparel



PRIMAGARD® Surgical 160 Mask

Level 3 | Latex Free | Plume Protector | Tieback

- Grey colour is the industry standard for laser masks.
- Fluid resistant to 21kPa (160mmHg), highest level of fluid protection based on AS4381:2015 standard.
- Suitable for surgical procedures in which a laser, harmonic scalpel or other smoke/vapor producing electro-surgical



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-6003	Plume Protector, Tieback, Grey, 4 Ply	50	6

techniques are being used and where there is a high risk of fluid splash.

- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 3 Barrier.

PRIMAGARD® Surgical 120 Mask

Level 2 | Latex Free | Tieback

- Fluid resistant to 16kPa (120mmHg).
- Suitable for procedures where there is a moderate risk of splash.
- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 2 Barrier.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PM4-305	Tieback, Blue, 3 Ply	50	6

PRIMAGARD® Procedural 120 Mask

Level 2 | Latex Free | Earloops

- Fluid resistant to 16kPa (120mmHg).
- Suitable for procedures where there is a moderate risk of splash.
- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 2 Barrier.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PM4-306	Earloop, Blue, 3 Ply	50	10

PRIMAGARD® Surgical 120 Mask

Level 2 | Latex Free | Anti-Fog | Tieback

- Anti-fog strip prevents moist air from rising up to fog glasses.
- Fluid resistant to 16kPa (120mmHg).
- Suitable for procedures where there is moderate risk of splash.
- Conforms to AS4381: 2015 & ASTM F2100-11.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PM4-311	Anti-Fog, Tieback, Blue, 3 Ply	50	6

- Level 2 Barrier.



Apparel

PRIMAGARD® Procedural 120 Mask

Level 2 | Latex Free | Anti-Fog | Earloops

- Anti-fog strip prevents moist air from rising up to fog glasses.
- Fluid resistant to 16kPa (120mmHg).
- Suitable for procedures where there is a moderate risk of splash.
- 12 Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 2 Barrier.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PM4-312	Anti-Fog, Earloop, Blue, 3 Ply	50	10

PRIMAGARD® Surgical 120 Mask

Level 2 | Latex Free | Anti-Fog | Tieback

- Anti-fog strip made of soft, open-celled foam, prevents moist air from rising up to fog glasses and eye protection.
- Fluid resistant to 16kPa (120mmHg).
- Suitable for surgical procedures where there is a moderate risk of fluid splash.
- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 2 Barrier.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-2492	Anti-Fog, Tieback, Denim Blue, 4 Ply	50	6

PRIMAGARD® Procedural 120 Mask

Level 2 | Latex Free | Anti-Fog | Anti-Glare | Visor | Earloop

- Attached visor has anti-glare and anti-fog properties, and wrap-around wings to protect the eyes from splashes from the front and sides.
- Fluid resistant to 16kPa (120mmHg).
- Suitable for procedures with moderate risk of splash.
- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 2 Barrier.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-1592	Anti-Glare, Anti-Fog, Visor, Earloop, Denim Blue, 4 Ply	25	4

Apparel



PRIMAGARD® Surgical 120 Mask

Level 2 | Latex Free | Anti-Fog | Anti-Glare | Visor | Tieback

- Attached visor has anti-glare and anti-fog properties, and wrap-around wings to protect the eyes from splashes from the front and sides.
- Fluid resistant to 16kPa (120mmHg).
- Suitable for surgical procedures where there is a moderate risk of fluid splash.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-2592	Anti-Glare, Anti-Fog, Visor, Tieback, Denim Blue, 4 Ply	25	4

- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 2 Barrier.

PRIMAGARD® PRIMASAURUS™ Standard Child mask

Non Rated | Latex Free | Earloop

- Small size designed for children.
- Used for respiratory etiquette to contain droplets from children who are symptomatic.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-1680	Earloop, Dinosaur Print, 3 Ply	50	10

PRIMAGARD® Standard 80 Mask

Level 1 | Latex Free | Tieback

- Very breathable.
- Fluid resistant to 11kPa (80mmHg).
- Suitable for surgical procedures where there is a minimal risk of fluid splash.
- Conforms to AS4381: 2015 & ASTM F2100-11.
- Level 1 Barrier.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PG4-2321	Fluid Resistant, Tieback, Blue, 3 Ply	50	6

Clear Visor

Latex Free



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
CVS001	33cm x 22cm	50	1
CVS002	33cm x 25cm	50	1

Apparel



Mask Summary Chart

REORDER	DESCRIPTION	LEVEL	PCS/UNIT	UNITS/CTN	TIE-BACK	EAR-LOOP	VISOR	ANTI-FOG	ANTI-GLARE	Delta P, AVG	FLUID RESISTANCE
PG4-1073	PRIMAGARD® Procedural 160 Mask	3	50	10		•		•		3.0	21kPa/160mmHg
PG4-1563	PRIMAGARD® Procedural 160 Mask	3	25	4		•	•	•	•	3.0	21kPa/160mmHg
PG4-4053	PRIMAGARD® Surgical 160 Mask	3	50	6	•			•		3.0	21kPa/160mmHg
PG4-3033	PRIMASOFT® Surgical Softbill Mask 160, Irritant Free	3	50	6	•					3.2	21kPa/160mmHg
PG4-2083	PRIMASAURUS™ Adult Dinosaur Print Mask 160	3	50	6	•					3.2	21kPa/160mmHg
PG4-4073	PRIMAGARD® Surgical Mask 160	3	50	6	•			•		3.0	21kPa/160mmHg
PG4-5073	PRIMAGARD® Surgical Mask 160, Visor	3	25	4	•		•	•	•	3.0	21kPa/160mmHg
PG4-6003	PRIMAGARD® Surgical Mask 160, Plume Protector	3	50	6	•					3.0	21kPa/160mmHg
PM4-305	PRIMAGARD® Surgical Mask 120	2	50	6	•					2.3	16kPa/120mmHg
PM4-306	PRIMAGARD® Procedural Mask 120	2	50	10		•				2.3	16kPa/120mmHg
PM4-311	PRIMAGARD® Surgical Mask 120	2	50	6	•			•		2.3	16kPa/120mmHg
PM4-312	PRIMAGARD® Procedural Mask 120	2	50	10		•		•		2.3	16kPa/120mmHg
PG4-2492	PRIMAGARD® Surgical Denim Blue Mask 120	2	50	6	•			•		2.4	16kPa/120mmHg
PG4-1592	PRIMAGARD® Procedural Mask 120, Visor	2	25	4		•	•	•	•	2.4	16kPa/120mmHg
PG4-2592	PRIMAGARD® Surgical Mask 120, Visor	2	25	4	•		•	•	•	2.4	16kPa/120mmHg
PG4-1680	PRIMASAURUS™ Child Dinosaur Print Face Mask		50	10		•				≤ 3.2	-

Dispensers

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
DIS001	Mask Dispenser - PM4-305, PM4-311, PG4-2321, PG4-2492	1	6
DIS002	Mask Dispenser - PG4-2083, PG4-4073, PG4-3033, PG4-6003	1	6
DIS003	Mask Dispenser - PG4-5073, PG4-2592, PG4-1592	1	10
DIS004	Mask Dispenser - PM4-306, PM4-312	1	10
PM4-BH01	Mask Holder - PM4-305, PM4-311, PG4-2321, PG4-2492	1	6
DISP-AP2	Apron Dispenser	1	6
DIS100	Thumb Loop Gown Wall Mount	1	100
RCD001	Round Cap Dispenser	1	5
DISP-AB	Acrylic Bandage Dispenser	1	4
VBD001	Emesis Bags Dispenser	1	12
VED001	VeraSoft® Handypack Dispenser 3mm/clear	1	48

All Dispensers are made available for use with Sentry Medical Products only




Miscellaneous



VBag® Emesis Bags

Sentry Medical Emesis Bags have a firm polypropylene twist and lock rim to provide easy control, and help prevent spillage after use. They isolate splash back and have a graduated indicative measurement system to 1500ml, allowing measurement of contents in the bag. On completion of use, the bag is easily disposed of with no mess. Simply twist the bag and lock into the keyhole slot on the rim.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
VB002	Red Ring Emesis Bag	50	10
VB011	Sickeze® 3pcs (6pks)	10	1
VBD001	Acrylic Emesis Bag Dispenser	1	12


Rims are colour coded red for infection control. Convenient dispenser also available. Measurement is approximate to the bag.



Pill Crusher & Pill Crusher Pouches

Sentry Pill Crusher is used to quickly crush tablets into easy-to-consume powder form in the corresponding pouch.

The pouch is formed with rib and beveled for ease of total powder removal.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PCS001	Pill Crusher	1	6
PCP001	Pill Crusher Pouch (Ribbed Beveled Seam)	1000	8

Miscellaneous



VeraSoft® Bath Wipes

Bathing System

VeraSoft® Bath Wipes are a convenient rinse free bathing system, used to reduce the risk of infection through control of cross contamination. These disposable and pre-moistened wash cloths are pH neutral, may be heated in the microwave, and have a gentle no rinse hydrating cleanser.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
VER001	20cm x 20cm Regular 8pcs	24	1
VER002	20cm x 20cm Antibacterial 8pcs	24	1

VeraSoft® Moist Wipes

Refresher Towel

VeraSoft® Moist Wipes are a refreshing and convenient disposable wipe, that are infused with aloe vera. Also ideal for use as an infant wipe.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
VER010	32cm x 22.5cm 64pcs	12	1
VER011	18cm x 20cm 25pcs	36	1
VED001	Acrylic 25 Pack Dispenser 3mm	1	48

VeraSoft® Multipurpose Towels

Dry

VeraSoft® Multipurpose Towels are low linting, highly absorbent and extremely strong and durable. Ideal for cleaning instruments and lining trays prior to sterilization. Available in two sizes, including a handy top shelf dispenser pack.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
VER020	30cm x 35cm 100pcs	12	1
VER021	30cm x 60cm 100pcs	6	1

VeraSoft® Shampoo Cap

Moist

VeraSoft® Rinse Free Shampoo Caps are the convenient way for cleansing hair without the hassle of rinsing out conventional shampoos and conditioners. Ideal for patients who have limited mobility.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
VER030	Shampoo Cap	30	1

Simple three step process:

1. Apply the shampoo cap
2. Massage hair through the cap until hair appears saturated
3. Remove cap and towel dry hair

Leaves hair feeling clean and smelling fresh.

Miscellaneous



Cotton Huck Towels

Non Sterile | Blue

Sentry Medical Blue Cotton Huck Towels are constructed from 100% cotton. Ideal for cleaning or hand drying. The towels are strong and durable, and can be used wet or dry.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
BCR001	Cotton Huck 45cm x 60cm	240	1

Bed Sheets

Draw Sheets | Stretcher Sheets

Two different Bed Sheets are available, Draw Sheets and Stretcher Sheets.

Sentry Medical Draw Sheets are constructed from impervious blue polyethylene, and are treated with anti-static solution during manufacturing.

The Basic Stretcher Sheet is white, made from polypropylene.

The Blue Stretcher Sheet is constructed from a non-woven, spunbond polyethylene fabric and has elastic on the 2 shorter edges.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
DS001	Plastic Draw Sheet 915mm x 1520mm	300	1
SSH001	Stretcher Sheet 1010mm x 2380mm	100	1
SSH002	Stretcher Sheet, Blue 1010mm x 2380mm	100	1

Pillow Sleeve

Non Woven

Sentry Medical Pillow Sleeves are light blue in colour, and constructed from a non-woven, spunbond polyethylene fabric.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PS001	Pillow Sleeve	100	4

Miscellaneous



Shroud White

Non Sterile | Latex Free

A polypropylene antistatic shroud covers or protects. Typically used in reference around burials. It is latex free.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
SHROUDWHT	Shroud White	20	1

Bluey Underpads

The Bluey Underpad has a waterproof backing and soft dry layering cover, with super absorbent material to lock away wetness.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
SCUP5	Bluey Underpad 5ply 42cm x 59cm	250	1

Body Bag With Handles

Latex Free | Envelope Zipper Style

Sentry Medical Body Bags with Handles are strong and durable during transportation and handling. They are highly resistant to leakage and contain a two way zipper on the top of the bag for easy access.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
BBA002	Standard 2.4m x 1.15m White	10	1
BBA003	Bariatric 1.4m x 2.6m White	3	1
BBA006	Paediatric 1.15m x 1.5m White	10	1
BBA007	Standard 2.4m x 1.15m Blue	10	1

Miscellaneous



Polythene Oversleeves White

Polythene oversleeves are ideal for a variety of medical applications.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
05189	Polythene Oversleeves White	1000	1

Patient Belongings Bag

White Plastic

Sentry Medical Patient Belongings Bags are constructed from a durable polyethylene. They can be clearly labelled with a patient name and room number. A draw string is provided at the opening of the bag for ease of closing and carrying.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
PB001	440mm x 610mm	300	1

Arm Glove

armRx® arm gloves are perfect for protecting plaster casts, dressings and “permanent” bandages from moisture and water penetration. Double seal, hook and loop closures make application and removal simple for carers or patients alike.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
122115	Arm Glove Standard (in Dispenser)	30	12

Miscellaneous



Foam Prep Sponge

Non Sterile

Sentry Medical Foam Prep Sponges are used for cleaning the incision site. Antiseptic solution is sprayed onto the site and wiped off with the foam sponge.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
FPS001	6cm x 5cm x 2.5cm	1000	1

Irrigation Bulb Syringe

Non Sterile

The Sentry Medical Irrigation Bulb Syringe is an easy and convenient way to irrigate patient wound sites. The device consists of a hard plastic syringe with a soft rubber bulb that is squeezed to commence irrigation. It can also be used in theatre.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
BSY001	60ml	50	1

Plastic Forceps

Non Sterile | Disposable

Plastic Forceps are intended for grasping and holding objects. The separated tip enables secure grasping while the interlocking teeth allow easy gripping on slippery or flimsy materials.




REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
FOR001	Blue Plastic Forceps BULK	1000	1

Disposable Dressing Tray Type 3

Non Sterile | Latex Free

Disposable Dressing Tray Type 3 is intended to hold various instruments and dressings needed to dress wounds.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
DDT003	Standard Type 3 Tray BULK	1000	1

« Indent Item

Miscellaneous



Tongue Depressor

Non Sterile | Latex Free

Sentry Medical Tongue Depressors are constructed from 100% non coloured plantation pine timber, and finished with a smooth polished edge.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
TDNS001	Wooden Tongue Depressor	100	50
TDNS002	Wooden Tongue Depressor Single Pack	50	20

NoPress® Eye Pressure Shield

Non Sterile | Latex Free

No Press® is designed to protect the anaesthetised patient's eyes from accidental pressure related injury due to positioning, draping and operational techniques or procedures.

The single piece shield coupled with thick medical grade foam distributes pressure evenly around the bony orbital margins. Biocompatible adhesive holds the shield firmly in place and integrated aeration holes minimise condensation and ensures visualisation of eyes.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
IVG-03	NoPress® Eye Pressure Shield	20	1

First Aid



Resuscitation Face Shield

Disposable

Sentry Medical Resuscitation Face Shields assist in the application of mouth to mouth resuscitation in emergency and First Aid situations. These shields have an oval shaped mouth piece incorporating a one way valve to allow air to be transferred into the patient without cross contamination.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
RES001	Face Shield	50	12
RES002	CPR Mask	30	1

Instant Ice Pack

An easy way to create instant ice to alleviate pain and reduce swelling from minor injuries.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ICP001	Instant Ice Pack Large 12cm x 24cm	24	1
ICP003	Instant Ice Pack Small 8.8cm x 16cm	40	1

Paper Tape

Non Sterile | Latex Free | Individually Wrapped

A low allergy paper tape that provides general purpose medical securement. Available in a variety of widths, with and without dispenser.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
HYP001	1.25cm x 9.1m	24	30
HYP002	2.5cm x 9.1m	12	30
HYP003	5cm x 9.1m	6	80
HYP004	7.5cm x 9.1m	4	30
HYP100	2.5cm x 9.1m with Dispenser	12	30
HYP101	1.25cm x 9.1m with Dispenser	24	30

First Aid



Alcohol Swabs

Non-woven swab moistened with 70% isopropyl alcohol ideal for cleansing skin e.g. prior to injection.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ALS001	Alcohol Swabs 3cm x 3cm	200	100

Povidone Iodine Antiseptic Wipes

A non-woven swab carrying 10% povidone iodine for skin and wound cleansing.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
POV001	Povidone Iodine Antiseptic Wipes 3cm x 3cm	200	50

Gloves Nitrile

Non Sterile | Latex Free | Nitrile Rubber | Powder Free

Protect yourself from contaminated items with non-sterile, latex free gloves.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
GNL200	Nitrile Powder Free Gloves Large (1 pair)	600	1
GNL210	Gloves Nitrile (5 pairs)	150	1

First Aid



Aqafy® Eyewash & Wound Irrigation

Sterile | Normal Saline

Aqafy® is a sterile solution of 0.9% sodium chloride, presented in an 'easy squeeze' plastic ampoule.

This solution is ideal

- for irrigating eyes to relieve discomfort
- to flush foreign substances from eyes
- to remove loose materials from wounds



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
EWS001	Aqafy® 15ml Saline Eye Wash Solution	360	1
EWS002	Aqafy® 30ml Saline Eye Wash Solution	180	1

Stainless Steel Forceps

Single use

Sentry Medical Forceps are forged from surgical grade stainless steel, ensuring exceptional strength and resistance to corrosion. Grooves on the forcep ends provide additional grip.

Ideal for First Aid kits.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
SSF001	Stainless Steel Forceps 125mm	100	1

Stainless Steel Sharp/Blunt Tipped Scissors

Single use

Sentry Medical Sharp/Blunt Tipped Scissors are forged from surgical grade stainless steel, ensuring sharpness and exceptional resistance to corrosion. These scissors have one tip sharp, one tip blunt for compatibility to different uses.

Ideal for First Aid uses.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
SSS001	Stainless Steel Scissors 125mm	50	16

First Aid



Splinter Probes

Hygienically Sealed

Sentry Medical Splinter Probes come in a packet of 5, and are packaged for individual use. Ideal for First Aid kits.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
SPR001	Disposable Splinter Probe 100pcs	1	50

Safety Pins

Safety pins, a helpful addition to every first aid kit.

Useful for securing bandages, slings, etc.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
SP003	Safety Pins Assorted Sizes 400 Dozens	1	1

Thermal Shock Blanket

Sentry Medical Thermal Shock Blanket offers effective protection against excessive heat and cold. Highly recommended as an integral part of any survival equipment, and in medical emergencies for shock treatment.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
THB001	130cm x 185cm	100	1

Blue Metal Detectable Dressing Strips

Sterile | Latex Free | Waterproof

The dressing strips are waterproof, extra wide, flexible, conformable and adhere firmly. The low adherent pad cushions the wound for comfort during wear.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
7007155	25mm x 72mm, 100 Strips	6	12
RBVM50	25mm x 72mm, 50 Strips	6	12

First Aid



First Aid Bag

Sentry Medical First Aid Bags are a black fold over bag with a velcro fastening. The bags have 'First Aid Kit' and the first aid symbol screen printed on one side. These bags measure 190mm x 318mm when opened.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
FABS005	Fold Over Black	1	1

Amputated Parts Bags

Transport amputated parts in one of these bags to ensure the part stays as clean as possible and gets to its destination. 3 sizes to carry small sized body parts.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
AMB001	Amputated Parts Bags 3pcs assorted size	450	1

Waste Bags



All Sentry Medical Clinical Waste Bags are constructed from a variety of low linear density polyethylene, high density polyethylene and ethyl-butylene-polyethylene. These materials contain organic colouring, provide strength at quoted thicknesses. These bags, produce no toxic by-products on incineration. A variety of sizes are available to fit small to large bins.

Sentry Medical also provides a full range of sizes in garbage bags and locker bags. All products are supplied in easy-to-use dispenser packs for convenient storage.

Clinical Waste Bag

Provided in two colours - yellow to visibly indicate cross-infection risk and purple to indicate cytotoxic risk - these bags are ideal for disposal of waste with need to identify particular hazards.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
CW003	56.5cm x 1m 55L 50um	200	1
CW004	51cm x 660mm 27L 20um	1000	1
CW010	80cm x 1.2m 100L 150um	50	1
CW017	Cytotoxic Purple 60cm x 84cm 70um	200	1

Garbage Bag

Made from HDPE for strength and security, whether the garbage you are managing is lightweight or heavy in nature or needs odour control, there will be a Sentry Medical Garbage Bag that will suit the purpose.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
GB002	75cm x 90cm 17um 73L Blue Lemon Scent	500	1
GB010	76cm x 107cm 25um 80L Natural	250	1
GB011	76cm x 95cm 25um 82L Black Heavy Duty	300	1

Locker Bag

Self-Adhesive

The Sentry Medical Locker Bag is the ideal size for management of daily waste at patients bedside. Self-adhesive strip makes attachment to locker or tray table easy and convenient.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
LB001	27.5cm x 20.5cm	1000	6

Advanced Wound Care



osmocel® Hydroporous Foam

Sterile | Non-Adhesive | Latex Free

osmocel® Hydroporous Foam is an ideal dressing to manage wounds with different needs and moisture levels

1. as a primary dressing on lightly exuding wounds, or
2. as a secondary dressing when other wound contact layers are needed.

The fine pores retain moisture at the wound surface, while excess exudate travels into upper pores of the foam to be eliminated as moisture vapour, increasing dressing wear time. Its insulation properties and low density are additional features for enhanced wound healing.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
AWF050	5cm x 5cm	10	10
AWF051	6cm x 8.5cm	5	10
AWF052	10cm x 10cm	10	10
AWF053	10cm x 20cm	5	10
AWF054	15cm x 15cm	5	10
AWF055	20cm x 20cm	«	5
AWF180	Vari-Site 12cm x 19cm	5	10

« Indent Item

osmocel® Si Silicone Hydroporous Foam

Sterile | Latex Free

osmocel® Si Silicone Hydroporous Foam is an all-in-one dressing that protects both skin and wound tissue.

Maintenance of a moist wound surface, absorption of excess exudate, provision of thermal insulation and incorporating gentle adherence, this dressing provides the benefits of modern wound healing and encourages less frequent dressing changes.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
AWS002	Bordered 5cm x 12.5cm	5	10
AWS006	Bordered 15cm x 20cm	5	10
AWS007	Bordered 10cm x 20cm	5	10
AWS008	Bordered 10cm x 25cm	«	5
AWS009	Bordered 10cm x 30cm	«	5
AWS050	Non-Bordered 5cm x 5cm	10	10
AWS051	Non-Bordered 7.5cm x 7.5cm	10	10
AWS052	Non-Bordered 10cm x 10cm	10	10
AWS053	Non-Bordered 15cm x 15cm	«	5
AWS054	Non-Bordered 15cm x 20cm	«	5
AWS100	Petal Bordered 7.5cm x 7.5cm	10	10
AWS101	Petal Bordered 10cm x 10cm	10	10
AWS102	Petal Bordered 15cm x 15cm	5	10
AWS150	Sacral Bordered 18cm x 18cm	5	10
AWS151	Sacral Bordered 23cm x 23cm	«	5
AWS160	Sacral Butterfly Bordered 18cm x 18cm	«	5
AWS161	Sacral Butterfly Bordered 23cm x 23cm	5	10
AWS180	Vari-Site 12cm x 19cm	5	10

« Indent Item

Advanced Wound Care



osmocol® Hydrocolloid Dressing

Sterile | Latex Free

osmocol® Hydrocolloid Dressings provide optimal healing conditions while allowing you to monitor healing progress – reducing the need for unnecessary dressing changes.

Post-operative and minor wounds benefit most from the undisturbed healing osmocol® supports, as it is transparent and occlusive.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
AWH001	Thin 5cm x 5cm	10	10
AWH002	Thin 5cm x 10cm	10	10
AWH003	Thin 5cm x 20cm	10	10
AWH004	Thin 10cm x 10cm	10	10
AWH005	Thin 15cm x 15cm	5	10
AWH006	Thin 20cm x 20cm	«	5
AWH050	Standard 10cm x 10cm	10	10
AWH051	Standard 15cm x 15cm	5	10
AWH052	Standard 20cm x 20cm	«	5
AWH075	Bordered 10cm x 10cm	10	10
AWH100	Bordered 14cm x 14cm	«	5
AWH150	Sacral Bordered 15cm x 18cm	«	5

« Indent Item

osmonate® Calcium Alginate Dressing

Sterile | Latex Free

osmonate® Calcium Alginate Dressing is a naturally derived gelling fibre dressing that assists with exudate and haemostasis management.

As the dressing gels, it conforms to contours of the wound bed, ensuring a moist wound environment is sustained, whilst its easy removal ensures minimal disturbance to new or fragile tissue. With 100% calcium ion availability, the dressing may assist in the control of minor bleeding.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
AWC001	Sheet 5cm x 5cm	10	10
AWC002	Sheet 10cm x 10cm	10	10
AWC003	Sheet 10cm x 20cm	5	10
AWC004	Rope 30cm	5	10

Advanced Wound Care



IMPERVIA® Neosorb Superabsorbent Dressing

Sterile | Latex Free

IMPERVIA® Neosorb Superabsorbent Dressing quickly absorbs and binds large amounts of wound exudate into its 'heart' of Superabsorbent polymers. The outer "envelope" is

- sealed on all sides preventing sideways leakage
- larger than absorbent core allowing outward, not upward, swelling
- is ultrasonically bonded (adhesive-free construction) increasing conformability and minimising risk of possible skin reactions.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
NWC040	10cm x 10cm	20	20
NWC041	15cm x 15cm	20	20
NWC042	15cm x 20cm	20	20
NWC043	20cm x 20cm	20	20
NWC044	20cm x 30cm	20	10
NWC045	20cm x 45cm	10	10

Zorflex® Wound Contact Layer

Sterile | Latex Free

Zorflex® Long Fibre Activated Carbon Cloth provides an excellent alternative to traditional antimicrobial dressings that may contain chemical agents such as silver, iodine or PHMB.

Zorflex® offers **three actions** that promote wound healing and patient wellbeing

1. Chemical free **antimicrobial** action
2. **Electrical conduction** for accelerated tissue formation
3. **Reduction** of MMPs, pain and odour.

Used as a primary wound contact layer, the secondary dressing of choice should be determined by the level of exudate.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ZF05	5cm x 5cm	10	14
ZF10	10cm x 10cm	10	10
ZF15	15cm x 25cm	10	5
ZF20	10cm x 20cm	10	10

Advanced Wound Care



BURNSHIELD® Hydrogel

Bottle | Sachet | Spray | Tube | Latex Free

BURNSHIELD® Hydrogel is a liquid treatment that provides pain relief to minor burns, and may minimise risk of infection. The **BURNSHIELD®** Gel is designed to cool and protect burns or other heat stress injuries, and prevents burn progression and further tissue damage.

Ideal for

- sunburn
- superficial/partial thickness burns
- remoistening of **BURNSHIELD®** dressings



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
550004	3.5ml Sachets	1200	1
550065	Blott Blister Sachets Pack (5)	« 100	1
550009	25ml Tube	72	1
550000	50ml Bottle	72	1
550016	Hydrogel Tube & Dispenser (12)	« 12	1
550081	50ml Hydrogel Spray	48	1
550003	125ml Hydrogel Spray	48	1

« Indent Item



BURNSHIELD® Dressings

Latex Free

BURNSHIELD® Dressings are a soft comfortable foam layer, saturated with **BURNSHIELD®** Hydrogel, specifically designed to cover and protect burnt areas.

The dressings cool, sooth and moisturise burnt areas and may provide protection against infection.

Ideal for **First Aid** use on any level of burn.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
901111	10 x 10cm	250	1
550006	10 x 10cm 15pc Dispenser Pack	« 10	1
900906	20 x 20cm	50	1
901001	20 x 45cm	30	1
901007	60 X 40cm	12	1
550013	Mini Burn Kit	100	1
901000	Face Mask 20cm x 45cm	30	1
901135	Face Mask 60cm x 40cm	120	1

« Indent Item



Compression Therapy Range



lymphlex® REDUCE - Two Layer Reusable Compression Therapy Kit

Non Sterile | Latex Free | No Clips

lymphlex® REDUCE Short Stretch Compression Therapy Kit contains all the components required to apply safe and efficacious compression therapy for the management of Venous Leg Ulcers. With padding and short stretch compression bandages that are washable up to 50 times, and individually replenishable single use items, lymphlex® REDUCE is both economically and environmentally responsible.

Kit Components:

1 Set (x2)

- 1 x lymphlex® AIR Reusable Foam Padding Bandage 8cm x 2.5m (LOR021)
- 1 x lymphlex® AIR Reusable Foam Padding Bandage 10cm x 2.5m (LOR022)
- 1 x lymphlex® Short Stretch Compression Bandage 8cm x 5m (LOR003)
- 1 x lymphlex® Short Stretch Compression Bandage 10cm x 5m (LOR004)

Shared Components

- 1 x lymtube® Tubular Overlay 10cm x 20m (LOR043)
- 1 x CoWrap® lite Lightweight Cohesive Bandage 8cm x 20m (LOR055)
- 1 x lymfix® Fixation Tape 2.5cm x 5m (LOR050)
- 1 x Laundry Bag

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
LOR100	Two Layer Reusable Compression Therapy Kit	1	4



lymphlex® AIR Reusable Foam Padding Bandage

Non Sterile | Latex Free | Reusable | No Clips

lymphlex® AIR is a fine pored, soft-to-the-touch foam bandage for use as reusable padding under compression therapy bandages. Independently tested to be washable up to 50 times.

lymphlex® AIR is an individually replenishable component of the lymphlex® REDUCE Compression Therapy Kit.

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
LOR021	8cm x 2.5m	1	24
LOR022	10cm x 2.5m	1	24



Compression Therapy Range



lymphlex® Short Stretch Compression Bandage

Non Sterile | Latex Free | No Clips

lymphlex® Short Stretch Compression Bandage is constructed from 100% cotton. Designed to provide a sub-bandage pressure of 30-40mm Hg when applied at full stretch, with 1/2 to 2/3 overlaps.

Ideal for use in phlebological conditions, lymphoedema management and soft tissue injury support.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
LOR003	8cm x 5m	1	300
LOR004	10cm x 5m	1	240

lymphlex® Short Stretch Compression Bandage is an individually replenishable component of the lymphlex® REDUCE Compression Therapy Kit. Independently tested to be washable up to 50 times.

lymtube® Tubular Overlay

Non Sterile | Latex Free | 100% Bamboo

lymtube® Tubular Overlay is knitted from OEKO-TEX® STANDARD 100 certified bamboo fibre, a more environmentally responsible fibre than cotton. It's smooth finish helps it glide on easily and its soft texture makes it more comfortable to wear than cotton.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
LOR042	8cm x 20m	1	60
LOR043	10cm x 20m	1	48
LOR044	12cm x 20m	1	48

lymtube®, the skin protection layer of the lymphlex® REDUCE Compression Therapy Kit, is an individually replenishable component. It may also be used as a lightweight, tubular retention bandage.



CoWrap® lite - Lightweight Cohesive Bandage

Non Sterile | Latex Free

CoWrap® lite is a lightweight, open weave, cohesive bandage designed as damage-free securement for lymphlex® AIR foam bandages.

It may also be used as a gentle retention device especially where latex is undesirable.

CoWrap® lite Lightweight Cohesive Bandage is an individually replenishable component of the lymphlex® REDUCE Compression Therapy Kit.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
LOR055	8cm x 20m	1	80

Compression Therapy Range



lymfix® Fixation Tape

Non Sterile | Latex Free

lymfix® Fixation Tape is constructed from plain cotton weave fabric with a zinc oxide, latex free adhesive.

lymfix® is designed to be used to secure lymphlex® compression bandages, and to support the “brown on brown” recommendation.

May also be used as a securing tape where latex is undesirable.

lymfix® Fixation Tape is an individually replenishable component of the lymphlex® REDUCE Compression Therapy Kit.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
LOR050	2.5cm x 5m	1	30

lymphlex® TLC - Two Layer Cohesive Short Stretch (Inelastic) Compression Kit

Non Sterile | Latex

lymphlex® TLC Short Stretch Cohesive Kit combines the skin friendly lymphlex® AIR+ padding layer and lymphlex® SU short stretch cohesive compression bandages to create a safe and easy to apply compression therapy kit for the short term management of Venous Leg Ulcers. It is the perfect compression therapy option where laundering is problematic, or outside of local infection control policies or guidelines.

REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
LOS100	Two Layer Cohesive Short Stretch (Inelastic) Compression Kit	1	20

Kit Components:

- 1 x lymphlex® AIR+ Padding Bandage 10cm x 4m (LOS022)
- 1 x lymphlex® SU Short Stretch (Inelastic) Cohesive Compression Bandage 8cm x 6m (LOS003)
- 1 x lymphlex® SU Short Stretch (Inelastic) Cohesive Compression Bandage 10cm x 6m (LOS004)



lymphlex® AIR+

Non Sterile | Latex Free

lymphlex® AIR+ is a foam padding bandage with an integrated skin protection layer. Designed to be skin friendly, this padding layer is ideal for single use under compression therapy bandages such as lymphlex® SU.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
LOS021	8cm x 4.5m	1	18
LOS022	10cm x 4m	1	24

Compression Therapy Range



lymphlex® SU - Short Stretch (Inelastic) Cohesive Compression Bandage

Non Sterile | Latex Free

lymphlex® SU Short Stretch (Inelastic) Cohesive Compression bandages are designed to provide support and management in the treatment of ulcers caused by vascular diseases such as venous insufficiency. They can also be used as maintenance therapy in the management of lymphoedema, once decongestion has been achieved.

Available in two widths, lymphlex® SU bandages provide less bulky, more comfortable compression and support for all anatomical zones.

REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
LOS003	8cm x 6m	1	300
LOS004	10cm x 6m	1	240



Surgical Consumables



ThreadX® XRD_φ Abdominal Sponges

Sterile | No Tapes

ThreadX® Abdominal Sponges are constructed from 100% cotton layered gauze, making them highly absorbent. A black x-ray detectable thread is heat bonded to the gauze.

Specifically designed for all open surgical procedures.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ABX002	22.5cm x 22.5cm 6 Ply 5pcs	e	60
ABX004	30cm x 30cm 6 Ply 5pcs	e	50
ABX006	45cm x 45cm 6 Ply 5pcs	« e	25
ABX008	30cm x 50cm 6 Ply 5pcs	e	30

φXRD - X-ray Detectable
 eDEHP Free
 «Indent Item

ThreadX® XRD_φ Abdominal Sponges

Non Sterile | Folded and Taped

ThreadX® Abdominal Sponges are constructed from 100% cotton layered gauze, making them highly absorbent. A black x-ray detectable thread is heat bonded to the gauze.

Specifically designed for all open surgical procedures.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ABS030	30cm x 30cm 6 Ply 5pcs BULK	« e	80
ABS031	30cm x 50cm 6 Ply 5pcs BULK	« e	32

φXRD - X-ray Detectable
 eDEHP Free
 «Indent Item

ThreadX® HYBRID XRD_φ Abdominal Sponges

Sterile | No Tapes

ThreadX® HYBRID Abdominal Sponges have almost double the absorbency of 100% cotton sponges. Consisting of a cotton gauze outer layer with an inner 'core' of non woven fabric. The black x-ray detectable thread is heat bonded to the gauze.

Ideal for all open surgical procedures, especially where extra absorbency is required.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ABX112	22.5cm x 22.5cm 5 Ply 5pcs	e	60
ABX114	30cm x 30cm 5 Ply 5pcs	e	50
ABX118	30cm x 50cm 5 Ply 5pcs	e	30

φXRD - X-ray Detectable
 eDEHP Free

ThreadX® HYBRID XRD_φ Abdominal Sponges

Non Sterile | Cotton and Non Woven

ThreadX® HYBRID XRD Abdominal Sponges are constructed from non-woven inner layers with an all-cotton gauze outer cover. A black x-ray detectable thread is heat bonded to the gauze.

These sponges have greater absorbency than 100% cotton. Designed for all open surgical procedures especially where improved fluid control is advantageous.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
ABX101	30cm x 30cm 5 Ply 5pcs	e	80
ABX103	30cm x 50cm 5 Ply 5pcs	e	32

φXRD - X-ray Detectable
 eDEHP Free

Surgical Consumables



ThreadX® XRD_φ Standard Gauze Swabs

Sterile

ThreadX® Standard Gauze Swabs are 100% cotton, with a conventional weave count, and a black x-ray detectable thread heat bonded to the gauze.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
XRD3252	10cm x 10cm 8 Ply 5pcs	e 250	1

Specifically designed for all open surgical procedures. French folded to reduce fraying and linting.

φXRD - X-ray Detectable

eDEHP Free

ThreadX® XRD_φ Optimum Gauze Swabs

Non Sterile | Sterile

ThreadX® Optimum Gauze Swabs are 100% cotton, with a high weave count making it highly absorbent. A black x-ray detectable thread is heat bonded to the gauze.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
XRD3253	10cm x 10cm 8 Ply Sterile 5pcs	e 250	1
XRD020	10cm x 10cm 8 Ply 5pcs	e 600	1

Specifically designed for all open surgical procedures. French folded to reduce fraying and linting.

φXRD - X-ray Detectable

eDEHP Free

ThreadX® XRD_φ Peanut Swab

Sterile | Adhesive Plastic Tray

ThreadX® Peanut Swabs are constructed from 100% cotton and are highly absorbent with a black x-ray detectable thread heat bonded to the gauze. Supplied in a plastic tray of 5 with an adhesive bottom for ease of use.



REORDER	DESCRIPTION	PCS/ UNIT	UNITS/ CTN
XPE001	Plastic Tray 5pcs	e 50	1

φXRD - X-ray Detectable

eDEHP Free

Surgical Consumables



ThreadX® XRD_φ Packing Gauze

Sterile

ThreadX® Packing Gauze is constructed from 100% cotton. Designed for open surgical procedures, an x-ray detectable thread is heat bonded to the gauze.

- Ideal for packing
- deep wounds,
 - undermining, and
 - tunnels



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
XRP001	1.25cm x 1m	e	100 1
XRP006	2.5cm x 1m	e	100 1
XRP007	2.5cm x 2m	e	100 1
XRP011	5cm x 1m	e	100 1
XRP012	5cm x 2m	e	100 1
XRP016	10cm x 1m	e	100 1
XRP017	10cm x 2m	e	100 1
XRV002	10cm x 2.5m	« e	100 1
XRV003	10cm x 3m	e	100 1

φXRD - X-ray Detectable

eDEHP Free

«Indent Item

ThreadX® XRD_φ Dental Pressure Pack

Sterile

ThreadX® Dental Pressure Packs are constructed from 100% cotton. A black x-ray detectable thread is heat bonded to the gauze. The dental pressure packs have a cotton tail for ease of removal and are highly absorbent.

Ideal for dental procedures.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
XDPO01	Small with Cotton Tail	e	100 1
XDPO02	Large with Cotton Tail	e	100 1
XDPO30	Large Ball with Cotton Tapes 2pcs	e	100 1

φXRD - X-ray Detectable


eDEHP Free

ThreadX® XRD_φ Throat Pack with Tail

Sterile

ThreadX® Throat Packs are constructed from 100% cotton. A green x-ray detectable thread is heat bonded to the gauze. The Throat Pack has a black tail for ease of visibility and removal.

Designed for open surgical procedures, these packs are highly absorbent.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
XRT001	5cm x 1m	e	100 1
XRT002	5cm x 2m	e	100 1

φXRD - X-ray Detectable

eDEHP Free

Surgical Consumables



ThreadX® XRD_φ Non-Woven Swabs

Sterile | Latex Free

ThreadX® Non-Woven swabs are constructed from a polyester and rayon blend, with a black x-ray detectable thread heat bonded to the material.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
XNW3150	10cm x 10cm, 4ply, 5pcs « _e	200	1

Specifically designed for all open surgical procedures. Highly absorbent and low linting.

φXRD - X-ray Detectable

_eDEHP Free

«Indent Item

ThreadX® XRD_φ Non-Woven Balls

Sterile | Latex Free

ThreadX® Non-Woven Balls are constructed from a polyester and rayon blend, with a black x-ray detectable thread heat bonded to the material.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
XNB020	XRD Non-Woven Balls 5pcs _e	600	1

Specifically designed for all open surgical procedures. Highly absorbent and low linting.

φXRD - X-ray Detectable

_eDEHP Free


Minor Procedure Drapes with Pigeonhole®

PE and Non-woven composite

Minor procedure drapes with Pigeonhole® are designed to provide a sterile field whilst providing access to the surgical site.

Ideal to

- isolate the area of surgical incision
- provide a protective barrier
- shield instruments from contamination during operative procedures.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
MDA001	Drape Adhesive 50cm x 60cm Non Sterile	500	1
MPD001	Drape Non Adhesive 50cm x 60cm Non Sterile	1000	1
MPD020	SENTURIAN® Minor Drape 50cm x 60cm Sterile «	50	1

«Indent Item

Surgical Consumables



vertex® BASIC Standard Surgical Gown

Sterile | AAMI Level 2

vertex® BASIC surgical gown is non reinforced and has a barrier protection of Level 2 (AAMI PB70:2012).

Ideal for use where there may be light exposure to blood or body fluids.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
GNS101	Small 1pk	20	1
GNS103	Medium 1pk	20	1
GNS105	Large 1pk	20	1
GNS107	X-Large 1pk	20	1
GNS109	XX-Large 1pk	20	1

vertex® ULTIMATE Surgical Gown

Sterile | AAMI Level 4

vertex® ULTIMATE surgical gown is made from BVB, a material that provides cool, breathable comfort combined with additional fluid repellancy, has a barrier protection of Level 4 (AAMI PB70:2012).

Ideal for use where there may be moderate to heavy exposure to blood or body fluids and/or need to wear for an extended period of time.



REORDER	DESCRIPTION	PCS/UNIT	UNITS/CTN
GNS141	Small 1pk	20	1
GNS142	Small 2 pk	20	1
GNS143	Medium 1pk	20	1
GNS144	Medium 2pk	20	1
GNS145	Large 1pk	20	1
GNS146	Large 2pk	10	1
GNS147	XL 1pk	20	1
GNS148	XL 2pk	10	1
GNS149	XXL 1pk	20	1
GNS150	XXL 2 pk	10	1

Community Support

Sentry Medical is dedicated to providing community support to numerous charitable funds both locally and internationally through the donation of product and funds. Below are some of the charitable organisations Sentry support.



Open Heart International provides specialist medical services in developing countries where men, women and children cannot access the necessary healthcare they so desperately require.

Inspired to action following a visit to Tonga, Open Heart International was established by three men as volunteers, who all worked for Sydney Adventist Hospital at the time.

Open Heart International coordinates teams of humanitarian health professionals and provides free life-transforming treatments including surgery, currently in Burns surgery, Cardiac, Orthopaedic, Ophthalmology and Women's Health. Projects have been undertaken throughout the Pacific, Asia and Africa.

Many of these patients would not have survived let alone be leading the happy and productive lives they currently are in their communities around the world.

Sentry Medical is pleased to continue to support Open Heart International through the donation of a number of products, contributing to the inspiring work they do each year. Sentry Medical donates a number of products from our extensive range to their work.

For more information please visit www.ohi.org.au



YWAM Medical Ships – Australia (YWAM MSA) is actively developing communities by addressing the health care and training needs in Papua New Guinea (PNG) alongside the priorities and vision of the PNG National Health Plan, PNG National Department of Education, and DFAT's commitment to the millennium Development Goals.


YWAM MSA is implementing programs with its Medical Ship and land-based teams in rural PNG communities in association with key stakeholders and partners. All of YWAM MSA staff, including senior project managers, are full-time volunteers.

YWAM MSA facilitates outreaches with its Medical Ship and land-based teams to rural and remote areas of PNG. The clinics include: primary health care, dentistry, immunisations, antenatal checks, malaria and haemoglobin testing and treatment, optometry, and ophthalmology. Training programs for rural healthcare workers and various health promotion initiatives for school children and the general community are incorporated into the clinics.

These clinics are primarily made up of medical professionals who volunteer their time, knowledge, and expertise for two weeks or more on board the YWAM Medical Ship.

Sentry Medical is pleased to support YWAM Medical Ships Australia through the donation of a number of products, contributing to the wonderful work they do each year in Papua New Guinea. Sentry Medical donates a number of products from our extensive range to do their work.

For more information about YWAM MSA or to volunteer with them, please visit www.ywamships.org.au



Jeans for Genes Day is a national day when people wear their jeans to raise much-needed funds for research into birth defects and diseases such as cancer, epilepsy, and a range of genetic disorders.

On the first Friday in August each year workplaces, schools and streets become a sea of denim in a united stance against childhood disease.

Every dollar raised on the day helps scientists at Children's Medical Research Institute discover treatments and cures, to give every child the opportunity to live a long and healthy life.

For more information please visit www.jeansforgenes.org.au

Community Support



Created in 2013, Restore One has a mission: To break poverty one community at a time, through House, Education, Employment and Health initiatives.

Below are some of Restore One's current projects.

PROJECTS:

The Restore One Cafe and Training Centre for Vulnerable Women Street 123, Number 23, Phnom Penh, Cambodia

If you are a poor, uneducated young person, a life of dignity, safety and security is hard to find. Many end up in forced labour or slavery as the only way to pay their debts and make a living. This cycle of poverty can affect generations of people to come, unless action is taken. The Restore One Cafe serves the best burgers in Phnom Penh, the capital of Cambodia! But our purpose is far greater than that. The cafe is a training ground for young people at risk. A 12month traineeship at the cafe equips them for a successful career in the booming hospitality and tourism industry. We are proud to be running a first-rate cafe that not only satisfies our customers, but changes lives in the process.

PROJECTS:

The Restore One Primary School Project Kampong Thom Province, Cambodia

In the West, we take quality schooling for granted. In Cambodia, education is the difference between a life of prosperity and a life-sentence of abject poverty, exploitation, and even slavery. We are developing a village school in the Kampong Thom Province that has experienced a number of challenges over the years.

These include:

- Malnourished students
- Cramped classroom conditions
- A general lack of resources
- Teachers who show up as little as two days a week, because they are unable to afford transport.

PROJECTS:

Houses For Families in Need Kampong Thom Province, Cambodia

Our most successful project so far has been the funding and construction of houses for poor rural families. Cambodia's monsoon season brings brutal rainstorms and flooding, so protection from the elements is an especially serious issue. In consultation with the local people, we have developed a design that is both simple and highly effective. The houses are made with a cost-effective yet reliable combination of bamboo and timber. They are also raised on stilts to avoid flooding. The houses are constructed as a collaboration between Restore One volunteers and local labour, to create an added benefit of employment and income for local tradespeople.

PROJECTS:

Micro-Business and Life Skills Training Jewellery Enterprise - Kampong Thom Province, Cambodia Textiles Enterprise - Phnom Penh, Cambodia

At Restore One, we are establishing two micro-business enterprises:

- A cost-effective jewellery business run by women from the village where we work in Kampong Thom Province.
- A textiles design and sewing business run by young mothers in a slum area of Phnom Penh.

The real expenditure of these projects is not the materials, but the life-skills training and education that we intend to deliver for these disadvantaged young women.

Our goal is to provide, for each project, a teacher to mentor the women in the following key areas:

- Health and hygiene principles
- Parenting skills
- Nutritional principles
- Dealing with family problems such as domestic violence.

Education in these areas will assist our young women in combating the cycle of poverty and disadvantage for themselves and their children.

For more information about Restore One or to volunteer with them, please visit www.restoreone.org.au