



Highly Absorbent



Comfortable



First Aid

Sentry Wound Dressings are constructed of a thick, highly absorbent cotton pad attached to a cotton gauze bandage. The dressings are intended to soak exudate away from the wound bed whilst protecting the open wound.

## HIGHLY ABSORBENT

Large absorption capacity and holds exudate away from the wound bed

## PROTECTS AGAINST POSSIBLE INFECTION

Large dressing sizes effectively cover wound and protect against debris and foreign matter



## Suitable for:

First aid

Hospitals and medical

All general wound dressings

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
WD001	Small	Sterile	12	48
WD002	Medium	Sterile	12	22
WD003	Large	Sterile	12	18

