



Low Linting



Available in 3 Grades



Dual Coloured

Sterifit Sterilisation Wraps are made from SMS bonded polypropylene. They are suitable for steam, ethylene oxide, and plasma sterilisation processes. Dual colouring allows for easy detection of strike-through, tears and punctures.

### DUAL COLOUR (LIGHT BLUE/BLEU)


For fast and efficient detection of tears

### STRONG AND DURABLE

Provides a high resistance to tearing and abrasion

### SUITABLE FOR ALL STERILISATION METHODS

Can be used in steam, ethylene oxide, or plasma sterilisation processes



## Suitable for:

### STANDARD

- Light to medium tray
- Linen packs
- Medium devices & instruments
- Basin sets


### HEAVY DUTY

- Large trays
- Instrument sets
- Medical devices
- Multiple basin sets

### ULTRA HEAVY DUTY

- Heavy surgical sets such as orthopaedic
- Heavy instrument sets

\*Reorder details on back



Reorder	Description	Size	Weight	Colour	Sterility	Pcs/Unit	Units/Ctn
<b>STANDARD</b>							
STS001<<	Sterifit - Standard	30cm x 30cm	43gsm	Light Blue/Mid Blue	Non-Sterile	1000	1
STS002	Sterifit - Standard	43cm x 43cm	43gsm	Light Blue/Mid Blue	Non-Sterile	500	1
STS003	Sterifit - Standard	61cm x 61cm	43gsm	Light Blue/Mid Blue	Non-Sterile	500	1
STS004	Sterifit - Standard	71cm x 71cm	43gsm	Light Blue/Mid Blue	Non-Sterile	300	1
STS005	Sterifit - Standard	91.5cm x 91.5cm	43gsm	Light Blue/Mid Blue	Non-Sterile	150	1
STS006	Sterifit - Standard	100cm x 106cm	43gsm	Light Blue/Mid Blue	Non-Sterile	150	1
STS007	Sterifit - Standard	100cm x 140cm	43gsm	Light Blue/Mid Blue	Non-Sterile	100	1
STS012<<	Sterifit - Standard	100cm x 150cm	43gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STS008	Sterifit - Standard	121.5cm x 121.5cm	43gsm	Light Blue/Mid Blue	Non-Sterile	100	1
STS009<<	Sterifit - Standard	130cm x 152.5cm	43gsm	Light Blue/Mid Blue	Non-Sterile	100	1
STS010<<	Sterifit - Standard	137cm x 183cm	43gsm	Light Blue/Mid Blue	Non-Sterile	100	1
STS011<<	Sterifit - Standard	213.5cm x 213.5cm	43gsm	Light Blue/Mid Blue	Non-Sterile	50	1
<b>HEAVY DUTY</b>							
STH001	Sterifit - Heavy Duty	71cm x 71cm	61gsm	Light Blue/Mid Blue	Non-Sterile	200	1
STH002	Sterifit - Heavy Duty	91.5cm x 91.5cm	61gsm	Light Blue/Mid Blue	Non-Sterile	150	1
STH003	Sterifit - Heavy Duty	106cm x 120cm	61gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STH004<<	Sterifit - Heavy Duty	106cm x 140cm	61gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STH005<<	Sterifit - Heavy Duty	121.5cm x 121.5cm	61gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STH006<<	Sterifit - Heavy Duty	130cm x 152.5cm	61gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STH007<<	Sterifit - Heavy Duty	137cm x 183cm	61gsm	Light Blue/Mid Blue	Non-Sterile	50	1
<b>ULTRA HEAVY DUTY</b>							
STU001<<	Sterifit - Ultra Heavy Duty	106cm x 120cm	71gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STU002<<	Sterifit - Ultra Heavy Duty	106cm x 140cm	71gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STU003<<	Sterifit - Ultra Heavy Duty	106cm x 160cm	71gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STU004<<	Sterifit - Ultra Heavy Duty	121.5cm x 121.5cm	71gsm	Light Blue/Mid Blue	Non-Sterile	50	1
STU006<<	Sterifit - Ultra Heavy Duty	137cm x 137cm	71gsm	Light Blue/Mid Blue	Non-Sterile	40	1
STU005<<	Sterifit - Ultra Heavy Duty	137cm x 183cm	71gsm	Light Blue/Mid Blue	Non-Sterile	40	1

<< Indent Item