

vertex[®] BASIC

VERTEX[®] BASIC STANDARD SURGICAL GOWN PACK



vertex[®] BASIC surgical gown is non reinforced and has a barrier protection of Level 2 (AAMI PB70:2012). Complies with ISO16603, ISO16604 and ISO6529.

PRODUCT FEATURES

- o Highly breathable and lightweight
- o Low linting
- o Anti-static
- o Microbial barrier

SUITABLE FOR

- o Light exposure to blood or body fluids
- o Operating theatres



REORDER	SIZE	PACK	FEATURE	PCS/UNIT	UNITS/CTN
GNS101	Small	1pk	AAMI Level 2	20	1
GNS103	Medium	1pk	AAMI Level 2	20	1
GNS105	Large	1pk	AAMI Level 2	20	1
GNS107	X-Large	1pk	AAMI Level 2	20	1
GNS109	XX-Large	1pk	AAMI Level 2	20	1

All product information is correct at date of publication

