

# GAMMEX<sup>®</sup> Latex Micro

S U R G I C A L   G L O V E S

## Latex, powder-free surgical gloves offering exceptional tactile sensitivity and dexterity

- **Heightened dexterity and tactile sensitivity:** Made 15% thinner than Ansell's standard latex surgical gloves for greater surgical finesse, these sterile GAMMEX<sup>®</sup> gloves deliver exceptional precision when performing delicate procedures
- **Elevated ergonomic comfort:** Ansell's unique anatomic glove formers deliver surgical gloves designed for enhanced comfort and movement to help reduce hand fatigue
- **Effortless donning and doffing:** A unique internal coating made possible by Ansell's donning aid technology enables easy donning and doffing under dry or damp hand conditions, be it the first don or during intra-operative procedures



**CAUTION:** Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

### Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Damp donning
- Prevention of eye fatigue

### FEATURED TECHNOLOGIES



## GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Latex Micro

SURGICAL GLOVES

## PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Brown
Thickness	Micro
Cuff Style	Straight with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Sodium di-n-butylthiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC), Zinc 2-mercaptobenzothiazole (ZMBT)

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330052055	330052060	330052065	330052070	330052075	330052080	330052085	330052090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu

**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com