

GAMMEX[®] Latex DermaShield[™]

Surgical Gloves

Latex, powder-free surgical glove you can rely on for everyday performance

Featuring Ansell's unique formers and proprietary formulation, this reliable latex surgical glove enables you to securely grip instruments while its smooth upper surface allows for easy suture manipulation.

KEY FEATURES AND BENEFITS:

- Smooth with micro-textured finish for secure handling of instruments under dry and damp conditions
- Made using Ansell's unique anatomic formers to enhance comfort and fit to help reduce hand fatigue
- The DERMASHIELD[™] Technology internal coating facilitates dry and damp donning as well as intra-operative changes
- Additional cleaning processes ensure low residues of both latex proteins and chemical accelerators to minimize the risk of developing latex Type I and chemical Type IV allergies and skin sensitivities

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- Brief, low-risk examination procedures
- Extra comfort during long, sustained procedures
- Sterile procedures outside of the O.R.
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES

 **A.R.T.**
Pinhole Barrier Technology

 **DERMASHIELD[™]**
Donning Aid Technology



GAMMEX[®] Innovation that Protects

Driven by innovation, GAMMEX[®] surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Latex DermaShield[™]

Surgical Gloves

PRODUCT INFORMATION

| | |
|-----------------------------------------------|---------------------------------------------------------|
| Material | Natural Rubber Latex |
| Powder Content | Powder-Free |
| Color | Natural |
| Thickness | Standard |
| Cuff Style | Straight with SUREFIT [™] Technology |
| External Glove Surface | Smooth with Micro-Textured Finish |
| Internal Glove Surface | Polymer Coated with DERMASHIELD [™] Technology |
| Allergy Prevention | None |
| Grip Level | Moderate |
| Double Gloving Recommendation | Outer glove |
| Tested For Use With Chemotherapy Drugs | No, not tested in accordance with ASTM D6978 |

ORDERING INFORMATION

| Size | 5.5 | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 |
|----------------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Ansell Product Code | 330063055 | 330063060 | 330063065 | 330063070 | 330063075 | 330063080 | 330063085 | 330063090 |

PACKAGING AND STORAGE

| | |
|-----------------------------|---------------------------------------------------------------------------------------------------------|
| Packaging | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case |
| Shelf Life | 3 years |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition |

STANDARDS AND CERTIFICATIONS

| | |
|------------------------------------------|----------------------------------------------------------|
| Product Standards | AS/NZS 4179, EN 420:2003 + A1:2009, ISO 10282, JIS T9107 |
| Quality / Environmental Standards | EN556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001 |

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 1-844-494-7854
 Email: insidesalesus@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cstac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com